

# Plus Size KOAlign OA Knee Wrap

SKU: KNB242

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



## Instructions For Use

**Step 1.** Unfasten the straps and open the brace. Locate the hinge and center the condyle pad on the side of the knee.

**Step 2.** Fasten the Velcro tabs in numbered order.

**Step 3.** Feed the bottom strap through the ring and fasten, ensuring the strap is above the calf muscle to prevent sliding. Fasten the thigh strap.

**Step 4.** Using the adjustment tool, locate the hex opening just above the hinge. Turn clockwise to relieve lateral pain or counterclockwise to relieve medial pain.



Still need help?  
**Perfect Fit Guarantee!**

Visit our website or scan the QR code to find out more.

### Care Instructions

- Spray with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

