

KOAlign II Adjustable OA ROM Knee Brace

SKU: KNB247

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Detach all buckles and set each strap to its loosest setting for first time application.

Step 2. While seated with the knee at a 90° angle, place condyle pad on the outside center of the knee directly on the skin or a knee sleeve.

Step 3. Ensure strap 1 is above the calf and adjust strap length for a secure fit. Strap 1 placement is key to preventing sliding.

Step 4. Fasten straps 2, 3 and 4 and adjust strap length for a secure fit.

Step 5. Using hex tool, insert into the slot above the hinge and turn counterclockwise for medial (inner) pain and clockwise for lateral (outer) pain.

**Test by walking. Additional adjustments can be made until the maximum pain relief is achieved.*

Step 6. Range of motion stops - use only if you need flexion / extension control. Use the screwdriver to loosen the screws. Insert desired stop into hinge. Tighten screw back into place.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Spray with a 50% Alcohol/Water Mixture
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach.  

