

Dorsal Plantar Fasciitis Night Splint

SKU: NTS264

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Detach all padding straps. Place foot into the brace with knee flexed to maximize the angle of dorsiflexion.

Step 2. Secure the forefoot strap snugly with the Velcro closure.

Step 3. Position the shell of the splint against the shin and wrap the calf strap around the lower leg to secure the brace. It is not recommended that you walk in this brace.



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse thoroughly and Air Dry.
- Do NOT dry clean, iron, or bleach.  

