

Universal Thumb Stabilizer Splint, Spica and Wrist Brace

SKU: WRB115

Prefer watching videos?

We've got you covered!

Visit our YouTube @BraceDirect or scan the QR code



Instructions For Use

Step 1. Open brace and position hand in thumb brace. The palmer stay should line up below the middle crease of the hand. The thumb stay should follow the length of the thumb and arm.

Step 2. Wrap the thumb strap around the thumb and secure to hook.

Step 3. Wrap the lacing straps around the wrist and loop, adjusting as necessary.

**Knuckles should be free bending.*

1



2



3



Still need help?
Perfect Fit Guarantee!

Visit our website or scan the QR code to find out more.

Care Instructions

- Hand Wash with a 50% Alcohol/Water Mixture
- Rinse thoroughly and Air Dry.
- Do NOT dry clean, iron, or bleach.

