

A lightweight walker boot providing support and comfort to patients, aiding rehabilitation capabilities.

The BeneCare Benekidz Walker is a lightweight walker boot providing support and comfort for children. The hook and loop closure provides easy adjustment and the fixed stirrups immobilises the ankle securely and comfortably. The rocker bottom sole not only allows for easy ambulation but also is impact resistant and shock absorbing increasing protection for the patient. The boot is extremely durable and its ultra-low profile design eliminates height difference and facilitates natural gait. The padded soft foam also ensures the lower leg is constantly comfortable.



Features & Benefits	Indications For Use
Rigid metal frame provides protection for active children	Replaces the need for heavy casts
Shock absorbing sole to improve comfort & protects the users from slips	Stable fracture of fibula, calcaneus and malleolus
Rocker bottom sole allows easier ambulation	Post ligament, soft tissue and tendon surgery
Lightweight design for kids to help them keep active	Post-surgery management
Hook & Loop closure provides easy adjustment	Toe walking
Ultra low profile design eliminates height difference and facilitates natural gait	
Padded soft foam keeps lower leg comfortable	
Excellent durability	



Benecare Medical
Branson Court, Branson Street,
Manchester, M40 7FP, United Kingdom
Email: info@benecaremedical.com
Tel: +44(0)161 273 6789
www.benecaremedical.com



Application Guidelines



1. Remove the liner from the boot and apply to the limb, ensure it is a perfect fit.

Secure on both the leg and the foot using the hook and loop fastening system.



2. Take the outer boot and position the stirrups in line with the malleolus (ankle bones).

When it is in the correct position, remove the pull up tabs.



3. Secure the stirrup to the liner using the adjustable straps, first on the foot. Followed by the three removable straps positioned evenly along the leg. To re-apply at home, leave the boot complete. Open the straps and the liner and position the limb inside, secure as above



Product Codes, Sizes & Colours

NHS Code	BeneCare Code	Size	UK Size	EU size
	BKWS	Small	5 - 7 (15cm)	22 - 24
	BKWM	Medium	7 - 12 (19cm)	24 - 30
	BKWL	Large	12 - 1 (21cm)	30 - 32.5

If you wish to speak to our expert Clinical Specialist Pat Radcliffe regarding any application techniques or material characteristics call +44 (0) 7817 800524 or email pat.radcliffe@benecaremedical.com

Buy Online.... www.benecaredirect.com

Benecare Medical
 Branson Court, Branson Street,
 Manchester, M40 7FP, United Kingdom
 Email: info@benecaremedical.com
 Tel: +44(0)161 273 6789
www.benecaremedical.com

B
 BENE CARE

