User Guide for Posture Corrector (Basic & Posture+)

Intended use:

You can start wearing for 15 minutes/day for the first 3 days and slightly increase the usage time after that. You will feel the posture brace pulling your shoulder back. This will function also to remember you to keep the proper posture, even when you are not wearing it.

For better results, we recommend you wear it for a minimum of two hours/day for a longer period. Please, consult the table below, where you can see the program by our specialist. Please don't wear the posture brace too tight; you need to allow your spine to adapt for about seven days.

Days	Time	Where are you now
1 to 3	15min / jour	1 2 3
4 to 7	30min / jour	4 5 6 7
8 to 11	50min / jour	8 9 10 11
12 to 15	70min / jour	12 13 14 15
16 to 22	100min / jour	16 17 18 19 20 21 22
23 to 31	120min / jour	23 24 25 26 27 28 29 30 31

Precautions:

When using any brace that restricts movement, controls structural positioning or limits muscular motion, the wearer must monitor the effect the brace has on their body to avoid risk of injury. Do not apply the posture corrector so tightly that it limits your blood supply or normal neurological function of the arms. If any problems arise, such as (but not limited to) numbness, swelling, discoloration, discomfort, pain, sensation or anything else unusual, stop using the posture corrector and consult a physician immediately.

Wash body parts covered by the posture corrector and the posture corrector itself, regularly in order to avoid skin irritations.

In case you suffer from back pain, shoulder pain or neck pain or had a recent surgery in this area, please talk to your doctor before using the posture corrector.

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 Zhongshan Chukuang Textile Co., Ltd
Baishawan Second Industrial Zone, East District,
Zhongshan City, China

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Shenzhen Youbest Testing Technology Co., Ltd.
3/F, Building 4, Zhenbao Industrial Park, Industrial West Road,
Longhua District, Shenzhen, Guangdong, P.R. China
0755-6111 7096
www.youbest-lab.com

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Anodyne GmbH