Clean Water System Installation

(Model #49-031 - with public water selection toggle)

- **Step 1-** Survey project and parts required. Parts can vary depending on circumstances.
- **Step 2-** Ensure main supply of air and water to delivery unit is off.
- Step 3- Locate the "main" supply of air going to the delivery unit (not the air "signal" to the master on/off). Splice the yellow tubing into this tubing, using the 3-port "tee" connector.
- **Step 4-** Locate the water supply to the delivery unit. At the delivery unit, disconnect this tubing from barb and attach the solid blue tubing onto this barb. Attach the water supply tubing disconnected (from unit) to the clear blue tubing.
- **Step 5-** Test for leaks and proper function consistent with the notes below.

Note:

- Ensure this system does not exceed 35 psi (see gauge on water system).
- When filling bottle- set selector switch to "bottled water," turn off switch on the water system and slowly unscrew bottle to relieve pressure before removing.
- Never over fill bottle, fill to approximately 2" from the bottle's top.





7310 Oxford Street • St. Louis Park, MN • 55426
Phone (800) 331-7993 • Fax (888) 729-1016
www.amerdental.com