
Clean Water System Installation

(Model #49-19- with public water selection toggle)

- **Step 1** Ensure main supply of air and water to delivery unit is off.
- Step 2- Locate the main supply of air going to the delivery unit. This will be a line that will have air in it only when the unit is on. Splice into this tubing using the 3-port "T" connector provided and connect the tubing marked "air inlet" to this "T".
- Step 3- Locate the water supply tubing to the delivery unit. At the junction box, disconnect this tubing from the barbed fitting and attach the tubing marked "fluid outlet" onto this barb. At the delivery system, disconnect the main water supply line and attach the tubing marked "fluid inlet" in its place.
- **Step 4** Turn main supply of air and water back on. Test for leaks and proper function consistent with the notes below.

NOTE:

- Ensure this system does not exceed 35 psi (see gauge on water system). You can also verify correct pressure by hitting the water button on your air water syringe. With the bottle set to 35psi, the water stream should be approximately 3 1/2' long.
- When filling bottle- set selector switch to "bottled water", turn off switch on the water system and slowly unscrew bottle to relieve pressure before removing.
- Never over fill bottle, fill to approximattely 2" from the bottle's top.





7310 Oxford Street • St. Louis Park, MN • 55426 Phone (800) 331-7993 • Fax (888) 729-1016 www.amerdental.com