INSTRUCTIONS FOR USE

Gingival Hemostatic Retraction Paste

A00661revG0418

Rx Only • For Safety Data Sheet (SDS) go to www.parkell.com

Stock No. S180 – Green Syringe Kit (7 Syringes + Tips)
 Stock No. S182 – Green Syringe Value Pack (25 Syringes + Tips)
 Stock No. S183 – Green Unit Dose Capsules (30)
 Stock No. S187 – Blu Unit Dose Capsules (30)
 Stock No. S186 – Blu Unit Dose Capsules (30) & Dispensing Gun (1)
 Stock No. S190 – Blu Syringe Kit (7 Syringes + Tips)
 Stock No. S191 – Blu Syringe Value Pack (25 Syringes + Tips)
 Stock No. S181 – Replacement Luer-Lock Tips (40)

STORE THIS PRODUCT IN A COOL, DRY PLACE.

Refrigeration will prolong shelf life, but avoid freezing. Bring to room temperature before use.

PRODUCT DESCRIPTION

Parkell's Dryz® and Dryz® Blu Gingival Hemostatic Retraction Pastes are dental products capable of temporarily controlling gingival bleeding, displacing marginal gingiva and drying the gingival sulcus around a tooth. The pastes allow the dental clinician to perform operative procedures without interference from blood, moisture or hyperplastic soft tissue.

Original Dryz Green and Dryz Blu are supplied in several ways:

1. In 0.85 gram syringes with hand plungers, or
2. In 0.24 gram single-patient use, unit dose capsules, designed for use with a Centrix®-compatible dispensing gun.

To simplify these instructions, both Original Green "Dryz" and "Dryz-Blu" are referred to here as "Dryz".

INDICATIONS FOR USE

Dryz Gingival Hemostatic Retraction Paste is indicated for temporary retraction and hemostasis of the marginal gingival tissues during clinical procedures such as taking subgingival and perigingival dental impressions, seating of fixed prostheses, and restorative dentistry.

CONTRAINDICATIONS FOR USE

Dryz Gingival Hemostatic Retraction Paste should not be swallowed, and is contraindicated for use on or by persons who are allergic to or sensitive to aluminum chloride, silicates or similar products. Symptoms may include dermatitis or inflammatory reactions. If this occurs, discontinue use and obtain medical care immediately. The product should not be used in the presence of, and is not a treatment for, active periodontal disease.

PRECAUTIONS FOR USE

- Dryz should be used in gingival sulci above the level of the junctional epithelium.
- The dispensing tip of the syringe or the tip of the unit dose capsule should be held at the mouth of the sulcus, without actually entering into it. The paste should be expressed slowly into the sulcus, so as not to inject the material under pressure into the periodontal tissues.
- After the material has been in the sulcus for approximately 2 minutes, it must be completely removed. The sulcus should be thoroughly rinsed with water and inspected for residual material so none is left behind on the gingiva or the tooth surface.

PRODUCT WARNINGS AND CLINICAL INFORMATION

- Store Dryz in a cool dark place prior to use, sealed in the foil package.
- After use, Dryz syringes may be 1) discarded or 2) resealed with the original cap tightly screwed back on, and placed into the foil package, to avoid drying out the material.
- Syringe dispensing needle tips are single-patient use only, and should be discarded after use.
- Unit dose capsules are single-patient use only, and should be discarded after use.
- Avoid sunlight or nearby heat sources. Shelf life will be reduced if not stored properly.

- Dryz Hemostatic Paste, like many agents used in dentistry to control bleeding, contains Aluminum Chloride. This may cause corrosion on some plastic or metal surfaces (even stainless steel), if left in extended contact. To protect dental instruments, always remember to rinse off Dryz from their surface immediately after each use.
- If product accidentally contacts skin or eyes, wash with copious amounts of water. Obtain medical care immediately if necessary.
- Do not ingest this product, and keep it away from children.

GENERAL INSTRUCTIONS FOR USE

1. After tooth preparation is complete, rinse the sulcus with water. Gently dry with compressed air at an angle to the open sulcus, not directly into it.
2. Open the resealable foil package by tearing off the top where indicated. Save the bag for later storage. Remove the syringe or a unit dose capsule from the package.
3. If dispensing from syringe, remember that Dryz is a thick material being pushed through a thin narrow tip. To ensure that the syringe needle tip does not come off during placement, follow these steps before use.
 - a. Remove the sealing cap from the syringe. Save the cap, as it may be used to reseat the syringe.
 - b. Depress the plunger to initiate the flow of a small amount of Dryz from the syringe orifice. Wipe off the excess.
 - c. Using the needle tip cover as a wrench, twist a needle tip onto the syringe until it is tight. If the tip is not fully locked on the syringe, the pressure may cause it to come off when the plunger is squeezed.
 - d. Remove the needle tip cover and carefully bend the needle to the desired angle using any suitable method. If the tip kinks, discard it and use another one.
 - e. Initiate the flow through the bent needle tip by slowly expressing a small amount of Dryz onto a pad.
4. If using a unit dose capsule, insert it properly into a Centrix-compatible gun, remove the cap and "start" the flow by extruding a small amount onto a pad first.
5. Express Dryz around the tooth into the mouth of the gingival sulcus, until it is slightly overfilled with material. Do not force the tip of either the needle or the capsule into the tissues. Allow the material to remain undisturbed in the sulcus for 2 minutes.
6. To increase gingival retraction, pre-made cotton compression caps (or standard cotton rolls cut into short sections) may be placed over the tooth, followed by 2 minutes of biting pressure by the patient's opposing teeth.
7. Dryz remains pliable, and does not fully dry out in the sulcus. After 2 minutes, rinse it away with water and air and remove any remaining material using a non-cutting hand instrument or cotton pledget. Rinse and dry again. The sulcus is now ready for the impression or restorative procedure.
8. After use, straighten the needle tip and re-cap it with the needle tip cover. Use the cover as a wrench to remove the tip from the syringe. Dispose of the covered tip or the unit dose capsule in accordance with applicable safety and environmental regulations. Needle tips and unit dose capsules are for single-use on one patient only. Each time Dryz is used on a different patient, apply a new needle tip or unit dose capsule.
9. If any Dryz remains in the syringe, you may reuse it by sealing it with the original plastic cap (not the needle tip), and placing it back

into the foil bag. The recapped syringe should be cleaned and disinfected with an EPA-registered low (HIV/HBV claim) to intermediate level (tuberculocidal claim) hospital disinfectant, as per the disinfectant manufacturer's instructions. Consult the Centers for Disease Control website at CDC.gov, referencing the most recent version of the "Guidelines for Infection Control in Dental Health-Care Settings". Store the syringe in the foil bag for use on the same patient in the future.

10. Alternately, the syringe may be discarded after a single use, in accordance with applicable regulations.

11. If you have any questions, call 1-800-243-7446 or e-mail info@parkell.com.

WARRANTY AND TERMS OF USE:

Parkell will replace defective material. This warranty is in lieu of all warranties of merchantability, fitness for purpose or other warranties, express or implied. Parkell does not accept liability for any loss or damage, direct, consequential or otherwise, arising out of the use of or the inability to use the product herein described. Before using, the user shall determine the suitability of the product for its intended use and the user assumes all risk and liability whatsoever in connection herewith. For full Warranty and Terms of Use information, please see www.parkell.com.

Safety Data Sheets (SDS) are available at www.parkell.com.
Parkell's Quality System is certified to ISO 13485.

EC REP European Authorized Representative
(Not a dealer/distributor):
Directa AB - P.O. Box 723, Finvids väg 8
SE-194 27 Upplands Väsby - Sweden

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