PULPDENT® Silane Bond Enhancer

PULPDENT SILANE is a single-component material which is used to promote increased bond strength between organic resins and porcelain. To achieve desired results, the etched porcelain surface must be completely dry and free of residual moisture before applying silane.

CONTROLLED DISPENSING

Syringe dispensing with Pulpdent's unique red dropper tips eliminates waste and provides controlled dispensing of silane.

Note: Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient. For additional information, refer to: <http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/DentalProducts/ucm404472.htm>

DIRECTIONS FOR USE

1. Carefully follow instructions for Pulpdent Porcelain Etch Gel, and etch porcelain surface to be treated. Rinse and dry.
2. Apply Pulpdent Dry-Rite™ drying agent to etched surface to assure that it is completely dry. Residual moisture will adversely effect bonding.
3. Place 22-gauge red dropper tip on silane syringe and apply Pulpdent Silane Bond Enhancer to treated porcelain surface. Allow to dry. A gentle flow of uncontaminated, oil-free and moisture-free air may be applied to accelerate drying. Cap syringe immediately after use. To check air line for contamination, blow air onto the surface of dental mirror. If streaks of water or droplets of oil appear, do not use the air line for this procedure.
4. Maintain a clean, dry field. Apply unfilled resin bonding agent and proceed with your preferred restorative composite or resin cement.

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PULPDENT®

12.0 Ecological Information

- 12.1 Ecotoxicity Follow good work practices and government regulations. Avoid release into environment.

13.0 Disposal Considerations

- 13.1 Regulations Follow all local and national government regulations in disposing material or contaminated packaging.

14.0 Transport Information

- 14.1 UN Number 1170
14.2 Technical name Ethyl alcohol
14.3 Packing group II
14.4 IATA class 3

15.0 Regulatory Information

- 15.1 EU Class IIa medical device under MDD 93/42/EEC.
15.2 US FDA Class II medical device.
15.3 Health Canada Class III medical device.

16.0 Other information

- 16.1 List of relevant R phrases
R11: Highly flammable
R36/ 37/38: Irritating to eyes, respiratory system and skin.
R66: Repeated exposure may cause skin dryness or cracking.
- 16.2 Hazard Statements
H225: Highly flammable liquid and vapor. Category 2.
H319: Causes serious eye irritation. Category 2.
H335: Specific Target Organ Toxicity (STOT), single exposure, respiratory tract, Category 3: May cause respiratory irritation.
H315: Causes skin irritation. Category 2.
EUH066: Repeated exposure may cause skin dryness or cracking.
- 16.3 Precautionary Statements
P210: Keep away from heat, sparks, open flame, hot surfaces. No smoking.
P403+P233: Store in a well-ventilated place. Keep container tightly closed.
P261: Avoid breathing fumes.
P280: Wear protective gloves/ clothing and eye protection.
P304+P340: If inhaled, remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305+P351+P338: If in eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353: If on skin or hair, remove contaminated clothing. Rinse skin with water.
P370+P378: In case of fire, use dry chemical, alcohol foam, or carbon dioxide for extinction.
- 16.4 Restrictions on use Pulpdent dental materials are to be sold to and used by dental professionals only.
- 16.5 Further information The information presented herein is believed to be factual as it has been derived from the works of persons believed to be qualified experts. However, nothing contained in this information is to be taken as a warranty or representation for which Pulpdent Corporation bears legal responsibility. The user should review any recommendations in the specific context of the intended use to determine whether they are appropriate.
- 16.6 Sources of key data National Institute for Occupational Safety (NIOSH)
US Occupational Safety and Health Administration (OSHA)
Eur-Lex European Union Law: Regulation (EC) No. 1272/2008 (CLP) and Regulation (EC) No. 1907/2006 (REACH).
Guidance on the compilation of safety data sheets. Version 1.1; December 2011. European Chemicals Agency
This Safety Data Sheet has been revised to meet the requirements of the GHS SDS format, Regulations (EC) No. 1272/2008 (CLP) and (EC) No. 1907/2006 (REACH). Specifically, Sections 2.1, 2.2, 3.2, 16.2, 16.3 have been deleted or revised.
- 16.7 Information which has been added, deleted or revised.

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Pulpdent Corporation		Revision Date: January 1, 2020
Safety Data Sheet		
Trade Name:	SILANE BONDING AGENT	
1.0 Commercial Product Name and Supplier		
1.1 Commercial product name/designation	Silane Bonding Agent	
1.2 Application / Use	Dental material for use by dental professionals.	
1.2.2 SIC	851 Human health activity	
1.2.3 Use Category	55	
1.3 Manufacturer, Importer		
1.3.1 Manufacturer	Pulpdent Corporation, 80 Oakland Street, P.O. Box 780, Watertown, MA 02472 USA Telephone: 1 617 926-6666 / Fax: 1 617 926-6262 Email: Pulpdent@pulpdent.com	
1.4 Emergency Telephone Number	1-800-535-5053 (24 Hour / USA)	
1.5 Authorized European Representative	Advena Ltd. Tower Business Centre Warwick CV34 6WE UK 2nd Floor, Tower Street Swatara, BKR 4013 Malta	

2.0 Hazards Identification			
2.1 Classification			
2.1.1 Classification according to Regulation (EC) No. 1272/2008 [CLP]	Hazard Class	Hazard Category	Hazard Statement
	Flammable liquid	2	H225
	Eye irritation	2	H319
	STOT SE	3	H335
	Skin irritation	2	H315; EUH066
2.1.2 Classification according to Directive 67/548/EEC (See SECTION 16 for full text of risk phrases)	Flammable (F) R11		
	Irritant (Xi) R36 / 37 / 38 - 66		
2.2 GHS Label Elements			
Hazard Pictograms 			
Signal Word: DANGER			
Restricted to use by dental professional only			
Hazard Statements:			
H225: Highly flammable liquid and vapor. Category 2.			
H319: Causes serious eye irritation. Category 2.			
H335: Specific Target Organ Toxicity (STOT), single exposure, respiratory tract, Category 3: May cause respiratory irritation.			
H315: Causes skin irritation. Category 2.			
EUH066: Repeated exposure may cause skin dryness or cracking.			
Precautionary Statements:			
P210: Keep away from heat, sparks, open flame, hot surfaces. No smoking.			
P403+P233: Store in a well-ventilated place. Keep container tightly closed.			
P261: Avoid breathing fumes.			
P280: Wear protective gloves/ clothing and eye protection.			
P304+P340: If inhaled, remove victim to fresh air and keep at rest in a position comfortable for breathing.			
P305+P351+P338: If in eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.			
P303+P361+P353: If on skin or hair, remove contaminated clothing. Rinse skin with water.			
P370+P378: In case of fire, use dry chemical, alcohol foam, or carbon dioxide for extinction.			

3.0 Composition				
3.1 Chemical characterization of the preparation:	Denatured ethyl alcohol preparation.			
3.2 Hazardous ingredients	CAS Number	Name of the Ingredient	Concentration	Classification per Regulation
	64-17-5	Ethyl alcohol	92%	Flammable (F); Irritant (Xi). R11- 36/ 37/38-66
	67-64-1	Acetone (denaturant)	7%	Flammable (F); Irritant (Xi). R11- 36/ 37/38-66

4.0 First Aid Measures		
4.1 General Information	May cause irritation of eyes or skin on contact. May cause irritation of respiratory tract if inhaled. Exposure to ethanol >1000 ppm may cause headache, drowsiness, lassitude, appetite loss. Show this safety data sheet to medical personnel. Get medical attention in case of uncertainty.	
4.2 Eye Contact	Keep eyelids apart, flush with running water for 15+ minutes. Get medical attention.	
4.3 Skin Contact	Remove contaminated clothing. Immediately wash with soap, running water. Use hand cream. Get medical attention if irritation persists.	
4.4 Ingestion	Rinse mouth with water. Do not induce vomiting. Get immediate medical attention. Never give anything by mouth to an unconscious person.	
4.5 Inhalation	Move to fresh air. If necessary, administer oxygen and/or artificial respiration; seek medical attention.	
4.6 Precautions for first responders	Ventilate the area. Wear safety glasses and gloves.	
4.7 Information for physicians	Symptoms Contact with material may cause irritation or redness in eyes or on skin. Inhalation may cause irritation of respiratory tract.	

Hazards	Exposure to ethanol >1000 ppm may cause headache, drowsiness, lassitude, appetite loss. Persons with chronic respiratory or skin disease are at increased risk with prolonged and/or repeated contact. Same as above under First Aid	
Treatment		
5.0 Fire Fighting Measures		
5.1 Suitable extinguishing media	Use dry chemical, alcohol foam, or carbon dioxide. Water may be ineffective, but should be used to keep fire-exposed containers cool.	
5.2 Extinguishing media to avoid	Water may be ineffective, but will keep fire-exposed containers cool.	
5.3 Special exposure	Carbon monoxide, carbon dioxide from incomplete combustion hazards in a fire	
5.4 Special protective equipment for fire-fighters	Self-contained breathing apparatus.	
6.0 Accidental Release Measures		
6.1 Personal precautions.	Wear chemical splash goggles and gloves.	
6.2 Environmental precautions	Avoid releasing large quantities into environment.	
6.3 Method for clean up	For small quantities: Ventilate area. Wear safety glasses, lab coat, gloves. Wipe up with absorbent material (paper or cloth towels). Rinse area of spill with water. Place all material in closed container away from heat, sparks, sun and oxidizers.	
7.0 Handling and Storage		
7.1 Handling	For use by dental professionals only. Remove applicator tip and recap immediately after use. Keep material tightly capped in original container. Do not use in presence of ignition sources. Take same precautions when container is emptied, as residual product is hazardous.	
7.2 Storage	Remove applicator tip after use. Keep tightly capped in original container. Store at cool room temperature in a well ventilated area. Avoid extreme temperatures (>27°C/80°F, <5°C/40°F), sparks, direct sunlight, oxidizing agents. Vapor may form flammable mixtures with air.	
7.3 Specific uses	Dental bond enhancer	
8.0 Exposure Controls / Personal Protection		
8.1 Exposure limit values	Ethanol: 1000 ppm. Acetone: 750 ppm	
8.2 Exposure controls		
8.2.1 Occupational exposure controls	No special equipment or ventilation required under normal conditions of use in dental practice For large quantities/prolonged exposure, use enclosure, local ventilation, dilution to reduce concentration below TLV.	
8.2.1.1 Respiratory protection	No special requirements under normal conditions of use. Good general ventilation is sufficient to control any airborne vapors.	
8.2.1.2 Hand protection	No special requirements other than surgical gloves.	
8.2.1.3 Eye protection	No special requirements other than safety glasses.	
8.2.1.4 Skin protection	No special requirements other than good hygiene and safety practices.	
8.2.1.5 Other controls	Emergency eye wash fountain should be close by and maintained. Wash hands after use. Do not eat, drink or smoke.	
8.2.2 Environmental exposure controls	Follow all government regulations.	
9.0 Physical and Chemical Properties		
9.1 Appearance / Color		
9.1.1 Color	Clear liquid	
9.1.2 Odor	Characteristic, sweet, ethanol odor	
9.2 Important health, safety and environmental information		
9.2.1 pH	Not applicable	
9.2.2 Boiling point	Boiling Point: 173°F / 78.3°C	
9.2.3 Flash point	43°F / 6°C (Tag closed cup)	
9.2.4 Ignition temperature	423°C	
9.2.5 Explosive properties	LEL: 3.3; UEL: 19	
9.2.6 Odor threshold	159 ppm	
9.2.7 Vapor pressure	44.6 mm Hg / 59 mbar / Id: E	
9.2.8 Specific gravity	0.795	
9.2.9 Solubility in water	Slightly	
9.2.10 Partition coefficient	Not determined	
9.2.11 Viscosity	Not determined	
9.2.12 Vapor density	1.59	
9.2.13 Evaporation rate	Not determined	
10.0 Stability and reactivity		
10.1 Conditions to avoid	Heat, sparks, open flame, any ignition source.	
10.2 Materials to avoid	Acetyl chloride and a wide range of oxidizing agents.	
10.3 Hazardous decomposition products	Carbon monoxide, carbon dioxide from incomplete combustion.	
10.4 Further information	Stable if stored and used as directed.	
11.0 Toxicological information		
11.1 Acute toxicity	Minimal health hazard under normal conditions of use in dental practice Ethanol: LD ₅₀ in young rats: 10.6 g/kg orally; LD ₅₀ in old rats: 7.06 g/kg orally. Acetone: LD ₅₀ in rats : 10.7 ml/kg orally	
11.2 Irritation and corrosiveness	May cause irritation/redness of eyes or skin on contact. May cause irritation of respiratory tract if inhaled.	
11.3 Sensitization	Not applicable.	
11.4 Sub-acute, sub-chronic and prolonged toxicity	No chronic health hazard under normal conditions of use.	
11.5 Carcinogenicity, Mutagenicity, Reproductive Toxicity	Not a carcinogen under normal conditions of use. Large quantities of ethanol, ingested over time, may be carcinogenic or a cause of Fetal Alcohol Syndrome. IARC has reported a relationship between habitual drinking of significant quantities of alcoholic beverages and cancer of oral cavity, pharynx, esophagus, liver.	
11.6 Empirical data	None available	
11.7 Clinical experience	0 and similar bond enhancers have been used in dentistry for decades with no reported problems.	