Betriebsanleitung
Operating instructions
Mode de fonctionnement
Istruzioni sul
funzionamento
Instrucciones de servicio
Upute za rad

Deutsch	04 - 15
English	16 - 27
Français	28 - 40
Italiano	41 - 52
Español	53 - 65
Hrvatski	54 - 63

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OPERATING INSTRUCTIONS CERAMILL THERM

Translation of the original operating instructions

- 1 Mains control switch
- 2 Indicator lamp heating
- 3 Digital display / actual temperature
- 4 Digital display / target temperature
- 5 Menu key
- 6 Programme start key
- 7 Selection keys
- 8 Status indicator LED programme sequence
- 9 Status indicator LED control system

PRODUCT DESCRIPTION

Ceramill Therm is a high-temperature chamber kiln for final sintering dental oxide ceramics, especially Ceramill ZI blanks. The sintering furnace exhibits 3 temperature programmes preset by the manufacturer, which are tailored to the Ceramill ZI framework material (Y-TZP ZrO₂, pre-sintered). An additional programme location can be freely programmed. The stackable sintering bowl filled with sintering balls (178311) serves to sinter the Ceramill ZI frameworks (178360).

UNPACKING AND CHECKING

Unpack the Ceramill Therm Sintering Furnace carefully. Two people must carry it when transporting it. Never lift up by the door! Remove all packaging material from the door and chamber. After the furnace has been delivered and unpacked, check that all the individual parts are complete. Check also for any possible damage in transit and the integrity of the three heating elements. If any damage should have arisen during transit, please notify the supplier immediately.

SETTING UP THE DEVICE/ASSEMBLY

Place the furnace in a well-ventilated room, not close to other heat sources and on a fireproof

base. See to it that there is sufficient free space around the furnace and that the ventilation slits are not covered.

CAUTION! Do not connect the power cable supplied with the device before the heating elements have been fitted.

No fumes arise when sintering Ceramill ZI frameworks, i.e. the furnace does not have to be placed under an extraction hood. Erect the furnace in such a way that the mains switch or mains plug are within easy reach so that the furnace can be quickly disconnected from the power supply in an emergency. The mains supply must have an earth (green and yellow cable).

CAUTION! The manufacture accepts no liability for any damage arising from improper installation.

COMMISSIONING


Detailed training to commission the furnace will take place during the Ceramill courses. In cases of doubt regarding the starting up of the furnace, please consult an authorised service provider/ electrician!

Before starting it up, check that the voltage specified on the label on the back of the casing corresponds to the mains voltage available.

Commissioning the furnace requires the following parts that have been delivered separately:

- _ 3 heating elements
- _ 6 bolts with nuts width across 10
- _ 12 washers
- _ 2 conductor lines
- _ Ceramic plate for furnace base
- _ Hexagon wrench key for rear wall screws

Fitting the heating elements

Ensure that the Ceramill Therm is not connected to the power supply. Remove the back of the casing. Remove the earth cable from the inside of the back wall.

Carefully insert the heating elements into the pipes provided for that purpose in the furnace chamber.

CAUTION! Proceed with utmost caution (risk of breakage)!


Align the contact wings of the heating elements vertically. Ensure that the heating element has a loose and stress-free position.

Electrical connection of the heating elements

Then screw the conductor lines loose with the aid of the bolts and washers.

Electrical wiring diagram (series connection)

If necessary, carefully tighten the screw connections by means of a spanner (size across flats 10). Steadying the screw by hand prevents the heating elements being burdened unrestricted.

Insert the earth cable on the inside of the back wall and screw the back of the casing shut again. Connect the Ceramill Therm to the power cord in the power point (16 A).

Inserting the base plate

Carefully insert the ceramic plate supplied into the furnace chamber with the smooth side up.

OPERATION/HANDLING

After connection to the mains power, the network control switch **1** will glow continuously.

Only after the control switch has been activated is the Ceramill Therm connected to the current control circuit; in doing so, the controller automatically conducts a lamp test (lamp test see general instructions).

Feeding the Ceramill Therm with zirconium frameworks

See the training documents relating to the Ceramill System and the instruction manual for Ceramill ZI zirconium blanks.


PROGRAMME OVERVIEW

The Ceramill Therm exhibits 3 programme locations for sintering Ceramill ZI frameworks of different indications as well as a freely programmable programme location for sintering dental oxide ceramics of other manufacturers (follow the instructions of the manufacturer!).

Programme description

_ Programme 1


Standard programme for Ceramill ZI.

Universal sintering programme for all common indications.

_ Programme 2


Fast sintering programme for Ceramill ZI single units.

Fast sintering programme for single units only. The furnace may be charged with max. one sintering bowl; otherwise, full sintering of the objects cannot be achieved.

_ Programme 3


Sintering programme with extended cooling period for Ceramill ZI.

Sintering programme with hold stage at 900°C and defined cooling rate. Especially for large and massive frameworks, where the connected residual blank is sintered together with the framework.

_ Programme 4

Freely programmable parameters.

4 successive segments can be assigned to the free programme location no. 4. Each segment contains one heating rate, temperature hold stage and, if defined, one cooling rate. The maximum adjustable sintering temperature is 1530°C; the heating rate is limited to 12°C/min. As a default, the programme is pre-set to 1500°C, 2 hours of holding time and a heating rate of 12°C/min (not specific to any material)

NOTE: If the maximum temperature of 1530°C is set, a maximum heating rate of 8°C/min is principally recommended to prevent the aging of the heating elements. Permanent use of the furnace at 1530°C leads to premature aging and vitrification of the heating elements!

PROGRAMMING AND HANDLING THE CONTROLLER OF THE FURNACE

Basic function of keys

Menu selection key:

- _ Selects the furnace parameters
- _ Saves altered furnace parameters

Start key:

- _ By pressing it for 2 seconds, the programme is started
- _ By pressing it for 2 seconds, the programme is halted
- _ By pressing it for approx. 5 seconds, the running or halted programme is stopped

- _ Selects parameters upwards in the programming mode (OP.P)

- _ By pressing the Menu selection and Start keys at the same time for more than 3 seconds, programming is enabled

Arrow keys:

- _ For individual adjustment of the selected parameter

- _ By pressing the Menu and Down keys at the same time, a lamp test is conducted (see general instructions)

- _ By pressing any key, the test stops

- _ By pressing the Start and Down keys at the same time, the respectively lowest parameter is reached

- _ By pressing the Start and Up keys at the same time, the respectively highest parameter is reached

SELECTING AND STARTING A PROGRAMME

After the furnace has been switched on

Standby „ready“

Call up the menu

Select the programme

Confirm the programme

Wait for standby (approx. 10 sec)

Start the programme (press the key for approx. 2 seconds)

Programme is running, actual and target temperatures are displayed.

- _ The furnace can only be started from the standby "rdy" mode!
- _ The selected programme needs to be confirmed by activating the Menu key ; otherwise, it will not be accepted!
- _ Even after the mains switch  has been activated, a selected programme is still saved and can be started directly from the standby mode.
- _ The controller automatically returns to the standby mode after 10 seconds.
- _ After starting the programme, the status indicator LED programm  lights permanently until the programme is terminated.

_ After starting the programme, the status indicator LED control system ❶ and the heating control lamp light ❷ synchronously until the holding time expires. If the furnace has reached the predefined target temperature, both LEDs may temporarily go out, as no additional heating is required due to the permanent temperature rise in the firing chamber.

Display of remaining time

_ By pressing the Start key ❶ for approx. 2 seconds during a running programme, the programme is halted. By pressing the key once again for 2 seconds, the programme continues. The programme can be stopped by pressing the Start key ❶ for 5 seconds.

Use the sintering pliers for removal (other sintering pliers may be used as well)

CAUTION! When removing the bowl, do not knock against the heating elements with the pliers!

A safety switch interrupts the heating process when the furnace door is opened! The sintering program will continue after the door has been closed again.

CAUTION! Do not open the furnace during the entire program sequence (frameworks, heating elements and insulation will receive a thermal shock)!

Halting the programme

_ By activating the Menu key ❷ during a running programme, segment information and the remaining time can be displayed. In doing so, the cooling period is not taken into account, if no cooling rate has been defined.

Display of current segment

DISPLAY OF PROGRAMME PARAMETERS

After switching on the furnace

Standby "ready"

Holding temperature 1 (°C). Here: 1450°C

Heating rate 1 (°C/min). Here: 8.0°C/min

Call up the menu and select the programme

Starting point of the programme "OFF": no fixed initial temperature (default setting)

Holding time 1 (hh.mm). Here: 2.00 hrs

End point of the programme (no additional segments)

Programme repetition. Here: 0 repetitions (default setting)

End temperature at which the furnace may be opened. Here: 200°C (default setting)

_ By using the Menu key , the respective programme parameters may be viewed in one direction only; they cannot be changed.

_ At the end of the parameter sequence, the standby mode is automatically displayed again.

PROGRAMMING THE CONTROLLER (PROGRAMME LOCATION 4)

Setting the controller from user mode (OP.b.) to programming mode (OP.P)

After the furnace has been switched on

Standby „ready“

Hold both keys simultaneously for approx. 3 seconds.

 - Press the key until OP.P appears.

Confirm the mode (password is required).

Call up the menu and select programme 4.

Password entry. Password: "4"

Call up the desired parameters
(see Chap. „Display of programme parameters“)

Confirm password.

Change the parameters.

Controller switches automatically to standby mode.

Press the Menu key until standby mode appears.

Hold both keys simultaneously for approx. 3 sec.

 - Press the key until O.P.U appears.

Confirm the mode.

User mode: parameters can no longer be changed (locked).

_ Hold the respective arrow keys  to reach the desired parameter value faster. As a result, the values are changed according to the decimal numbering scheme (after 10 counting units, the respective preceding value is changed).

MAINTENANCE/CARE

_ Wipe the upper surface of the furnace with a moist cloth when necessary. In no event must water enter the inside of the chamber! Do not clean with organic solvents.

_ Regular maintenance/calibration of the heating elements and insulation is not necessary. A visual check of the heating elements before starting the sintering program is recommended!

Exchanging the heating elements

_ If the indicator lamp of the heating  shines continuously after the start-up of the sintering program, this means that the heating element has failed.

_ If one heating element fails, always replace all 3 elements! If a heating element should fail again, the heating elements already used can be fitted again.

General instructions

_ During the installation of the heating elements in the furnace, a small amount of refractory material may adhere to the heating elements; this is harmless and will vitrify during the first few sintering processes.

_ Lamp test: in case not all of the two LEDs  and  flash up, the lamps are defective; please contact the customer service.

_ If many coloured objects are sintered, it is recommended to occasionally run a sintering programme with empty furnace.

_ Each furnace runs through one programme

cycle for test purposes before leaving the factory. This causes discolouration of the firing chamber and the panels.

- After 1-2 years, the heating elements may slightly adhere to the brick box. Should they need to be exchanged after this period, release the heating elements prior to pulling them out to prevent damage of the refractory material.

GENERAL MESSAGES IN THE DIGITAL DISPLAY

Automatic adjustment of the target and actual temperature values.

Cooling phase after sintering process (undefined cooling rate).

End of programme. The furnace may be opened.

Error message thermal element (contact the customer service).

Wrong password. Re-entry required.

Technical data:

L/B/H:	654 x 382 x 525 mm
Capacity of firing chamber:	1,1l
Power input:	2 kW
max. Temperatur:	1530 °C
Weight:	35 kg
E- fuse:	16 A
Power supply:	220-230 V 50/60 Hz

Packing list:

Order No.:	
178300	Ceramill Therm
	3 heating elements made of silicon carbide in separate protective packaging
	ceramic plate for furnace base
	mains cable
	sintering bowl and sintering balls
	Hexagon wrench key (4mm)

Accessories/Components:

Order No.:	
178360	Stackable sintering bowl (in stacked condition for use with Ceramill Therm 178350 only)
178311	Sintering balls 200 g
178361	Pliers for stackable sintering bowl
178313	Thermal element
178362	Heating element silicon carbide CT II, 1 piece)

**WARNING AND SAFETY INSTRUCTIONS:**

Ceramill Therm was developed for use in the dental industry. Commissioning and operation must be performed precisely according to the instructions for use. The manufacturer accepts no liability for damage caused by inappropriate use or improper handling. Furthermore, the operator is obliged to check the machine before using it. The operator himself is responsible for checking that the machine is suitable and fit for use for the designated purposes, insofar as these purposes are not specified in the instructions for use.

All repairs and/or servicing jobs, apart from those described under Maintenance/care, may only be carried out by the manufacturer or authorised service providers.

- _ **Before carrying out any repair or maintenance work on the furnace, always ensure that the furnace is not connected to the power supply!**
- _ **Wear protective equipment when carrying out repair work on the insulation.**
- _ **Handle the heating elements with great care (risk of breakage)!**

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Amann Girrbach AG
Amann Girrbach GmbH