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Working instructions

**Combilabor Induction Heated
Vacuum Pressure Casting Machine**

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These working instructions apply to::			
Order No.	Type	Features	Date
660 04 331	Heracast iQ	integrated water cooling and vacuum pump	01/2001

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Caution ! High frequency !

Warning for persons with heart pacemakers !

Read working instructions prior to use !

**We cannot accept any guarantee claims or assume
 liability if the machine is used for other purposes as
 stated or for damage resulting from non-compliance
 with these working instructions!**

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Confirmation of the manufacturer

We herewith confirm that the Heracast iQ is interference-screened in accordance with the regulations of VDE 0871 B, paper 163/1984 Vfg. 1046 .

The machine has been registered with the German Telekom AG. They reserve the right to verify production conforms with the regulations.

Hanau, January 2001
 Heraeus Kulzer GmbH & Co. KG,
 Grüner Weg 11, D-63450 Hanau

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1. Information on safe operation

Symbols

Symbol	Explanation
	Safety-relevant chapters and sections in these working instructions have been marked with this symbol.
	Information within the working instructions on the optimum use of the unit.
	Hot surface – risk of getting burned
	Caution! High frequency! Warning for persons with heart pacemakers

General

Please check the unit for transport damage and, if necessary, report the damage to the forwarder within 24 hours after receiving the unit.

EC Statement of Conformity

Herewith we, **Heraeus Kulzer GmbH & Co. KG, Grüner Weg 11, 63450 Hanau**, confirm that the following unit due to its intended use and the version marketed by us corresponds to the relevant basic safety and health requirements of the EC guideline. This statement will become invalid in case of any modification of the unit that is not coordinated with us.

Designation of the machine: Induction casting machine
Heracast iQ

Type of machine: Heracast iQ

Serial number.: 01 0005 ff.

Relevant EC guidelines:

Guideline on machines 89/392/EWG, appendix II A Elektromagnetic Compatibility 89/336/EWG

The unit fulfills the safety requirements of

- DIN VDE 0750 T1 05.82
- DIN VDE 0411 T1 06.86
- DIN 1988 T4 12.88

Working instructions

Using these working instructions and based on the work to be performed by the operator, easily understandable, written instructions have to be prepared and published in the language of the operators concerning the work at and with this machine.

Unit book

We recommend to keep a unit book. All tests as well as all essential works (e.g. repair work, modifications) must be documented in this book..

Safety information

With these laboratory units the safety concerning the protection of persons, the environment and the material to be processed mainly depends on the behaviour of the persons operating the unit.

Prior to operation read the working instructions carefully, adhere to the information provided in order to avoid errors and damage, in particular damage to the health.

In addition to the information in these working instructions, relevant national laws and guidelines must be observed for setting up and operating this unit (technical connection requirements of the electrical supply companies, etc.).

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The unit must not be used for heating up food.

The unit must not be used for works which include the release of combustible gases and vapors which burn with air or which may form a dangerous, explosive mixture.

The unit is not suitable for thermal treatment of hazardous dusts and fibres either.

The metal surface around the casting chamber heats up during continuous operation. Do not touch this surface.

 Due to the waste gases released during working, adequate ventilation must be provided.

 Power cable and plug must be checked for damage prior to operation. If any damage exists, the unit must not be connected to the mains.

Works at the electrical equipment of the unit must only be performed by the **authorized Heraeus service** and in the safe condition (voltage cleared).

Only permissible original spare parts must be used. The use of different parts holds unknown risks and must be avoided at any rate.

Proper function and safety of the unit are only guaranteed if the required tests, maintenance and repair work have been performed by **Heraeus service agents** or by adequately trained personnel.

Heraeus Kulzer GmbH & Co. KG **will not accept any liability** for damage to the unit resulting from inexpert repair which has not been performed by the **Heraeus service agents** or by adequately trained personnel or if no original spare parts or accessory parts have been used during the exchange of these parts.

2. Use in accordance with specifications

The induction casting machine Heracast iQ is a laboratory unit for casting of precious metal dental alloys and CoCrMo alloys for dental applications with a liquidus temperature of up to 1500°C except for titanium.

Precision castings with meltable alloy quantities of 130 g in the graphite crucible and 100 g in the ceramic crucible are possible. Alloy quantity for CoCrMo and non-precious metal bonding alloys in the ceramic crucible: up to 60 g

 We recommend the exclusive use of Original Heraeus Kulzer crucibles which are especially matched for this type of application. The use of other material leads to the exclusion of guarantee claims in case of damage to the unit or the molten material. In view of the various causes for faulty castings, we do not grant any guarantee for such cases.

Normally, the unit is suitable to be set up and operated in the following fields:
Commercial and industrial laboratories, schools, universities, hospitals, etc.

The unit has been designed for continuous operation.

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Working rules

Personal protective equipment such as hand, face and body protection must be worn; jewelry must be taken off prior to working.

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The unit must not be:

- used to heat up food.
- used for drying processes or thermal treatment which may lead to the release of combustible gases and vapors which burn with air or which may form a dangerous, explosive mixture.
- The unit is not suitable for thermal treatment of hazardous dusts and fibres either.

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3. Description

Features

- Safe and material-appropriate melting and casting under vacuum
- All types of precious metal alloys and almost all types of non-precious alloys (excluding titanium alloys) as well as aluminium and beryllium containing alloys can be melted and cast (→ hazardous to health).
- Temperature control via electronic power control and time
- Long lifetime of ceramic and graphite crucibles

☞ In order to avoid errors and problems during working we would like to ask you to read these working instructions as well as the golden booklets

“Precision-fit dental casting” and “A guide to processing errors”

carefully and store them at a place ready to hand. Further „Golden booklets“ are available on request. Tel. Customer Service, see last page.

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4. Scope of delivery/Original equipment

Scope of delivery

- 1 x Unit, working instructions, guarantee card, test certificate
- 1 x Filling hose with connector and container
- 1 x Muffle holder
- 1 x Pressure hose 13 mm, incl. 2 hose clamps and a socket 13 mm
- 2 x Door signs Caution! High frequency Germ./Engl.

Original equipment

- 1 x Small alloy spoon
- 1 x Box cont. small melting powder pellets
- 1 x Spruing aid – plexiglass half-pressing
- 1 pack 6 pcs ceramic crucibles for CL-IG/IM/IQ/195
- 1 pack 6 pcs ceramic crucibles NPM-crucibles for CL-IG/IM/IQ/195 (not for USA)
- 1 x Crucible box CL-IG/IM/IQ
- 1 Pack 10 pcs graphite inserts for CL-IG/IM/IQ/195
- 3 x each Casting ring (X3, X6, X9)
- 1 x each Cone former (X3, X6, X9)

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Photo 2) Front panel with control keys and display

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No. Description

1. Front panel with control elements
2. Viewing glass - casting chamber
3. Front limit switch - chamber lock
4. Locking lever – casting chamber (with carrying handle) and bow
5. Viewing glass for oil level of the pump
6. Holder for special tool muffle lift
7. Display
8. Adjustment knob
9. Keyboard
10. Main switch
11. Port for modem
12. Mains connection
13. Compressed air connection
14. Graduated viewing glass for cooling water level
15. Connection to fill/empty the cooling water tank
16. Filter cover air outlet
17. Rear limit switch
18. Viewing glass - casting chamber
19. Ceramic crucible
20. Graphite insert
21. Shielding plate
22. Mould support
23. Chamber gasket
24. Mould
25. Adjustable screw for size of mould
26. Clamping device, guide and snap-in locking device for fixation of the mould
27. Filter cover - air supply

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6. Location and installation

Transport

Carefully transport the unit horizontally to prevent the pump oil from leaking and damaging the unit. Packings and units must not be stacked. Shocks must be avoided!

For dimensions and weight refer to section "Technical data".

Unpacking

Remove straps. If required, screw in the carrying handles at the corners of the unit.

Set-up

Location: Table with load bearing capacity of at least 80 kg.

Table area : (w x h x d) 100 x 50 x 60 cm

The casting machine must be placed on a solid, skid-proof surface (laboratory desks, racks) so that a horizontal safe position is ensured. The surrounding temperature may amount to 40°C. Air inlet and outlet openings in the housing of the unit (rear and lower surface) must not be covered or blocked. Minimum distance to the wall: 100 mm.

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

Important information

- The rear side of the unit must be at least 10 cm away from the wall to avoid blocking the air outlet of the built-in fan.
- Do not place flammable material underneath the casting chamber.
- To transport the unit, leave the carrying handle under the chamber. Do not lift the unit at the chamber! Risk of damaging the unit!
- There must not be any objects within the turning range (horizontal and vertical) of the chamber.
- During operation - turning movement of the chamber 90° downwards!
- During turning movement do not touch/operate/hold chamber and lock.
- Do not lean on the chamber.
- Risk of damage and injury in case of non-adherence!

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Minimum distances to adjacent areas or pieces of equipment:

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Mains connection

- Mains: 200 - 250 V (AC), 1 P/N/PE, 50/60 Hz, according to VDE requirements and those of your local electricity board.
- Mains fuse: Separate safety fuse 16 A inert or safety cut-out fuse C 16 A
- Mains connection: The unit should not be connected via a fault current circuit breaker. If the use of a fault-current circuit breaker is specified by the local electricity board, type 30 mA should be used.

Compressed air supply

Note: The compressed air must be clean and dry!

- Line pressure min. 4 bar, max. 7 bar! Higher line pressure (even short-term) can result in damage to the inner valves! To avoid this, an optional pressure reducer filter combination must be used. Order No.: 646000743. In case of moist compressed air, this pressure reducer-filter combination must also be connected in series. The pressure reducer includes accessories for unit mounting or wall mounting.
- In case of moist compressed air a water separator/filter, Order No. 64600743 must be interconnected.
- In case of moist compressed air a water extractor must be interconnected Best. Nr. 64600743
- Line cross-section (inner) min. 10 mm
- Rapid pressurization (< 2 sec.) is essential for the mould filling behaviour. Non-compliance can result in faulty castings. The use of a separate compressed air tank (< 10 mm) in the direct vicinity of the casting machine is recommended for small line cross-sections. Best.Nr.: 64601059.

Note: Any guarantee claims shall be excluded in case of malfunctions or damage resulting from inadequate compressed air supply!

Filling the cooling system

For safety reasons the inner cooling system is only prefilled with a small quantity of antifreezing and anti-corrosion agent. To fill the system completely, the enclosed bottle must be used! For this purpose open the cap of the bottle and fill the container with tap water (do not use distilled water). Connect the hose to the connection piece at the rear side of the unit. Hold the bottle well above the unit. Press the plastic bottle slightly and shortly. Due to the difference in height the unit will be filled automatically; if required, make a hole in the container! Wait until the bottle is empty and remove the hose from the fill-in connection piece. Screw the lid off the bottle and fill the bottle with water again. Connect the bottle to the filling hose and continuously check the filling level of the reservoir (14) on the right side during filling up. If the water level has reached the value "MAX.", stop filling and remove the filling hose. For this purpose activate the unlocking lever above the connection with a tool, if required.

Important: Do not fill over the "MAX." marking. Risk of damaging!

After initial operation, the cooling water level must be checked again and the unit must be filled up with pure tap water. For maintenance work it has turned out to be helpful to reuse the **empty plastic bottles** to empty the cooling water tank; please do not discard these bottles.

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Room ventilation

The room in which the unit is operated must have sufficient technical ventilation. The unit must not be operated in recesses that can not be ventilated. If several units are to be placed in one room, special ventilation measures may be required (e.g. zone ventilation).

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7. Operation

The following pages are to provide basic information and hints which are essential for successful and error-free working. Please observe the order of the working steps.

Mostly, casting errors can not be attributed to the machine. In such cases please refer to the golden booklet: "A Guide to Processing Errors".

Putting into operation

1. Turn on the compressed air supply for the unit

2. Main switch "On" (I)

The test cycle starts after switching on the unit. The vacuum pump evacuates the chamber. Then the chamber is turned automatically and pressurization is started. Prior to turning back the chamber into the stand-by position, the chamber is deaerated automatically. In case of a defect or an error, a respective message is displayed. Once the test run has been successfully completed, the unit is ready for operation.

Information

In case of questions concerning error messages and the elimination of these, please contact our technical service, Telephone (0 61 81) 35 - 58 80.

Information

If an error can not be eliminated, please contact the responsible service agent. The addresses can be found on the last page of these instructions.

Unauthorized opening of the machine include unknown risks and are not permissible.

In case of questions concerning application-technical problems, please contact our Application-technical department, Telephone (0 61 81) 35 - 58 89.

Putting out of operation

- Remove crucible and mould from the chamber.
- Let the unit cool down for approx. 3 min. before switching it off (then turn off the cooling water supply).
- Switch off the unit; unplug the unit if it is not to be operated for longer periods.
- Remove contamination from the crucible and the chamber.

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Errors and causes


If a malfunction occurs during the test run or during operation, the error is shown in the display:

Error resp. warning message	Warning	Error
UDC too high	> 210 V	> 220 V
IDC too high	> 16,5 A	> 17A
Inner temperature too high	> 50°C (or < 3°C)	> 60°C (or < 3°C)
Oscillator temperature too high	> 70°C (or < 3°C)	> 75°C (or < 3°C)
Temperature of cooling liquid too high	> 50°C (or < 3°C)	> 65°C (or < 3°C)
PDC too high	> 2050 W	> 2100 W
15 V-supply too high	> 15.70 V	> 15.99 V
15 supply too low	< 14.30 V	< 14.00 V
Deviation of set-actual power	If -after 3 s -deviation > 5% Conditions: • if graphite crucible • if chamber in stand-by position	If - after 3 s- deviation > 5% Conditions: • if graphite crucible • if chamber in stand-by position • if unit in Auto-Test mode
Vacuum not built up completely	Vacuum > 500mbar after 15s	Vacuum > 500 mbar after 25s or if in "automatic test mode" Vacuum > 50 mbar after 45s
Check compressed air	Additional message in case of vacuum error	Additional message in case of vacuum error
Pressure too high	If pressure > 3.5 bar after 4s and after finishing the run-up phase	If pressure > 3.7 bar after 4s and after finishing the run-up phase
Pressure too low	If pressure < 2.9 bar after 4s and after finishing the run-up phase	If pressure < 2.7 bar after 4s and after finishing the run-up phase
Chamber movement		Chamber has not reached target position after 4s; Deviation >15 digits from set value
Chamber open		Chamber open Conditions: • Vacuum • Pressure • Chamber movement if required, check chamber limit switch and cable run
Check crucible	if - for graphite crucible - after 3s UDC > 180V	
Pressure build-up too slow	Pressure < 2.9bar after 4s	
Error during reading of EEPROM		No reply of EEPROM
EEPROM checksum error		SRS checksum error during reading of EEPROM
EEPROM data error		Incorrect data during reading the EEPROM
Error ADC	If no reply of A/D converter after more than 20 sec	
Temperature power pack 1	Signal from power pack 1	Signal from power pack 1
Temperature power pack 2	Signal from power pack 2	Signal from power pack 2
Error buffer overrun	Number of warnings exceeds the number of those which can be displayed	Number of warnings exceeds the number of those which can be displayed
Energy control		Automatic deactivation after 240000 Ws = 240 kW (e.g. 2 KW *2min)
Check pressure sensor	If during Auto-Test or automatic test mode a surrounding pressure of < 700 mbar or >1200 mbar is measured	
Auto-Test unsuccessful	If tried to change to casting menu and Auto-Test unsuccessful	
Memory test	Error during check of program memory	
Chamber position		If chamber position outside permissible range Chamber adjustment range tolerance ± 15 digits
Watchdog		Unauthorized abortion of program
Water flow		Step 0 autom. after 15 min. without HF-function, error message deactivated Step 1 < 200 ml Step 2 < 400 ml Step 3 < 400 ml Error display if error lasts at least 4 sec
Pressure build-up too slow		If pressure > 0.2 bar 4s after deaerating
Leakage		if at 50 mbar vacuum drops at least 30mbar in 5 seconds or if during pressurization a drop of at least 0.3 bar in 5 seconds

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8. Working with the Heracast IQ

After switching on and successful completion of the Auto-Test, the user is automatically directed to the main menu. In this menu additional options are available.

8.1 Main menu – casting (diagram)

Description

The crucibles “graphite” or “ceramic” can be selected with the adjustment knob. For this purpose turn the knob to the left or the right until the marking “*” marks the desired value. Confirm your selection by slightly pressing the adjustment knob (action point).

Proceed identically in the menu “Enter quantity/liquidus temperature”. The values can be increased or lowered by turning the adjustment knob to the left or right.

All preset data are displayed in the “Casting menu”. To adjust a different vacuum value, also turn the knob to the left or the right (without pressing).

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Start the premelt now by pressing the key under the menu item "VS". Press "Stop" key after finishing the premelt.


The main melt can always be carried out automatically or manually.

- Automatic main melt / automatic turning chamber is automatically turned (after a given or calculated delay time)***
- Manual main melt / automatic turning chamber is turned automatically (decision of the operator)

The decision for automatic or manual operation mode depends on various factors and is finally made by the operator depending on the alloy that is used.

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Description

The settings in the Setup menu allow individual configuration of the user surface with country- and user-specific conditions and are mainly performed using the adjustment knob.


The individual menus can be selected one after the other by turning the adjustment knob.

The current menu item is marked by "*". To select the options of a menu item, the adjustment knob must be pressed slightly.

Values can now be adjusted according to your wishes. If you want to end the settings of a menu item, press the adjustment knob again. By turning the adjustment knob to the left or right, you can select another menu item.

The Setup menu can be quit in two ways:

- **Save**
All your Setup changes will be saved and run
- **Cancel**
The Setup changes you have made will be deleted and the unit maintains its previous settings.

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Description

The password is also entered with the adjustment knob (press after each individual number). The user password is „00000“. It has been prestored and needs only to be confirmed by pressing the “OK” key. Current menus are marked with “*“. Scrolling through the individual options is possible by turning the knob to the left or the right. Press the knob to select an option. Rights concerning your password can be taken from the table above.

Note

⚠ Customers do not have access to levels 2 and 3 since the settings performed in these levels might affect the proper function of the unit so that damage in case of maladjusting can not be excluded.

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9. Teleservice, configuration of modem

9.1 Teleservice

The casting machine Heracast iQ features a serial port and a system processor. In conjunction with a modem these two features allow for Teleservice. It is possible to log in each Heracast iQ casting machine all over the world if this casting machine is also connected to the telephone network via a modem.

The Teleservice offers the advantage of a fast remote diagnosis or a simple transfer of a software update.

Depending on the telephone network, an analog modem or an ISDN modem is required for Teleservice. The modem is connected to the serial port of the casting machine. The modem is connected with the telephone network.

Authorized service agents may contact the connected machine, check the machine and prepare a diagnosis.

9.2 Configuration of modem

The menu item "Configuration of modem" checks the compatibility of the connected modem. If the modem is compatible, "OK" is displayed in the menu.

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10. Melting and casting

The technical procedure is described in the following. For alloy-specific application information, see Alloys, section 11.

10.1 General

Information brochures on casting according to the Heraeus system:

- Precision-fit dental casting
- A Guide to processing errors

Heraeus spruing system

Information can be found in section 12 of these working instructions.

10.2 Suitable investment materials

No graphite-containing investment materials must be used. The graphite content can result in degassing or in damage to the alloy. We recommend to use our graphite-free, phosphate-bonded Heravest investment materials.

10.3 Premelting

In order to obtain uniform casting conditions, all alloys are pre-melted. The mould is placed into the casting chamber after pre-melting the alloy.

Important notice:

The alloy must be observed continuously during pre-melting. Generally, each melt may only be observed through the blue glass because of the high luminous intensity. Risk of getting blinded! The pre-melting process must only be interrupted when all the alloy has melted! The alloy has a spherical shape; there are no protruding edges of the molten metal. Splitting up of the oxide layer is without any relevance during melting in the graphite crucible. When melting large quantities (more than 50 g), small individual portions must be pre-melted. The metal should only be melted so that filling in of the next portion is possible. Only in the last pre-melt all the alloy is melted.

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Exception: CoCrMo alloys (see section 11.3)

10.4 Casting

After premelting, the mould is rapidly placed in the casting chamber and locked in place (ideally approx. 20 – 40 seconds); the chamber is closed and the melting process is started by pressing the key “HS auto.“ or “HS man.“

Important: The interruption between pre- and main melt should not exceed one minute to prevent the melt and the mould from cooling down too much.

Each crucible should only be used for one alloy to avoid mixing of alloys. Use a felt pen to mark the crucibles according to the intended use.

Graphite insert

☞ When using graphite inserts, the alloy is heated under continuous observation by the user up to the set point of time. After the set time is reached, the chamber is turned automatically or manually.

If - in case of manual turning - casting is not started immediately after the set time is reached, the alloy may be damaged through overheating.

- If graphite inserts lose some of their height or the upper border becomes thin and brittle, they must be exchanged → minimum weight of graphite insert approx. 7 g.
- If a lot of melting powder has accumulated in the graphite insert, no melting powder should be added for the next casting.
- Prior to each casting blow out graphite insert (cleanliness!)

Ceramic crucible

When using ceramic crucibles, ceramic bonding resp. palladium-based alloys are heated until the oxide layer splits open. When casting CoCrMo and NPM alloys, melting is continued until the shadow disappears (see picture phase 5, section 13) and, depending on the type of alloy, casting is performed with or without delay time.

- The sprues and the casting buttons must be cut prior to melting. When filling the crucible, make sure that the alloy is placed in the crucible as far down as possible to obtain optimum and even melting power.
- When melting CoCrMo alloys, the correct position of the cylinders must be ensured.

11. Alloys

11.1 Gold casting and high gold content ceramic bonding alloys Universal and silver-palladium alloys

Processing: Alloy quantities: 5 g to 130 g
Gold casting alloys: Use ceramic crucible with graphite and melt with melting powder pellet (small!)
Ceramic bonding alloys: Use ceramic crucible with graphite insert! Melt without melting powder pellet!

In the Heracast iQ the melting powder pellets are used in exactly the reverse order as in the CL-G, CL-G 77, G 94, G 97 und G 2002 since the alloys are melted under vacuum which leads to reduced formation of protective gas.

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11.2 Reduced precious metal ceramic bonding and palladium-based alloys

Processing: Use ceramic crucible without melting powder pellet.

Alloy quantity: 15 g to 100 g (Min. quantity of 20 g for reduced precious metal content bonding alloys)

11.3 CoCrMo and non-precious metal alloys

Processing: Use NPM ceramic crucible without melting powder pellet

Important: The NPM ceramic crucibles are exclusively suitable for casting NPM alloys and feature a longer lifetime than normal ceramic crucibles. If NPM crucibles are used for precious metal alloys, the silicon might be damaged.

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
Caution!


- Under no circumstances must the machine be left unattended!
- The melt must always be observed!
- The viewing glass must be used for each melting process!
- Risk of getting blinded due to the high luminous intensity of the melt.
- No hot melting crucible must remain in the coil after working.

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• Heraenium NF:

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Preparations:

- Place ceramic crucible into the coil
- Fill alloy into the crucible
- Close chamber

Premelting (phase 1, 2, 3)

- Press key “VS”
- Vacuum between 250 mbar
- Observe heating process of the alloy (phase 1, 2, 3)
- when phase 3 sets in press “Stop” key

Insert the casting mould

- Open chamber

- Position the casting mould

- Close chamber

Main melt

- Press key “HS”

- Then select between “HS auto.” or “HS man.”

- In “HS auto.” mode the delay time expires and the chamber is turned automatically.

- In “HS man.” mode observe melting process of the alloy (phase 4, 5)

! Only in “HS man.” mode !**• Heraenium CE:**


Immediately turn casting chamber after reaching phase 6 (disappearance of the shadow) - key “Schwenken” (turn)

• Heraenium EH:


3 seconds after phase 6 has been reached, turn chamber key “Schwenken” (turn)

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



14. Maintenance

Proper function and operational reliability of the unit are only ensured if the required tests, maintenance and service work are performed by Heraeus Kulzer service agents or personnel instructed by Heraeus Kulzer.

Heraeus Kulzer GmbH & Co. KG will not accept any liability for damage - in particular personal injury - resulting from improper maintenance or repair work which has not been performed by Heraeus service agents or adequately trained personnel or if no original spare parts or accessory parts have been used during the exchange of parts.

We recommend to conclude a maintenance contract with our service agents; an offer can be requested there (see service agents, section 19).

Maintenance and care

Note:

Prior to maintenance and service work the unit must be switched off and unplugged!

Strict adherence to the maintenance intervals is recommended to avoid faulty castings and damage to the unit.

Maintenance work types:

a) daily

- Check water level
- Check water level and remove contaminations (e.g. milky consistency)
- Clean the inside and the outside of the viewing glass using a soft cloth.
- Remove alloy residues, etc. from the rubber seal of the casting chamber (vacuum cleaner or blow out with compressed air).
- In case of moist compressed air: empty the water extractor of the optional pressure reducer through the release valve on the bottom of the casting machine.

b) weekly (of after 100 castings)

- Clean the inside of the casting chamber, the mould holder and the rubber seal when the unit has cooled down.
- Check function and proper activatability of the chamber lock. Adjustment and cleaning may only be carried out by adequately trained personnel.

c) annually (or after 10.000 castings)

Note:

The following work includes interventions into the unit and may only be performed by adequately trained service personnel!

- Change oil of the vacuum pump.
- Clean oil mist filter of the pump.
- Check casting chamber lock and support, check welding seams.
- Check circulation of cooling water and flow control instrument (manual test program). If required, change cooling water.
- Change filter for vacuum and compressed air in the valve box
- Check pressure and vacuum hoses.
- Check activatability of the turning mechanism, wear and firm seat of the gears
- Performance check (manual test run) with graphite crucible (2000 W).
- Clean filter mats (rear and bottom surface) (16/27), if required, exchange them

d) after 30.000 casting or after 5 years

- Exchange locking handle and chamber lock.

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Tests

Works at the electrical equipment of the unit must only be performed by the authorized Heraeus service agents or adequately trained personnel and in the safe condition (voltage cleared). Only permissible original spare parts must be used.

Maintenance

Permissible parts and accessories:
Proper function and operational reliability of the unit are only ensured if permissible original spare parts are used. The use of other parts includes unknown risks and must be avoided at any rate.

Disposal

The unit is designed to be operated for 10 years. As far as the disposal of spare parts or the unit is concerned, please contact Heraeus Kulzer at Hanau, Service Dept., addresses see last page.

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
15. Repair / Maintenance work

The following work includes interventions into the unit and must only be performed by adequately trained service personnel!

15.1 Exchanging the pump oil

The unit must be switched off and unplugged before filling in or filling up oil! Open the unit.

- Place paper towel underneath the pump head to absorb leaking oil.
- Loosen oil discharge screw (B) using an Allan key.
- Hold cup under the oil discharge opening, empty the pump.
- Tighten the oil discharge screw.
- Loosen the oil filling screw (A) using an Allan key.
- Fill in the required quantity (max. 60 ml), tighten the screw.
- Check oil level in the viewing glass (C).

15.2 Exchanging the pump filter

The unit must be switched off and unplugged to exchange the pump filter (D)! Open the unit.


- Carefully remove the pump filter by using pliers (turn to the left).
- Place O-ring over the thread, screw in new filter (clockwise), do not overturn!

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
15.3 Exchanging the filter in the valve box

The unit must be switched off and unplugged to exchange the filter in the valve box!

Disconnect compressed air supply! Open the unit.

- Carefully open the unit using an Allan key.
- Change prefilter element and main filter element

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15.4 Exchanging the cooling water

- Connect the hose of the bottle to the respective socket (13) at the rear of the casting machine.
- Remove the adhesive strip from the deairing hole of the bottle.
- Press the plastic bottle once and seal the small hole with the adhesive strip again.
- Hold the bottle clearly lower than the casting machine, the tank is emptied automatically.
- For filling, see section 6 "Filling the cooling system".

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
16. Technical data

	Value	Unit
Dimensions (w x h x d)	660 x 500 x 550	mm
Required table area (w x h x d)	1000 x 500 x 560	mm
Mains connection:	200 - 250	V
Rated frequency:	50/60	Hz
Rated power:	2,6	KVA
Generator power:	2	KW
Compressed air supply:	4 - 7	bar
Minimum cross-section of compressed air supply (inner):	10	mm
Weight:	65	kg
Vacuum:	< 50	mbar
Noise level	< 70	dB (A)
Protective measures - protection class:	1	
Protection type:	IP 10	

Fuse protection: For the connection to the mains the electrotechnical rules and the technical regulations of the local Electricity Board have to be observed. Cut-out fuse 16 A or automatic circuit breaker C16 A.

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



Alloy		Alloy Type		Alloy Type	
Hera PF	Graphite + MP	Gold casting / Universal	890°C		
Hera SG	Graphite + MP	Gold casting / Universal	895°C		
Bio Maingold SG, Hera GG	Graphite + MP	Gold casting / Universal	920°C		
Maingold SG, Bio-Maingold IT	Graphite + MP	Gold casting / Universal	930°C		
Maingold G, MP	Graphite + MP	Gold casting / Universal	960°C		
Maingold GV	Graphite + MP	Gold casting / Universal	970°C		
Mainbond A, Maingold OG, Bio-Maingold TK	Graphite + MP	Gold casting / Universal	990°C		
Alba SG	Graphite + MP	Gold casting / Universal	1000°C		
Mainbond EH	Graphite + MP	Gold casting / Universal	1010°C		
Bio-Maingold I, Bio-Heranorm	Graphite + MP	Gold casting / Universal	1035°C		
Maingold I, Keramikgold PKF	Graphite + MP	Gold casting / Universal	1050°C		
Maingold W	Graphite + MP	Gold casting / Universal	1060°C		
Hera KF, Mainbond KF	Graphite + MP	Gold casting / Universal	1070°C		
Bio-Herador GG	Graphite	High gold content / silver palladium	1110°C		
Herador GG	Graphite	High gold content / silver palladium	1125°C		
Bio-Herador SG/N	Graphite	High gold content / silver palladium	1130°C		
Herador C	Graphite	High gold content / silver palladium	1135°C		
Herador MP	Graphite	High gold content / silver palladium	1140°C		
Alba KF	Graphite	High gold content / silver palladium	1165°C		
Herador S/SG	Graphite	High gold content / silver palladium	1150°C		
Herador PF	Graphite	High gold content / silver palladium	1160°C		
Herador G, H	Graphite	High gold content / silver palladium	1200°C		
Heradent E	Graphite	High gold content / silver palladium	1235°C		
Heradent NH	Graphite	High gold content / silver palladium	1260°C		
Albaloy	Ceramic	Red. prec. met. cont. / palladium-based		6 sec after splitting up	
Albabond / E / EH / U / GF / A / B / C	Ceramic	Red. prec. met. cont. / palladium-based		6 sec after splitting up	
Herabond N	Ceramic	Red. prec. met. cont. / palladium-based		6 sec after splitting up	
Heralay G / U, Herador P	Ceramic	Red. prec. met. cont. / palladium-based		6 sec after splitting up	
Heranium CE	Ceramic	CoCrMo / NPM bonding		0 sec after shadow	
Heranium EH	Ceramic	CoCrMo / NPM bonding		3 sec after shadow	
Heranium NF (Pre-melt end after alloy has completely melted)	Ceramic	CoCrMo / NPM bonding		6 sec after shadow	

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	High gold content bonding alloys	Casting Pd-based alloys Gold-Pd alloys	NPM / CoCrMo	
<ul style="list-style-type: none"> • Check compressed air supply, if required, open • Switch on main switch • Auto-Test • Unit is ready for casting 	<ul style="list-style-type: none"> • Select menu graphite crucible and gold casting • Enter liquidus temperature and quantity • Place ceramic crucible with graphite insert in the coil • Adjust muffle size and height • Add metal in small pieces into the crucible • Close chamber • Press key "VS" (pre-melt) • Check vacuum and change, if required (recommendation approx. 100 mbar) • Observe melt, interrupt pre-melt with Stop key after the alloy has completely melted • Open casting chamber and add melting powder pellet • Place muffle in chamber and lock in place, close chamber, • Press key "HS auto" • When the time has expired, the chamber is turned • After approx: 60 sec the chamber is turned back, open chamber and remove casting mould 	<ul style="list-style-type: none"> • Select menu graphite crucible and high gold content bonding alloy • Enter liquid. temp. and quantity • Place ceramic crucible with graphite insert in the coil • Adjust muffle size and height • Add metal in small pieces into the crucible • Close chamber • Press key "VS" (pre-melt) • Check vacuum and change, if required (recommendation approx. 100 mbar) • Observe melt, interrupt pre-melt with Stop key after the alloy has completely melted • Open casting chamber (no melting powder pellet) • Place muffle in chamber and lock in place, close chamber, • Press key "HS auto" • When the time has expired, the chamber is turned • After approx. 60 sec the chamber is turned back, open chamber and remove casting mould 	<ul style="list-style-type: none"> • Select menu ceramic crucible and Pd-based alloy or Gold-Pd alloy • Enter delay time 6 s • Place ceramic crucible in chamber • Adjust muffle size and height • Add metal in small pieces into the crucible • Close chamber • Press key "VS" (pre-melt) • Check vacuum and change, if required (recommendation approx. 100 mbar) • Observe melt, after splitting up of the oxide layer interrupt pre-melt with "Stop" key • Open casting chamber (no melting pellet) • Place muffle in chamber and lock in place • Press key "HS" • After splitting up of the oxide layer press key "HS auto." • When the time has expired, the chamber is turned • After approx. 60 sec the chamber is turned back, open chamber and remove casting mould 	<ul style="list-style-type: none"> • Select menu ceramic crucible and NpM • Enter delay time depending on the alloy • Place ceramic crucible in chamber • Adjust muffle size and height • Ensure horizontal, parallel, placing of metal cylinders • Close chamber • Press key "VS" (pre-melt) • Check vacuum and change, if required (recommendation approx. 250 mbar) • Observe melt, depending on alloy type interrupt pre-melt with "Stop" key • Open casting chamber (no melting pellet) • Place muffle in chamber and lock in place? • Press key "HS" • After shadow appears, press key "HS auto." • When the time has expired, the chamber is turned • After approx. 60 sec the chamber is turned back, open chamber and remove casting mould
<p>Putting unit out of operation</p>				
<ul style="list-style-type: none"> • Remove crucible after casting • the unit can be switched off after 5 minutes. • If required, turn off compressed air supply 				

159
120
160

112
7
113
115
118

136
135
101
102
103

118
146
156
161
162
104
157
105

58 74 75 73 72 78 79 80 61 163 60 62 64 66 1
63 68 84 143 65 136 135 90
152 138 123 147 122 124 88 144 106 107 111 112 114 83 80 79 91
136 135 80 79 92 94 95 93 100 97 99
154 131 130 132 148 149 150