SafeBasics• True Fit Thin™



TECHNICAL DATA SHEET

FEATURES & BENEFITS

- The unrivaled strength of nitrile delivers enhanced barrier protection and exceptional fit, while providing the tactile sensitivity necessary for the most complex procedures.
- Designed to fit like a glove should, while providing the tactile sensitivity needed for the most complex procedures.
- o Textured fingertips for a secure grip, even when wet.

RECOMMENDED USE

 Non-sterile, single use, medical examination gloves are intended to prevent transmission of a wide variety of diseases to both patients and health care personnel. Specifically, gloves are intended to act as a barrier between the healthcare provider and patients to prevent contamination between patient and caregiver.



True Fit Thin

ISO 13485 certified manufacturer. Class 2 Medical Device. Health Canada registered medical device manufacturer.

STORAGE CONDITIONS

These gloves should be disposed of after each use. Store in a cool, dry place. Avoid excessive heat (40 °C or 104 °F).

Recommended storage temperature: 10-30°C. Shield open box from direct sun and fluorescent light. Shelf Life: 3 years.

PRODUCT DETAILS

Product Name: Medicom SafeBasics® True Fit Thin™

Powder-Free Nitrile Medical

Examination Gloves

Product Type: Single-use, non-sterile, ambidextrous

Material: Nitrile
Colour: Blue
Quantity per Box: 300

Specs: 3 mil, 3.5 g (Blue), 2.5 mil, 3.2 g (Pink)

Country of Origin: Malaysia

GLOVE SPECIFICATIONS									
	CHARAC	TERISTIC	ACCEPTABLE LIMIT						
Dimension & Tolerance			Extra-Small (XS)	60 - 80					
			Small (S)	70 - 90					
			Unisize (S-M)	75 - 95					
		Width (mm)	Medium (M)	85 - 105					
			Large (L)	100 - 120					
			Extra-Large (XL)	110 - 130					
			XX-Large (XXL)	120 - 140					
	Length (mm)		XS & S ≥ 220 Others: ≥ 230						
	Thickness (mm)	Palm	≥ 0,05						
		Finger	≥ 0,05						
Physical Requirement	Before Ageing	Tensile Strength	≥ 14 MPa						
		Ultimate Elongation	≥ 500%						
	A.C. A .	Tensile Strength	≥ 14 MPa						
	After Ageing	Ultimate Elongation	≥ 400%						
Powder Residue	Powder Free Gloves		≤ 2,0 mg/glove						
		Powdered Gloves	≤ 10,0 mg/dm2						







CHEMOTHERAPY CHART						
CHEMOTHERAPY DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME					
Carboplatin / Carboplatine, 10,0 mg/ml	>240 minutes					
Carmustine, 3,3 mg/ml	11.3 minutes					
Cisplatin (BCNU), 1,0 mg/ml	>240 minutes					
Cyclophosphamide (Cytoxan), 20,0 mg/ml	>240 minutes					
Dacarbazine (DTIC), 10,0 mg/ml	>240 minutes					
Doxorubicin Hydrochloride / Chlorhydrate de doxorubicine, 2,0 mg/ml	>240 minutes					
Etoposide (Toposar) / Étoposide (Toposar), 20,0 mg/ml	>240 minutes					
Fluorouracil / Fluorouracile, 50,0 mg/ml	>240 minutes					
Ifosfamide, 50,0 mg/ml	>240 minutes					
Methotrexate / Méthotrexate, 25,0 mg/ml	>240 minutes					
Mitomycin C / Mitomycine C, 0,5 mg/ml	>240 minutes					
Mitoxantrone, 2,0 mg/ml	>240 minutes					
Paclitaxel (Taxol), 6,0 mg/ml	>240 minutes					
Thiotepa / Thiotépa, 10,0 mg/ml	7.4 minutes					
Vincristine Sulfate / Sulfate de vincristine, 1,0 mg/ml	>240 minutes					

^{*}The aforementioned chemicals have been tested with these gloves.

WARNING: DO NOT USE WITH CARMUSTINE AND THIOTEPA.

CAUTION: Testing showed an average breakthrough time of 11.3 minutes with Carmustine and 7.4 minutes with Thiotepa 10.0 mg/ml.

Gloves used for protection against chemotherapy drug exposure should be selected specifically for the type(s) of chemicals used. Users should review drug labeling or material safety data sheets for the drugs being used to determine an adequate level of protection.

The gloves have been tested according to ASTM D6978. Tested for use with chemotherapy drugs.

Store in a cool dry place. Avoid excessive heat (40 °C/104 °F). Shield open box from direct sun and fluorescent lighting.



Appendix 1 (List of items)											
Item Code	Description	Packaging Information	Size	Color	Treatment	Marketed in USA as	Marketed in Canada as				
1185	Medicom SafeBasics "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 8 Boxes/Case, 2400 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Non-Sterile	Medical Device Class 1	Medical Device Class 2				
1186	Medicom SafeBasics "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 8 Boxes/Case, 2400 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large	Pink	Non-Sterile	Medical Device Class 1	Medical Device Class 2				

AMD Medicom Inc. | 2555 Chemin de l'Aviation, Montreal, Quebec, Canada H9P 2Z2

APP1 04.02-04EN_5M Date: 11-Jan-2023 Page 1 / 1