

Redcap™

Emergency first response for chemical splashes.



* Excludes product and packaging.

aura³: act fast with the right kit.

Our New Aura3 chemical splash kit is better for the workplace and better for the planet, providing quick and easy access to Redcap eyewash buffered phosphate solution for neutralising acid and alkali splashes to the skin and eyes.. Aura3 represents a first for the first aid industry, the first 100% recycled and fully recyclable first aid kit box on the market.

The fresh new box design feature colour coded, embossed, and screen-printed lids. The new graphics stand out from a distance, enabling quick identification in an emergency, especially in a busy workplace. Compact and convenient, the integral bracket allows easy placement and fixing.

After eye or skin contact with acid or alkalis, it is only seconds before the pH value is significantly altered, and the corrosive damage may lead to permanent injury. Acid and alkalis must be quickly neutralised to a pH value equal to the normal level in the eye or on the skin. The use of ordinary water or 0.9 % sodium chloride will only provide a limited degree of effective dilution with an increased flushing time.

The neutralising effect of Redcap, reduces flushing time and volume of fluid needed to dilute acid and alkali splashes*, therefore reducing patient discomfort and secondary risk to the first aider.



Integral bracket



Single clip



Embossed icon



Laser etching, no labels



Includes eye bath



and eye pad dressings

Better for the workplace, better for the planet.



+44 (0) 08456 448808 | sales@reliancemedical.co.uk | www.reliancemedical.co.uk

RELIANCE MEDICAL
'Leading First Aid'



REV210323

Easy to use.



Before approaching a casualty that has chemical burns, you must first protect yourself by wearing the appropriate PPE.

If needed, call the emergency services.

Identify the chemical.



CHEMICALS ON THE SKIN

Carefully remove contaminated clothing. Drench the area with fluid for at least 20 minutes.

CHEMICALS IN THE EYES

Place the casualty in a comfortable position with their head inclined, so the affected eye is lower..

Irrigate the whole eye, including under the eyelid. Buffered Phosphate Solution offers a better chance of saving an eye than standard Sodium Chloride (Saline), Eye Wash. As well as diluting and flushing out the eye, the Phosphate has a neutralising effect, this helps to restore PH levels of the eye back to normal.

Contaminant:

ACIDS, ALKALIS

- **Action:** To quickly neutralise the pH value and prevent corrosive damage to the eyes and skin.
- **Directions:** Rinse with Redcap™ Buffer solution. Continue rinsing with saline eye wash or water until medical help arrives.



Quicker neutralisation, especially on the face, can relieve the distress of feeling 'underwater' from flushing, the first aider is also better protected from becoming a secondary victim whilst in contact with the contaminant.

The Aura3 Chemical Splash kit including Redcap Buffered Phosphate Solution is perfectly suited to enable first responders to achieve this goal.

This kit contains an easy-to-squeeze bottle, allowing the eye to be irrigated easily and efficiently, eye pods and eye pads. The box can be easily wall mounted with the integrated bracket.



5970 EYE WASH AND SKIN FLUSH KIT
in Medium Red Aura3 Box

SIZE: 25.5cmH x 28cmW x 10cmD
UNIT: 1
MOQ: 10

CONTENTS

5970

DOUBLE

Dressings, Eye Pad	2
Redcap™, Bottles 500ml	2
Eye Wash, Pods 20ml	2
Mirror	1

**It is recommended to continue flushing with water or saline eye wash until medical help arrives.*



When you see this symbol with our products, you are guaranteed we have reduced plastic in the mould, used 100% recycled materials, 100% recyclable product.



+44 (0) 08456 448808 | sales@reliancemedical.co.uk | www.reliancemedical.co.uk

RELIANCE MEDICAL
'Leading First Aid'

