Technical	
Data	
Sheet	

Product specification

1. Product info.	Name: COVID-19 Antigen Detection Kit-Nasal Swab Model: 1test/kit,5tests/kit, 25tests/kit
2. Description	This product is used to qualitatively detect the novel coronavirus (SARS-CoV-2) antigen in human anterior nasal swabs .
3. Features	Our products have the characteristics of high sensitivity, good specificity and convenient operation.
4. Intended use	This product is used to qualitatively detect the novel coronavirus (SARS-CoV-2) antigen in human anterior nasal swabs . It can be used for rapid investigation of suspected novel coronavirus cases, and can also be used as a reconfirmation method for nucleic acid detection in discharged cases. A positive test result requires further confirmation, and a negative test result cannot rule out the possibility of infection. This kit is for home use by laymen in a non-laboratory setting (such as person's home or certain non-traditional sites such as offices, sporting events, airports, schools etc.). It is recommended to conduct a comprehensive analysis of the condition based on the patient's clinical manifestations and other laboratory tests.

Technical information

	Component name	Material
	Test card	nitrocellulose membrane
1. List of Materials	Extraction tube&Tube Cap	PE+PP
	Sample collection swab	РР
2. Latex free	Yes	
3. Characteristic	Sensitivity	97.1%
	Specificity	99.2%
	Accuracy(Total coincidence rate)	98.4%
	Minimum detection limit	When the virus content is greater

	than 400TCID ₅₀ /ml, the positive detection rate is greater than 95%. When the virus content is less than 200TCID ₅₀ /ml, the positive detection rate is less than 95%, so the minimum detection limit of this product is 400TCID ₅₀ /ml.
Precision	Three consecutive batches of reagents were tested for precision. Different batches of reagents were used to test the same negative sample 10 times in succession, and the results were all negative. Different batches of reagents were used to test the same positive sample 10 times in succession, and the results were all positive.
HOOK effect	When the virus content in the sample to be tested reaches 4.0*10 ⁵ TCID ₅₀ /ml, the test result still does not show the HOOK effect.

The test card of this product was tested with virus samples within the following concentration range, and the results showed that there was no cross-reaction.

No.	Name	Concent ration	No.	Name	Concentr ation
1	HCOV-HKU1	10⁵TCID₅ ₀/ml	17	H7N9	10 ⁵ TCID ₅₀ /ml
2	Staphylococcus aureus	10 ⁶ TCID₅ ₀/ml	18	H5N1	10⁵TCID ₅₀ /ml
3	Group A streptococci	10 ⁶ TCID₅ ₀/ml	19	Epstein-Barr virus	10 ⁵ TCID ₅₀ /ml
4	Measles virus	10⁵TCID₅ ₀/ml	20	Enterovirus CA16	10⁵TCID ₅₀ /ml
5	Mumps virus	10⁵TCID₅ ₀/ml	21	Rhinovirus	10 ⁵ TCID ₅₀ /ml
6	Adenovirus type 3	10⁵TCID₅ ₀/ml	22	Respiratory syncytial virus	10 ⁵ TCID ₅₀ /ml
7	Mycoplasmal pneumonia	10 ⁶ TCID₅ ₀/ml	23	Streptococcus pneumoniae	10 ⁶ TCID ₅₀ /ml
8	Paraimfluenzavi rus,type2	10⁵TCID₅ ₀/ml	24	Candida albicans	10 ⁶ TCID ₅₀ /ml
9	Human metapneumovir us	10⁵TCID₅ ₀/ml	25	Chlamydia pneumoniae	10 ⁶ TCID ₅₀ /ml
10	Human coronavirus OC43	10⁵TCID₅ ₀/ml	26	Bordetella pertussis	10 ⁶ TCID ₅₀ /ml
11	Human coronavirus 229E	10⁵TCID₅ ₀/ml	27	Pneumocystis jiroveci	10 ⁶ TCID ₅₀ /ml
12	Bordetella parapertusis	10 ⁶ TCID₅ ₀/ml	28	Mycobacterium tuberculosis	10 ⁶ TCID ₅₀ /ml

4. Cross-reactivity

-	13	Influenza B	10 ⁵ TCID ₅	29	Legionella	10 ⁶ TCID ₅₀
_		Victoria STRAIN	₀/ml		pneumophila	/ml
_	14	Influenza B	10 ⁵ TCID ₅	30	MERS-	10 ⁵ TCID ₅₀
_		YSTRAIN	₀/ml	50	coronavirus	/ml
-	15	Influenza A	10 ⁵ TCID ₅	31	Human	10 ⁵ TCID ₅₀
		H1N1 2009	₀/ml	51	coronavirus NL63	/ml
	16	Influenza A	10 ⁵ TCID ₅			
_		H3N2	₀/ml			

The test results do not be interfered with the substance at the following concentration:

	No.	Name	Conce ntratio n	No.	Name	Conc entrat ion
_	1	Whole Blood	4%	9	Oseltamivir	5mg/ ml
-	2	Ibuprofen	1mg/ ml	10	Naphazoline Hydrochloride Nasal Drops	15%
_	3	tetracycline	3µg/m I	11	Mucin	0.5%
_	4	chloramphenicol	3µg/m I	12	Compound Benzoin Gel	1.5m g/ml
_	5	Erythromycin	3µg/m I	13	Cromolyn glycate	15%
_	6	Tobramycin	5%	14	Deoxyepinephrine hydrochloride	15%
_	7	menthol	15%	15	Afrin	15%
_	8	Mupirocin	10mg/ ml	16	Fluticasone propionate spray	15%

6. Sterilization metho	d	Non-sterile	
7. Shelf Life		24 months	
8. Packaging specification	6.1 Sales unit	1、5、25	Per kit

REV	12.12.2021	01 Date
-----	------------	---------