



Pulse Oximeter INV-430J/PE

Precautions

- Do not attempt to repair the Pulse Oximeter by yourself. Only certified engineers should maintain and repair it.
- Change the contact position between the Oximeter probe and the finger periodically if you are monitoring your SpO₂ levels and pulse rate for more than 2 hours.
- This product is not designed to be used by newborn babies.
- Seek for medical care if the measured value goes beyond the normal range and you are sure that the instrument is not malfunctioning.
- The pulse oximeter uses infrared light (invisible to your eyes) to measure your SpO₂ levels. Hence, please do not stare at the light-emitting components of the Oximeter, as that could cause harm and/or potentially blind your eyes.
- For details about clinical limitations and contraindications, please carefully consult relevant medical literatures.

The following factors may affect the accuracy of measurement:

- The Oximeter is used in an environment involving high-frequency devices, such as high-frequency electric knives and CT apparatuses.
- The probe of the Oximeter is placed on the same arm that a blood pressure cuff arterial duct or intravenous injection.
- The user suffers from hypotension, severe vascular atrophy, severe anemia, or low oxygen.
- The user is in sudden cardiac arrest or shock state.
- The user is wearing nail polish or artificial nails.

Warning: Do not use the Oximeter in an environment with any flammable gases, flammable anesthetic, or other flammable substances.

Warning: Keep unit and lanyard away from children as the included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard.

Warning: Do not throw the batteries into fire, as that could cause an explosion.

Warning: Do not attempt to charge the included batteries, as that could cause leakage, fire disaster, or even explosion. Dispose the used batteries in accordance to the local laws and regulations.

Warning: Do not use the Oximeter in an MRI or CT environment.

Warning: Do not operate the Oximeter if it is wet. Avoid moving the oximeter from a cold to a hot and humid environment.

Warning: Install the batteries properly before powering on the Oximeter for normal use. Please remove the batteries if you are not planning to use the Oximeter for a long time.

Warning: Close the battery cover when the instrument is in use.

Symbol Conventions


Symbol	Description
	BF-type application part
	Caution: Please see this manual.
%SpO ₂	Symbol of oxygen saturation
bpmPR	Symbol of pulse rate
	Temperature limitation
	When end users abandon this product, they must send the product to the collection place for recycling.
	Manufacturer information

Overview

SpO₂ stands for peripheral capillary oxygen saturation. Oxygen saturation is defined as the ratio of oxyhemoglobin (HbO₂) to the total concentration of hemoglobin (i.e. Oxyhemoglobin + reduced hemoglobin) present in the blood. It is an important physiological parameter involved in respiration and circulation. The Pulse Oximeter feature herein is small, portable, non-invasive and easy to use. The user only needs to insert a finger into the chamber to measure his/her SpO₂ and pulse rate.

Working Principles, Expected Usage, and Applicable Scope

Oxygenated blood absorbs light at 660nm (red light), whereas deoxygenated blood absorbs light preferentially at 905nm (infra-red). A pulse oximeter works by passing a beam of red and infrared light through a pulsating capillary bed and then measure the amount of red and infrared light emerging from the tissues via a sensor. To improve accuracy, the NV-430J/PE uses a proprietary algorithm to collect data from pulsatile arterial blood and ignores local noise from the tissues. The relative absorption of light by oxyhemoglobin (HbO) and deoxyhemoglobin is then processed according to the Beer-Lambert's law and a qualitative measurement of the users' oxyhemoglobin status i.e. Oxygen saturation level (SpO₂) derived.



This device is intended for non-medical use in healthy people to monitor their pulse and blood oxygen levels for sports and/or aviation only. Do not use for continuous monitoring.

Limitations

Spot check pulse oximeter takes 1-2 reading every second and takes an average of 3-6 readings before displaying the result. Hence, you will have to wait at least 3-6 seconds before a result will be displayed on the monitor. If the pulse oximeter fails to detect blood flow for the first few seconds, then the results will be delayed accordingly.

The pulse oximeter does not work for people with naturally small fingers or low peripheral blood flow. This is because the pulse oximeter would not have enough blood to take a measurement. Please note that prolonged use of a pulse oximeter, hypotension, vasoconstriction, hypothermia and certain medications can lead to low or restricted blood flow. Nail polish can also impede the pulse oximeter from taking a reading. Finally, it is generally accepted that the saturation percentage is unreliable on the steep part (around 60 mm Hg) of the oxyhemoglobin dissociation curve.

Innovative features of the INV-430J/PE

The INV-430J/PE attempts to ameliorate some of the above problems by 1) Increasing the sensitivity of the sensor so that it can measure the user's SpO₂ and Pulse Rate (PR) at blood perfusion Index (PI) as low as 0.2%.

- Comes with anti-motion technology so that the pulse oximeters can still function even when the finger moves.
- Using a self-adjusting smart spring system to fit the finger snugly but not so tight as to impede blood flow that might lead to a low PI.
- Using soft, hypoallergenic medical grade silicon in the finger chamber so that the pulse oximeter will not cause discomfort.
- Using solid ABS plastic to block ambient light from reaching the sensor so that the signal to background ratio and the accuracy of the measurement is dramatically increased.
- Having a plethysmograph that allows the user to see and gauge the amount of blood detected by the pulse oximeter. This will allow the user to know if the finger is properly inserted into the finger chamber and if the pulse oximeter is taking the measurements properly (See below for details)


Operation Guide



Switch on the power supply of the pulse Oximeter

Stick one finger completely into the finger chamber of the Oximeter. The fingernail should be facing upward. Release the clip and press the power key to power on the Pulse Oximeter. The INV-430J/PE pulse oximeter



comes with anti-motion technology so that the pulse oximeter can function even if there is finger movement. However, for best results, try to keep your finger still during measurement. It is also not advisable to use this instrument during sports activities as movement may lead to inaccuracies. Once the reading stabilizes, read the measured values of oxygen saturation and pulse rate on the OLED screen. The Oximeter will automatically shut down 10 seconds after you remove your finger.


 **If you do not insert your finger completely into the chamber the measurement will be inaccurate.**

 **Replace the batteries when the batteries are low and the symbol  flickers on the screen.**

Schematic Diagram of Display

The following figure shows the information display on the OLED screen of the Oximeter in normal detection state:

- 1 Symbol of oxygen saturation
- 2 Measured value of oxygen saturation
- 3 Bar graph
- 4 Plethysmograph
- 5 Measured value of pulse rate
- 6 PI
- 7 Symbol of pulse rate
- 8 Battery power indication



Power-On Key/Function Key Operations

PRESS (<0.5 sec) the power-on/function key to turn on the pulse oximeter. Once it is turned on, **HOLD** (>0.5 sec) the power-on/function key to enter into the menu interface. **PRESSING** (<0.5 sec) on the power-on/function key will allow you to scroll through options while **HOLDING** (>0.5 sec) it will select an item. "*" indicates the option you are currently at (Please see figure below).

Alarm Sound Setting

To turn the alarm on/off, **PRESS** the power-on key/function key to move "*" to **Alm**. **HOLD** the power-on key/function key to switch the option from on to off or vice versa. When **Alm** is set to on and the measured values of the blood oxygen saturation (SpO₂) and/or pulse rate (PR) go beyond the stated upper or lower limit (See Alarm Range Setting to set parameters for Alarm), the alarm will go off. When **Alm** is set to off, the alarm is deactivated. The alarm will not go off regardless of your SpO₂ and/or PR values.

Brightness Setting

To change brightness, **PRESS** the power-on key/function key to move "*" to **brightness**. **HOLD** the power-on key/function key to set the brightness. There are 5 brightness settings. 5 is the brightest. The factory default setting is 4. Please note that a brighter setting will drain the

battery faster.


V2.01.E	
Alm setup	*
Alm	on
Beep	off
Demo	off
Restore	ok
Brightness	4
Exit	

Interface 1

V2.01.E	
Sounds Setup	*
SpO ₂ Alm Hi	100
SpO ₂ Alm Lo	94
PR Alm Hi	130
PR Alm Lo	50
+/-	+
Exit	

Interface 2

Plethysmograph



← Plethysmograph

The INV-430J/PE comes with two invaluable features, the plethysmograph and the perfusion index, that help the user determine if the readings are accurate and reliable. The plethysmograph indicates the amount of blood flow detected by the pulse oximeter and the perfusion index (PI) indicates the strength of your pulse. Each wave in the plethysmograph (see above) corresponds to a heartbeat and the wave amplitude corresponds to the amount of blood detected by the pulse oximeter flowing through the blood vessels. The pulse oximeter is optimized when the height of the wave amplitude is consistent throughout as shown in the figure above. That is when you should take the reading. The pulse oximeter can function with a PI reading as low as 0.2. However, if your PI is below 0.2%, this means that your blood perfusion is too low for a reliable read. Please warm your hands to increase blood flow and retake your measurement to make sure that your finger is positioned right under the LED lights. Most of the time, a higher PI will give you a more reliable reading.

Alarm Range Setting

To change the alarm range settings, **HOLD** the power-on/function key to enter Menu Interface 1. With the "*" next to **Alm Setup**, **HOLD** the power-on key/function key to enter in Menu Interface 2. "*" should be at **Sounds Setup**. Press the power-on key/function key to move "*" to the option you desire to change. **SpO₂ Alm Hi** and **SpO₂ Alm Lo** refer to the upper and lower limit of the SpO₂ levels that you wish the alarm to go off at respectively. **PR Alm Hi** and **PR Alm Lo** refer to the upper and lower limit of the Pulse Rate levels that you wish the alarm to go off at respectively.

While the "*" symbol stays on the +/- option, hold the function key to set the option to + or -. In + mode, select the corresponding upper or lower limit option and hold the function key to increase the upper or lower limit; in - mode, hold the function key to decrease the upper or lower limit.

Move "*" to the **Exit** option, and hold the function key to return to the monitoring interface.

About the Display

PRESS (<0.5 sec) the function key in measurement mode to change the display mode. The six different modes are shown in the figure to the right. Choose the mode that is the most convenient for you.



Cleaning

To clean the instrument, power off and remove the batteries first. Then clean the outer surface of the instrument (including the LED screen) using a piece of dry soft cloth dipped with 75% medical alcohol. **Do not immerse the unit in alcohol.**

Caution: Do not use any strong dissolving agent such as acetone.

Caution: Do not rub the body of the instrument using materials such as steel wire balls or polished metal objects.


Caution: Ensure that there is no washing liquid on the surface of the instrument.

Caution: Do not allow liquid to flow into the instrument during cleaning.

Caution: Do not immerse any part of the instrument into any liquid.

Disinfection

Before measurement with the instrument, wipe the rubber finger pad inside the chamber using a piece of dry soft cloth dipped with 75% medical alcohol. Clean the finger to be measured using the medical alcohol for disinfection purposes before and after use.

 Do not disinfect the instrument by using high-temperature/high-pressure disinfecting gas.

Maintenance

- Remove the batteries from the battery slot and properly store them if you do not plan to use the Oximeter for a long period of time.
- Replace the batteries if they are low on power.
- Clean the Oximeter and the fingertip before use to ensure accurate reading.
- Store the Oximeter between 14 and 122°F (-10 to +50°C) and at humidity levels no greater than 93%.
- Periodically check the Oximeter for damage.
- Avoid using the Oximeter in an environment with flammable gases or using it in an environment where the temperature or humidity is

excessively high or low.

- Check the accuracy of the oxygen saturation and pulse rate readings by using an appropriate calibration apparatus.

Technical Specifications

1. Dimensions: 62.0 mm (Width) × 37.0 mm (Depth) × 32.0 mm (Height)
Weight: 42.5 g (including the height of the two AAA dry batteries)
2. Peak wavelength range of the light emitted from the probe: red light 660 nm ± 3; infrared light 905 nm ± 5.
3. Maximum optical output power of the probe: 1.2 mW for infrared light (905 nm).
4. Normal working condition

Working Temperature	5°C to 40°C (41°F to 104°F)
Relative Humidity	15% to 80%, non-condensing
Atmospheric Pressure	70 kPa to 106 kPa
Rated Voltage	DC 3.0 V

5. Default values and conditions of alert

Parameter	Value
Oxygen saturation	Upper limit: 100 Lower limit: 94
Pulse rate	Upper limit: 130 Lower limit: 50
Alert condition	When the alert switch is on and the actual measured value goes beyond the preset alert parameter range, the Oximeter gives an alert sound.

6. Technical parameters

Parameter	Value
Display range	Oxygen saturation 35% to 100%
	Pulse rate 25 bpm to 250 bpm
Resolution	Oxygen saturation 1%
	Pulse rate 1 bpm
Measurement precision	Oxygen saturation ±2% (70% to 100%) No requirement (≤ 69%)
	Pulse rate ±2 bpm
Alert range	Oxygen saturation Upper limit: 50% to 100% Lower limit: 50% to 100%
	Pulse rate Upper limit: 25 bpm to 250 bpm Lower limit: 25 bpm to 250 bpm
Alert error	Oxygen saturation ± 1% of the preset value
	Pulse rate The greater of ±10% of the preset value and ±5 bpm
PI	Weak PI Min. 0.2%

Safety Type


Anti-electric-shock type: internal power supply device

Anti-electric-shock degree: BF-type application part

Running mode: continuous working

Waterproof grade: IP22

6.4 Storage and Transportation



Innovogroups

Phone: +1-858-888-9781

cs@innovogroups.com

<http://www.innovogroups.com>