HYAL*SYSTEM

AUTO CROSS-LINKED HYALURONIC ACID 2.0%

ACP

- EN | Package insert for the physician
IT | Foglietto illustrativo per il medico
ES | Prospecto para el médico
DE | Gebrauchsinformation für den Arzt
FR | Notice d'information destinée au médecin
UA | Листок-вкладиш в упаковці для лікар
EL | Ένθετο συσκευασίας για τον ιατρό
BG | Вложка на опаковката за лекаря
PL | Ulotka dla lekarza
EE | Pakendi infoleht: teave arstile
RO | Prospectul produsului pentru medic
LV | Lietošanas instrukcija ārstam
LT | Medikui skirtas informacinis lapelis

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Product Information

HYAL SYSTEM ACP is a sterile, highly viscous gel obtained by condensation of hyaluronic acid, a main component of human connective tissue and of epithelial and mesothelial tissues. Thanks to its prolonged residence in the application site, it slowly releases hyaluronic acid that, due to its natural hydrating properties, improves and maintains turgidity and skin elasticity for a longer time with respect to hyaluronic acid-based preparations.

Composition

Principal component: ACP (auto-cross-linked hyaluronic acid) 2.0%

Other components: sodium chloride, water for injection.

Indications

HYAL SYSTEM ACP is indicated for a medium lasting corrective action on minor soft tissue defects such as facial fine lines. It improves skin turgidity and elasticity, giving a rejuvenating effect to the skin. It is also indicated for the management of hypertrophic scars (pathological or not) resulting from burns, surgical wounds and skin trauma.

Instructions for use

1. Bring the product to room temperature.
 - a. **Treatment of soft tissue defects**
2. Remove the protective tip cap of the syringe.
3. Insert one of the needles provided in the box. It is recommended to institute personalized treatments to fulfill patient's needs. It is suggested to use the 4 mm needle to treat areas such as the frontal, peri-ocular, peribuccal and the mandibular supero-inferior line, while it is preferable to use the 13 mm needle to treat other areas such as the zygomatic antero-lateral region and the whole sottozygomatic area.
4. Inject intradermally the gel, superficially, at mid level or deeply, according to the aimed results. Different injection techniques could be used at physician's discretion and in different facial areas. It is generally recommended a linear technique or multiple injections in the middle or deep dermis. For wide areas of the face, a linear crossed technique is suggested. Fine lines can be treated with superficial injections, using linear and/or micropomphi injection techniques.
5. After injection, gently massage the treated area to improve product distribution in the tissue.
 - a. **Treatment of hypertrophic scars**
2. Remove the protective tip cap of the syringe.
3. Insert the 13 mm needles provided in the box.
4. Infiltrate the scar with the product in the intra lesion dermal planes.
5. Infiltrate the same amount of product in the subdermal level.
6. Repeat the procedure until complete scar infiltration, using a maximum of two ml of product for 25 cm² of scar area, up to obtain a transitory ischemia indicated by the transition to a white color of the lesion.

Contraindications

HYAL SYSTEM ACP is contraindicated in patients with a known hypersensitivity to hyaluronan or its derivatives.

As per any other intra-dermal injection, avoid application in anatomical sites with infection or inflammatory processes.

Undesirable effects

The treated area may show reactions commonly related to the injection, such as tumefaction, mild pain and erythema, that spontaneously resolve within few hours after application. In rare instances, small ecchymosis may be induced; they spontaneously resolve within few days.

Interactions

HYAL SYSTEM ACP can be used with topical anaesthetics containing lidocaine.

HYAL SYSTEM ACP must not be used concomitantly with disinfectants containing quaternary ammonium salts because hyaluronic acid can precipitate in their presence.

Storage conditions

Keep refrigerated (+2 / +8 °C)

Do not freeze.

Warnings and Precautions

- Do not overcorrect the skin defects.
- Administer in aseptic conditions. As any other percutaneous treatment, application of HYAL SYSTEM ACP may imply a risk of infection.
- The device is supplied in pre-filled single use syringes, packaged in a protective blister.
- The syringe is single use. Discard any unused portion.
- The needles are single use. Discard after use and do not sterilize again.
- The needles may be damaged if used for multiple injections on the same patient, in that case it is suggested to replace them with similar needles.
- If this product is reprocessed and/or reused, Fidia farmaceutici cannot guarantee functionality, material structure, or cleanliness or sterility of the product. Reuse could lead to illness, infection and/or serious injury to the patient or user.
- Dispose according to applicable national regulations.
- Keep out of reach of children.
- Do not use if the package and, in particular, the syringe and needle blisters are damaged.
- Do not utilize the product after the expiration date reported on the packaging
- No clinical data are available on possible interactions with other aesthetic correction treatments (implant of biomaterials, lifting, botulinum toxin injections, laser, chemical peeling) done in the previous three months.
- No clinical data are available on possible interactions with permanent fillers.
- No clinical data are available on the use of HYAL SYSTEM ACP in patients affected by autoimmune conditions (LES, scleroderma, etc.) neoplasia or with ongoing immunosuppressor therapy; depending on the nature of the disease, the physician shall therefore choose the administration and the associate treatment on an individual case-by-case basis.
- No clinical data are available on the use of HYAL SYSTEM ACP in women who are pregnant or in supposed pregnancy or during lactation.

How it is supplied

Cardboard box containing one single-use 1 ml prefilled syringe with backstop, one sterile 30 G, 13 mm needle and one sterile 30 G, 4 mm needle (TSK Laboratories, Japan - EC Rep.: EMERGO Europe - CE 0123).

The box is also provided with a sticker for the physician reporting name of the product, lot number and expiry date.

The contents of the syringe are sterilized using steam.

The needles are sterilized using gamma radiation.

This leaflet was last revised in April 2019

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Fidia farmaceutici S.p.A.