

# Dikma QuEChERS

## 1. Dikma QuEChERS

# QuEChERS

**Quick, Easy, Cheap, Effective, Rugged and Safe  
Method for Determining Pesticide Residues**



# Fruits and Vegetables



# Milk and Honey



# Animal Tissues



# Tobacco



# What Makes Traditional Multiresidue Methods Inefficient?

- **Complicated**
- **Laborious**
- **Time consuming**
- **Require high amount of solvents**
- **Expensive**



# Ideal Multiresidue Method

- Fast and easy to perform
- A minimum amount of chemicals
- Good selectivity to avoid complicated cleanup procedures
- Cover sufficiently broad spectrum of analytes



# Benefits of QuEChERS Method

- **High recoveries**
- **Accurate results**
- **High sample throughput**
- **Low non-chlorinated solvent usage**
- **Reduce reagent costs and staff exposure to hazardous solvents**
- **Reduce glassware usage & labor costs**
- **Broad utility & ease of use**

# QuEChERS

- **Primarily for analysis of pesticides from hydrated matrices (80-95% water content)**
- **Multi-residue pesticide analysis of food and environmental samples can be problematic due to the wide range of chemical properties encountered with pesticide residues**
- **The complex sample matrix may contain abundant quantities of chlorophyll, lipids, sterols and other components that can interfere with good sample analysis**
- **Use of the QuEChERS method reduces these problems**

# QuEChERS Method

- **Consists of a liquid-liquid micro extraction**
- **Partitioning and extraction of polar analytes is aided by MgSO<sub>4</sub>**
- **The preferred solvent is acetonitrile**
- **Acetonitrile provides extraction of the broadest range of organic compounds without co-extraction of large amount of lipophilic material and is highly compatible with GC-MS and LC-MS applications showing the fewest interferences**
- **Followed by dispersive SPE sample clean-up to remove unwanted matrix materials**

# Modifications to the Original Method

- **Some modifications to the original QuEChERS method have been introduced to ensure efficient extraction of pH dependent compounds and to minimize degradation of base and acid labile pesticides**
- **Buffering with citrate salts has been introduced in the micro extraction to adjust the pH to a compromise value of 5 to 5.5, where most acid and base labile pesticides are sufficiently stabilized. To improve stability of base-labile compounds in the sample extracts, a small amount of formic acid is added to the final extract after cleanup**
- **Acidic pesticides are directly analyzed from the raw extract before PSA cleanup since they would be adsorbed and not released by the sorbent**
- **In another modification introduced by Schenck, a graphitized carbon black (GCB) PSA cartridge is used to remove plant pigments without the loss of planar compounds**

# Four QuEChERS Methods

- The original QuEChERS method introduced in 2003, uses NaCl to enhance extraction (reduce polar interferences)
- Dispersive AOAC 2007.01 uses 1% acetic acid in acetonitrile, sodium acetate as a buffer replacing sodium chloride
- The dual phase column – this method variation introduces the use of PSA and GCB cartridges to remove high levels of chlorophyll and plant sterols in the final extract without the loss of planar pesticides (polar aromatics) using a 3:1 acetone:toluene solvent mixture
- The European version (EN 15662) is similar to the AOAC method, except the extraction uses NaCl, sodium citrate dihydrate and disodium citrate sesquihydrate instead of NaOAc. *Acetic acid is not used*

# Original QuEChERS

10 g sample + 10 mL MeCN

add internal standard

add 4 g MgSO<sub>4</sub> + 1 g NaCl  
shake vigorously for 1 min

centrifuge for 5 min

1 mL of the upper layer  
+ 25 mg PSA + 150 mg MgSO<sub>4</sub>  
mix for 30 s  
centrifuge for 1 min

extraction

dispersive SPE  
clean-up

# Buffered QuEChERS (AOAC 2007.01)

15 g sample + 15 mL 1% HOAc in MeCN

add internal standard

add 6 g MgSO<sub>4</sub> + 1.5 g NaOAc  
shake vigorously for 1 min

centrifuge for 5 min

1 mL of the upper layer  
+ 50 mg PSA + 150 mg MgSO<sub>4</sub>  
(+ 50 mg C<sub>18</sub>)  
mix for 30 s  
centrifuge for 1 min

extraction

dispersive SPE  
clean-up



# Buffered QuEChERS (EN 15662)

extraction

10 g sample + 10 mL MeCN

add internal standard

add 4 g  $\text{MgSO}_4$  + 1 g NaCl  
+ 1 g  $\text{Na}_3\text{Citr}\cdot 2\text{H}_2\text{O}$  + 0.5 g  $\text{Na}_2\text{HCitr}\cdot 1.5\text{H}_2\text{O}$   
(+ 0.6 mL 5N NaOH for lemons, limes etc.)  
shake vigorously for 1 min

centrifuge for 5 min

dispersive SPE  
clean-up

1 mL of the upper layer  
+ 25 mg PSA + 150 mg  $\text{MgSO}_4$   
(+ 2.5 or 7.5 mg GCB for matrices with a high  
content of carotenoids and chlorophyll)  
mix for 30 s, centrifuge for 1 min

# The Dual Phase Variation

- **Matrix plant pigments often interfere with analysis**
- **To reduce these interferences, graphitized carbon can be added to the dispersive solid-phase cleanup tubes**
- **Carbon however may result in a loss of planar (polar aromatic) pesticides**
- **Cleanup of plant pigments with minimum loss of planar pesticides can be accomplished by using a dual-phase cartridge containing PSA and GCB**

## Dual Phase Cartridge Clean-up Procedure

1. Pre-rinse cartridge with 5 mL of toluene
2. Add an aliquot of the supernatant to cartridge
3. Start collection
4. Elute with 6-12 mL 3:1 acetone:toluene
5. Concentrate for GC-MS analysis or concentrate to dryness and reconstitute in mobile phase for LC analysis



## **Other Variations** – additional cleanup options

- **C18 – removes lipids**
- **Graphitized carbon black (GCB) – removes pigments (and planar compounds)**
- **ChloroFiltr<sup>®</sup> - removes chlorophyll**



# QuEChERS

## BREAKING IT DOWN



# Acetonitrile

- **Extracts the greatest number of analytes with the least number of interferences**
- **Can be used on GC and LC**

# Acetic Acid

- **1% acetic acid in acetonitrile, when combined with sodium acetate, prevents base sensitive analytes from breaking down during extraction**
- **Compromises PSA**
- **Work best with LC-MS-MS analyses**
- **Tailing issues with GC**

## **Sodium Citrate Dibasic Sesquihydrate and Sodium Citrate Tribasic Dihydrate**

- **Buffers extraction to prevent break down of pH sensitive analytes**
- **Does not require acetic acid**
- **Does not compromise PSA (primary secondary amine)**



# Anhydrous Magnesium Sulfate

- In extraction, acid in partitioning and improves recoveries of polar analytes
- In clean-up, works as a desiccant



# Sodium Chloride

- **Help to reduce polar co-extractables**
- **Can reduce recovery of polar analytes**

# PSA

**Used in clean-up to reduce levels of organic acids, sterols, some sugars and lipids**



# Endcapped C18

**Used in clean-up to aid PSA in removal of hydrophobic interferants such as lipids**



# Aminopropyl

- **Similar properties to PSA, but lower exchange capacity**
- **Less likely to damage base sensitive analytes**



# Graphitized Carbon Black (GCB)

- **Used in clean-up to remove pigments and polyphenols**
- **GCB binds planar analytes such as acephate, bromophos, carbendazim, chlorthiophos, cyprodinil, ethoprop, fonofos, leptophos, methamidophos, pyrimethanil, and thiabendazole and will remove these planar analytes**
- **More GCB = Lower recovery of planar analytes**

# Advantages of Buffered Methods

- **Buffered extract to protect base sensitive analytes such as folpet, dichlorofluanid, chlorothalonil, pymetrozine, dicofol, captan, tolyfluanid**
- **Extract protected**
  - LC – formic acid (PSA is basic)
  - GC – toluene and magnesium sulfate (prevent thermal breakdown)

# Disadvantage of AOAC Method

- **Acetic acid – PSA removes organic acids**
- **Less clean-up**
- **Tailing issues on GC**





# Cartridge or Dispersive Clean-up?

- **Cartridge (not QuEChERS)**
  - Cleaner extracts
  - Takes longer
  - Uses more solvent
  - Requires a manifold and accessories
- **Dispersive (QuEChERS)**
  - Quick
  - Easy
  - Not as clean as a cartridge
  - Requires a centrifuge

# PUTTING IT ALL TOGETHER



1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48 49 50 51 52 53 54 55 56 57 58 59 60 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75 76 77 78 79 80 81 82 83 84 85 86 87 88 89 90 91 92 93 94 95 96 97 98 99 100


# The QuEChERS Method for Pesticide Residues



1) Shake sample with solvent and salts



2) Centrifuge for 1 min



3) Mix a portion with a sorbent



5) Analyze Pesticides



4) Centrifuge for 1 min



**Add hydrated vegetation  
to salts & solvent**




**Shake**



Courtesy of Dr. Frank Schenck, USFDA

# Centrifuge






← MeCN Extract

← Plant Matrix

← Aqueous


# Recoveries in the QuEChERS method

229 pesticides analyzed by GC-MS and LC-MS/MS



Courtesy of Dr. Steven Lehotay, USDA-ARS

## Recoveries of 15 Pesticides in Different Matrices




**No differences were found vs. matrix for individual pesticides or concentration**




# Repeatability in the QuEChERS method


229 pesticides analyzed by GC-MS and LC-MS/MS



# Why PSA Cleanup – GC/FPD asparagus extract



# Matrix Components Removed by SPE Sorbents



# Pesticide Analytes

**GC amenable pesticides are capitalized. LC/MS/MS pesticides are not capitalized. Analytes that can be analyzed under either instrument are underlined**

<u>acephate</u> *	acetamiprid	Acrinathrin	aldicarb	aldicarb sulfone
aldicarb sulfoxide	Aldrin	azaconazole	azamethiphos	<u>azinphos-methyl</u>
<u>azoxystrobin</u>	Bifenthrin	<u>bitertanol</u>	Bromopropylate	<u>bromuconazole</u>
Bupirimate	<u>buprofezin</u>	butocarboxim	butocarboxim sulfone	butocarboxim sulfoxide
Cadusafos	<u>carbaryl</u>	carbendazim	<u>carbofuran</u>	3-hydroxy-carbofuran
chlorbromuron	( $\alpha$ -, $\gamma$ -)Chlordane	( $\alpha$ -, $\beta$ - Chlorfenvinphos	Chlorpropham	Chlorpyrifos
Chlorpyrifos-methyl	Chlorthaldimethyl	Chlorothalonil*	Chlozolate	clofentezine
Coumaphos	cycloxydim*	( $\lambda$ -)Cyhalothrin	cymoxanil	Cypermethrin
<u>cyproconazole</u>	<u>cyprodinil</u>	(2,4'-4,4'-)DDE	(2,4'-4,4'-)DDT	Deltamethrin
demeton	demeton-O-sulfoxide	demeton-S-methyl	demeton-S-methyl sulfone	desmedipham
Diazinon	<u>dichlofluanid</u> *	Dichlorobenzophenone	<u>dichlorvos</u>	diclobutrazole
Dicloran	dicrotophos	Dieldrin	<u>Diethofencarb</u>	<u>difenoconazole</u>
Diflufenican	<u>dimethoate</u>	dimethomorph	<u>diniconazole</u>	Diphenyl
Diphenylamine	<u>disulfoton</u>	<u>disulfoton sulfone</u>	diuron	<u>dmsa</u>
<u>dmst</u>	dodemorph	$\alpha$ - Endosulfan	$\beta$ -Endosulfan	Endosulfan sulfate

# Pesticide Analytes

**GC amenable pesticides are capitalized. LC/MS/MS pesticides are not capitalized. Analytes that can be analyzed under either instrument are underlined**

ethiofencarb sulfoxide	Ethion	ethirimol	<u>Ethoprophos</u>	<u>etofenprox</u>
Etridiazole	Famoxadone	<u>fenamiphos</u>	<u>fenamiphos sulfone</u>	<u>Fenarimol</u>
Fenazaquin	fenbuconazole	<u>fenhexamid</u> *	Fenithrothion	<u>fenoxycarb</u>
Fenpiclonil	Fenpropathrin	Fenpropidine	<u>fenpropimorph</u>	<u>fenpyroximate</u>
<u>Fenthion</u>	<u>fenthion sulfoxide</u>	Fenvalerate	florasulam*	Flucythrinate I & II
Fludioxonil	flufenacet	Flufenconazole	<u>flusilazole</u>	Flutolanil
Fluvalinate	Fonophos	fosthiazate	Furalaxyl	furathiocarb
<u>furmecyclox</u>	Heptachlor	Heptachlor epoxide	Heptenophos	Hexachlorobenzene
<u>hexaconazole</u>	hexythiazox	imazalil	imidacloprid	Iprodione
iprovalicarb	isoprothiolane	isoxathion	<u>kresoxim-methyl</u>	Lindane
linuron	<u>Malathion</u>	<u>malathion oxon</u>	Mecarbam	<u>mephosfolan</u>
Mepronil	Metalaxyl	metconazole	<u>methamidophos</u> *	Methidathion
<u>methiocarb</u>	methiocarb sulfone*	methiocarb sulfoxide	methomyl	methomyl-oxime
metobromuron	metoxuron	Mepanipyrim	Mevinphos	<u>monocrotophos</u>
monolinuron	<u>myclobutanil</u>	nuarimol	Ofurace	<u>omethoate</u>
<u>oxadixyl</u>	oxamyl	oxamyl-oxime	oxydemeton-methyl	paclobutrazole

# Pesticide Analytes

**GC amenable pesticides are capitalized. LC/MS/MS pesticides are not capitalized. Analytes that can be analyzed under either instrument are underlined**

Parathion	Parathion-methyl	<u>penconazole</u>	<u>pencycuron</u>	<i>cis</i> - Permethrin
<i>trans</i> -Permethrin	phenmedipham	<i>o</i> -Phenylphenol	<u>Phorate</u>	<u>phorate sulfone</u>
Phosalone	Phosmet	Phosmet-oxon	phosphamidon	Phthalimide
<u>picoxystrobin</u>	Piperonyl butoxide	<u>pirimicarb</u>	<u>pirimicarb-desmethyl</u>	Pirimiphos-methyl
prochloraz	Procymidone	<u>profenofos</u>	Prometryn	Propargite
Propham	<u>propiconazole</u>	<u>propoxur</u>	Propyzamide	Prothiofos
pymetrozine*	Pyrazophos	<u>pyridaben</u>	<u>pyridaphenthion</u>	<u>pyrifenox</u>
<u>pyrimethanil</u>	Pyriproxyfen	Quinalphos	Quinoxifen	Quintozene
sethoxydim*	spinosad	<u>spiroxamine</u>	<u>tebuconazole</u>	tebufenozide
<u>Tebufenpyrad</u>	<u>tetraconazole</u>	Tetradifon	Tetrahydrophthalimide	Terbufos
Terbufos sulfone	thiabendazole	thiacloprid	thiamethoxam	thiodicarb
thiofanox	thiofanox sulfone	thiofanox sulfoxide	thiometon	thiometon sulfone
thiometon sulfoxide	thiophanate-methyl	Tolclofos-methyl	<u>tolyfluanid</u> *	<u>triadimefon</u>
<u>triadimenol</u>	Triazophos	trichlorfon	tricyclazole	tridemorph
<u>trifloxystrobin</u>	<u>trifluminazole</u>	Trifluralin	<u>Triphenylphosphate</u>	vamidothion
vamidothion sulfone	vamidothion sulfoxide	Vinclozolin		

# Key Steps

- 1. Hydrate samples to 80% or higher**
- 2. Add solvent prior to adding salts**
- 3. Use internal standard**
- 4. Use matrix matched calibration standards**
- 5. Use buffer when necessary**
- 6. Use analyte protectants (formic acid or toluene)**
- 7. Beware of GCB**



Thank You!

