

Ultra-low input linked-read library method generates highly accurate and economical long-range information from short reads for broad applications

Tom Chen¹, Peter Chang¹, Long Pham¹, Donna Wu¹, Gloria Mo¹, Yu Xia¹, Chris Boles², Yong Wang¹

1. Universal Sequencing Technology, Carlsbad, CA 92011

2. Sage Science, Inc., Beverly, MA 01915

email: tchen@universalsequencing.com



Abstract

Transposase Enzyme Linked Long-read Sequencing (TELL-Seq™) is a simple and scalable NGS library technology to generate barcode linked-reads for genome scale haplotype phasing and de novo sequencing. In a PCR tube, under a standard NGS laboratory setting and without the need for any expensive protocol-specific instrument, TELL-Seq Whole Genome Sequencing (WGS) Library Prep Kit will generate an Illumina sequencing library in 3 hours from as low as 0.1ng DNA input. Multiple samples can be easily processed in parallel as needed. The library protocol can be adjusted and used for various sized genomes ranging from bacteria to human. The TELL-Seq library will replace both fragment library and mate-pair library altogether and become a new standard library method for WGS.


Genome Assembly

Sample	<i>Escherichia coli</i> MG1655	<i>Escherichia coli</i> DH10B	<i>Arabidopsis thaliana</i> 3	<i>Homo sapiens</i> NA12878
gDNA input (ng)	0.5	0.1	3	5
Average gDNA length (Kb)	39	50	18	50
Assembler	TuringAssembler		SuperNova	
Global/local k-mer sizes	105/69	105/31	N/A	N/A
Genome fraction (%)	99.9	99.3	97.4	93.7
Duplication ratio	1.009	1.012	1.061	1.075
Largest alignment	4,637,716	3,459,910	5,613,314	23,573,913
Total aligned length	4,678,679	4,705,787	119,712,720	2,771,334,486
NA50	4,637,716	3,459,910	1,272,225	4,302,918
# misassemblies	0	1	246	1,987
# mismatches per 100 kbp	5.37	12.31	14.17	114.07
# indels per 100 kbp	0.39	0.24	5.96	25.61
# N's per 100 kbp	0.0	0.0	4,499.1	6419.8
# contigs (>= 5000 bp)	3	2	382	13,584
# contigs (>= 10000 bp)	3	2	200	7,508
Largest contig	4,640,256	4,657,354	10,292,118	109,183,970
Total length (>= 1000 bp)	4,725,322	4,822,052	126,210,917	2,979,445,840
N50	4,640,256	4,657,354	4,677,832	31,462,027
GC (%)	50.75	50.67	36.29	40.91


Whole Exome Phasing




How does TELL-Seq Work?



Whole Genome Phasing




Targeted Full-Length Gene Diplotyping



TELL-Seq WGS Library Workflow is Simple



Structure Variation Detection



Technical Indices	TELL-Seq™
Instrument cost	Zero
Library prep cost	Low
Input DNA	0.1 ng to 5 ng
Processing time	3 hours
Barcode capacity	>2 billion
Genome size	≥ 0.1 Mb

References

Chen, Z et al. 2020. Ultralow-input single-tube linked-read library method enables short-read second-generation sequencing systems to routinely generate highly accurate and economical long-range sequencing information. *Genome Research* 30: 898-909.