SILVER CALCIUM ALGINATE DRESSINGS & ROPES

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

Silver Calcium Alginate Dressings and Ropes are composed of calcium alginate, sodium alginate, and sodium silver phosphate. Vitále® STRONG™ Calcium Alginate Dressings are made of absorbent fibers produced by a special process. The dressing absorbs exudate and forms a gel-like covering over the wound. This maintains a moist wound healing environment and keeps the dressing from adhering to the wound. Vitále® STRONG™ Silver Calcium Alginate Dressings have been tested in vitro and found to be effective against bacteria such as: Staphylococcus aureus, Escherichia coli, Enterococcus faecalis, Pseudomonas aeruginosa, Klebsiella pneumoniae, and Staphylococcus epidermidis.

Vitále® STRONG™ Calcium Alginate Dressings with Antibacterial Silver are sterilized by exposure to Gamma Radiation, and must not be re-sterilized.

INDICATIONS

May be used as a primary wound dressing for moderate to heavily exuding partial to full thickness wounds, including: postoperative wounds, trauma wounds, leg ulcers, pressure ulcers, diabetic ulcers, and graft and donor sites.

Preservative agent	Sodium Silver Phosphate		
Active against pathogens	Gram positive and Gram negative bacteria		
Spectrum of activity	GRAM NEGATIVE BACTERIA		
	Bacterial	Escherichia coli ATCC 8739	Klebsiella pneumoniae ATCC 4352
			Pseudomonas aeruginosa ATCC 9027
	24 hour reduction %	>99.9%	>99.9%
	GRAM POSITIVE BACTERIA		
Spectrum of activity	Bacterial	Staphylococcus aureus ATCC 6538	Staphylococcus epidermidis ATCC 12228
			Enterococcus faecalis ATCC 19615
	24 hour reduction %	>99.9%	>99.9%
Concentration on the device	Each dressing contains 0.3% (w/w) silver ion		

CAUTION

- This product is single use. Do not use if pouch is damaged or open.
- Vitále® STRONG™ Silver Calcium Alginate Dressings and Ropes should be changed daily to avoid colonization of the dressing. The dressing should not be used on dry wounds.

- Frequent or prolonged use of this product may result in permanent discoloration of skin.
- Silver Calcium Alginate Dressing with STRONG Fiber Blend should not be used on patients with a known sensitivity to silver.
- Not intended for use on third degree burns.

IN VITRO STUDY RESULTS

The device is non-irritating, non-sensitizing, non-cytotoxic when tested in vitro. The device has shown antibacterial preservative effect up to 7 days against both Gram positive and Gram negative bacteria.

SIDE EFFECTS & INTERACTIONS

Silver Calcium Alginate Dressing with STRONG Fiber Blend has been tested for in-vitro biocompatibility, including cytotoxicity, sensitization, skin irritation and systemic toxicity. All tests were performed in accordance with the ISO 10993 Standard Series for Biological Evaluation of Medical Devices. The results of these studies indicate that Vitále® STRONG™ Silver Calcium Alginate Dressings are safe for their intended use.

STORAGE

Vitále® STRONG™ Silver Calcium Alginate Dressings and Ropes must be stored in a cool & dry area. The products are sterile unless the package is opened or damaged.

APPLICATION

- Clean the wound with sterile physiological saline using existing clinical protocol. Dry the surrounding skin gently.
- Choose the right size dressing to cover the wound area.
- Cover and secure the Vitále® STRONG™ Silver Calcium Alginate Dressing with a secondary dressing.
- Dressing change frequency will depend on patient condition and the level of exudates. While exudates are moderate to heavy, dressing changes should be made in 24 hours.
- Gently remove the secondary dressing. If the wound appears dry, saturate the dressing with sterile saline prior to removal.

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