

ANKLE BRACE

QUICK START GUIDE

OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Vive Ankle Brace as easy as possible.



To see all FAQ's in one place visit vhealth.link/rxO

HOW TO USE

Step 1:

Insert your foot into the brace so the Vive logo is over the front of your shin.



Step 2:

Fasten the top straps and adjust to your comfort level.



Step 3:

Attach provided ankle extender straps if needed. Extenders can be cut to desired length for a custom fit.



Step 4:

Wrap the bottom straps around your ankle and fasten them securely.



vive®