

vive®

BOXER SPLINT

Owner's Manual - SUP2053BLK



vivehealth.com

OVERVIEW

The Vive Boxer Splint provides support for breaks, fractures, or sprains of the 4th and 5th metacarpals (ring and pinky fingers). The splint support contains two aluminum splint bars to encase and immobilize the injured fingers and help facilitate the healing process.



We are constantly answering questions and recording helpful videos to make using your Vive Boxer Splint as easy as possible. Check out the included link and QR code to help you through the process.



To see all FAQ's in one place
visit vhealth.link/00i

APPLICATION

1. Insert your hand through the wrist and palm straps.



2. Position the bottom splint bar under the pinky and ring fingers that you wish to immobilize, and wrap the top splint bar over the same fingers.



3. Tighten and secure the palm strap.



4. Tighten and secure the wrist strap.



5. Tighten and secure the finger strap.



6. Adjust the straps as needed to reach the desired fit.

BEND THE SPLINT FOR CUSTOM FIT (ASSISTANCE REQUIRED)

The splint can be bent if a different angle is desired. To bend:

1. Grip the splint at each end and bend slightly.



2. Refit the splint using the steps above.
3. Repeat until the desired splint angle and fit is achieved.

CARE

The splint support comes with two aluminum splint bars installed in the sleeves.



Remove the splint bars prior to washing the fabric splint support. Hand wash with cold water and mild soap. Rinse and then hang dry.



GOT MORE QUESTIONS?

Check out our list of Frequently Asked Questions at vhealth.link/00i for helpful answers.



And if that doesn't answer your question, our customer service team would love to help! Feel free to connect with them by phone, e-mail, or chat on our website.



service@vivehealth.com



1-800-487-3808



vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109