

CRISPR/Cas9 nuclease



Catalog No.: 9000

Expressed in: *E. coli*

Contents:

CRISPR/Cas9 is provided at a concentration of 10 pmol/μL (10 μM)

Background:

CRISPR/Cas9 is a *Streptococcus pyogenes* Type II Cas9 endonuclease that can be programmed by a small RNA to site-specifically digest DNA. The specificity of cleavage is determined by the sequence of a single guide RNA (sgRNA) that contains bases that are complimentary to the target site. The simplicity of the CRISPR/Cas9 endonuclease system enables precise *in vitro* DNA digestion.

Application Notes:

CRISPR/Cas9 can be programmed to cleave DNA precisely by the sequence of the sgRNA that is loaded onto the protein to cleave DNA virtually anywhere. This enables site-specific *in vitro* engineering of genomic, plasmid, or PCR amplified DNA.

*These products are intended for research use only, not for diagnostic use. The safety and efficacy of these products in diagnostic or other clinical uses has not been established.

Protein Details:

CRISPR/Cas9 is expressed in *E. coli* fused to a solubility-enhancing protein. The molecular weight of this fusion protein is 171 kDaltons.

Shipping and Storage:

CRISPR/Cas9 is stored at -20 °C in 50% glycerol, 50 mM Tris-HCl, 50 mM KCl, 1 mM DTT, 1 mM EDTA, pH 7.5. *Can be supplied in a glycerol-free buffer as a custom order.*

CRISPR/Cas9 is shipped on dry or blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided.

Quality Control:

- CRISPR/Cas9 activity: >90% cleavage of 1 nM, 4 kb PCR product in 60 minutes at 37 °C using 40 nM targeted gRNA and 1 pmol CRISPR/Cas9 in 50 μL (20 nM).
- Purity: >95% as determined by SDS-PAGE analysis
- <0.05 ng contaminating host DNA per *

***In vitro* digestion of DNA with CRISPR/Cas9 nuclease**

Product Overview:

CRISPR/Cas9 is a recombinant form of *Streptococcus pyogenes* Type II Cas9 nuclease that site specifically digests DNA using the complementarity sequence of a small RNA bound to the protein to guide the nuclease to its target.

A protocol is provided for *in vitro* digestion of double stranded DNA with CRISPR/Cas9 and a single guide RNA.

Useful Information:

The DNA, guide RNA and nuclease-free water for this protocol are not supplied with the CRISPR/Cas9.

Information for the design of a Cas9 nuclease guide RNA containing a target sequence can be viewed at the Addgene website.

10X Reaction Buffer is used in the protocol below. 10X Reaction Buffer is 500 mM NaCl, 100 mM Tris-HCl, 100 mM MgCl₂, 10 mM DTT, pH 7.9. This buffer can be substituted with 10X TA Buffer. 10X TA Buffer is 330 mM Tris acetate, 660 mM KCl, and 100 mM magnesium acetate, pH 7.5.

General Protocol:

1. Assemble the reaction components at room temperature in a microfuge tube.

2. Mix the components at 37 °C for 1 hour.
3. Stop the heating at 70 min.

Component	Volume
Nuclease-free water	20 µL
10X Reaction Buffer	3 µL
300 nM guide RNA	3 µL (30 nM final)
1 µM CRISPR/Cas9	1 µL (~30 nM final)
Pre-incubate the above mixture at RT for 10 min before adding DNA	
30 nM DNA	3 µL (30 nM final)
Total reaction volume	30 µL

reaction and incubate hour.
reaction by °C for 15