Bedienungsanleitung / Manual

AXION STIM-PRO X9+

Digitaler 4-Kanal

Nerven- und Muskelstimulator

INHALTSVERZEICHNIS / CONTENT

2 – 33	Deutsche Bedienungsanleitung	DE
34 – 65	Instruction manual in English	EN
66 – 97	Manuel d'utilisation Français	FR
98 – 129	Manual de Instrucciones Español	ES
130 – 161	Istruzioni per l'uso Italiano	IT
162 – 193	Instrukcja Obsługi Polski	PL

EINLEITUNG

DE

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INTRODUCTION

EN

The primary functions of the STIM-PRO X9+ are pain relief and pain treatment with the TENS mode and electrical muscle stimulation with the EMS mode. The device is used by the patient himself directly.

Explanation of TENS

There is nothing “magic” about **T**ranscutaneous **E**lectrical **N**erve **S**timulation (TENS). TENS is intended to be used to relieve pain. The TENS unit sends comfortable impulses through the skin that stimulate the nerve (or nerves) in the treatment area. In many cases, this stimulation will greatly reduce or eliminate the pain sensation the patient feels. Pain relief varies by individual patient, mode selected for therapy, and the type of pain. In many patients, the reduction or elimination of pain lasts longer than the actual period of stimulation (sometimes as much as three to four times longer). In others, pain is only modified while stimulation actually occurs. You may discuss this with your physician or therapist.

How TENS works

Transcutaneous Electrical Nerve Stimulation is a non-invasive, drugfree method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however, in most patients it is effective in reducing or eliminating the pain.

Explanation of EMS

Electrical **M**uscle **S**timulation is an accepted and proven way of treating muscular injuries. EMS has also been

used successfully for many years in the field of high-performance sport for muscle regeneration and muscle development.

How EMS works

The brain sends impulses to a muscle which causes the muscle to work. EMS generates artificial stimulation of the muscle by the use of electrical pulses that stimulate the muscles located beneath the skin. The muscle is unable to distinguish between an impulse from the brain or one that comes from outside. The muscle reacts with a contraction and is trained passively. Then when the pulse ceases, the muscle relaxes and the cycle starts over again (stimulation, contraction and relaxation).

Indications

People use TENS to relieve pain for several different types of illnesses and conditions. They use it most often to treat muscle, joint, or bone problems that occur with illnesses such as osteoarthritis or fibromyalgia, or for conditions such as low back pain, neck pain, tendinitis, or bursitis. People have also used TENS to treat sudden (acute) pain, such as labor pain, and long-lasting (chronic) pain, such as cancer pain.

EMS is widely used in hospitals and sports clinics for the treatment of muscular injuries and for re-education of paralyzed muscles, to prevent atrophy in affected muscles and improve muscle tone and blood circulation. It works for:

Relaxation of muscle spasms, prevention or retardation of disuse atrophy, increasing local blood circulation, muscle re-education, immediate post-surgical stimulation of calf muscles to prevent venous thrombosis, maintaining or increasing range of motion.

IMPORTANT SAFETY INFORMATION

Read instruction manual before operation. Be sure to comply with all “CAUTIONS” and “WARNINGS” in the manual. Failure to follow instructions can cause harm to user or device.

EN

SIDE EFFECTS

TENS

Side effects of TENS machines are generally mild, even with extended use. If any effects are especially bothersome, a physician can evaluate the continued use of the device.

Skin irritation and minor burns are potential side effects. If this occurs, discontinue use and consult your physician.

EMS

Side effects from EMS are generally very mild. Patients may experience tenderness for a day or so afterward where muscles are already sore from inflammation. Shortterm redness and skin irritation may also occur. Those who use implanted medical devices and pregnant women should avoid this therapy.

Skin irritation and minor burns are potential side effects. If this occurs, discontinue use and consult your physician.

CAUTIONS & CONTRAINDICATIONS

TENS

- ▶ Only use original accessories from the manufacturer for your device.
- ▶ Do not use this device for undiagnosed pain symptoms until consulting a physician.
- ▶ Patients with an implanted electronic device, such as a cardiac pacemaker, implanted defibrillator, or any other metallic or electronic device should not undergo TENS treatment without first consulting a doctor.
- ▶ Patients with heart disease, epilepsy, cancer, fever or any other health condition should not undergo TENS treatment without first consulting a physician.
- ▶ Stimulation delivered by this device may be sufficient to cause electrocution. Electrical current of this magnitude must not flow through the thorax or across the chest because it may cause a cardiac arrhythmia. **Do not place electrodes above the heart!**
- ▶ Do not place electrodes on the front of the throat as spasm of the Laryngeal and Pharyngeal muscle may occur. Stimulation over the carotid sinus (neck region) may close the airways, make breathing difficult, and may have adverse effects on the heart rhythm or blood pressure.
- ▶ Do not place electrodes on your head or at any sites that may cause the electrical current to flow transcerebrally (through the head).
- ▶ This device should not be used while driving, operating machinery, close to water, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.

- ▶ Always turn off the device before applying or removing TENS accessories to the skin.
- ▶ Isolated cases of skin irritation may occur at the site of electrode placement following long term application. If this occurs, discontinue use and consult your physician, if the irritation does not disappear.
- ▶ If TENS therapy becomes ineffective or unpleasant, stimulation should be discontinued until its use is re-evaluated by a physician.
- ▶ Keep this device out of the reach of children.
- ▶ The units have no AP/APG protection. Do not use it in the presence of explosive atmosphere and flammable mixture.
- ▶ Only apply electrodes to healthy, uninjured skin.

EMS

- ▶ Only use original accessories from the manufacturer for your device.
- ▶ Safety of powered muscle stimulators for use during pregnancy has not been established.
- ▶ Caution should be used for patients with suspected or diagnosed heart problems.
- ▶ Caution should be used for patients with suspected or diagnosed epilepsy.
- ▶ Caution should be used in the presence of the following.
 - When there is a tendency to hemorrhage following acute trauma or fracture.
 - Following recent surgical procedures when muscle contraction may disrupt the healing process.
 - Over the menstruating or pregnant uterus; and
 - Over areas of the skin which lack normal sensation.

- ▶ Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium for sensitive skin, or alternate electrode placement.
- ▶ Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.
- ▶ Powered muscle stimulators should be kept out of the reach of children.
- ▶ Powered muscle stimulators should be used only with the leads and electrodes recommended for use by the manufacturer.
- ▶ Portable powered muscle stimulators should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.

WARNINGS

- ▶ The unit should be used only under the continued supervision of a physician.
- ▶ The unit devices have no curative value.
- ▶ The unit is a symptomatic treatment and as such suppresses the sensation of pain, which would otherwise serve as a protective mechanism.
- ▶ The unit is not effective for pain of central origin. (This includes headache)
- ▶ Electronic monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when electric stimulation is in use.
- ▶ Caution should be used in applying electric stimulation to patients suspected of having heart disease. Further clinical data is needed to show there are no adverse results.
- ▶ Electrodes should not be placed over the eyes, in the mouth, or internally.
- ▶ Do not use while sleeping.
- ▶ Do not use during pregnancy unless directed by your physician.
- ▶ The long-term effects of chronic electrical stimulation are unknown.
- ▶ Stimulation should not be applied over the carotic sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
- ▶ Stimulation should not be applied transcerebrally (through the head).
- ▶ Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.
- ▶ Stimulation should not be applied over, or in proximity to, cancerous lesions.

Steps to select a program (P)

Select P mode with


In program mode you can choose between 37 different programs for various areas of application. These are fully pre-programmed - just set the intensity for the connected channels and start the application!

EN

The program starts immediately when you set the intensity!

Most programs run in continuous mode, recognizable by a "C" shown in the lower right corner of the display. Some programs also have a pre-defined duration. The device ends the application automatically when the time expires. Please use these programs for the pre-defined duration only and do not extend or shorten the duration of the application.

To set a program, press the „MODE“ button until you reach P mode, recognizable by a P shown in the bottom-left corner of the display.

Press  or  to select the desired programs (P1 – P37). You find an overview over all available programs on pages 42 to 46.

When you chose the right program, please set the intensity for the connected channels (1-4).

Increase intensity 

Decrease intensity 

In P mode, the parameters are set by default and can not be changed!

Program mode (P)

All parameters are set by default. The selected program is automatically activated for all channels.

TENS – STIMULATION

EN

The programs have proven to be helpful in most cases for the suggested applications. However, they might be used for the treatment of other injuries or discomforts as well. Every patient responds individually to stimulation current.

TENS – Programs P1 - P16

P01	Classic TENS-mode
	To be applied in most existing types of pain (acute and chronic). Gate-Control-Effect.
P02	Low frequent TENS (Burst)
	Release of endorphins. Application for example in case of radiating pain on arms/legs/feet, muscular deeper-lying pain, lumbago-ischialgia syndrome, knee arthrosis
P03	Modulated TENS (massaging effect)
	Against habituation effect when applied for longer periods. Can be applied for most types of pain (acute and chronic).
P04	Alternating frequencies 15/2Hz
	Application for example in case of stiff joints, neck pains, lumbago, menstrual problems.
P05	Alternating frequencies 80/2 Hz
	Perfect for a variety of different pain conditions.
P06	The special program for nausea, feeling of sickness, travel sickness, seasickness
	The electrodes should be placed above the acupuncture point C6.
P07	Migraine and sensible parts of the body
	Due to the reduced duration of pulse these applications are to be specifically used in cases of tension headache, face pain, neck pain, herpes zoster, migraine.

TENS – Programs P1 - P16

P08	Chronic pain
	Application in case of pains in the neck, hip, shoulders, elbows, rheumatic pains, back pain (lumbago), menstrual problems and knee arthrosis.
P09	Rheumatic pain
	Pains due to rheumatoid arthritis, radiating pain on arms/legs/feet, through modulation suitable for long-term use.
P 10	Sweep-Mode 2 – 125 Hz
	In cases of therapy-resistant complaints, stimulation both on sensory and motor basis, against habituation effect.
P 11	Frequency change (2 – 100 Hz) and pulse duration modulation
	Applicable for most types of acute and chronic pain. With anti-habituation effect.
P 12	Frequency modulation
	Application in case of pain in the back, knee, shoulders, legs, elbows and neck, rheumatic pains.
P 13	Prof. Han - stimulation
	Muscular spasms, allodynia, duration: 30 minutes.
P 14	High-frequent TENS
	Acute pains, tensions, Gate-Control-Effect.
P 15	Combi-TENS-Mode
	4-Phases TENS. In cases of persisting pains. Consists of high and low frequent TENS as well as frequency/pulse duration/intensity modulation. Duration: 40 minutes.
P 16	Frequency-Wobbel-mode
	Pains in the wrist and fingers, long-term application. Duration: 40 minutes.

EMS – STIMULATION

Please note that the respective programs consist of several phases. Programs 17 – 28 contain a warming-up phase of approx. 60 seconds, followed by the stimulation itself which may also consist of various phases.

EN

EMS – Programs P17 - P30

P17	EMS – Strength (beginners) for small muscle groups
P18	EMS – Strength (beginners) for large muscle groups
P19	EMS – Strength (fully trained persons) for small muscle groups
P20	EMS – Strength (fully trained persons) for large muscle groups
P21	EMS – Muscle mass (beginners) for small muscle groups
P22	EMS – Muscle mass (beginners) for large muscle groups
P23	EMS – Muscle mass (fully trained persons) for small muscle groups
P24	EMS – Muscle mass (fully trained persons) for large muscle groups
P25	EMS – Endurance (beginners) for small muscle groups
P26	EMS – Endurance (beginners) for large muscle groups
P27	EMS – Endurance (fully trained persons) for small muscle groups
P28	EMS – Endurance (fully trained persons) for large muscle groups

P29	EMS – Relax I
	muscle regeneration, muscle relaxation – moderate intensity
P30	EMS – Relax II
	muscle regeneration, muscle relaxation – moderate intensity

MASSAGE - Programs P31 - P32

P31	Massage I
	consists of three different massaging phases
P32	Massage II
	consists of three different massaging phases

EN

What are small and large muscle groups?

The muscles can be divided into small and large muscle groups. Here you will find some examples for a better understanding:

Large muscle groups:

- ▶ thigh musculature
- ▶ buttock musculature
- ▶ abdominal musculature
- ▶ back musculature
- ▶ chest musculature

Small muscle groups:

- ▶ calf musculature
- ▶ biceps and triceps
- ▶ shoulder musculature
- ▶ forearms
- ▶ fingers and toes

STIMULATION OF THE PELVIC FLOOR MUSCLES

EN

Programs P33 - P37

Application only possible in conjunction with a vaginal or rectal probe (not part of the scope of delivery).

Noticeable but pleasant contractions of the pelvic floor musculature. Please alternate the programs in use to prevent habituation effects. Always use contact gel to ensure optimal stimulation.

Session duration: 20 minutes

Application: once or several times a day as necessary

Note: Please read the safety warnings carefully before using pelvic floor training.

INDICATIONS

A Probe electrode is the accessory of incontinence stimulator intended for the use of pelvic floor muscle stimulation. It sends gentle electrical impulses generated by the stimulator via the electrodes on the surface of each probe. By exercising the muscle group passively, it can help relieve male or female incontinence and support treatment.

Electrical stimulation via the pelvic nerves is a recognized treatment alternative for urinary incontinence. When treating stress incontinence, the aim of the electrical stimulation is to improve the function of the pelvic floor muscles. For urge incontinence, the aim is to inhibit involuntary bladder contractions by stimulating the nerves in the pelvic floor.

SAFETY WARNINGS & CONTRAINDICATIONS (pelvic floor training)

EN

- ▶ Read the instructions carefully before using the device.
- ▶ Keep this device out of the reach of children.
- ▶ Do not place electrodes in the area of the carotid arteries!
- ▶ Do not place electrodes on the front of the throat!
- ▶ Be careful with stimulation over metal implants.
- ▶ Do not place electrodes above your heart.
- ▶ Do not use electrodes on open wounds / injuries.
- ▶ Do not use if you have fever.
- ▶ Do not use if you have arrhythmia.
- ▶ Do not use if you have epilepsy.
- ▶ Never use the device when you are operating machines.
- ▶ Never use near explosive or flammable gases.
- ▶ Always TURN OFF the device before placing / removing electrodes.
- ▶ Do not use on overflow incontinence
- ▶ Do not use for cervical cancer
- ▶ Do not use if there are metal implants in the area of application
- ▶ Do not use in case of severe local inflammation
- ▶ Do not use when using a metal containing Intrauterine pessary (contraceptive device)
- ▶ Do not use in case of severe sensory disturbances in the area of application.
- ▶ Do not use for higher level haemorrhoid disease (when using a rectal probe)

- ▶ For users of pacemakers and during pregnancy, the device must only be used under medical supervision (!).
- ▶ In case of overflow incontinence, the unit should not be used.
- ▶ The device must not be used in case of cervical cancer, fistulas, fistulas, or the inversion of uterus / vagina.
- ▶ The device must not be used in case of infections in the vaginal or rectal area.
- ▶ For hygiene reasons, the probe may only be used by one person.
- ▶ The unit must not be used if the pelvic floor is completely denervated.

EN

PERFORMING PELVIC FLOOR EXERCISES

EN

1. Connecting the wires

Unpack the vaginal / rectal probe, clean it under running water and dry it. Connect the wires to the device and attach it to the probe.

2. Inserting the probe

Begin by emptying your bladder. The exercises should be carried out when lying down with your knees slightly bent. Apply contact gel to the probe and insert it slowly and carefully.

3. Program selection, performing the stimulation

Turn on the device. Select the desired program (P33 - P37) and start the stimulation by increasing the intensity by pressing the buttons \triangle ∇ of the respective channel.

The program starts immediately when you set the intensity!

4. Finishing the stimulation

Turn off the device and slowly remove the probe. **Attention: To prevent broken cables please grab the probe by the shaft and do not pull on the cables.** Use mild soap and water to clean the probe.

OPERATING THE DEVICE

Power On/Off button and intensity control


EN

The device can be switched on and off by pressing the On/Off button. The intensity can be reduced or increased by pressing the intensity buttons. The intensity can be set per channel from 0 to 100mA in steps of 1mA.

There are 4 LED lights at the top of the device. They shine yellow to indicate activity of the corresponding channel.

Mode control

There are 5 TENS modes available (B, N, M, SD1, SD2) and 3 EMS modes (C, S, A) as well as the P mode (programs 1-37). The modes can be selected by pressing the "MODE" button. If a TENS mode is selected, the LCD display shows "TENS". If an EMS mode is selected, "EMS" is displayed.

Set control

By pressing the  control, you may enter the setting you intend to adjust. You may start to set the value by pressing the  or  controls when the value is flashing.

Increment control

When pressing this button, the parameter will increase. This button controls the increase of settings.

Decrement control

When pressing this button, the parameter will decrease. This button controls the decrease of parameter.

Dual timer

The unit has a timer of 1-60 minutes and Continuous mode. It can be adjusted by pressing the  and the  or  controls. The application will end when the time is up.

Displaying time (L) for channels 1+2 / (R) for channels 3+4

Low battery indicator

If this symbol is shown, the battery is low and the batteries need to be replaced immediately.

How to set personal TENS parameters


The settings can be adjusted according to the following steps.

1. Turn on the device

After the electrodes are placed firmly on skin and the lead wires are plugged in the socket of device, press the  button.

2. Select mode

Select a mode by pressing the "MODE" control. The mode you select will be shown on the display. There are 5 modes available:


(B) Burst, (N) Normal, (M) Modulation, SD1 and SD2.

When a TENS-mode is selected, the LCD-Display will show "TENS". After a mode is selected, press  control to enter next setting. Then press the  or  control to change the settings. The settings are stored in the unit as soon as selected.


3. Set pulse width

The pulse width can be set from 30 μ s to 400 μ s. Press (SET) control to enter this setting. Then press the (▲) or (▼) control to change the settings. If no instructions regarding the pulse width are given in therapy, set the control to 70 – 120 μ s.

4. Set pulse rate

The pulse rate is adjustable from 2Hz to 150Hz. Press (SET) control to enter this setting. Then press the "Increment" or "Decrement" control to change the settings. If no instructions regarding the pulse width are given in therapy, set the control to 70 – 120 μ s.

5. Set timer




There are two adjustable timers available.

- (L) controls the application duration for the channels CH 1&2
- (R) controls the application duration for the channels CH 3&4

The settings are adjustable from 1 to 60 minutes. With the setting (C) the machine runs in continuous mode. Press (SET) control to enter the settings. Press (▲) or (▼) control to adjust the setting. Press the (▲) control when the Timer shows 60 minutes to set the timer to continuous stimulation. Both Timers can be adjusted in the same way.

(L) Timer CH1 & CH2 (R) Timer CH3 & CH4

Continuous

EN

6. Adjusting the intensity


The intensity can be adjusted from 0 – 100. Press or control to adjust the setting of the respective channels 1-4.

How to set personal EMS parameters

The settings can be adjusted according to the following steps:

1. Turn on the device

After the electrodes are placed firmly on skin and the lead wires are plugged in the socket of device, turn the on/off control.

2. Select mode


Select a mode by pressing the “Mode” control.

The mode you selected will show up on the top of liquid crystal display. There are 3 modes available: **C (Constant)**, **S (Synchronous)** and **A (Asynchronous)**.

When an EMS mode is selected, the LCD display will show “EMS”. After a mode is selected, press control to enter next setting. Then press the or control to change the settings. The settings are set as soon as they are selected.


3. Set ramp time

The ramp time controls the time of output current that increase from 0 to the setting level, and from the setting value back to 0. The ramp time is adjustable from 1 to 8 seconds.

4. Set contraction time (ON-Time)


The ON time controls the contraction time. By pressing the **SET** control, the contraction time is adjustable from 0 – 90 seconds. This time includes the ramp time twice, once for the ramp up and once for the ramp down.

Example: 10s ON time & 3s ramp time = 3s ramp up, 4s full intensity, 3s ramp down

5. Set pause time (OFF-Time)


The Off Time controls the time of relaxation between two cycles. By pressing the **SET** control, the relaxation time can be adjusted. Both channels' stimulation is cycled on and off by the contraction and relaxation settings. The range is adjustable from 0 to 90 seconds. In Alternate mode, the OFF time should be equal or more than the ON time.

(OFF TIME ≥ ON TIME)


6. Set pulse width

The pulse width is adjustable from 30 to 400 μ S. Press **SET** control to enter this setting. Then press the **▲** or **▼** control to change the settings. If no instructions regarding the pulse width are given in therapy, set the control to 70 – 120 μ s.

7. Set Pulse Rate

The pulse rate is adjustable from 1Hz to 150Hz. Press **(SET)** control to enter this setting. Then press the **(▲)** or **(▼)** control to change the settings. If no instructions regarding the pulse rate are given in therapy, set the control to 70 – 120 Hz.

EN

8. Set timer


There are two adjustable timers available.

(L) controls the application duration for the channels CH 1&2


(R) controls the application duration for the channels CH 3&4

The settings are adjustable from 1 to 60 minutes. With the setting **(C)** the machine runs in continuous mode. Press **(SET)** control to enter the settings. Press **(▲)** or **(▼)** control to adjust the setting. Press the **(▲)** control when the timer shows 60 minutes to set the timer to continuous stimulation. Both timers can be adjusted in the same way.


(L) Timer CH1 & CH2

(R) Timer CH3 & CH4

Continuous


9. Adjusting the intensity

The intensity can be adjusted from 0 – 100. Press **(▲)** or **(▼)** control to adjust the setting on the respective channels 1-4.

CONSTRUCTION OF THE DEVICE


Front

EN

Back

EN



PATIENT COMPLIANCE METER

This unit can store 60 sets of operation records. Total treatment time up to 999 hours can be stored.

EN



Check & delete individual records

Press “Mode” control and turn on the power simultaneously. The LCD will show the number of records and operation time. Press the  or  button to check each record. To delete a record, press  control for 3 seconds.

Check & delete accumulative record

At the individual records menu, press “Mode” control to switch to accumulative record menu. Press the “SET” control first, then press the “Mode” control simultaneously for 3 seconds and all of the records will be deleted followed by a beeper sound.


BATTERY RECHARGING

When the low battery indicator is shown, the batteries need to be replaced. Remove the batteries if the device is not likely to be used for some time.

EN

- ▶ **Make sure the power is off.**
- ▶ **Open the battery compartment (push the cover downwards)**
- ▶ **Remove the batteries**
- ▶ **Insert the new batteries. Please check that the polarity is correct!**
- ▶ **Close the battery compartment**

MAINTENANCE, TRANSPORTATION AND STORAGE OF THE DEVICE

EN

- ▶ Non-flammable cleaning solution (containing 70% of alcohol) is suitable for cleaning the device.
- ▶ Stains and spots can be removed with a cleaning agent.
- ▶ Do not submerge the device in liquids or expose it to large amounts of water.
- ▶ Return the device to the carrying box with sponge foam to ensure that the unit is well-protected before transportation.
- ▶ If the device is not to be used for a long period of time, remove the batteries and put it back into the carrying box and keep it in a cool, dry place.
- ▶ The packed TENS device should be stored and transported under the temperature range of -20°C ~ $+60^{\circ}\text{C}$, relative humidity 20% ~ 95%, atmosphere pressure 500 hPa ~ 1060 hPa.

SAFETY CHECKS

For safety reasons, review the following checklist once a week.

- ▶ Check the device for deformation of the housing or damage to the output sockets.
- ▶ Make sure that the descriptions and labels are not damaged.
- ▶ Check that the LED is on when the device is turned on.
- ▶ Check the cables and electrodes for damage.

- ▶ The device must be subjected to safety checks and maintenance by authorized technicians before use and each re-use, but at least every 24 months.
- ▶ Please consult your distributor if there are any problems with the device and accessories.
- ▶ The manual must always be carried with the device.

MALFUNCTIONS

Should any malfunctions occur while using your device, check the following points. If none of these can solve the problem, please contact the customer service.

- ▶ Check if the controls are correctly adjusted. Check if the values suggested by the physician have been adjusted.
- ▶ Check if the cables are correctly connected to the device.
- ▶ The display does not turn on? Replace the batteries.
- ▶ Check the cables for possible damage. The device may only be used when subjected to a load.
- ▶ There is no stimulation current? Check if the electrodes are correctly connected to the same channel and if the intensity is high enough.

CONFORMITY TO SAFETY STANDARDS

The STIM-PRO X9+ device is in compliance with the EN60601-1:2006/A1:2013 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance & EN60601-1-2:2015 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances - Requirements and tests

TECHNICAL DESCRIPTION

01	Channel	4 channels – with adjustable intensity
02	Intensity	Adjustable, 0 – 100 mA (in steps of 1 mA) with 500 ohm load on each channel
03	Pulse amplitude	Asymmetrical, Bi-phasic square pulse
04	Output voltage	0 – 50 V
05	Power source	4pcs 1.5V AA-batteries
06	Size	13.8cm (L) x 7.8cm (B) x 2.8cm (H)
07	Weight	276 grams (batteries included)
08	Pulse rate	Adjustable, from 1 – 150 Hz, 1 Hz/step
09	Pulse width	Adjustable, from 30 – 400 microseconds, 10 μ /step
10	ON time	Adjustable, from 2 – 90 seconds, 1 Sec./step
11	OFF time	Adjustable, from 0 – 90 seconds, 1 Sec./step
12	Ramp time	Adjustable, from 1 – 8 seconds, 1 Sec./step
13	Mode	5 TENS modes: B (Burst) N (Normal), M (Modulation), SD1 (40% modulation), SD2 (70% modulation) – 3 EMS modes: C (Constant), 5 (Synchronous), A (Alternate), P (Program-Mode) – 37 TENS programs (set by default)
14	TENS Burst mode B	Burst rate: Adjustable, 0.5 – 5 Hz – Pulse Width adjustable, 30 – 400 μ s – Frequency fixed = 100 Hz
15	TENS Normal mode N	The pulse rate and pulse width are adjustable. It generates continuous stimulation based on the setting value
16	TENS Modulation mode M	Modulation mode is a combination of pulse rate and pulse width modulation. The pulse rate and width are automatically varied in a cycle pattern. The pulse width is decreased by 50% from its original setting in 0.5 second, then the pulse rate is decreased by 50% from its original setting in 0.5 second. Total cycle time is 1second. In this mode, pulse rate (1 – 150 Hz) and pulse width (30 – 400 μ s) are fully adjustable.

17	TENS SD 1 mode	The SD1 (Strength-Duration) mode consists of automatic modulation intensity and pulse width in 40% range. The intensity is always increasing while the pulse width is decreasing and vice versa. The intensity is decreased by 40% while the pulse width is increased by 40% in 5 seconds. In the next 5 seconds, the intensity is increased by 40% while the pulse width is decreased by 40%. Total cycle time is 10 seconds. Pulse rate (1 – 150 Hz) and pulse width (30 – 400 μ s) are fully adjustable.
18	TENS SD 2 mode	The SD2 (Strength Duration) mode consists of automatic modulation intensity and pulse width in 70% range. The intensity is always increasing while the pulse width is decreasing and vice versa. The intensity is decreased by 70% while the pulse width is increased by 70% in 4 seconds. In the next 4 seconds, the intensity is increased by 70% while the pulse width is decreased by 70%. Total cycle time is 10 seconds. Pulse rate (1 – 150 Hz) and pulse width (30 – 400 μ s) are fully adjustable.
19	EMS C Constant mode	Constant stimulation based on setting value. Only pulse width and pulse rate are adjustable in this mode.
20	EMS S Synchronous mode	Stimulation of all channels occurs synchronously. The "ON" time including "Ramp Up" and "Ramp Down" time. Therefore, the setting of ON Time should be no less than two times of the "Ramp" time in this mode.
21	EMS A Asynchronous mode	Alternating stimulation. Used for the training of both agonist and antagonist. First, channels 1&3 run simultaneously, then channels 2&4. The duration of one cycle is the duration of the ON time, including the ramp time. Please note that the ON time can not be shorter than twice the ramp time.

21	TENS/EMS P program mode	See section Program mode (P) – pages 42 – 46
22	Patient compliance meter	This unit can store 60 sets of operation records. The maximum recordable time is 999 hours.
23	Low battery indicator	A low battery indicator will show up on the LCD when the battery is low.
24	Operating condition	Temperature: 0°C - 40°C Relative humidity: 30 % - 75 % Air pressure: 700hPa - 1060hPa
25	Notice	All technical values include a tolerance of +/- 5%

SCOPE OF DELIVERY

Each STIM-PRO X9+ is fully equipped with the following standard equipment.

- ▶ **1 piece** **TENS EMS combo device**
- ▶ **8 pieces** **Self-adhesive electrodes**
- ▶ **4 pieces** **1.5V - AA batteries**
- ▶ **1 piece** **Storage case**
- ▶ **4 pieces** **Electrode cables**
- ▶ **1 piece** **Instructions manual**

WARRANTY

EN

All TENS models carry a warranty of 24 months from the date of delivery. The warranty applies to the stimulator only and covers both parts and labor relating thereto. The warranty does not apply to damage resulting from improper handling, the failure to follow the operating instructions, loss or dropping.

Manufacturer:
axion GmbH
Mollenbachstr. 13
71229 Leonberg, Germany
www.axion.shop/en

1.  Advice for operation
2.  Degree of electrical protection Type BF
3.  Do not insert the plug into the power supply socket of 230V
4.  IP protection class
5.  Do not dispose of via residual waste
6.  Increment
7.  Decrement
8.  Read the instruction manual
9.  Direct current (DC)
10.  Manufacturer
11.  Serial number

Manual Version: 1.2

Release Date of the Manual: December 2018

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