bob

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Nature is what takes place
around us since the beginning of time,
even before mankind.

The ancient Egyptians, observing an Ibis
(which uses its long beak to draw water
and introduce it into the intestine),
understood the natural need
of internal purity.

Nature knows itself
and it is always perfect.

*Looking inwards and feeling perfect
is a right, and it is natural...*

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INDEX

PACKAGE CONTENTS	27
WHAT IS BOB	28
WHY IT IS USED	29
CONTRAINDICATIONS	30
WARNINGS	31
INSTALLATION OF SPECIAL FILTER TO THE FAUCET	33
HOW TO USE BOB	35
FREQUENCY AND DURATION OF CLEANSING	38
SUGGESTIONS FOR RECTAL CLEANSING	39
MAINTENANCE AND STORAGE	41
ADDY UNIVERSAL ADAPTOR	42
ADDY WARNINGS	43
ADDY INSTRUCTIONS	44
SPARE PARTS AND ACCESSORIES	46
DISCLAIMER	47

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PACKAGE CONTENTS

- **Special Aerator Filter (A1, A2) bag, containing:**
 - Special filter holder with internal thread (X)
 - Special filter holder with external thread (Y)
 - Aerator Filter(K)
 - Gasket (W - with rectangular section) for Special filter holder with internal thread
 - Gasket (Z - with “step” section) for Special filter holder with external thread
- **Attachment (B) with 3-mt. / 10’ spiral tube (C), fitting (D), handle (F) and end fitting (G).**
- **Rigid aluminum tubes (H1-H2) with premounted nozzle (I).**
- **Bag containing n. 20 mini cannulas (L).**
- **Key for removing filter holder (M).**
- **Universal adapter Addy (O)* with fastening key (P)***
- **Bob carrying bag.**
- **Multi-language user manual.**

The letters in the brackets refer to the general diagram. Item “E” not included, can be purchased separately as an accessory.

Included only in the **Bob Travel packages.*

bob

WHAT IS BOB

***Bob** is a medical device for anal and rectal cleansing. **Bob** connects directly to your faucet thanks to its special filter and valve for the automatic control of water pressure (patented system).*

bob

WHY IT IS USED

Bob is used for deep anal cleaning and to help remove fecal matter residue present on the walls of the rectum-sigma. Bob can be used in cases of constipation, bloating, flatulence, intestinal discomfort in general and for proper rectal hygiene.

CONTRAINDICATIONS

***Bob** should not be used in cases of: Enteritis proctitis (Inflammation of the small and large intestines and of the anus); Acute abdominal pathologies; Pregnancy; Bleeding hemorrhoids and anal fissures; Perforation of the rectum; Rectum cancer; Hematochezia; Melena; Recent colon, rectum or anus surgery; Heart/respiratory failure; Epilepsy and psychosomatic syndromes; Diarrhea; Dismetabolic Syndrome; Arterio-venous disorders; Diverticular disease in its active phase; Nausea or vomiting; Serious state of dehydration; Prostate surgery.*

In addition to the pathologies listed above, always consult your physician before using this medical device.

WARNINGS

- Before using **Bob**, check the integrity of the device and its components, as well as the correct assembly of all parts;
- If you experience any discomfort after using the device, stop using of **Bob** immediately and seek medical advice from your doctor;
- **Bob** should be used with warm water, never to exceed 98.6°F. Ideal water temperature should be maintained between 29°C and 36°C (84°F – 96°F). After connecting **Bob** to the faucet and without pressing the handle lever (F), use your hand to check the temperature of the water coming out of the valve for the automatic control of water flow (B - Fig. 10). Keep in mind that the sensitivity to heat or cold of the rectal channel is considerably greater than that of your hand.
- If at any point while using **Bob** if you feel that the water temperature may be too hot or too cold, immediately stop the water flow by releasing the handle lever (F), and adjust the temperature from the faucet. NOTE: some water heaters may take a while to start; to solve this issue, simply open the hot water from another faucet (e.g. bath tub) to kick start the water heater.
- The rectal cannulas are disposable; every rectal cannula must be used one time only, even if the same person uses it.

- *Water used for the rectal cleansing must be clean (normal drinkable water supplied by the domestic water network). If you believe that the water coming out of your tap is not of good quality or is not drinkable, we recommend using a **Bob** medical filter with double membrane of 1.2 microns (E), which can be purchased separately.*
- *Colon cleansing should always be performed with light to moderate water flow, especially during the first washes. Water flow and pressure can be adjusted directly from the faucet. In case of diverticulosis in particular, if advised by your doctor and while keeping in mind all other counterindications, the water flow must be adjusted to a minimum.*
- ***Bob** is not intended for use and it should not be used by anyone under the age of 16 without adult supervision and the authorization of your doctor. Do not use the device on someone who is sleeping or unconscious, an infant, or a patient with altered mental status or decreased skin sensation.*

INSTALLATION OF SPECIAL FILTER TO THE FAUCET

*NOTE: The installation of the special filter holder on the faucet must be done only the first time you use **Bob**: the exclusive **Bob aerator filter does not need to be removed and replaced after each use** because it allows normal use of the tap (fig.5), providing water aeration and filtration. (The special filter holder **Bob** connects to all faucets with removable filter holder according CEN-EN 246 standards.)*

*To connect **Bob** special aerator filter (A1/A2) to your faucet (sink, bidet, or shower hose**):*

1 - Unscrew the aerator filter holder from your faucet, making sure the gasket is also removed (fig. 1). If necessary, use the key (Q).

2a - If your faucet has an internal thread (fig. 2a), insert the aerator filter (K) with rounded side up inside the filter holder with external thread (Y) and place the gasket with step section (Z) on the aerator filter (fig. 2a).

2b - If your faucet has an external thread (fig. 2b), insert the aerator filter (K) with rounded side up inside the filter holder with internal thread (X) and insert the gasket with rectangular section (W) (fig. 2b).

3- Thread the new filter holder to the tap (fig. 3a/3b) until it stops using the connector (B - fi g. 4) with dry hands, turning counterclockwise.

to follow>

bob

*** To connect **Bob** to the shower hose, use the Multi-function Shower Joint (accessory - can be purchased separately. See the “Spare Parts and Accessories” below for more info).*

***Bob** can be connected to the tap also via the universal adapter **Addy** (O) without replacing the aerator filter, ideal while travelling (**Addy** is included only with **Bob Travel**; you can purchase it separately for **Bob** - see “Spare Parts and Accessories”.)*

HOW TO USE BOB

1. Thread the rigid tube (H1) onto the rigid curved tube (H2) and then onto the handle (F) (Fig. 15) so that the water flow is in line with the handle (fig. 16). Make sure the gray washers are mounted on the threads of the two rigid tubes.
2. If the water spout of the nozzle (I) is not aligned with the handle (F), as shown in Fig. 16a, slightly loosen the end fitting of the handle (G) and rotate the rigid tube (H1 + H2), to align the nozzle (I) to the handle (F) (Fig. 16b). Thread on again the handle's end fitting (G) to block the rotation of the rigid tube (H1 + H2).
3. For a deeper cleansing of the rectum, use the mini rectal cannulas by threading one on the nozzle (I).
4. With dry hands, insert the quick connector (B) to the special filter holder (A1 - A2)) or Addy (O) previously mounted on your faucet, and rotate counterclockwise (Fig. 17a - 17b) to secure Bob to the faucet
5. Turn on the water and, without pressing the handle lever (F), regulate the desired temperature. We recommend warm water (29°C/84°F - 36°C/96°F approx). **Never use water at a temperature higher than 36 ° C/96.8°F.** Use your hand to check the temperature of the water coming out of the control valve (B - Fig. 18)
6. Direct the nozzle in the toilet or sink, and press the handle lever handle (F) to check that water is flowing from the nozzle and

bob

that the temperature is stable. Release the handle lever to stop the water flow.

7. Sit on the toilet, position the rigid tubes (H1 + H2) of Bob between your legs so that the water flow is in line with the sphincter (Fig. 19). In order to keep the nozzle closer to the sphincter, lightly press the handle lever down. If you are using the rectal cannulas, gently insert the tip into the anus.

Stop using if the tip is hard to insert.

Forcing the cannula into the rectum can cause serious injury. In case of bleeding or pain, get immediate medical care.

8. While remaining seated, gradually open the water flow to start the cleansing. Always start the water flow slowly when you let water in the rectum. (CAUTION: check “Frequency and duration” and “Suggestions for rectal cleansing” for additional instructions).

9. Each cleanse consists of these 3 phases:

a) Filling (flow of the water in to the rectum);

Stop filling with water immediately if you feel full, even if this is the first filling and performed only for a couple of seconds;

b) Holding (water is held inside the rectum until the urge to evacuate is strong, usually 1 to 5 minutes);

to follow>

bob

*c) Emptying (evacuation); during evacuation, the rectal cannula may be left in or taken out of the anal duct; if **no urge is felt after 5 minutes try to empty the bowel; call a doctor immediately is no water comes out after 30 minutes** as this might be a sign of dehydration.*

10. After the cleansing has ended, remain seated and close the water flow by releasing the handle lever (F) and closing the faucet.

11. Make sure you expel all water from the colon. If all water is not expelled, subsequent light leakage or slight discomfort may occur.

12. Disconnect the quick connector (B) from the faucet filter holder and wash the rigid tube with the nozzle by unscrewing them from the handle, cleaning them with care and rinsing thoroughly (Fig. 22). Discard the rectal cannula.

*13. Store **Bob** in its case (see “maintenance and storage”).*

For a correct execution of the cleanse, it is recommended to alternate phases of filling with phases of holding, each of the duration of 10-15 seconds;

- *For the first cleansings, repeat the steps of filling and retention at least twice. For the subsequent washings, repeat these steps for a maximum of five times. In any case, immediately stop the water injection if you feel a sensation of bloating (even if the water injection occurs only once and for a few seconds).*
- *When the last phase of water intake is complete, try to retain the water inside (without any particular efforts and, in any case, no longer than 4-5 min.) before evacuating.*
- *The total amount of water used during the cleanse is subjective and it depends on your personal body structure and condition. Normally, for each cleanse, you should not insert more than 1,5 lt / 1.5 qt of water overall (1,5 lt / 1.5 qt passes through the pipe on average in about 40 seconds with the tap and water flow fully open and at a normal water pressure of about 2.5-3 bar);*
- *The total duration of the cleansing is variable and subjective. It depends, not only on the time needed for the water injection, but also on the time it takes for one to retain the water inside the rectum, and the time required for the evacuation which may vary depending on the individual.*
- *It is recommended to do no more than one wash per week for the first month and no more than 3 washes per week afterwards, unless instructed otherwise by your doctor.*

SUGGESTIONS FOR RECTAL CLEANSING

- *To facilitate evacuation, take a deep breath, contract the stomach and massage the area above the pubic bone and the left side of the stomach for about 1-2 minutes. The evacuation will occur naturally when the intestine is ready and the stool is softened.*
- *The sensitivity to water pressure and temperatures are subjective and therefore the user must evaluate and adjust them, although never to exceed the parameters set in these instructions.*
- *During the first uses, it is recommended that **Bob** be used with a light to moderate flow of water. In the case of excessive water flow, reduce it via the flow regulator;*
- *After each cleansing, it is advisable to take milk enzymes in order to reintegrate intestinal bacterial flora.*
- *In several cases, it may occur that the cleansing effect may be particularly effective and that the stimulus of evacuation may occur several times over the following 2-3 hours. Therefore, it is recommended, especially during the first cleansings, that a restroom be readily available for use during these following 2-3 hours. If this effect should occur, for successive cleansings, the quantity of water used should be reduced.*

to follow>

bob

- *Inside the quick connector there is a valve that automatically checks the water flow - patented system - which will automatically drain excess water directly into the basin/ bidet. You may hear a sound similar to a whistle when the water flow is open: this is normal and it means that the valve is working properly. To eliminate this sound simply decrease the water flow from the faucet.*

MAINTENANCE AND STORAGE

- *After each cleansing, blow air into the fitting (B) to expel any remaining water from the pipes in order to keep **Bob** in a clean and hygienic condition (Fig. 23).*
- *The rigid pipes (H1 and H2) should be detached from the handle and washed after each use before storing the product. It is also advisable to clean all the other parts (Fig. 22);*
- *In order to avoid damaging the rubber or plastic parts, do not clean the product with aggressive detergents;*
- *Before using the product, connect **Bob** to the faucet and let hot water flow through the device and out of the nozzle for at least 30 seconds;*
- *Between uses it is recommended that the product be kept well clean and dry inside the case.*

bob

ADDY UNIVERSAL ADAPTOR

*Included only with **Bob - Travel** version. It can be purchased separately as an accessory for the **Bob** (basic version.)*

***Addy** (O) allows you to connect **Bob** to the tap in a quick and easy way, without replacing the aerator filter holder (A1 / A2). It is ideal for use while travelling or in case of non-standard tap threads. **Addy** can be connected to all the taps with water filter holder with an outside diameter from 16mm to 25mm and minimum length of 5 mm (fig.16).*

ADDY WARNINGS

- *When using Addy, always open the water flow gradually and up to a moderate flow max. For a correct use of Bob it is sufficient to open the faucet until the water begins to flow from the valve located in the connector.*
- *Do not touch or move Addy while the faucet is on.*
- *For the first few moments, there could be some minor water leakage from the upper part of Addy but this should quickly stop. If the leakage persists or becomes excessive, repeat the process of connecting Addy to the faucet.*
- *In case of leakage from your faucet's aerator filter while Addy is connected, check that the your filter is tighten well on the faucet.*
- *If Addy is not used for a long time, the screws may be harder to move: without disconnecting Addy from the faucet, try to unscrew and tighten the screws several times in order to facilitate their movement again, without force. Should the screws not move, spray some lubricant in the screw holes where the fastening key would be inserted.*
- *Excessive pressure during screwing/ unscrewing and/or repeated use of Addy on the same faucet may cause slight scratches on the faucet's filter holder (for the extended use of Bob on the same faucet, use special aerator filters).*
- *If Addy is left connected to the faucet, check that the hold is strong before each use, making sure that the 3 screws are well tightened.*

ADDY INSTRUCTIONS

1. *Open the three screws so the tips are not showing inside the Addy opening, and check that the gasket is properly in place, without visible deformation (fig. 17).*

2. *Holding Addy with the gasket and opening facing upwards (fig. 18), push it up against the faucet until the gasket completely adheres to the filter holder, making sure that:*

- *one of the 3 screws is aligned with the faucet and the two remaining screws are to the side, and none of the screws is under the neck of the faucet (Fig. 19);*

- *Addy follows the inclination of the faucet's filter holder, allowing the gasket to completely adhere and ensure tightness (fig. 20);*

- *Addy is as close as possible to the faucet's filter holder (Fig. 21).*

Special Cases:

- *Some filter holders have a slight indentation: make sure that the tips of the screws do not lean against the indentation in order to prevent Addy from disconnecting.*

- *In case of faucets that are slightly flared/rounded (ex. kitchen faucets), mount the screws on the straight side of the filter and never on the flared parts to prevent Addy from disconnecting from the faucet (Fig. 22).*

3. *Holding Addy to the faucet with the aid of your thumb (fig. 23), use the key (P) to tighten the three screws to a secure hold on the faucet, following the direction indicated by the arrow located on the cylindrical housing of the screws (fig. 24).*

bob

- 4. Once you tighten one of the screws all the way, double check the other two screws for tightness, never forcing them (fig. 24). While tightening the screws, be careful not to damage the inner gasket so that the hold is not compromised (Fig. 17). You are now ready to connect your device.*
- 5. When you are finished, disconnect your device and unscrew one of the three screws to remove Addy from the faucet (Fig. 13). Do not force the screws once they've reached the end of the travel.*
- 6. Before storing Addy, dry it and leave the screws slightly unscrewed, making sure however that the screw tips do not deform the gasket (Fig. 17).*

*You can purchase the complete **Bob** product line and accessories on our website:
www.bobispop.com*

SPARE PARTS AND ACCESSORIES

Mini rectal cannulas: to install on the nozzle for a deeper cleansing. Available in packs of 30 units.

Rectal cannulas: available in refill packages of 30 units

Double filter membranes of 1.2 microns: for filtering water from bacteria and suspended impurities larger than 1.2 microns.

Universal adapter Addy: allows the use of Bob without removing your filter holder from the faucet. Its use is ideal for traveling, in hotels or in the case of faucets with non-standard threading.

Multifunction shower joint: to connect Bob to the shower hose.

Metal filter holder for faucet: available with standard internal and external threads.

Adapters: for faucets with non-standard threads.

Extension: 3-meter / 10 feet spiral tube to connect Bob to the faucet when the distance from the faucet to the toilet is longer than 3 meters / 10 feet covered by the provided tube.

DISCLAIMER

The manufacturer, distributors and retailers refuse any liability for any damages that may, directly or indirectly, occur to people or property from non-compliance to all the requirements specified in this user manual, its installation, use, counter indications and maintenance.

***Bob** products are not and do not substitute any equipment that is used in hydro-colon therapy centers.*

The company reserves the right to make partial or complete modifications to the product, which may slightly differ from the descriptions in the present manual.

For any clarification or information please visit www.bobispop.com or contact customer care via email at info@bobispop.com (or supportusa@bobispop.com only for USA and Canada).

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Int. Pat. Pending

 Water Powered srl - Via Rep. Marinare, 128,
80147, Napoli NA, Italy

 Water Powered srl uses only solar power