User's Manual

# KetoSens® BT

Blood  $\beta$ -Ketone Monitoring System

- Easy Data Transfer by Bluetooth Connection
- 2.5" Large, Backlit Display
- Strip Ejector for Easy Ejection of Used Test Strips

For single user / home use only

PKTAA-0000005 REV0 2021/01

Ketosis Range Card on the last page of the English User's Manual

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# KetoSens® BT

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## Welcome to the KetoSens® BT Blood $\beta$ -Ketone Monitoring System

Thank you for choosing the KetoSens® BT Blood  $\beta$ -Ketone Monitoring System. The system provides you with rapid and convenient blood  $\beta$ -Ketone *in vitro* (i.e., outside the body) diagnostic monitoring. You'll receive an accurate result in eight seconds with a small 0.5  $\mu$ L blood sample.

### Bluetooth® wireless technology

Bluetooth® wireless technology is used by some smartphones and many other devices. KetoSens® BT meter uses Bluetooth® wireless technology to pair\* and to send your  $\beta$ -Ketone results to a smartphone. The KetoSens® BT meter is designed to work with the SmartLog™ mobile app. When using the KetoSens® BT System, we suggest you pair your KetoSens® BT meter with your smartphone and track your results.

Manufacturer: i-SENS, Inc.  
43, Banpo-daero 28-gil, Seocho-gu, Seoul 06646, Korea  
www.i-sens.com

Contact Customer Service:

**1-800-429-5001**

9 am–9 pm EST, Monday–Saturday

At all other times or in case of emergency, please contact your healthcare professional or emergency medical response.

- No part of this document may be reproduced in any form or by any means without the prior written consent of i-SENS.
- The information in this manual is correct at the time of printing. However, i-SENS reserves the right to make any necessary changes at any time without notice as our policy is one of continuous improvement.

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## **WARNING:**

Turn off the Bluetooth® feature in areas where the use of wireless devices is restricted, such as hospitals, some healthcare professional offices, and airplanes.

## **Trademarks**

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by i-SENS, Inc. is under license. All other trademarks and trade names are those of their respective owners.

## **FCC Radiation Exposure Statement:**

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with minimum distance 8 inches (20 cm) between the radiator and your body.

Any changes or modifications not expressly approved by the manufacturer could void the user's authority to operate this equipment.

### **Note**

The SmartLog™ mobile app may not be compatible with all smartphones.

# Table of Contents

## Information

Important Information: Read This First! .....	6
Specifications .....	10
KetoSens BT Blood $\beta$ -Ketone Monitoring System .....	12
Inserting or Replacing the Batteries .....	13
KetoSens Blood $\beta$ -Ketone Test Strip .....	15
KetoSens BT Blood $\beta$ -Ketone Meter .....	17
KetoSens BT Blood $\beta$ -Ketone Meter Display .....	18

## Preparation

Setting Up Your System.....	20
Adjusting the Date and Time.....	20
Setting Up Bluetooth® .....	21
Adjusting the Date and Time.....	23
Setting the Sound On/OFF .....	25
Turning on the Strip Expiration Date Indicator .....	26
Setting the Strip Expiration Date Indicator .....	27
Checking the System.....	28
Control Solution Testing.....	29
Comparing the Control Solution Test Results.....	32

## Testing

Using the Lancing Device.....	33
Preparing the Lancing Device.....	34
Preparing the Meter and Test Strip .....	36
Applying Blood Sample .....	36
Discarding Used Lancets.....	40
HI Message.....	41

## Additional Functions

Meter Memory .....	42
Viewing Test Results Stored in the Meter's Memory .....	42
Setting the Alarm Function .....	44
Setting the Time Alarms (alarm 1–3) .....	44
Caring for Your System.....	46

## Maintenance

Understanding Error Messages .....	52
General Troubleshooting.....	54
Performance Characteristics .....	55
EMC Table .....	57
Warranty Information .....	58

## Important Information: Read This First!

To receive safe and optimum system benefits, please read the entire manual contents before using the system.

### Intended use:

The KetoSens BT Blood  $\beta$ -Ketone Monitoring System is intended to be used for the quantitative measurement of  $\beta$ -Ketone (beta-hydroxybutyrate) level in capillary whole blood drawn from the fingertip. The KetoSens BT Blood  $\beta$ -Ketone Monitoring Systems are for self-testing outside the body (for *in vitro* diagnostic use) in the home as an aid to monitor the effectiveness of diabetes control program. The system is not intended for use in the diagnosis of or screening for diabetes mellitus and is not intended for use on neonates. This system is intended to be used by a single person and should not be shared.

The KetoSens Blood  $\beta$ -Ketone Test Strips work with the KetoSens BT Blood  $\beta$ -Ketone Meter to quantitatively measure Blood  $\beta$ -Ketone in capillary whole blood samples drawn from the fingertip.

The KetoSens Control Solutions are intended for use with the KetoSens BT Blood  $\beta$ -Ketone meter and KetoSens Blood  $\beta$ -Ketone test strip to check that the meter and the test strip are working together properly and that the test is performing correctly.

- Please use this device only for the intended use described in this User's Manual.
- Please follow the suggested cleaning and disinfection procedures described in this User's Manual.
- KetoSens Blood  $\beta$ -Ketone Test Strips are intended for self-testing only. Follow the instructions provided by your healthcare professional on the safe disposal of used test strips.

## Important Safety Information:

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members. Do not use on multiple patients.
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed pre-cleaning and disinfection.
- “FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication” (2010) <http://wayback.archive-it.org/7993/20170111013014/http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
- “CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens” (2010) <http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

## Limitations

- Do not use for diagnosis of or screening for diabetes mellitus.
- The system should not be used to test neonates.
- Do not test samples other than fresh capillary whole blood obtained from the fingertip.
- Do not use the system at altitudes above 10,000 feet (3,000 meters).
- Do not use when hematocrit is outside the acceptable hematocrit range for testing of 20 % to 55 %.
- Severe dehydration (excessive water loss) may cause inaccurate results. If you believe you are suffering from severe dehydration, consult your healthcare professional immediately.
- For *in vitro* diagnostic use only.
- Critically ill patients should not be tested with this device.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock.
- Incorrect result may occur in individuals who are dehydrated.
- The meter and lancing device are for single patient use. Do not share these items with anyone, including other family members! Do not use on multiple patients!
- Do not reuse; each test strip is for single use only.
- Do not use when humidity is higher than 90 % and lower than 10 %, as extremes in humidity may affect results.
- For single patient use only.
- For over-the-counter use.

## Important Information

- The KetoSens BT Blood  $\beta$ -Ketone Monitoring System is intended for self-testing outside the body (*in vitro* diagnostic use).
- $\beta$ -Ketone in blood samples reacts with the chemical in the test strip to produce a small electrical current. The KetoSens BT meter detects the electrical current which reflects the amount of  $\beta$ -Ketone in the blood sample.
- The KetoSens BT Blood  $\beta$ -Ketone Meter is designed to minimize code related errors in monitoring by using the no-coding function.
- The KetoSens BT Blood  $\beta$ -Ketone Meter should be used only with the KetoSens Test Strips.

*If you need assistance, please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. At all other times or in case of emergency, contact your healthcare professional or emergency medical response.*

This device is not intended for use in healthcare or assisted-use settings such as hospitals, physician offices or long-term care facilities because it has not been cleared by the FDA for use in these settings, including for routine assisted testing or as part of glycemic control procedures.

Use of this device on multiple patients may lead to the transmission of Human Immunodeficiency Virus (HIV), Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), or other bloodborne pathogens.

# Specifications

## Product specifications

Measurement range	0.1–8.0 mmol/L
Sample size	Minimum 0.5 $\mu$ L
Test time	8 seconds
Sample type	Fresh capillary whole blood (fingertip)
Assay method	Electrochemical Amperometry
Battery life	1,000 tests
Power	Two 3.0 V lithium batteries (disposable, type CR2032)
Memory	1,000 test results
Size	106 x 58 x 17 (mm)
Weight	72.6 g (with batteries)
Bluetooth®	<ul style="list-style-type: none"><li>• Frequency range: technology 2.4–2.4835 GHz</li><li>• Operating range distance: maximum 32 feet (10 meters) unobstructed</li><li>• Operating channels: 40 channels</li><li>• Security encryption: 128-bit AES (Advanced encryption standard)</li></ul>

## Operating ranges

Temperature	50–104 °F (10–40 °C)
Relative humidity	10–90 %
Hematocrit	20–55 %

## Storage/Transport conditions

Temperature	Meter (with batteries)	32–122 °F (0–50 °C)
	Test strip	39–86 °F (4–30 °C)
	Control solution	46–86 °F (8–30 °C)
Relative humidity	Test strip	20–80 %

# KetoSens BT Blood $\beta$ -Ketone Monitoring System

**KetoSens BT Blood  $\beta$ -Ketone Monitoring System includes the following items:**

- KetoSens BT Blood  $\beta$ -Ketone Meter
- User's Manual
- Quick Reference Guide
- Bluetooth® Pairing Quick Guide
- Batteries (2)

**KetoSens BT Blood  $\beta$ -Ketone Monitoring System may include the following items:**

- KetoSens Blood  $\beta$ -Ketone Test Strips
- Lancets
- Lancing Device
- Logbook
- Carrying Case

The KetoSens Blood  $\beta$ -Ketone Control Solution Level A and B are sold separately and can be purchased by contacting Customer Service.

- Check all the components after opening the KetoSens BT blood  $\beta$ -Ketone monitoring system package. The exact contents are listed on the main box.
- The cable for data management software can be ordered separately. Please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.


## Inserting or Replacing the Batteries

The KetoSens BT meter uses two 3.0 V lithium batteries. Before using the meter, check the battery compartment and insert batteries if empty.

When the  symbol appears on the display while the meter is in use, the batteries should be replaced as soon as possible. The test results may not be saved if the batteries run out completely.



### Step 1

Make sure the meter is turned off. Push the cover in the direction of the arrow to open the battery compartment.


### Step 2

Remove the old batteries one by one by lifting with the index finger and pulling it out with your thumb and index finger as shown in the figure on the right. Insert two new batteries with the + side facing up and make sure the batteries are inserted firmly.

### Step 3

Place the cover on the battery compartment. Push down until you hear the tab click into place.


### Note

Removing the meter batteries will not affect your stored results. However, you may need to reset your meter settings. See page 20.

## KetoSens Blood $\beta$ -Ketone Test Strip

The KetoSens BT Blood  $\beta$ -Ketone Monitoring System measures blood  $\beta$ -Ketone quickly and accurately. It automatically absorbs the small blood sample applied to the narrow edge of the strip.

### Warning:

- The KetoSens test strip should be used only with fresh capillary whole blood samples drawn from the fingertip.
- Do not reuse test strips.
- Do not use test strips past the expiration date.
- Test strips in new, unopened foil packet can be used up until the expiration date printed on the foil packet if the test strips are used and stored according to its storage and handling methods.

- Store foil packet in a cool and dry place between 39–86 °F and 20–80 % Relative Humidity.
- Keep the foil packet of test strips away from direct sunlight or heat.
- Store unused test strips in their original foil packet to avoid damage or contamination.
- Open the foil packet to take out a test strip for testing and use the test strip immediately. Do not use the test strip if the foil is damaged or torn.
- Avoid getting any liquid or moisture in the test strip vial or box. This can affect the test strips and cause inaccurate test results.
- Do not apply samples other than capillary whole blood or control solution to the test strip.
- Handle test strips only with clean and dry hands.
- Do not bend, cut, or alter test strips in any way.
- For detailed storage and usage information, refer to the KetoSens test strip package insert.
- Keep the meter and testing supplies away from young children. The test strips and foil packet may cause choking hazards.

#### **Note**

You can get more safety information at FDA Public Health Notification (<http://wayback.archive-it.org/7993/20170111013014/http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>) or at CDC Clinical Reminder (<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>).

# KetoSens BT Blood $\beta$ -Ketone Meter

## Data Port

Used to transfer data from the meter to a computer with a cable

## ◀, ▶ Button

Turns the meter on, selects or changes information

## Test Strip Port

Insert test strip here

## Display

Shows results, messages

## S Button

Turns the meter on/off, confirms menu selections, and changes information


## Test Strip Ejector

Slide down to discard the used strips

## Note

- The unit of measurement is fixed and it cannot be changed by the user.
- The cable for data management software can be ordered separately. Please contact Customer Service: 1-800-429-5001 (9 am–9 pm EST, Monday–Saturday).

# KetoSens BT Blood $\beta$ -Ketone Meter Display

<b>1 Memory recall mode</b>	appears when test results stored in the memory are displayed
<b>2 Bluetooth® symbol</b>	appears when Bluetooth® is turned on
<b>3 Mute symbol</b>	appears only when the sound is set to OFF
<b>4 Control Solution flag</b>	appears when the control solution test results are saved or displayed
<b>5 Blood insertion symbol</b>	indicates meter is ready for the application of a drop of blood or control solution
<b>6 alarm</b>	appears when the time alarm has been set
<b>7 mmol/L</b>	unit for measuring blood $\beta$ -Ketone
<b>8 Battery symbol</b>	indicates meter batteries are running low and need to be replaced

**9 KETONE**

appears when the KetoSens strip is inserted, and when the test results are saved or displayed

**10 Test results**

**11 Month/Day/Hour/Minute**

 **Note**

It is recommended to check if the display screen on the meter matches the illustration above every time the meter turns on. Do not use the meter if the display screen does not exactly match the illustration as the meter may show incorrect results.

## Setting Up Your System

Press and hold the **S** button for 3 seconds to enter the SET mode. After all settings are finished, press and hold the S button for 3 seconds to turn off the meter.

Press the ◀ or ▶ button to change values. Press and hold the ◀ or ▶ button to scroll faster.

### Adjusting the Date and Time

#### Step 1 Entering SET Mode

Press and hold the **S** button for 3 seconds to enter SET mode. After all the segments flash across the screen, 'SET' will show up.

Press the ◀ or ▶ button to select 'YES' and press the **S** button to go to the next step.

#### Note

Follow steps 2–3 to pair your meter and smartphone. Pairing allows the meter to communicate wirelessly with your smartphone. Ensure that devices are within the maximum Bluetooth® range (10 meters). Before pairing your meter and smartphone, download and install the SmartLog™ mobile app on your smartphone.

## Setting Up Bluetooth®

### Step 2 Bluetooth® Pairing

#### Note

For details about the Bluetooth® pairing process, refer to the Bluetooth® Pairing Quick Guide enclosed in the KetoSens BT product box.

① **Meter:** Your KetoSens BT meter will have Bluetooth® feature turned on. If you want to turn off Bluetooth®, press the ◀ or ▶ button and **S** to confirm your selection. The meter will go to step 3, Setting the Year.

② **Smartphone:** If you **want to pair** (connect) your meter to your smartphone, launch the SmartLog™ mobile app and find the **Accessories** menu on your smartphone.

③ **Smartphone:** Select the KetoSens BT model from the meters list, and then tap **Bluetooth Register > Next**.

④ **Meter:** With the meter turned off, press and hold the ▶ button for three seconds to enter pairing mode.

The Bluetooth® symbol and **bT** will appear and **YES** will blink at the bottom of the screen.

⑤ **Meter:** Press the **S** button to select **YES** and set the Bluetooth® pairing type setting. The pairing type number **2** and **tYPE** (type) will appear, and the Bluetooth® symbol will blink at the top.

⑥ **Smartphone/Meter:** Follow the on-screen instructions on the SmartLog™ mobile app to search your meter, and complete the pairing process.

⑦ **Meter:** When the pairing process is complete, the meter will display **End** and it will be automatically turned off after three seconds.

#### **Note**

Some content or menus may differ depending on your smartphone's operating system or SmartLog™ version.

#### **Note**

The SmartLog™ mobile app may not be compatible with all smartphones.

Some smartphones, especially those that are not tested or approved by i-SENS, may be incompatible with your meter. Visit [www.i-sens.com/smartlog](http://www.i-sens.com/smartlog) for more information about supported smartphones. You can also scan the QR code on the back cover of this user manual.

## Adjusting the Date and Time

### Note

- If your meter and smartphone are paired, the meter communicates with the smartphone, automatically synchronizing the current date and time. In this case, continue to press the **S** button to skip the date and time setting and go on to step 9.
- To set the date and time manually, press and hold the **S** button to enter the SET mode, and continue to press the **S** button to go to step 3.

### Step 3 Setting the Year

Press the ◀ or ▶ button until the correct year appears. When the present year appears, press the **S** button to confirm your selection and to go to the next step.

### Step 4 Setting the Month

A number indicating the month will blink on the bottom left corner of the screen. Press the ◀ or ▶ button until the correct month appears. Press the **S** button to confirm your selection and to go to the next step.

### Step 5 Setting the Date

Press the ◀ or ▶ button until the screen displays the correct date. Press the **S** button to confirm the date and to go to the next step.

### Step 6 Setting the Time Format

The meter can be set in the AM/PM 12-hour or the 24-hour clock format.

Press the ◀ or ▶ button to select a format. The AM•PM symbol is not displayed in the 24-hour format. After selecting the format, press the **S** button to go to the next step.

### Step 7 Setting the Hour

Press the ◀ or ▶ button until the correct hour appears.

After the hour is set, press the **S** button to go to the next step.

## Step 8 Setting the Minute

Press the ◀ or ▶ button until the correct minute appears. After setting the minute, press the **S** button to go to the next step.

## Setting the Sound On/OFF

### Step 9

On pressing the ◀ or ▶ button the screen will display 'On' or 'OFF'. Press the **S** button to confirm the selection.

The meter will beep in the following instances if the sound is set to On.

- When you press a button to turn on the meter
- When the test strip is inserted in the meter
- When the blood sample is absorbed into the test strip and the test starts
- When the test result is displayed
- When it is time for a pre-set blood  $\beta$ -Ketone test

If the sound is set to OFF, none of the sound functions will work. After setting the sound, press the **S** button to go to the next step.

### Note


The  symbol is displayed only when the sound is set to OFF.

## Turning on the Strip Expiration Date Indicator

### Step 10

This setting allows you to turn the strip expiration date indicator on or off. This setting turns the function on or off only. See page 27 to set the strip expiration date.

When 'EP' appears on the screen, press the ◀ or ▶ button. The screen will display 'On' or 'OFF'. Press the **S** button to confirm the setting. If you do not want to set the indicator, press the **S** button while the screen displays 'OFF'.

### Note

If the pre-set expiration date is reached, the meter will display 'EP' when the test strip is inserted and when the test result is displayed right after the test.

If the expiration date is set to October of 2023, the meter will display 'EP' at the beginning of October, 2023.

# Setting the Strip Expiration Date Indicator

## Step 1 Entering the Expiration Date Setting

Press and hold the ◀ and ▶ buttons at the same time for 3 seconds to enter the expiration date settings. After all segments flash across the screen, 'EP' will show up.


### Note

The strip expiration date is printed on the test strip box.

## Step 2 Setting the Year

A number indicating the year will blink in the left corner of the screen. Press the ◀ or ▶ button until the correct year appears.

Press the **S** button to confirm the year and set the month.

## Step 3 Setting the Month

A number indicating the month will blink at the bottom of the screen. Press the ◀ or ▶ button until the correct month appears.

Press and hold the **S** button for 3 seconds to confirm the month and turn off the meter.

## Checking the System

You may check your meter and test strips using the KetoSens Control Solutions (level A and B).

The KetoSens Control Solution contains a known amount of  $\beta$ -Ketone and is used to check that the meter and the test strips are working properly.

KetoSens Control Solution ranges are printed on the KetoSens test strip box. Compare the result displayed on the meter to the KetoSens Control Solution range printed on the test strip box.

Before using a new meter or a new box of test strips, you should conduct a control solution test following the procedure on page 29–31 with two different levels of solutions (Control A and B).

### Note

- Use only the KetoSens Control Solution (available for purchase separately).
- Check the expiration date printed on the bottle. When you first open a control solution bottle, record the discard date (date opened plus three (3) months) in the space provided on the label.
- Make sure your meter, test strips, and control solution are at room temperature before testing. Control solution tests must be done at room temperature, 68–77 °F (20–25 °C).
- Before using the control solution, shake the bottle, discard the first 1 or 2 drops and wipe the tip clean.
- Close the control solution bottle tightly and store at a temperature between 46–86 °F (8–30 °C).

## You should perform a control solution test:

- When you want to practice the test procedure using the control solution instead of blood,
- When using the meter for the first time,
- Whenever you open a new box of test strips,
- If the meter or test strips do not function properly,
- If your symptoms are inconsistent with the blood  $\beta$ -Ketone test results and you feel that the meter or test strips are not working properly,
- If you drop or damage the meter.

## Control Solution Testing

### Step 1

Insert a test strip into the meter's test strip port, with the contact bars facing upwards. Gently push the test strip into the port until the meter beeps. Do not force the test strip into the meter.

The ▲ symbol will show up.


### Step 2

Press and hold the ► button for 3 seconds till the 📱 symbol appears on the display.

### Step 3

Shake the KetoSens Control Solution bottle well before each test. Remove the cap and squeeze the bottle to discard the first drop. Then wipe the tip with a clean tissue or cloth. Dispense a drop of control solution onto a clean non-absorbent surface. It helps to squeeze a drop onto the top of the cap as shown. After the ▲ symbol appears on the display, apply the solution to the narrow edge of the test strip until the meter beeps. Make sure the confirmation window fills completely.


### Note

The meter may switch off if the control solution is not applied within 2 minutes of the ▲ symbol appearing on the screen. If the meter turns off, remove the strip, reinsert, and start from step 2.

### Step 4

A test result will appear after the meter counts down from 8 to 1.

The result is stored in the meter's memory but it is not included in the averages.

## Step 5

Compare the result displayed on the meter to the KetoSens Control Solution range printed on the test strip box.

The result should fall within that range.

Follow the instructions provided by your healthcare professional on the safe disposal of used test strips.

### Control Solution Range

Control A: X.X–X.X mmol/L

Control B: X.X–X.X mmol/L

## Caution

The range printed on the test strip box is for the KetoSens Control Solution only. It has nothing to do with your blood  $\beta$ -Ketone level.

## Note

The KetoSens Control Solution can be purchased separately. Contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.

## Comparing the Control Solution Test Results

The test result of each control solution should be within the range printed on the label of the test strip box. Repeat the control solution test if the test result falls outside of this range. Out of range results may occur in the following situations:


Situations	Do This
<ul style="list-style-type: none"><li>• When the control solution bottle was not shaken well,</li><li>• When the meter, test strip, or the control solution were exposed to high or low temperatures,</li><li>• When the first drop of the control solution was not discarded or the tip of the bottle was not wiped clean,</li><li>• When the meter is not functioning properly.</li></ul>	Repeat the control solution test by referring to the <b>Note</b> on page 28.
<ul style="list-style-type: none"><li>• When the control solution is past the expiration date printed on the bottle,</li><li>• When the control solution is past its discard date (the date the bottle was opened plus three (3) months),</li><li>• When the control solution is contaminated.</li></ul>	Discard the used control solution and repeat the test using a new bottle of control solution.

If results continue to fall outside the range printed on the test strip box, the KetoSens test strip and KetoSens BT meter may not be working properly. Do not use your system and contact Customer Service: 1-800-429-5001 (9 am–9 pm EST, Monday–Saturday). *At all other times, contact your healthcare professional.*

## Using the Lancing Device

You will need a lancing device in order to collect a blood sample. You may use the lancing device contained in the KetoSens BT Blood  $\beta$ -Ketone Monitoring System or any other medically approved lancing device.

- The lancing device is for use by a single user only and should not be shared with anyone including other family members.
- Use a soft cloth or tissue to wipe the lancing device. If necessary, a small amount of alcohol on a soft cloth or tissue may be used.

### **⚠ Caution: To reduce the chance of infection**

- All components that come into contact with blood samples should be considered to be biohazards capable of transmitting viral diseases, even after disinfection.
- Do not use a lancet that has been used by another person and never share your meter and lancing device with anyone.
- Wash your hands thoroughly with soap and water after handling your meter, lancing device, or test strips.
- Keep your meter and lancing device clean. See the pre-cleaning and disinfection section on page 46.


### Note

Repeated puncturing at the same sample site may cause pain or skin calluses (thick hard skin). Choose a different site each time you test.

## Preparing the Lancing Device


### Step 1

Wash hands and sample site with soap and warm water. Rinse and dry thoroughly.


### Step 2

Unscrew and remove the adjustable tip.



### Step 3

Insert a new disposable lancet firmly into the lancet holder. Twist off the protective cover of the lancet and set it aside, then replace the adjustable tip. Keep the protective cover to replace on top of the used lancet after testing.

## Step 4


Turn the adjustable tip until it is aligned with the load confirmation window and release button as shown in the diagram.

## Step 5

The lancing device has six puncture depth settings, numbered 0 through 5 (0 for a shallow puncture, 5 for a deeper puncture).

Choose a depth by rotating the top portion of the adjustable tip until the desired number aligns with the arrow.


## Note

0 = least penetration of lancet into the skin.

5 = most penetration of lancet into the skin.

## Step 6

To cock the lancing device, hold the body of lancing device in one hand and pull the sliding barrel with the other hand. The device is loaded when you feel a click and the load confirmation window turns red.

## Note

The skin depth to retrieve samples will vary for various people at different sample sites. The lancing device's adjustable tip allows the best depth of skin penetration for an adequate sample size.

## Preparing the Meter and Test Strip

### Step 7

Insert a test strip with the contact bars facing upwards into the meter's test strip port. Gently insert the test strip until the meter beeps. Do not force the test strip into the meter.


The  symbol will appear on the screen.

## Applying Blood Sample

### Step 8


Obtain a blood sample using the lancing device. Place the device against the pad of the finger. Press the release button. Remove the device from the finger. Wait a few seconds for a blood drop to form. A minimum volume of 0.5 microliter is needed to fill the confirmation window (actual size of 0.5  $\mu\text{L}$  : ●).

## Step 9

After the  symbol appears on the screen, apply the blood sample to the narrow end of the test strip till the meter beeps. If the confirmation window is not filled in time because of abnormal viscosity (thickness and stickiness) or insufficient volume, the Er4 message may appear.

It is recommended to place the test strip vertically into the blood sample site as shown below.

## Caution

Do not allow any foreign substances, such as dirt, blood, or water, enter into the meter. The meter may be damaged or may malfunction. Follow the warning information provided below to prevent possible damage to the meter.

- Do not apply the blood sample directly to the test strip port.
- Do not apply the blood sample to the test strip while holding the meter in a way that the tip of the test strip faces upwards. The blood sample may run down the surface of the test strip and flow into the test strip port.
- Do not store your meter in unsanitary or contaminated sites.


### Note

The meter may switch off if the blood sample is not applied within 2 minutes of the  symbol appearing on the screen. If the meter turns off, remove the strip, reinsert, and start from step 8.

### Step 10

At this time, the display segments will rotate clockwise while the blood is going in.

The test result will appear after the meter counts down from 8 to 1. The result will be automatically stored in the meter's memory. If the test strip is removed after the test result is displayed, the meter will automatically turn off after 3 seconds. Discard used test strips safely in a disposable container.

## Note

To transmit  $\beta$ -Ketone data using the Bluetooth® feature,

- The Bluetooth® feature on the meter must be turned on,
- The meter and a smartphone must be paired,
- The SmartLog™ mobile app must be launched.

The meter will transmit data in the following cases,


- When the strip is ejected after measuring,
- When the meter is turned on (only when untransmitted data exists).

## Discarding Used Lancets



### Step 1

Unscrew the lancing device tip.

### Step 2

Push the lancet ejector forward with the thumb to dispose of the used lancet. Follow the instructions provided by your healthcare professional on the safe disposal of used lancets.

### Note


- Check for damages before using the lancet. If they have been damaged, please discard it and use other lancet.
- The lancet is very sharp. Please keep away from children.
- Keep the lancets in a cool and dry place.
- The lancet is for single use only. Never share or reuse a lancet. Always dispose of lancets properly.

# HI Message


## HI Message

The meter displays results between 0.1–8.0 mmol/L. 'HI' appears when the blood  $\beta$ -Ketone level is greater than 8.0 mmol/L. If 'HI' is displayed again upon retesting, please contact your healthcare professional immediately.

## Lo Message

'Lo' appears when a test result is less than 0.1 mmol/L.

## Note

If you continue to get HI results and think the meter may not be functioning properly, contact Customer Service: **1-800-429-5001 (9 am–9 pm EST, Monday–Saturday)**. At all other times or in case of emergency, contact your healthcare professional or emergency medical response.

# Meter Memory

The KetoSens BT meter can save up to 1,000  $\beta$ -Ketone test results with time and date. If the memory is full, the oldest test result will be deleted and the latest test result will be stored.

## Viewing Test Results Stored in the Meter's Memory

### Step 1

Press any button to turn the meter on.  
After all the segments flash across the screen, the latest test result will be displayed with the current date and time.


### Note

If there are no results saved in memory, only 'KETONE' will be displayed on the screen when you turn on the meter. Refer to the image on the right.

## Step 2

Use the ► button to scroll through the stored test results, starting from the most recent and ending with the oldest. Press the ◀ button to return to the result seen previously.

After checking the stored test results, press the **S** button to turn off the meter.

### Note

The control solution test results saved with the  symbol will be displayed with the  symbol when you review the stored test results.

# Setting the Alarm Function

Three types of alarms can be set in the KetoSens BT meter. The alarms ring for 15 seconds and can be silenced by pressing any button or by inserting a test strip.

## Setting the Time Alarms (alarm 1–3)

### Step 1

Without inserting a test strip, press and hold the ◀ and **S** buttons simultaneously for 3 seconds to enter the time alarm setting. 'alarm 1' will be displayed while 'OFF' blinks on the screen.

### Step 2

On pressing the ▶ button, 'alarm 1' is set and 'On' is displayed on the screen. Press the ▶ button again to cancel 'alarm 1'. 'OFF' will blink on the screen.


### Step 3

Press the ◀ button to adjust the time of alarm 1. A number representing the time will blink on the screen. Press the ▶ button to set the time.


### Step 4

On pressing the ◀ button, the number indicating the minute will start blinking. Press the ▶ button to set the minute.

### Step 5

Press the **S** button to finish and to enter the 'alarm 2' mode.  
Repeat steps 2 to 5 to set the remaining time alarms (alarm 2 and 3).

### Step 6

Press and hold the **S** button for 3 seconds to finish and turn the meter off.

## Caring for Your System

- To minimize the risk of transmission of blood-borne pathogens, the pre-cleaning and disinfection procedure should be performed as recommended in the instructions below.
- Wash your hands thoroughly with soap and water after handling the meter, lancing device, or test strips.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be disinfected prior to use by the second person.

### Pre-cleaning and Disinfection:

Pre-cleaning procedure is needed to clean dirt as well as blood and other body fluids on the exterior of the meter and lancing device before performing the disinfection procedure.

The disinfection procedure is needed to prevent transmission of blood-borne pathogens.

- For the meter and lancing device, this pre-cleaning and disinfection procedure should be performed **once per week**.

#### Note

The life span of a KetoSens BT meter is 5 years. We recommend cleaning and disinfecting both the meter and lancing device at least once per week. We have validated a total of 260 cleaning and disinfecting cycles (260 Cleaning and 260 disinfection cycles) to represent weekly cleaning and disinfecting over the use life of your meter and lancing device.

1 pre-cleaning and 1 disinfection cycles per week \* 52 weeks per year \* 5 years = 260 pre-cleaning and 260 disinfection cycles

- We have validated **Clorox Healthcare Bleach Germicidal Wipes** with 0.55 % sodium hypochlorite as the active ingredient, for disinfecting the KetoSens BT meter and lancing device. It has been shown to be safe for use with the meter and lancing device. This disinfectant is available commercially in towelette form. In addition to KetoSens BT Blood  $\beta$ -Ketone Monitoring System instruction, please read the instructions provided by the manufacturer of Clorox Healthcare Bleach Germicidal Wipes before using it.


<b>Name</b>	<b>Clorox Healthcare® Bleach Germicidal Wipes</b>	
<b>Manufacturer</b>	<b>Clorox® Professional Products Company</b> [Phone] 1 800 537 1415 [Website] <a href="http://www.cloroxprofessional.com">www.cloroxprofessional.com</a>	
<b>EPA registration number</b>	67619-12	
<b>Active ingredients</b>	Sodium Hypochlorite: 0.55 %	

### Note

The disinfectant products can be purchased through online retailers (e.g. Amazon or Walmart.com) or by telephone at the manufacturer. To find out where to purchase the disinfectant product, please contact the manufacturer or visit their website as listed above.

## Pre-cleaning and Disinfection Procedures:

- ① Open the cap of the Clorox Healthcare Bleach Germicidal Wipes container and pull out 1 towelette and close the cap.

- ② Wipe the entire surface of the meter 3 times horizontally and 3 times vertically using one towelette to pre-clean blood and other body fluids.

Front


Back


Both sides

- ③ Dispose of the used towelette in a trash bin.


- ④ Pull out 1 new towelette and wipe the entire surface of the meter 3 times horizontally and 3 times vertically using a new towelette to remove bloodborne pathogens.

- ⑤ Dispose of the used towelette in a trash bin.

- ⑥ Allow exteriors to remain wet for 1 minute.
- ⑦ Repeat the same procedure for the lancing device (step ① to step ⑥).

- \* **After the pre-cleaning and disinfection procedure, the control solution should be tested to confirm that the meter works properly before using the meter. Control solution tests should be performed with two different levels of control solutions (Level A and B). Verify that the test results are within the range printed on the test strip box. See pages 29–31 for how to do a control solution test.**

#### Note

If any of the following deterioration signs appear after pre-cleaning or disinfecting, please stop using the system and contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. *At all other times, contact your healthcare professional.*

- When the inscriptions on the exterior of the meter (or lancing device) have been removed
- When the color of the meter (or lancing device) has changed or faded
- When cracks or roughness develop on the meter (or lancing device)
- When a part of the segment on the meter display does not flash
- When control solution test results do not fall within the stated range on the test strip box.

### **Caution**

- Do not use other cleaners or disinfectants because other chemicals have not been validated and may damage the meter.
- Do not get fluids inside the meter through the test strip port, data transmission port or battery compartment. Never immerse the meter or hold it under running water because this will damage the meter.

### **Caution: Storage and Handling**

- Do not expose the meter to direct sunlight, heat, or excessive humidity for an extended period of time.
- Do not let dirt, dust, blood, or water enter into the meter's test strip port.
- Do not drop the meter or subject it to strong shock.
- Do not try to fix or alter the meter in any way.
- Keep the meter in a cool and airy place.
- Store all meter components in the carrying case to prevent loss.
- Avoid getting any liquid or moisture in the test strip vial or box. This can affect the test strips and cause inaccurate test results.
- Do not apply samples other than capillary whole blood or control solution to the test strip.
- Use of the meter adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the meter and the other equipment should be observed to verify that they are operating normally.

- Use of accessories, transducers and cables other than those specified or provided by i-SENS could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF Communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the KetoSens BT meter, including cables specified by i-SENS.

 **Note**

You can get additional information or technical assistance by calling our Customer Service Center: 1-800-429- 5001 (9 am–9 pm EST, Monday–Saturday). *At all other times, contact your healthcare professional.*

## Understanding Error Messages

	<p>A used test strip was inserted.</p> <p>&gt; <i>Repeat the test with a new test strip.</i></p>
	<p>The blood or control solution sample was applied before the  symbol appeared.</p> <p>&gt; <i>Repeat the test with a new test strip and wait until the  symbol appears before applying the blood or control solution sample.</i></p>
	<p>The temperature during the test was above or below the operating range.</p> <p>&gt; <i>Move to an area where the temperature is within the operating range 50–104 °F (10–40 °C) and repeat the test after the meter and test strips have reached a temperature within the operating range.</i></p>
	<p>The blood sample is of insufficient volume or blood is too thick.</p> <p>&gt; <i>Repeat the test using a new test strip.</i></p>
	<p>This error message may appear when the wrong blood <math>\beta</math>-Ketone test strip is used instead of the KetoSens Blood <math>\beta</math>-Ketone test strip.</p> <p>&gt; <i>Repeat the test after inserting a KetoSens test strip.</i></p>

There is a problem with the meter.

- > *Do not use the meter. Contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. At all other times, contact your healthcare professional.*

There is a problem with Bluetooth® communication.

- > *Contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. At all other times, contact your healthcare professional.*

An electronic error occurred during the test.

- > *Repeat the test with a new test strip. If the error message persists, contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.*

 **Note**

If the error messages persist, contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. At all other times or in case of emergency, contact your healthcare professional or emergency medical response.

## General Troubleshooting

Problem	Troubleshooting
The display is blank even after inserting a test strip.	<ul style="list-style-type: none"><li>• Check whether the test strip is inserted with the contact bars facing up. Check if the strip has been inserted completely into the test strip port.</li><li>• Check if the appropriate test strip was used. Check whether the batteries are inserted with the + side facing up.</li><li>• Replace the batteries.</li></ul>
The test does not start even after applying the blood sample on the test strip.	<ul style="list-style-type: none"><li>• Check if the confirmation window is filled completely.</li><li>• Repeat the test after inserting a new test strip.</li></ul>
The test result does not match the way you feel.	<ul style="list-style-type: none"><li>• Repeat the test after inserting a new test strip.</li><li>• Check the expiration date of the test strip.</li><li>• Perform control solution test.</li></ul>

### Note

If the problem is not resolved, please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. At all other times or in case of emergency, contact your healthcare professional or emergency medical response.

## Performance Characteristics

The performance of KetoSens BT Blood  $\beta$ -Ketone Monitoring System has been evaluated in laboratory and in clinical tests.

**Accuracy:** The accuracy of the KetoSens BT Blood  $\beta$ -Ketone Monitoring System was assessed by comparing blood  $\beta$ -Ketone results obtained by patients with those obtained using a Randox Monaco, a laboratory instrument. The following results were obtained by diabetic patients at clinic centers.

### Capillary Blood Sample

System accuracy results for capillary blood ketone

Slope	0.994
Y-intercept	-0.007
Correlation coefficient (r)	0.970
Number of Subjects	185
Range tested	0.017–1.054 mmol/L

System accuracy results for capillary blood ketone

Within $\pm 0.15$ mmol/L	Within $\pm 0.225$ mmol/L	Within $\pm 0.3$ mmol/L
183/185 (98.9 %)	185/185 (100 %)	185/185 (100 %)

**Precision:** Precision studies were performed in a laboratory using the KetoSens Brand of Blood  $\beta$ -Ketone Monitoring Systems.

<i>Within Run Precision</i>		
Blood average	0.4 mmol/L	SD = 0.027 mmol/L
	1.1 mmol/L	SD = 0.056 mmol/L

Blood average	3.4 mmol/L	CV = 3.7 %
	5.2 mmol/L	CV = 3.5 %
	6.9 mmol/L	CV = 3.6 %

<i>Between Run Precision</i>		
Control average	0.6 mmol/L	SD = 0.049 mmol/L
	2.2 mmol/L	CV = 3.8 %
	4.0 mmol/L	CV = 4.3 %

This study shows that there could be a variation of up to 4.3 %.

### Reference:

1. A. Rewers, Current Controversies in Treatment and Prevention of Diabetic Ketoacidosis, *Advances in Pediatrics* 57 (2010) 247–267.
2. CLSI EP5-A3(Evaluation of Precision Performance of Quantitative Measurement Methods;Approved Guideline-Third Edition)
3. CLSI EP6-A (Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical A Approach; Approved)
4. CLSI EP9-A3(Measurement Procedure Comparison and Bias Estimation Using Patient Samples;Approved Guideline-Third Edition.
5. CLSI EP30-A Characterization and Qualification of Commutable Reference Materials for Laboratory Medicine; Approved Guideline.
6. CLSI EP7-A2, 2005 (Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition

## EMC Table

The following tables contain the Manufacturer's declaration and additional information required by IEC 60601-1-2:2014 (Fourth Edition).

Phenomenon	Basic EMC standard or test method	Port tested	Test level/ requirement
Radiated disturbance	CISPR 11: 2015 +A1: 2016	Enclosure	Group 1, Class B
Electrostatic Discharges (ESD)	IEC 61000-4-2:2008	Enclosure	±8 kV Contact ±2 kV, ±4 kV, ±8 kV, ±15 kV Air
Radiated RF Electromagnetic Fields	IEC 61000-4-3:2006 + A1:2007+A2:2010	Enclosure	3 V/m 80 MHz to 2.7 GHz 80%, 1 Hz AM RF Wireless Comm. (Refer to test report clause 1.15)
Electric Fast Transients and bursts	IEC 61000-4-4:2012	AC Mains	±2 kV AC, 100 kHz PRR
Surges	IEC 61000-4-5:2014	AC Mains	±0.5 kV, ±1 kV L1 to L2 (DM)
Conducted Disturbances, induced by RF fields	IEC 61000-4-6:2013	AC Mains	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz
Power-frequency Magnetic Field	IEC 61000-4-8:2009	Enclosure	30 A/m

## Warranty Information

---

### Manufacturer's Warranty

i-SENS, Inc. warrants that the KetoSens BT Meter shall be free of defects in material and workmanship in normal use for a period of five (5) years. The meter must have been subjected to normal use. The warranty does not cover improper handling, tampering, use, or service of the meter. Any claim must be made within the warranty period.

i-SENS, Inc. will, at its discretion, repair or replace a defective meter or meter part that is covered by this warranty. As a matter of warranty policy, i-SENS will not reimburse the consumer's purchase price.

### Obtaining Warranty Service

To obtain warranty service, you must return the defective meter or meter part along with proof of purchase to your nearest i-SENS sales or customer service representative.

## MEMO

## MEMO

## MEMO

## MEMO

# MEMO

## KetoSens® BT


### $\beta$ -Ketone Measurement for Type 1 Diabetic Users

<b>Lo-0.5 mmol/L</b> Re-test for blood ketones only if the blood glucose levels are higher than 120 mg/dL. Revert to routine blood glucose monitoring.	<b>1.5-3.0 mmol/L</b> Refer to medical staff for assessment for DKA*. Re-test blood glucose and ketone 2 hourly.
<b>0.5-1.5 mmol/L</b> Re-test for blood glucose and ketone 2 hourly. Report to medical staff if the ketone levels do not fall as the insulin dosage may need reviewing. Continue to monitor 4 hourly until the patient is back on a subcutaneous insulin regime. Resume testing for blood ketone level if blood glucose level is above 120 mg/dL.	<b>3.0 mmol/L and higher</b> Monitor blood ketone and glucose hourly and refer to the medical team for assessment for DKA*. It is expected that blood ketone level will reduce by 0.4 mmol/L per hour. If this is not being achieved, ask a doctor to review management.

## MEMO

**KetoSens® BT**

### Keto Diet Ketone Blood Level Ranges

\* May require medical attention.  
Please check with your medical provider.

References: <http://www.diabetes.co.uk/keto/measuring-ketosis-on-a-keto-diet.html>  
<http://litfl.com/ketoacidosis/>  
<http://www.nhs.uk/conditions/diabetic-ketoacidosis/>

Please check the back of this card for important information about ketone levels.

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The image shows the error code 'Er-5' displayed on a digital screen with a light gray background.

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The image shows the error code 'Er-6' displayed on a digital screen with a light gray background.

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The image shows the error code 'Er-7' displayed on a digital screen with a light gray background.

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The image shows the error code 'Er-8' displayed on a digital screen with a light gray background.

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### **Referencia:**

1. A. Rewers, Current Controversies in Treatment and Prevention of Diabetic Ketoacidosis, *Advances in Pediatrics* 57 (2010) 247–267.
2. CLSI EP5-A3(Evaluation of Precision Performance of Quantitative Measurement Methods;Approved Guideline-Third Edition)
3. CLSI EP6-A (Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical A Approach; Approved)

- 
4. CLSI EP9-A3(Measurement Procedure Comparison and Bias Estimation Using Patient Samples;Approved Guideline-Third Edition.
  5. CLSI EP30-A Characterization and Qualification of Commutable Reference Materials for Laboratory Medicine; Approved Guideline.
  6. CLSI EP7-A2, 2005 (Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition

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