# VETMATE™

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Auto Backlight LCD

PGUAA-0000064 REV0 2019/10

## Welcome to the VetMate™ Blood Glucose Monitoring System

Thank you for choosing the VetMate™ Blood Glucose Monitoring System. The system provides you with safe, fast, and convenient blood glucose *in vitro* (i.e., outside the body) diagnostic monitoring. You can obtain accurate results of your dog/cat in just 5 seconds with a small (0.4 µL) blood sample.

- No part of this document may be reproduced in any form or by any means without the prior written consent of i-SENS.
- The information in this manual is correct at the time of printing. However, i-SENS reserves the right to make any necessary changes at any time without notice as our policy is one of continuous improvement.

*Manufacturer:* i-SENS, Inc.  
43, Banpo-daero 28-gil, Seocho-gu,  
Seoul 06646, Korea  
[www.i-sens.com](http://www.i-sens.com)

Contact Customer Service:

**1-800-429-5001**

9 am–9 pm EST, Monday–Saturday  
At all other times or in case of emergency,  
contact your veterinarian/veterinary  
surgeon or emergency veterinary response.

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## Preparation

Setting Up the System .....	14
Adjusting the Date and Time .....	14
Setting the Sound On/OFF .....	16
Setting VET Mode (Animal Mode).....	17
Turning on the Hypoglycemia (HYPO) Indicator .....	18
Setting the 'Test Result Reset' (Deleting the saved test results).....	19
Checking the System.....	20
Control Solution Testing.....	21
Comparing the Control Solution Test Results.....	24

## Testing

Using the Lancing Device.....	25
Preparing the Lancing Device .....	26
Preparing the Meter and Test Strip.....	27
Applying Blood Sample .....	28
Discarding Used Lancets .....	31
HI and Lo Messages .....	32
Target Blood Glucose Ranges.....	33
Transferring Test Results.....	34

## Additional Functions

Meter Memory .....	35
Viewing Averages Stored in Memory .....	35
Viewing Test Results.....	36
Setting the Alarm Function.....	37
Setting the Time Alarms (Alarm 1–5) .....	37

## Maintenance

Understanding Error Messages .....	39
General Troubleshooting .....	42
Performance Characteristics.....	43
EMC Table .....	46
Warranty Information.....	48

## Important Information: Read This First!

For optimum safety and benefits, please read the entire manual contents before using the system.

### Intended Use:

VetMate Blood Glucose Monitoring System is used for use with VetMate Blood Glucose Meters to quantitatively measure glucose in capillary and venous whole blood samples drawn by veterinarians or veterinary surgeons for the management of diabetes in cats and dogs.

- The VetMate Blood Glucose Monitoring System is intended for use in monitoring glucose samples taken from the dog or cat.
- Glucose in blood samples reacts with the chemical in the test strip to produce a small electrical current. The VetMate Meter detects this electrical current and measures the amount of glucose in the blood sample.
- The VetMate Blood Glucose Meter is designed to minimize code related errors in monitoring by using the no-coding function.
- The VetMate Blood Glucose Meter should be used only with the VetMate Blood Glucose Test Strips.
- An abnormally high or low red blood cell count (hematocrit level over 65 % or below 15 %) may produce inaccurate results.

*If you need assistance, please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. At all other times or in case of emergency, contact your veterinarian/veterinary surgeon.*

## Specifications

### Product specifications

Measurement range	20–600 mg/dL
Sample size	Minimum 0.4 µL
Test time	5 seconds
Sample type	Fresh capillary whole blood, Fresh venous whole blood
Calibration	Plasma-equivalent
Assay method	Electrochemical
Battery life	3,000 tests
Power	Two 3.0 V lithium batteries (disposable, type CR2032)
Memory	500 test results
Size	92.4 x 51.4 x 21.3 (mm)
Weight	55 g (with batteries)

### Operating ranges

Temperature	41–113 °F (5–45 °C)
Relative humidity	10–90 %
Hematocrit	15–65 %

### Storage Conditions

Glucose meter (with batteries)	32–122 °F (0–50 °C)
Test strip	34–86 °F (1–30 °C)

## VetMate Blood Glucose Monitoring System

### VetMate BGM System includes the following items:

- \* VetMate Blood Glucose Meter
- \* User's Manual
- \* Quick Reference Guide
- \* Batteries

### VetMate BGM System may include the following items:

- \* VetMate Blood Glucose Test Strips
- \* Lancets
- \* Lancing Device
- \* Logbook
- \* Carrying Case


- The VetMate Blood Glucose Control Solution Level A and B are sold separately and can be purchased by contacting Customer Service.
- Check all the components after opening the VetMate Blood Glucose Monitoring System package. The exact contents are listed on the main box.
- The USB cable for data management software can be ordered separately. Please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. *At all other times, contact your veterinarian/veterinary surgeon.*

## Replacing the Batteries

The VetMate meter comes with two 3.0 V lithium batteries. Before using the meter, check the battery compartment and insert batteries if empty. When the **+ -** symbol appears on the display while the meter is in use, the batteries should be replaced as soon as possible. The test results may not be saved if the batteries run out.



### Step 1

Make sure the meter is turned off. Push the cover in the direction of the arrow to open the battery compartment.


### Step 2

Remove the used batteries one at a time. Slip your index finger under the battery to lift and pull out as shown. Insert two new batteries with the + side facing up and make sure the batteries are inserted firmly.

### Step 3

Place the cover on the battery compartment. Push it down until you hear the tab click into place.

**Note:** Removing the meter batteries will not affect the stored results. However, you may need to reset the meter settings. See pages 14–19.

## Caring for the System

Use a soft cloth or tissue to wipe the meter exterior. If necessary, dip the soft cloth or tissue in a small amount of alcohol.

Do not use organic solvents such as benzene, acetone, or any household and industrial cleaners that may cause irreparable damage to the meter.

### Caution:

- Do not expose the meter to direct sunlight, heat, or excessive humidity for an extended period of time.
- Do not let dirt, dust, blood, or water enter into the meter's test strip port.
- Do not drop the meter or submit it to strong shock.
- Do not try to fix or alter the meter in any way.
- Strong electromagnetic radiation may interfere with the proper operation of this device. Keep the device away from sources of strong electromagnetic radiation, especially when measuring your pet's blood glucose.
- Keep the meter in a cool and well ventilated place.
- Store all the meter components in the carrying case to prevent loss and help keep the meter clean.
- Use of the meter adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the meter and the other equipment should be observed to verify that they are operating normally.


- Use of accessories, transducers and cables other than those specified or provided by i-SENS could result in increased electromagnetic emissions or decreased electromagnetic immunity of the meter and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the meter, including cables specified by i-SENS. Otherwise, degradation of the performance of the meter could result.

### Disposal of the meter

If you need to throw your meter away, you should follow existing policies and procedures of your own country or region. For information about correct disposal, please contact your local council or authority. If you need assistance, contact Customer Service or visit [www.i-sens.com](http://www.i-sens.com).

## VetMate Blood Glucose Test Strip

The VetMate Blood Glucose Monitoring System measures blood glucose quickly and accurately. It automatically absorbs the small blood sample applied to the narrow edge of the strip.

### Warning!

- Apply fresh whole capillary blood to the test strip. Fresh venous whole blood may also be used only if drawn by veterinarian/veterinary surgeon.
- Fresh venous whole blood specimens containing the anticoagulants EDTA and Heparin is acceptable. Iodoacetate or fluoride/oxalate are not recommended.
- Do not reuse test strips.
- Do not use test strips past the expiration or discard date.
- Test strips in new, unopened vials and test strips in vials that have been opened can be used up until the expiration date printed on the test strip box and vial label if the test strips are used according to its storage and handling methods.
- Store test strips in a cool and dry place at a temperature between 34–86 °F (1–30 °C) and 10–90 % Relative Humidity.
- Keep test strips away from direct sunlight or heat and do not freeze.
- Store test strips only in their original vial.
- Close the vial tightly after taking out a test strip for testing and use the strip immediately.
- Handle test strips only with clean and dry hands.
- Do not bend, cut, or alter test strips in any way.
- For detailed storage and usage information, refer to the VetMate test strip package insert.

### Caution:

- Keep the meter and testing supplies away from young children or pets.
- Drying agents in the vial cap may be harmful if inhaled or swallowed and may cause skin or eye irritation.

## VetMate Blood Glucose Meter

### Data Port

Used to transfer data from the meter to a computer with a cable

### ▲ Button

Selects or changes information

### S Button

Turns the meter on/off and confirms menu selections

### Display

Shows results, messages

### Ejection Button

Slide down to discard the used test strips

### ▼ Button

Selects or changes information

### Test Strip Port

Insert test strip here

### Note:

- The cable for data management software can be ordered separately. Please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.
- The unit of measurement is fixed and it cannot be changed by the user.

## VetMate Blood Glucose Meter Display

### Display Segment Check

It is recommended to check if the display screen on the meter matches the illustration below every time the meter turns on. Do not use the meter if the display screen does not exactly match the illustration as the meter may show incorrect results.

### mem

Appears when test results stored in the memory are displayed

### alarm


Appears when the time alarm has been set

### Hypoglycemia symbol

Appears when the test result is below the pre-set hypoglycemic level

### Mute symbol

Appears only when the sound is set to OFF

### Battery symbol

Indicates meter battery is running low and needs to be replaced

### check

Appears when the control solution test results are saved or displayed

### Test results

### Blood insertion symbol

Indicates the meter is ready for the application of a drop of blood or control solution

### Animal symbol

Shows the set animal (dog/cat) mode, Can be changed in SET mode or immediately prior to blood insertion

## Setting Up the System

Press and hold the **S** button for 3 seconds to enter SET mode. After all settings are finished, press and hold the **S** button for 3 seconds to turn off the meter. Press the **▲** or **▼** button to reach the accurate value. Press and hold the **▲** or **▼** button to scroll faster.

### Adjusting the Date and Time


#### Step 1 Entering SET Mode

Press and hold the **S** button for 3 seconds. After all the segments flash across the screen, 'SET' will be displayed on the screen. Press the **S** button again to go to the next step.

#### Step 2 Setting the Year

Press the **▲** or **▼** button to adjust until the correct year appears. After setting the year, press the **S** button to confirm your selection and go to the next step.

#### Step 3 Setting the Month

A number indicating the month will blink on the bottom left corner of the screen. Press the **▲** or **▼** button until the correct month appears. Press the **S** button to confirm your selection and go to the next step.


#### Step 4 Setting the Date

Press the **▲** or **▼** button until the screen displays the correct date. Press the **S** button to confirm the date and go to the next step.


#### Step 5 Setting the Time Format

The meter can be set in the AM/PM 12-hour or the 24-hour clock format. Press the **▲** or **▼** button to select a format. The AM•PM symbol is not displayed in the 24-hour format. After selecting the format, press the **S** button to go to the next step.


## Step 6 Setting the Hour

Press the ▲ or ▼ button until the correct hour appears. After the hour is set, press the **S** button to go to the next step.

## Step 7 Setting the Minute

Press the ▲ or ▼ button until the correct minute appears. After setting the minute, press the **S** button to go to the next step.

## Setting the Sound On/OFF

### Step 8

On pressing the ▲ or ▼ button, the screen will display On or OFF. Press the **S** button to confirm the selection.

The meter will beep in the following instances if set to On:

- When the test strip is inserted in the meter,

- When the blood or control solution sample is absorbed into the test strip and the test starts,
- When the test result is displayed,
- When you press the **S** or ▲ button to turn on the meter,
- When it is time for a preset blood glucose test.



If the sound is set to OFF, none of the sound functions will work. After setting the sound, press the **S** button to go to the next step.

**Note:** The  symbol is displayed only when the sound is set to OFF.

## Setting VET Mode (Animal Mode)

### Step 9

In this setting, VET mode can be selected; dogs or cats. On pressing the ▲ or ▼ button, the screen will display the  or  symbol. Press the **S** button to confirm the setting.


## Turning on the Hypoglycemia (HYPO) Indicator

### Step 10

In this setting, the hypoglycemia (possible low blood sugar) level can be selected.

'HYPO' will be displayed with the  symbol. On pressing the  or  button, the screen will display 'On' or 'OFF'. Press the **S** button when 'On' appears to set the hypoglycemia level. Press the  or  button until the desired hypoglycemia level between 20 and 90 mg/dL appears. Press the **S** button to confirm the level and go to the next step.



If you do not want to set the indicator, press the **S** button while the screen displays 'OFF'.


**Caution:** Before setting the hypoglycemia level, ask your veterinarian/veterinary surgeon to help you decide what the hypoglycemia level is.

## Setting the 'Test Result Reset' (Deleting the saved test results)


### Step 11

In this setting, the test results saved in the meter can be deleted. Please note that if you select 'YES', the saved test results will be deleted and cannot be restored.

'dEL' and 'no1' will be displayed on the screen. Press  to switch between 'YES' or 'no'. Press  to select a desired dog (no 1-4) or cat (no 1-4) to delete the saved test results.

To delete all the saved test results, press the  button until the  and  symbols, and 'ALL' appear on the screen. Press the  button until 'YES' appears on the screen. When 'YES' blinks, press the **S** button. All the saved test results will be deleted from the meter. Refer to the process shown in the figure below.

If you do not want to delete the results, press the **S** button while the screen displays 'no'. Then, the screen will return to step 2. See page 14.

# Checking the System

You may check the meter and test strips using the VetMate Control Solution. The VetMate Control Solution contains a known amount of glucose and is used to check that the meter and the test strips are working properly. The test strip vials have VetMate Control Solution ranges printed on their labels. Compare the result displayed on the meter to the VetMate Control Solution range printed on the test strip vial. Before using a new meter or a new vial of test strips, you may conduct a control solution test following the procedure on pages 21–23.

## Notes:

- Use only the VetMate Control Solutions.
- Check the expiration date printed on the bottle. When you first open a control solution bottle, record the discard date (date opened plus 3 months) in the space provided on the label.
- Make sure the meter, test strips, and control solution are at room temperature before testing. Control Solution tests must be done at room temperature (68–77 °F/20–25 °C).
- Before using the control solution, shake the bottle. discard the first few drops and wipe the tip clean.
- Close the control solution vial tightly and store at a temperature between 46–86 °F (8–30 °C).

## You may do a control solution test:


- When you want to practice the test procedure using the control solution instead of blood,
- When using the meter for the first time,
- Whenever you open a new vial of test strips,
- If the meter or test strips do not function properly,
- If your pet's behaviors are inconsistent with the blood glucose test results and you feel that the meter or test strips are not working properly,
- If you drop or damage the meter.

**Caution:** If all the results you get on testing three times continuously are within the range printed on the test strip vial, the meter and test strips are working properly and you may use them for the blood glucose test.

## Control Solution Testing

### Step 1

Insert a test strip into the meter's test strip port, with the contact bars facing upwards. Gently push the test strip into the port until the meter beeps. Be careful not to bend the strip while pushing it in. The ▲ symbol will be displayed on the screen.


When pressing and holding the ▲ button for 3 seconds, the control solution testing will be ready.

## Step 2

Shake the VetMate Control Solution bottle before each test. Remove the cap and squeeze the bottle to discard the first drop. Then wipe the tip with a clean tissue or cloth. After the ▲ symbol appears on the display, apply the solution to the narrow edge of the test strip until the meter beeps. Make sure the confirmation window fills completely.

**Note:** The meter may switch off if the control solution sample is not applied within 2 minutes of the ▲ symbol appearing on the screen. If the meter turns off, remove the strip, reinsert, and start from step 1.

## Step 3

A test result will appear after the meter counts down from 5 to 1. When 'check' is displayed, the result is stored in the meter's memory but it is not included in the averages.

## Step 4

Compare the result displayed on the meter to the range printed on the test strip vial. The result should fall within the range. Used strips should be discarded safely in disposable containers.

Control Solution Range		GDH-FAD
Dog	Cat	
Control A: 139~209	140~210	mg/dL
Control B: 245~368	236~354	mg/dL

### Caution:

- If the results are not within the range on the test strip vial, the meter and test strip may not work properly. Then, stop using the meter and contact Customer Service: 1-800-429-5001 Mon-Sat, 9 am-9 pm EST.
- The range printed on the test strip vial is for the VetMate Control Solution only. It does not have any connection to the blood glucose level of your dog/cat.

## Using the Lancing Device

### Comparing the Control Solution Test Results

The test result of each control solution should be within the range printed on the label of the test strip vial. Repeat the control solution test if the test result falls outside of the range. Out of range results may occur in the following situations:


Situations	Do This
<ul style="list-style-type: none"><li>• When the control solution bottle was not shaken well,</li><li>• When the meter, test strip, or the control solution were exposed to high or low temperatures,</li><li>• When the first drop of the control solution was not discarded or the tip of the bottle was not wiped clean,</li></ul>	Repeat the control solution test by referring to the 'Notes' on page 20.
<ul style="list-style-type: none"><li>• When the control solution is past the expiration date printed on the bottle,</li><li>• When the control solution is past its discard date (the date the bottle was opened plus three (3) months),</li><li>• When the control solution is contaminated.</li></ul>	Discard the used control solution and repeat the test using a new bottle of control solution.

If results continue to fall outside the range printed on the test strip vial, the VetMate Test Strip and Meter may not be working properly. Do not use your system and contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. *At all other times, contact your veterinarian/veterinary surgeon.*

**Note:** The VetMate Control Solution can be purchased separately. Please contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. *At all other times, contact your veterinarian/veterinary surgeon.*

You will need a lancing device in order to collect a blood sample.

You may use the lancing device included in the VetMate Blood Glucose Monitoring System or any other medically approved lancing device.

- Use a soft cloth or tissue to wipe the lancing device. If necessary, a small amount of alcohol on a soft cloth or tissue may be used.


**Caution:** To avoid infection when drawing a sample, do not use a lancet more than once, and:

- Do not use a lancet that has been used by other dogs/cats.
- Always use a new sterile lancet.
- Keep the lancing device clean.

## Preparing the Lancing Device


### Step 1

Wash your hands and sample site with soap and warm water. Rinse and dry thoroughly.


### Step 2

Unscrew and remove the lancing device tip.


### Step 3

Firmly insert a new lancet into the lancet holder. Hold the lancet firmly. Gently twist to pull off protective disk. Save disk to recap lancet after use. Replace lancing device tip.


### Step 4

Turn the adjustable tip until it is aligned with the load confirmation window and release button as shown in the diagram.


### Step 5

The lancing device has six puncture depth settings (0 for a shallow puncture, 5 for a deeper puncture). Choose a depth by rotating the top portion of the adjustable tip until the desired number aligns with the arrow. A beginning setting of three (3) is recommended.

### Step 6


To cock the lancing device, hold the body of lancing device in one hand and pull the sliding barrel with the other hand. The device is loaded when you feel a click and the load confirmation window turns red.

## Preparing the Meter and Test Strip


### Step 7

Insert a test strip with the contact bars facing upwards into the meter's test strip port. Push the strip in gently until the meter beeps. Be careful not to bend the test strip. The  symbol will appear on the screen.

## Step 8

In this step, VET mode can be selected; dogs or cats. You can change the pet type between dogs and cats by pressing and holding the ▼ button for 3 seconds. Press the ▲ or ▼ button to select a number for your dog or cat. Test results of up to 4 different dogs and cats can be distinguished by selecting numbers between 1–4 in this step.

## Applying Blood Sample

### Step 9


Obtain a blood sample using the lancing device. Place the device against the test site: the marginal ear vein in both cats and dogs; paw pad in cats and dogs; leg callus, usually dogs; inner or outer lip, dogs only. Press the release button. Remove the device from the sample site. Wait a few seconds for a blood drop to form. A minimum volume of 0.4 microliter is needed to fill the confirmation window (actual size of 0.4  $\mu\text{L}$ : ●).

**Note:** Petroleum jelly may be applied prior to lancing the site. This can help the blood to form into a droplet and avoid spilling into the hair.

### Step 10

After the ▲ symbol appears on the screen, apply the blood sample to the narrow end of the test strip till the meter beeps. If the confirmation window is not filled in time because of abnormal viscosity (thickness and stickiness) or insufficient volume, the Er4 message may appear. It is recommended to place the test strip vertically into the blood sample site as shown below.

Good Sample

Insufficient Sample

**Note:** The meter may switch off if the blood sample is not applied within 2 minutes of the ▲ symbol appearing on the screen. If the meter turns off, remove the strip and reinsert it, and apply the blood sample after the ▲ symbol appears on the screen.

### Step 11

The test result will appear after the meter counts down from 5 to 1. The result will be automatically stored in the meter's memory.


Gently apply pressure to the test site with sterile gauze or cotton wool to help stop the blood flow.

### Step 12

Slide the ejection button to discard the used strips safely in disposable container.


If the test strip is removed after the test result is displayed, the meter will automatically switch off after 3 seconds.

## Discarding Used Lancets



### Step 1

Unscrew the lancing device tip.

### Step 2

Stick the lancet into the saved protective disk. Push the lancet ejector forward with the thumb and simultaneously pull out the sliding barrel to dispose of the used lancet in a proper biohazard container.


**Caution:** The lancet is for single use only. Never share or reuse a lancet. Always dispose of lancets properly.

## HI and Lo Messages

### HI Message


The meter displays results between 20–600 mg/dL. 'HI' appears when the blood glucose level is greater than 600 mg/dL and indicates severe hyperglycemia (much higher than normal glucose levels).

If 'HI' is displayed again upon retesting, please contact your veterinarian/veterinary surgeon immediately.

### Lo Message

'Lo' appears when a test result is less than 20 mg/dL and indicates severe hypoglycemia (very low glucose levels). If 'Lo' is displayed again upon retesting, please contact your veterinarian/veterinary surgeon immediately.

**Note:** If the error messages persist even though your pet does not have hyperglycemia or hypoglycemia, please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am –9 pm EST. At all other times or in case of emergency, contact your veterinarian/veterinary surgeon or emergency veterinary response.

## Target Blood Glucose Ranges

Consult a veterinarian/veterinary surgeon for the target glucose values for each animal.

Reminders	Target ranges from veterinarian/ veterinary surgeon
Time of day	

**Unexpected Values :** Low or high blood glucose readings can indicate a potentially serious medical condition. If your pet's result is unusually high or low, or do not match the way your pet behaves, repeat the test with a new test strip. If the reading is inconsistent with the symptoms or if the blood glucose result is less than 60 mg/dL for dogs and less than 70 mg/dL for cats<sup>1)</sup> or higher than 180 mg/dL for dogs and higher than 200 mg/dL for cats<sup>2)</sup>, contact your veterinarian/veterinary surgeon.

### Reference

1. Edward F, Nelson RW, Reusch C, Scott-Moncrieff JC Canine & Feline Endocrinology. 4th Ed. St. Louis (MO): Saunders; 2015.
2. Rucinsky R, et al. AAJA Diabetes Management Guidelines for Dogs and Cats. J Am Anim Hosp Assoc 2010;46:215-224.

## Transferring Test Results

Test results stored on VetMate meter can be transferred from the meter to a computer using VetLog software and cable. The meter screen displays 'Pc' when it is connected to the computer using the data cable.

You are responsible for properly securing and managing your PC. If you suspect an adverse cybersecurity event related to VetMate meter, contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.

## Meter Memory

The VetMate meter can save up to 500 glucose test results with time and date. If the memory is full, the oldest test result will be deleted and the latest test result will be stored. The meter calculates and displays the averages of total test results from the last 14 days.

### Viewing Averages Stored in Memory


#### Step 1

Press the ▲, ▼ or S button to turn the meter on. The current date and time will be displayed at the bottom of the screen for 2 seconds. A 14-day average value for your dog assigned to 'no 1' will be followed.

#### Step 2

Press the ▲ button to scroll through the 14-day average values of pets assigned to each number. Press the ▼ button to scroll back through the averages see previously.

## Setting the Alarm Function

Five time-set alarms (alarm 1–5) can be set in the VetMate meter:

### Viewing Test Results


#### Step 1

Press the **▲** or **S** button to turn the meter on. The current date and time will be displayed on the bottom of the screen for 2 seconds. A 14-day average value for your dog assigned to 'no 1' will be followed.

#### Step 2


Use the **▼** button to scroll through the test results, starting from the most recent and ending with the oldest. Press the **▲** button to return to the result seen previously. The test results will be displayed with the dates and time following selected numbers (no 1–no 4). After checking the stored test results, press the **S** button to turn off the meter.

### Setting the Time Alarms (Alarm 1–5)

#### Step 1

Without inserting a test strip, press the **▲** and **S** buttons simultaneously for 3 seconds to enter time alarm mode. 'alarm 1' will be displayed while 'OFF' blinks on the screen.

#### Step 2

On pressing the **▼** button, 'alarm 1' is set and 'On' is displayed on the screen. Press the **▼** button again to cancel 'alarm 1' and 'OFF' will blink on the screen.

#### Step 3

Press the **▲** button to adjust the time of 'alarm 1'. A number representing the time will blink on the screen. Press the **▼** button to set the time.

## Understanding Error Messages

### Step 4

On pressing the ▲ button, the number indicating the minute will start blinking. Press the ▼ button to set the minute.

### Step 5

Press the S button to finish and to go to 'alarm 2' mode. Repeat steps 2 to 5 to set the remaining time alarms (alarm 2–5).

### Step 6

Press the S button for 3 seconds to finish and turn the meter off.

Message	What It Means	What To Do
	A used test strip was inserted.	Repeat the test with a new test strip.
	The blood or control solution sample was applied before the ▲ symbol appeared.	Repeat the test with a new test strip and wait until the ▲ symbol appears before applying the blood or control solution sample.
	The temperature during the test was out of the operating range.	Move to an area where the temperature is within the operating range (41–113 °F/ 5–45 °C) and repeat the test after the meter and test strips have reached a temperature within the operating range.

Message	What It Means	What To Do
	The blood sample has abnormally high viscosity or insufficient volume	Repeat the test with a new test strip.
	A non-VetMate test strip was used. Test strip was not inserted properly.	Repeat the test with a VetMate test strip. Gently push the test strip into the test strip port with the contact bars facing upwards until the meter beeps.
	There is a problem with the meter.	Do not use the meter. Contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.

Message	What It Means	What To Do
	An electronic error occurred during the test.	Repeat the test with a new test strip. If the error message persists, contact Customer Service at 1-800-429-5001 Mon–Sat, 9 am–9 pm EST.

**Note:** If the error messages persist, contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. *At all other times or in case of emergency, contact your veterinarian/veterinary surgeon or emergency veterinary response.*

## General Troubleshooting

Problem	Troubleshooting
The display is blank even after inserting a test strip.	<ul style="list-style-type: none"> <li>• Check whether the test strip is inserted with the contact bars facing upwards. Check if the strip has been inserted completely into the test port.</li> <li>• Check whether the batteries are inserted with the + side facing upwards.</li> <li>• Replace the batteries.</li> </ul>
The test does not start even after applying the blood sample on the strip.	<ul style="list-style-type: none"> <li>• Check if the confirmation window is filled completely.</li> <li>• Repeat the test with a new test strip.</li> </ul>
The test result does not match your expectation.	<ul style="list-style-type: none"> <li>• Repeat the test with a new test strip.</li> <li>• Check the expiration or discard date of the test strip.</li> <li>• Perform control solution test.</li> </ul>

**Note:** If the problem is not resolved, please contact Customer Service: 1-800-429-5001 Mon–Sat, 9 am–9 pm EST. *At all other times or in case of emergency, contact your veterinarian/ veterinary surgeon or emergency veterinary response.*

## Performance Characteristics

VetMate Blood Glucose Monitoring Systems are calibrated to yield results equivalent to plasma glucose concentrations.

**Accuracy:** The accuracy of the VetMate Blood Glucose Monitoring Systems was assessed by comparing blood glucose results obtained by patients with those obtained using a YSI Model 2300 Glucose Analyzer, a laboratory instrument. The following results were obtained by diabetic patients at clinic centers.

### For dogs

Slope	1.0341
Y-Intercept	3.3029
Correlation coefficient (r)	0.9928
Number of measurements	882
Range tested	28.2–382 mg/dL

Accuracy results for glucose concentration < 100 mg/dL

Within ± 5 mg/dL	Within ± 10 mg/dL	Within ± 15 mg/dL
124/306 (40.5 %)	241/306 (78.8 %)	294/306 (96.1 %)

Accuracy results for glucose concentration ≥ 100 mg/dL

Within ± 5 %	Within ± 10 %	Within ± 15 %
346/576 (60.1 %)	521/576 (90.5 %)	575/576 (99.8 %)

Accuracy results for glucose concentration between 28.2 mg/dL and 382 mg/dL

Within 15 mg/dL or 15 %
869/882 (98.5 %)

### For cats

Slope	1.0533
Y-Intercept	4.9766
Correlation coefficient (r)	0.9935
Number of measurements	600
Range tested	28.9–496.3 mg/dL

Accuracy results for glucose concentration < 100 mg/dL

Within ± 5 mg/dL	Within ± 10 mg/dL	Within ± 15 mg/dL
114/210 (54.3 %)	189/210 (90 %)	209/210 (99.5 %)

Accuracy results for glucose concentration ≥ 100 mg/dL

Within ± 5 %	Within ± 10 %	Within ± 15 %
204/390 (52.3 %)	331/390 (84.9 %)	386/390 (99 %)

Accuracy results for glucose concentration between 28.9 mg/dL and 496.3 mg/dL

Within 15 mg/dL or 15 %
595/600 (99.2%)

**PRECISION:** Precision studies were performed in a laboratory using the VetMate BGM Systems.

### For dogs

<i>Within Run Precision</i>			
Interval 1. 30–50	36.9 mg/dL	SD = 2.3 mg/dL	–
Interval 2. 51–110	60.1 mg/dL	SD = 4.2 mg/dL	–
Interval 3. 111–150	123.7 mg/dL	SD = 3.3 mg/dL	CV = 2.6 %
Interval 4. 151–250	187.6 mg/dL	SD = 8.5 mg/dL	CV = 4.5 %
Interval 5. 251–400	286.9 mg/dL	SD = 11.2 mg/dL	CV = 3.9 %

### For cats

<i>Within Run Precision</i>			
Interval 1. 30–50	34.6 mg/dL	SD = 2.7 mg/dL	–
Interval 2. 51–110	64.9 mg/dL	SD = 2.6 mg/dL	–
Interval 3. 111–150	134.2 mg/dL	SD = 4.4 mg/dL	CV = 3.3 %
Interval 4. 151–250	191.1 mg/dL	SD = 8.6 mg/dL	CV = 4.5 %
Interval 5. 251–400	304.6 mg/dL	SD = 12.2 mg/dL	CV = 4.0 %

### For dogs and cats

<i>Between Run Precision</i>			
Interval 1. 30–50	42.2 mg/dL	SD = 2.3 mg/dL	–
Interval 2. 96–144	121.5 mg/dL	SD = 4.7 mg/dL	CV = 3.9 %
Interval 3. 280–420	357.8 mg/dL	SD = 9.8 mg/dL	CV = 2.7 %

## EMC Table

The following tables contain the Manufacturer's declaration and additional information required by IEC 60601-1-2:2014 (Fourth Edition).

Phenomenon	Basic EMC standard or test method	Port tested	Test level/requirement
Mains terminal disturbance voltage	CISPR 11:2015 +A1:2016 EN 55011:2016 +A1:2017	AC Main	Group1, Class B
Radiated disturbance	CISPR 11:2015 +A1:2016 EN 55011:2016 +A1:2017	Enclosure	Group1, Class B
Harmonic Current Emission	IEC 61000-3-2:2014 EN 61000-3-2:2014	AC Main	Class A
Voltage change, Voltage fluctuations and Flicker Emission	IEC 61000-3-3:2013 EN 61000-3-3:2013	AC Main	Pst: 1 Plt: 0.65 Tmax: 0.5 dmax: 4 % dc: 3.3 %
Electrostatic Discharge Immunity	IEC 61000-4-2:2008 EN 61000-4-2:2009	Enclosure	± 8 kV/Contact ± 2, ± 4, ± 8, ± 15 kV/Air
Radiated RF Electromagnetic Field Immunity	IEC 61000-4-3:2006 +A2:2010 EN 61000-4-3:2006 +A2:2010	Enclosure	10 V/m, 80 MHz - 2.7 GHz 80% AM at 1 kHz
Immunity to Proximity Fields from RF wireless Communications Equipment	IEC 61000-4-3:2006 +A2:2010 EN 61000-4-3:2006 +A2:2010	Enclosure	Table 9 in IEC 60601-1-2: 2014

Phenomenon	Basic EMC standard or test method	Port tested	Test level/requirement
Electrical Fast Transient/Burst Immunity	IEC 61000-4-4:2012 EN 61000-4-4:2012	AC Main	± 2 kV, 100 kHz repetition frequency
Surge Immunity	IEC 61000-4-5:2014 EN 61000-4-5:2014	AC Main	Line to Line ± 0.5 kV, ± 1 kV Line to Ground ± 0.5 kV, ± 1 kV, ± 2 kV
Immunity to Conducted Disturbances Induced by RF fields	IEC 61000-4-6:2013 EN 61000-4-6:2014	AC Main	3 V 0.15 MHz - 80 MHz 6 V in ISM and amateur radio bands Between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Power Frequency Magnetic Field Immunity	IEC 61000-4-8:2009 EN 61000-4-8:2010	Enclosure	30 A/m 50 Hz and 60 Hz
Voltage dips	IEC 61000-4-11:2004 EN 61000-4-11:2004	AC Main	0 % $U_T$ ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % $U_T$ ; 1 cycle and 70 % $U_T$ ; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11:2004 EN 61000-4-11:2004	AC Main	0 % $U_T$ ; 250/300 cycle

## Warranty Information

---

### Manufacturer's Warranty

i-SENS, Inc. warrants that the VetMate Meter shall be free of defects in material and workmanship in normal use for a period of five (5) years. The meter must have been subjected to normal use. The warranty does not cover improper handling, tampering, use, or service of the meter. Any claim must be made within the warranty period.

i-SENS will, at its discretion, repair or replace a defective meter or meter part that is covered by this warranty. As a matter of warranty policy, i-SENS will not reimburse the consumer's purchase price.

### Obtaining Warranty Service

To obtain warranty service, you must return the defective meter or meter part along with proof of purchase to your nearest i-SENS sales or customer service representative.

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