



**WWW.OSUNGUSA.COM**

**Instructions For Use  
FOR  
OCCLUSAL PLANE PLATE, ARP1**

**1. Purpose of use**

An instrument that records the positional relationship of the jaw joint of the upper jaw arch in the dentist and reproduces the relationship with the opening axis of the upper jaw on the occlusal machine.

**2. Preparation**

Before using this product, be familiar with the instructions and follow the instructions.

- 1) Significant deformation of the product's appearance, pinholes, attachment of foreign substances, etc. should be checked.
- 2) Follow CDC recommendation for sterilization: Minimum exposure periods for sterilization of wrapped healthcare supplies are 30 minutes at 121°C (250°F) in a gravity displacement sterilizer or 4 minutes at 132°C (270°F) in a pre-vacuum sterilizer.
- 3) The instrument should not come in contact with other instruments during sterilization and drying.

**3. Method of use**

- 1) Follow the method of use during dental procedures in the hospital.
- 2) Let the patient stand so that the patient's head is horizontal. Put the bite on your front teeth. Check the level and remove the product from the patient's mouth.

**4. Storage and keeping conditions**

- 1) After use, wash with a hand or washing machine using a cleaning agent and water to remove blood or foreign substances.
- 2) Store in a dry and cool place out of direct sunlight.
- 3) Check the appearance to check the wear condition or deterioration of the product.

**5. Cautions in use or precaution in use**

- 1) The sterilizer manufacturer's instructions should be followed in order to avoid staining when the instrument is sterilized. Keep the cover closed until the steam is completely removed.
- 2) Do not use the product other than the intended purpose, store at room temperature.

**PRODUCT WARRANTY:**

Replacement warranty is limited due to the degradation of materials at high temperatures after repeated autoclaving.

11/2/2023

15020 West Drive, Houston, TX 77053  
Email: Support@OsungUSA Tel: (844)-870-2781 Fax (847) 589-1298  
www.osungUSA.com