

**Disposable Respirator Mask** Technical Specification Sheet

KN95 GB2626 2006

Proudly developed and tested in collaboration with the Canadian EMS community.





# CANADIAN DESIGNED FOR THE FRONTLINE

The respirator mask is designed to help provide respiratory protection for the wearer. The respirator has a filter efficiency level of at least 95% against particulate aerosols free of oil. It is disposable and can fit a wide range of face sizes.

Always follow User Instructions and use in manners as indicated.

# **KEY FEATURES**

- Synthetic Blood Penetration:ASTM Level 3 Third party tested
- Non-sterile Flammability Rating Class I
- Adjustable Nose Clip
- +90% Fit Test Demographic Across Both Models
- Bacterial Filtration Efficiency +99% -Third party tested

# **MATERIAL COMPOSITION**

- Rubber Elastic Straps
- Aluminum Nose Clip
- Melt Blown Fabric (40%)
- Non Woven Fabric (60%)

## **INSTRUCTIONS FOR USE**



**DO NOT** use for gases or vapors (i.e. anesthetic gases such as isoflurane or vapors from sterilants such as glutaraldehyde.)

**DO NOT** use in any manner not indicated in the User Instructions.

### TIME USE LIMITATION

Respirator may be used until damaged, breathing becomes difficult, or contaminated with blood or body fluids. **Discard after every use.** Follow national, provincial, local, and facility infection control guidance and policies.

### **SHELF LIFE AND STORAGE**

2 years from the date of manufacture.

Use By date on box in YYYY-MM-DD format.

Store respirators in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals.

Store in temperatures between -4 °F (-20°C) and +86°F (+30°C) and not exceeding 80%RH.

### **WARNING**



This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting any disease or infection. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Follow all local regulations. Misuse may result in sickness or death.





25pcs Respirators per box | Carton Quantity 40pcs

### Additionally certified as:

CE Notified Body #: 0370

Conformity to: REGULATION (EU) 2016/425 PERSONAL PROTECTIVE EQUIPMENT

Performance Level: FFP2 NR

Standard: EN 149:2001 + A1:2009

Respiratory protective devices. Filtering half masks to protect against particles. Requirements, testing, marking.

FDK Supply Canada inc. is the sole Manufacturers Representative of Shenzen Fudakin Technology co., LTD