

BIORELEASE®

Wound & Burn Dressing

Description

BioRelease is a transparent, non-ionic hydrogel that maintains a moist environment and helps dissolve wound biofilms. The product is thermo-reversible; it is water-thin at room temperature and becomes a thick hydrogel at body temperature, forming a perfect mold for the wound bed or tunnel. BioRelease is non-cytotoxic.

Ingredients

Poloxamer, phenoxyethanol, octenidine, water, propanediol

Indications for Use

BioRelease is intended for the management of wounds such as stage I-IV pressure ulcers, partial and full-thickness wounds, diabetic foot and leg ulcers, post-surgical wounds, first and second-degree burns, and grafted and donor sites.

Contraindications

BioRelease should not be used if there is a hypersensitivity to any ingredients.

Warning and Precautions

BioRelease can be used during pregnancy as healthcare providers allow. Avoid contact with eyes. No specific drug interaction studies have been undertaken; therefore, use with other products should follow healthcare providers' instructions.


Do not use with saline solution. Use caution if using in combination with any of the following ingredients (Hypochlorous acid, carbomer, acrylates/C10-C30 alkyl acrylate crosspolymer, dentin, bovine serum albumin, copolymer M239144, sodium lauryl sulfate, heat-killed microbes, triethanolamine, sodium hypochlorite, EDTA, ethanol, and nystatin) as it may change the effectiveness of the antimicrobial properties.













Instructions For Use

BioRelease should be applied directly to the wound bed and then covered with the appropriate secondary dressing as indicated.

Storage Conditions

The product should be stored between 65F and 85F

 **Manufactured by**
McCord Research, Inc.
580 N Madison, Ave Suite 4 • North Liberty, IA 52317
www.mccordresearch.com

	Manufacturer		Keep Dry
	Catalog Number		Caution
	Batch Code		Temperature Limits
	Expiration Date		Don't Use If Package Is Damaged
	Prescription Required		Keep Away From Sunlight
	Medical Device		Consult Instructions