

Healing Multiple Open Wounds

PATIENT BACKGROUND

2017

- Split thickness graft right lower leg with later failure of grafts
- Congenital lymphedema

2018

- PVD (peripheral vascular disease) (HCC)
- Atherosclerosis of native artery of both lower extremities with intermittent claudication (HCC)
- Varicose veins of left lower extremity with ulcer of ankle with fat layer exposed (HCC)

2019

- Thigh shingles
- Bell's palsy

2020

- Open wound of right knee with artificial knee with complication, subsequent encounter

- Popliteal artery stenosis, right (HCC)
- Tibial artery stenosis, right (HCC)

2021

- Mitral valve replacement

2022

- Acute respiratory failure with hypoxemia (HCC)
- Acute on chronic systolic CHF
- Acute PE
- Multiple open Wounds in bilateral low ext
- Moderate protein-calorie malnutrition (HCC)
- Hypercholesterolemia
- S/P MVR (mitral valve replacement)
- UTI (urinary tract infection)
- Sepsis (HCC)
- Pacemaker implanted

TREATMENT PLAN

The patient was changed to:

BioRelease® Foaming Cleanser
AgFresh®

The wounds were cleaned and dressings were changed 3 times a week.

OUTCOME

Two wounds completely resolved. Two wounds have reduced 60-70% in size and healthy viable tissue can now be seen over the entire area of each wound.

Total product costs to date less than \$300.

8/28/22




8/28/22



9/29/22



9/29/22



11/8/22



2/22/23




PRODUCT INFORMATION


PART NUMBER	NAME	SIZE	DESCRIPTION	CASE UPC	CARTON UPC	CARTON QTY	CASE QTY
MCR71359	AgFresh®	30mL	Silver Wound and Burn Dressing	754684518572	614409567967	5 EA	5 CARTON
MCR71360	BioRelease®	1.7oz	Wound and Burn Dressing	754684518565	614409567882	1 EA	24 CARTON
MCR71361	BioRelease® Foamer	5oz	Sprayable Wound and Burn Dressing	754684518558	754684518541	1 EA	12 CARTON

McCord Wound Treatment Protocol


Wounded ECM



Synthetic ECM



Closed Wound



Wounds With Bioburden

1. Follow instructions for BioRelese® Foaming Cleanser.
2. Apply a thin coat (a nickel in thickness) of BioRelese® Hydrogel.
3. Cut the corner of the AgFresh® packet and squeeze a bead of the product down the center of the wound.
4. Cut the corner of the BioRelese® packet and squeeze a bead of the product over the AgFresh® Hydrogel.
5. Spread the two products evenly over the wound to fill. The AgFresh® will incorporate into the entire area by itself.
6. Cover with secondary dressing as indicated.

Wounds Without Bioburden

1. Debride wound as indicated.
2. Apply BioRelese® Foaming Cleanser into wound bed liberally and allow to soak in...in preparation for BioRelese® Hydrogel.
3. Cool BioRelese® Hydrogel. in a cup of ice water until it is water thin (about 65° F). BioRelese® Hydrogel will form a thick hydrogel with body contact. This will allow filling in of all crevices before gelling.
4. Squeeze bottle to release product into wound making sure to fill any tunnels.
5. Fill the entire wound bed with product. Spread, if required, to get uniform coverage.
6. Cover with secondary bandage as indicated.

BioRelese® Foaming Cleanser

1. Debride wound as indicated.
2. Spray BioRelese® Foaming Cleanser liberally into wound bed, and 2-3 inches on to wound margin, and allow to soak in. Dab excess product.
3. BioRelese® Foaming Cleanser is a no-rinse formula – The product is designed to remove debris below the wound bed in preparation for BioRelese®.
4. BioRelese® Foaming Cleanser is an antimicrobial product that effectively kills hard to reach pathogens found just below and at the wound's surface.