

**MCCORD**

# BioRelease<sup>®</sup> Foamer

## Wound & Burn Cleansing Dressing

### Description

BioRelease Foamer is a non-ionic, non-cytosine wound cleanser that dissolves biofilms. It can be applied directly to the wound or onto a secondary dressing to aid in moist wound management. If indicated, cover wound with BioRelease or secondary dressing as directed by the healthcare provider.

### Ingredients

Poloxamer, phenoxyethanol, ozenidine, water, propanediol

### Indications for Use

BioRelease Foamer is intended for the cleansing of wounds such as stage I-IV pressure ulcers, partial and full thickness wounds, diabetic foot and leg ulcers, post surgical wounds, first and second-degree burns, grafted and donor sites.

### Contraindications

BioRelease Foamer should not be used if there is a hypersensitivity to any of the ingredients

### Warning and Precautions

BioRelease Foamer can be used during pregnancy as allowed by healthcare provider. Avoid contact with eyes. No specific drug interaction studies have been undertaken; therefore use with other products should follow healthcare providers instructions. Do not use with saline solution. Use Caution if using in combination with any of the following ingredients (Hypochlorous acid, carbomer, acrylates/C10-C30 alkyl acrylate crosspolymer, dentin, bovine serum albumin, copolymer M239144, sodium lauryl sulfate, heat-killed microbes, triethanolamine, bark cork, sodium hypochlorite, EDTA, ethanol, and andnystatin) as it may change the effectiveness of the cleansing properties.

[https://www.researchgate.net/publication/304365219\\_Evidence\\_of\\_incompatibility\\_for\\_topical\\_anionic\\_agents\\_used\\_in\\_conjunction\\_with\\_chlorhexidine\\_gluconate\\_A\\_systematic\\_review](https://www.researchgate.net/publication/304365219_Evidence_of_incompatibility_for_topical_anionic_agents_used_in_conjunction_with_chlorhexidine_gluconate_A_systematic_review)

Reference note: This list is from the Chlorhexidine review, which has the same properties as our products.

### Instructions For Use

BioRelease Foamer should be applied directly to the wound and then covered with the appropriate dressing as recommended by healthcare provider.

Cleanse the entire wound as shown by healthcare provider.

### Storage Conditions

Product should be stored between 65F and 85F

Manufactured by  
McCord Research, Inc.  
580 N Madison  
Ave Suite 4  
North Liberty, IA 52317  
[www.mccordresearch.com](http://www.mccordresearch.com)



	<b>Manufacturer</b>		<b>Keep Dry</b>
	<b>Catalog Number</b>		<b>Caution</b>
	<b>Batch Code</b>		<b>Temperature Limits</b>
	<b>Expiration Date</b>		<b>Don't Use If Package Is Damaged</b>
	<b>Prescription Required</b>		<b>Keep Away From Sunlight</b>
	<b>Medical Device</b>		<b>Consult Instructions</b>