


Mecca Care Shower Seat 400x330mm



SKU



CHROME
NRCR0003CH



MATTE BLACK
NRCR0003MB



BRUSHED NICKEL
NRCR0003BN



GUN METAL
NRCR0003GM



BRUSHED GOLD
NRCR0003BG



BRUSHED BRONZE
NRCR0003BZ

Features

1. Ergonomic design folds neatly up against wall.
2. High quality, high strength materials and finishes.
3. Compact size suitable for all shower sizes.
4. Over 150kg capacity when installed into brick/masonry or multiple wall studs.
5. AS1428.1 Approved for ambulant bathrooms tested to 150kg.
6. Gun metal, brushed gold and brushed bronze colours are PVD coated for long term durability.

WARNING:

These seats are designed to be mounted onto wall studs or masonry walls. Installation is recommended to be completed by qualified trades person who complies with the AS 1428.1-2009 design for access and mobility standard.

It is recommended rail location is determined by a qualified occupational therapist or trained professional to the AS 1428.1-2009 standard. Correct installation is required to ensure safe, long term operation. Incorrect mounting can result in reduced capacity of the rail which can result in injury to user and damage to the wall.



MECCA CARE
DDA CARE 32MM

Installation Instruction

Cautions

- 1: Please confirm whether the wall is a load-bearing before drilling;
- 2: Before drilling, please confirm whether there is a water pipe at the hole position;
- 3: Check whether the product accessories are complete before installation;
- 4: Make sure the drill bit is 10mm in size.
- 5: Make sure the hole depth is 6.5cm
- 6: The recommended height from floor to top of seat is about 470-480mm.

| | |
|--|--|
| | <p>1. Ensure wall is suitable for installing the shower seat before you begin installation. Place wall fixing plate on wall in correct position and mark out holes.</p> <p>Tips: Recommend height from floor to top of seat is 470-480mm. Then recommend height from floor to wall fixing plate is 453-463mm, refer to left picture.</p> |
| | <p>2. Firstly, drill 6 holes on the wall based on the screw holes of our shower seat (with a drill bit of 10mm and a hole depth of 6.5cm). Apply silicon to the hole and use a hammer to knock the nylon bolt cap into the hole on the wall.</p> |
| | <p>3. Place the shower seat on the wall, align the screw holes with the wall holes, and then tighten the screws with a socket wrench.</p> |
| | <p>4. Cover the base cover tightly.</p> |
| <p>*Note: you can choose the style without slow down function, and the installation method is the same</p> | <p>5. After installation, fold up the seat when not to use it and unfold it when to use.</p> |

***Note: The drilling bit is 10mm, and the drilling depth is greater than or equal to 6.5cm**