



DEFENDER

FOOT DEFENDER FITTING TIPS




STEP 1

Deflate the air cell by pushing the silver button

Place patient in a tall diabetic sock that extends above mid calf

Remove AFO front closure

Place patient's heel firmly into the back of the boot




STEP 2

Apply front AFO closure, starting with the hook & loop near the toes, work your way up

Be sure that the top straps are not too tight, allowing for a 2 finger gap at the shin.

Note the boot is not intended to fit tight



STEP 3

We don't recommend patients with neuropathy or wounds use the pneumatic system for fitting; for all others, pump air (red button) 3-8 times for a snug fit at the heel.

Release the air before fitting or removing the boot. **Never exceed 8 pumps total** as it may trigger the safety valve on the pump and render it unusable.

This is a safety feature for patients with neuropathy. The pneumatic pump also can restrict the blood flow down to the foot so check with your medical professional before using.

✓ JUST RIGHT

Note the AFO angle to the tibia (shin)



✗ TOO TIGHT

Wrong angle to the tibia (shin)



Scan QRCode
YouTube Patient Instruction Link

We don't recommend patients with neuropathy or wounds use the pneumatic system for fitting.