Lansinoh®

Compact Single Electric Breast Pump

Instructions for Use

Important Safeguards

The Lansinoh® Compact Single Electric Breast Pump is a personal care item and for hygienic reasons should never be shared between mothers or resold from one mother to another. Either of these actions will void the warranty.

If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breastmilk using the Lansinoh® Compact Single Electric Breast Pump will not reduce or remove the risk of transmitting the virus to your baby through your breastmilk.

When using electrical products, especially when children are present, basic safety precautions should always be maintained, including the following:

Intended Use of Breast Pump

Lansinoh Compact Single Electric Pump is intended to express and collect the breastmilk of a nursing woman for the purpose of feeding collected breastmilk to a baby. Lansinoh Compact Single Electric Pump is intended for a single user.

Indications of Breast Pump

Lansinoh Compact Single Electric Pump enable mothersto provide breast milk for babies who cannot feed directly at the breast for a variety of reasons, such as latch-on problems, premature or unwell babies.

Contraindication of Breast Pump

There are no known contraindications.

Warning: To reduce the risk of burns, electrocution, fire, or injury to persons:

- This product should never be left unattended when plugged in.
- Close supervision is necessary when this product is used near children.
- Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep the cord away from heated surfaces.
- Never use while sleeping or drowsy.
- Never drop or insert any object into any opening or tubing.
- Do not use this product if it is broken or any parts are missing.
- Do not allow small children or pets to play with the pump unit, AC adapter, etc.
- Do not use outdoors with the AC adapter or operate where aerosol spray products are being used or where oxygen is being administered.

Packaging materials are not part of the product. For your child's safety, please remove and discard all packaging materials prior to use but read and retain the instructions for future reference.

Save these instructions and your receipt for future reference.

Read All Instructions Before Using

Danger: To reduce the risk of electrocution:

- Always unplug this product immediately after using.
- Do not use while bathing.
- Do not place or store this product where it can fall or be pulled into a tub or sink.
- Do not place or drop pump into water or other liquid.
- Do not reach for an electrical product that has fallen into water. Unplug immediately.

1. QUICK START GUIDE	2
2. PUMP	2
2.1 Breast Pump Parts and Accessories	2
2.2 Cleaning Instructions	4
2.3 Assembling the Breast Pump	4
2.4 Using the Breast Pump	4
2.5 Breast Flange Sizing	5
2.6 Pumping Tips	5
2.7 Troubleshooting	6
3. BREASTMILK STORAGE GUIDELINES	7
4. TECHNICAL DATA	8
4.1 Power Supply	8
4.2 Control System	8
4.3 Regulatory Compliance	8
5. CUSTOMER SERVICE	8
6. WARRANTY	8

Please do not return this pump to the retailer.

Lansinoh's Customer Service team is here to help. If you experience any issues with the product, please call our Customer Service Department at +86 400 821 0385, email customerservice@lansinoh.cn, visit our website www.lansinoh.cn, or contact your country distributor.

Health protection & Hygiene i.e. for health protection and hygiene reasons, breast pumps may not be returned or exchanged once the hygiene seal has been broken. If your product is faulty, please contact our Customer Service Department.

1. QUICK START GUIDE

DO NOT wash tubing.

Before first use, clean all parts that will come into contact with milk according to instructions on page 3, section 2.2.

Assembly

See detailed instructions on page 3, section 2.2. Proper assembly is essential to ensure proper function and suction.

Insert purple diaphragm so that feet align with the cutouts on the breast flange body

Attach clear cover and twist so tubing connector is at the back

Attach breast flange

Insert valve

Attach bottle

Attach tubing

Secure connector

Attach to electric pump unit

Only use the Lansinoh® Compact Single Electric Breast Pump with Lansinoh® Breastmilk Expression Set.

Getting Started

Hold Breastmilk Expression Set to your breast so your nipple will extend straight toward the nipple tunnel of the breast flange. Make sure the breast completely fills the ComfortFit® breast flange so that a vacuum is formed and no air escapes. Make sure the set is tilted slightly downward to allow milk to flow into the bottle.


Using Let-Down and Expression Phases


- Press the power button  to start the pump in Let-Down Phase . Let-Down Phase will run automatically for 2 minutes.
- To begin the Expression Phase sooner, press the Let-Down/Expression toggle button  to move to the Expression Phase. Press the Let-Down/Expression toggle button  again at any point to switch back to Let-Down Phase.
- Let-Down Phase is a faster pumping style that mimics the short, rapid sucking your baby uses to initiate milk flow.
- Expression Phase mimics the slower, deeper sucking of feeding, and features 5 adjustable suction levels that can be tailored for your comfort.
- The breast pump features LED lights beneath the purple keypad to indicate which setting you are on. Let-Down Phase is identified on the keypad with slowly dimming LED lights. Expression Phase is identified on the keypad with still LED lights.

2. PUMP

2.1 Breast Pump Parts and Accessories

Breast pump unit

Breastmilk Expression Set	
a	Diaphragm Cap
b	Diaphragm
c	ComfortFit® Breast Flange (25mm)
d	Breast Flange Body for Electric Pump
e	White Valve
f	Soz (160ml) Bottle
g	Sealing Cap
Breast Pump Parts and Accessories	
h	Silicone Tube (800mm)
i	Pump Connector
j	AC Adapter
Manual Breast Pump Parts	
l	Handle
m	Stem
n	Silicone Diaphragm
o	Breast Flange Body for Manual Pump

2.2 Cleaning Instructions

*Unless otherwise instructed by your health care provider, follow these instructions.

**DO NOT wash tubing or connector.

Cleaning of Breastmilk Expression Set Before First Use

1. Always wash and dry hands thoroughly with soap and water before handling the Breastmilk Expression Set.
2. Disassemble the Breastmilk Expression Set and place parts in a large pot.
3. Cover the parts in water making sure there is room for all of the parts to float freely in the pot.
4. Bring the water to the boil and allow the parts to a boil for 5 minutes.
5. After boiling remove the parts from the pot with tongs and place the pieces neatly on a clean paper towel or in a clean drying rack and allow them to air dry.
6. Allow the parts to dry completely before assembly and use. Partially dried parts may affect the operation of

the pump.

7. Once parts are dry they can be stored in a clean plastic tub with a lid or in a clean plastic zipper bag.

Cleaning of Breastmilk Expression Set Before Each Use

1. Always wash and dry hands thoroughly with soap and water before handling the Breastmilk Expression Set.
2. Before assembling the Breastmilk Expression Set, be sure that all of the parts are completely dry.
3. Visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. These may affect the functionality of the pump.
4. Should you have any questions or concerns about your parts please contact Lansinoh at **+86 400 821 0385**, Wechat account:lansinoh520, visit our website www.lansinoh.cn, for further guidance.

Cleaning of Breastmilk Expression Set After Each Use


1. Disassemble the Breastmilk Expression Set, setting aside tubing and connectors
2. Rinse all parts with cold water immediately after use. Place pump parts in a clean wash basin used only for washing infant feeding equipment. Do not place pump parts directly in the sink because germs in sinks or drains could contaminate the pump.
3. Hand wash all parts of the Breastmilk Expression Set using warm soapy water. A teat brush may be used to scrub the hard plastic pieces of the Breastmilk Expression Set, but should not be used on the white valve or diaphragm.
4. Rinse parts with running clean, warm water, or by submerging in fresh water in a separate basin that is only used for cleaning infant feeding equipment.
5. Place parts in a pot and cover with water making sure there is room for all of the parts to float in the pot.
6. Bring the water to a boil and allow to boil for 5 minutes.
7. After boiling remove the parts from the pot with tongs and place the pieces neatly on a clean paper towel or in a clean drying rack and allow them to air dry. Do not use a dish towel to rub or pat items dry because doing so may transfer germs to the items.
8. Allow the parts to dry completely before assembly and use. Partially dried parts may affect the operation of the pump.
9. Once parts are dry they can be stored in a clean plastic tub with a lid or in a clean plastic zipper bag.
10. Clean your wash basin and bottle and/or teat brush after each use.
11. Electric unit should be wiped down with a clean paper towel or soft cloth after use. The electric unit should never be put into water or other liquids for cleaning.

WARNING: Use caution when cleaning the small white valve. It is important for proper suction. Do not use excessive force and do not clean with a baby bottle brush or a nipple brush. Gently wash by hand in hot, soapy water and rinse under hot running water.

WARNING: It is highly recommended not to use any method other than the above mentioned cleaning methods.

To use steam or microwave sterilisers: Follow the manufacturers instructions. Ensure all openings are face down in the steriliser. Ensure you wash your hands before removing equipment from the steriliser.

2.3 Assembling the Breast Pump

2.4 Using the Breast Pump


The Lansinoh® Compact Single Electric Breast Pump is very easy to use, offering several settings that you can customise to maximise your milk production and comfort. There are 2 elements which you can customise whilst expressing, and the breast pump features LED lights beneath the purple keypad to indicate which setting you are on.

Let-Down and Expression Phases







When babies breastfeed, they start with quick sucking to stimulate your milk let-down and flow. Once your milk has let down, the baby draws the milk out with slower, deeper sucks. To simulate this natural rhythm, the The Lansinoh® Compact Single Electric Breast Pump offers two phases – the Let-Down Phase and the Expression Phase.

- Let Down-Phase is a faster pumping style which mimics the short, rapid sucking your baby uses to initiate milk flow
- Expression Phase mimics the slower, deeper sucking of feeding and features 5 adjustable suction levels that can be tailored for your comfort.

The switch between the 2 phases can be easily operated using the Let-Down/Expression toggle button. Let-Down Phase is identified on the keypad with slowly dimming LED lights. Expression Phase is identified on the keypad with still LED lights.

When the pump is turned on, it starts in Let Down-Phase and will automatically switch to the Expression Phase after 2 minutes. If your milk starts flowing before the 2 minute transition, you can switch to Expression Phase at any point by pressing the  button. You can always return to the Let-Down Phase when needed by pressing the  button again.

Operating the Breast Pump

1. Ensure that the breast pump has been properly cleaned. Wash your hands thoroughly and make sure your breasts are clean.
2. Ensure the breast pump and Breastmilk Expression Set have been correctly assembled and the silicone tube is securely attached to the pump via the tubing connector and to the diaphragm cap.
3. Hold the Breastmilk Expression Set to your breast with the centre of your nipple extending towards the tunnel of the ComfortFit® breast flange.
4. Make sure the breast completely fills the ComfortFit® breast flange so that a vacuum is formed and no air escapes.
5. The Breastmilk Expression Set and nipple tunnel must tilt slightly downward to allow milk to flow naturally toward the bottle.
6. Press and hold the On/Off button  for a few seconds to turn the pump on.
7. The breast pump will always start in Let-Down Phase Level 3 and will automatically switch to Expression Phase after 2 minutes. If your milk starts flowing before the 2 minute transition, you can switch to Expression Phase at any point by pressing the  button. You can always return to the Let-Down Phase when needed by pressing the  button again.
8. The Expression Phase has 5 adjustable suction levels. You can increase the suction level by pressing the  button. You can decrease the suction level by pressing the  button. A maximum of 5 LED lights will illuminate depending on the selected suction level.
- **Note: If at any time the suction becomes uncomfortable, reduce it by pressing the  button.**
9. When a pumping session is finished, turn the pump off with the On/Off Button . Ensure there is no suction on the nipple before removing the ComfortFit® breast flange from the breast. If necessary, break the suction

by inserting a finger between the breast flange and the breast. Once the suction is broken, the ComfortFit® breast flange can be removed safely.

- Carefully remove the bottle from the breast flange body and close the bottle with the sealing cap provided. Please follow the breastmilk storage guidelines as outlined in Section 3.

● **Note: If the pump unit is left on and unattended for more than 60 minutes, it will automatically power off.**

Power Options

The pump may be operated using the following power supply modes:


- **AC adapter**
Connect the AC Adapter to the power connector socket located on the On/Off button side of the pump. Plug the AC Adapter into an available electrical outlet. Electrical safety can only be achieved when the breast pump is connected to the AC Power Adapter supplied with the product (See “Technical Data” section).
- **Power bank (not included)**
Connect the power bank USB cable to the power connector socket located on the On/Off button side of the pump. Electrical safety can only be achieved when the breast pump is connected to the power bank specified power supply (See “Technical Data” section).

We recommend using your pump with the AC adapter provided.


The Lansinoh® Compact Single Electric Breast Pump can easily switch to Manual Breast Pump.

Note: Prior to assembling the Lansinoh® Manual Breast Pump, ensure that all parts have been cleaned.

1. Wash your hands thoroughly.
2. Insert the stem into the silicone diaphragm.

3. Place the silicone diaphragm and stem on top of the breast flange body (for manual breast pump). Make sure it fits securely around the rim by pressing down with your fingers to ensure a perfect seal. (These parts are designed to fit snugly.)


4. Gently secure the white valve onto the breast flange body from underneath, making sure white valve is flush with the breast flange body. It is not necessary to use excess force when attaching the white valve.
5. Align the groove in the handle with the breast flange body and gently press the handle down. Push until you hear it click.
6. Adjust the stem and handle so that the stem fits the first notch of the handle, known as the Let-Down Phase setting.

7. Attach the ComfortFit® breast flange to the breast flange body.
 - To switch to Expression Phase, gently press the “PUSH” part at the tip of the handle to set the handle to the second(bottom) notch of the stem. Continue to squeeze and release the handle to find the appropriate suction pressure.
 - To return to the Let Down Phase from the Expression Phase,(1) Pull the tip of the stem towards you,(2) gently squeeze the handle, and set it in the first(top) notch.

2.5 Breast Flange Sizing

When you are pumping only your nipple and a small amount of your areola should be pulled into the breast flange tunnel. If a large amount of areola is being pulled into the tunnel, this indicates that the breast flange is too large. You need a smaller size.

When you are pumping your nipple should move in and out of the tunnel freely. If you notice that your nipple is touching or rubbing along the sides of the tunnel, this indicates that the breast flange is too small. You need a larger size.

The breast cushion you received with the pump is 25mm. 21mm, 28mm, 30.5mm and 36mm ComfortFit® breast cushions are available for purchase. Simply visit www.lansinoh.cn, or call +86 400 821 0385 or contact your country distributor for guidance.

2.6 Pumping Tips

Breastmilk is produced on a supply and demand basis. Expressing your breastmilk regularly will stimulate your body and help you maintain your milk supply. If you have any breastfeeding or breast pump related questions, consult an International Board Certified Lactation Consultant (IBCLC), breastfeeding counselor, or other healthcare provider such as your midwife or health visitor.

Top Tips for Pumping Sessions

- If you and baby are apart, look at their pictures or videos, or have something with your baby's scent on it. These cues will help to trigger let-down.
- If baby is with you, skin-to-skin contact and scent will help to trigger let-down. You may pump one breast while nursing baby on the other. This allows you to maximize pumping effectiveness by using the let-down reflex stimulated by your baby's sucking.
- Give yourself sufficient time to pump. Stress can negatively impact milk production.
- Stay hydrated; having a glass of water before and during each session is important for milk production.
- Many mums find that they have a high milk supply in the morning after sleeping, so you may find it helpful to try a morning pumping session.
- Some mums can achieve more than one let-down during a pumping session. You may be able to achieve a second by returning to the Let-Down Phase by pressing the Let-Down/Expression toggle button.
- Track the date, time, and duration of your feeding and pumping sessions as well as the amount of milk expressed from each breast. Keeping these records can help you identify trends or patterns so you can establish your optimal feeding and pumping routine. If you're having any difficulties, it can be helpful to share this information with your healthcare professional.
- You can save time and prevent messy milk transfers or spills by pumping directly into Lansinoh® Breastmilk Storage Bags. Just ensure that your Lansinoh® Breastmilk Storage Bag has white and purple flowers and a white line above the logo. If your bag does not have this design, it is not compatible with pumping directly into the bag. NOTE: Pumping directly into bags requires two hands to safely remove milk storage bags. Visit www.lansinoh.cn to learn more.

2.7 Troubleshooting

Issue	Solution
The pump does not turn on	<p>If you experience an interruption in the normal operation of the electrical and/or mechanical functioning of the pump, simply do the following:</p> <ul style="list-style-type: none"> ● Check to confirm the AC adapter or power bank is securely connected to the pump. ● If you are using the AC adapter, unplug it from the electrical outlet, wait a few seconds, and then plug it back into the electrical outlet. ● If you are using a power bank, unplug it from the pump, wait a few seconds, and then plug it back in. <p>If you encounter further problems, please visit www.lansinoh.cn, call +86 400 821 0385, or contact your country distributor for guidance.</p>
I experience pain when I express	<p>If you feel any discomfort, you may need a different size breast flange. Visit www.lansinoh.cn, or call +86 400 821 0385, or contact your country distributor.</p>
I do not feel any suction or the suction is very low	<p>If it feels like the suction has decreased, try the following:</p> <ul style="list-style-type: none"> ● While pump is on, remove the tubing from the diaphragm cap and place a finger over the tubing connector. If you feel suction, the pump is operating properly, but the Breastmilk Expression Set may not be correctly assembled or the white valve may be damaged. If you don't feel any suction, try unplugging the pump, wait a few seconds, and plug it back into the electrical outlet. ● Make sure the Breastmilk Expression Set has been properly assembled. Check the assembly information in Section 2.3, page 4. ● Confirm that all connections (especially the tubing) are secure. Make sure the breast flange forms a complete seal around the breast. ● Inspect the white valve for a tear or pinhole. This component is critical to achieving proper suction. Replace the white valve with one of the spare valves provided with the product.

3. BREASTMILK STORAGE GUIDELINES

Breastmilk Storage and Thawing Guidelines for Healthy Term Babies

Check with your healthcare professional or an International Board Certified Lactation Consultant for specific storage instructions. When freezing breastmilk, store in the back of the freezer. Do not store breastmilk in the freezer door as the temperature can fluctuate when the door is opened and closed.

Before storing, date all breastmilk and use the oldest breastmilk first.

Milk Storage Guidelines, Full term infant			
WHERE	TEMPERATURE	TIME	ACTION
At Room Temperature	60-80°F (16-26°C)	4 hours (ideal) up to 6 hours (acceptable)	Contents should be covered and kept as cool as possible; covering the container with a damp towel may keep milk cooler
Insulated Cooler Bag	39-59°F (4 - 15°C)	24 hours	Keep ice packs in constant contact with milk containers; limit opening cooler bag
In a Refrigerator	<39°F (<4°C)	4 days (ideal) up to 5 days (acceptable)	Store milk in the back of the main body of the refrigerator
Freezer	0°F (-18°C)	3 months (ideal) up to 6 months (acceptable)	Store milk away from sides and toward the back of the freezer where temperature is most constant
Deep Freezer	-4°F (-20°C)	6 months (ideal) up to 12 months (acceptable)	

Source: La Leche League International 2017

These breastmilk storage guidelines are a recommendation. Contact your health professional or consult your country's national guidelines for further information.

Thawing Breastmilk

Thaw frozen milk in the refrigerator or place the containers under cool running water, gradually warming the water until the milk is thawed. Thawed milk can be kept in a refrigerator for up to 24 hours.

DO NOT refreeze milk once it is thawed.

Gently swirl thawed milk before feeding to mix the layers that have separated. Breastmilk is not homogenized and the "cream" will rise to the top. Breastmilk can vary in colour and consistency depending upon the time of the day that the milk was expressed, what you have eaten, and the age of the baby at the time of pumping. For further advice, please contact your health care professional.

DO NOT microwave breastmilk!

Microwaving alters the composition of the milk. Microwaving also may create "hot spots" in the milk that could burn the infant.

DO NOT thaw breastmilk in hot or boiling water.

4. TECHNICAL DATA

Power Supply

AC Adapter operation: 100V/240V wall plug type AC adapter, galvanically isolated with 5V DC output voltage.

Pump power dissipation: 10 watts maximum

Voltage regulation: Linear voltage regulator with low dropout voltage and short circuit protection.

Power bank operation: CE marked and 5V/2A power bank can be used.

Control System

The expected service life of the pump is 500 hours.

Control unit: Microcontroller with integral program memory

Program memory: 8K flash ROM, 1K RAM memory
There are 2 phases: Let Down and Expression selectable from keypad.

Adjustable suction level: 5 suction levels are selectable and controlled from keypad.

Suction vacuum level(Expression): Minimum 120 mmHg, maximum 270 mmHg

Regulatory Compliance

SYMBOLS	MEANING OF THE SYMBOLS
	Follow instructions for use
	Identifies the manufacturer
	Indicates the date of manufacture

SYMBOLS	MEANING OF THE SYMBOLS
	Type BF applied part
	Class II equipment
	Keep the device dry
	Keep the device away from sunlight
	Indicates the temperature limitation for operation, transport and storage
	Indicates the humidity limitation for operation, transport and storage
	Indicates the atmospheric pressure limitation for operation, transport and storage
IP22	Indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water
	This pump is made of various metal and plastics. Before disposal, the device is to be rendered unusable and it must be disposed of as sorted waste in accordance with local waste electrical and electronic equipment regulations. To avoid the harmful effects on the environment and public health align with your defined local return and collection rules for waste electrical and electronic equipment (incl. batteries).
	Indicates that the packaging contains products intended to come in contact with food
	Indicates to handle the fragile device with care

 **Warning: Children have STRANGLED in cords.**
Keep the cords out of the reach of children.

LANSINOH MAKES NO EXPRESS WARRANTY OF ANY KIND REGARDING THE PRODUCT OTHER THAN THOSE WARRANTIES EXPLICITLY SET FORTH HEREIN. ANY IMPLIED WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR ANY IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, TO THE EXTENT PERMITTED BY LAW, SHALL BE LIMITED IN DURATION TO A PERIOD OF 90 DAYS FROM THE DATE OF PURCHASE BY THE ORIGINAL PURCHASER.

In the event the product is found by Lansinoh® to not meet the above limited warranty, as purchaser's sole and exclusive remedy Lansinoh® will repair or, at the option of Lansinoh®, replace this product without charge for such replacement or parts of labour.

NOTE: The purchaser shall bear all expense for returning this product to Lansinoh®. This warranty does not apply to any product that has been subjected to misuse, abuse, or alteration. Only Lansinoh® replacement parts and accessories are recommended for use with Lansinoh® breastfeeding products. Improper or incorrectly performed maintenance or repair or use of non-Lansinoh® replacement parts or accessories may invalidate this warranty.

TO THE EXTENT PERMITTED BY LAW, LANSINOH SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, REPLACEMENT COSTS RESULTING FROM THE BREACH OF ANY WRITTEN OR IMPLIED WARRANTY.

If you wish to make a claim under this warranty, please telephone +86 400 821 0385, email customerservice@lansinoh.cn, or contact your country distributor. You will be given a Freepost address to return the product, if necessary, together with your original dated bill of sale or other dated proof of purchase and a brief statement of the problem.

Declaration for Electromagnetic Compatibility (EMC)

For information use only

FCC Statement

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Recommended separation distance between Portable and mobile RF communications equipment and the model @ 3Vrms

The Lansinoh® Compact Single Electric Breast Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Lansinoh® Compact Single Electric Breast Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Lansinoh® Compact Single Electric Breast Pump as recommended below, according to the maximum output power of communications equipment.

5. CUSTOMER CARE

Lansinoh is here to support you in your breastfeeding journey, so if you have any questions, please call our Customer Service Department at **+86 400 821 0385**, email customerservice@lansinoh.cn, visit our website www.lansinoh.cn, or contact your country distributor.

6. WARRANTY

Limited Consumer Product Warranty

This product is warranted to the original consumer purchased (the "purchaser") to be free from defects in material and workmanship which are not commercially acceptable for the period of one year from the date of purchase on the pump mechanism and the period of 90 days on all other product components. Warranty coverage terminates if you sell or otherwise transfer this product to another person. This warranty gives you specific legal rights and you may also have other rights, which vary by location.

Separation distance according to frequency of transmitter m
 Rated Maximum output power of transmitter W

Wireless Statement

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance of 3.5m away from the equipment.

Stacking

This device should not be used adjacent to or stacked with other equipment. And if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.

This product is a personal care product and therefore cannot be returned once opened. If you experience any problems with this product, please telephone our free customer service line at +86 400 821 0385 in China or contact your country distributor, or contact your country distributor. Our representative will be happy to assist you.

Serious events experienced from using this product should be reported immediately to Lansinoh Laboratories and to the local competent Authority.

**Manufacturer: LANSINOH LABORATORIES SAG.
 GERIC. TAS. SAN. TIC. LTD. STI.**

**Zafer SB Mah. Gunduz Sk. No:13 E.S.B. 35410
 Gaziemir – İzmir - Turkey**

**Lansinoh Laboratories UK Ltd | lansinoh.co.uk
 Ground Floor, West Mayesbrook House,
 Lawnswood Business Park,
 Leeds, LS16 6QY,
 UK**

Lansinoh® is a registered trademark of Lansinoh Laboratories, Inc. ©2021 Lansinoh Laboratories, Inc. All rights reserved.

Date IFU issued: 13 July 2021

Version: Second Version

Recommended separation distance between Portable and mobile RF communications equipment and the model @ 3Vrms			
The Lansinoh® Compact Single Electric Breast Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Lansinoh® Compact Single Electric Breast Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Lansinoh® Compact Single Electric Breast Pump as recommended below, according to the maximum output power of communications equipment.			
Separation distance according to frequency of transmitter m			
Rated Maximum output power of transmitter W	150 khz to 80 Mhz	80 Mhz to 800 Mhz	800 Mhz 2,5 Ghz
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,36
100	11,67	11,67	23,33

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All Rights Reserved
TAS, SAN, TIC, LTD, STI
Lansinoh Laboratories Saglik Gereçleri Tasarim San. Tic. Ltd. Sti.