argotherm 2

DE Betriebsanleitung	3 - 19
EN User Manual	20 - 36
FR Mode d'emploi	37 - 53
IT Istruzioni d'uso	54 - 71
ES Modo de empleo	72 - 88
PT Instruções de operação	89 - 106

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Tab. 4

Table of Contents

1	Explanation of Symbols	21
2	General Safety Instructions	22
3	Suitable Personnel	22
4	Machine Specifications	23
4.1	Argotherm 2	23
4.2	Program overview, shortcut keys .	26
5	Installation	27
5.1	Setting up	27
5.2	Protective gas connection	28
5.3	Compressed-air connection	29
5.4	Electrical connection	29
6	Application and Operation	30
6.1	Function test	31
6.2	Starting the sintering process	31
6.3	Finalizing the Sintering Program ..	32
7	Cleaning and Maintenance	33
7.1	Cleaning	33
7.2	Daily checks	33
7.3	Monthly maintenance	33
8	Malfunctions, Repairs and Warranty ..	34
8.1	Malfunctions	34
8.2	Repairs	34
8.3	Warranty	34
9	Environmental Protection	35
10	Technical Data and Accessories/Spare Parts	36

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

When there is a risk of burns from hot surfaces, the exclamation mark in the warning triangle is substituted by a heat-alert symbol.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
_	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

DANGER:

This product does not comply with the ATEX Directive.

- ▷ Do not use in explosive atmospheres!
- ▷ Do not operate with explosive gases or mixtures!

WARNING:

Fire hazard!

- ▷ In case of unexpected processes in the furnace (for example, strong smoke development or odour nuisance), switch the furnace off immediately and allow the furnace to cool down to room temperature.

CAUTION:

A defective machine can lead to malfunctions! On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

WARNING:

Operation with power sources, products, consumables, auxiliary aids, etc., which are subject to the Ordinance on Hazardous Substances or in any way have an impact on the health of the operating personnel, is not permitted.

NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

NOTE:

- Increased wear of door gasket, insulation, heating rods and furnace housing!
- ▷ Avoid opening the furnace when in hot condition (above 200 °C)!

3 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 Machine Specifications

4.1 Argotherm 2

4.1.1 Delivery Scope

- _ High-temperature furnace Argotherm 2
- _ Heating rod (4 x)
- _ insulating stone, furnace ceiling
- _ Exhaust air pipe
- _ Ceramic tube
- _ Floor plate
- _ Insulating stone, front
- _ Sintering compartment Argovent (2)
- _ Gas regulator/gauge
- _ Connection hose for gas supply
- _ Connection hose for compressed-air supply
- _ Mains cable
- _ Allen key

▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

4.1.2 Intended Use

Argootherm 2 is a high-temperature furnace for sintering milled CoCrMo units. The compact furnace is used as a table-top model.

This high-temperature furnace has been specially developed for sintering the Sintron material. Only the Sintron materials, whose properties and sintering temperatures are known, may be used. Any other or further-going use, such as the processing of other products than those intended for, as well as the handling of hazardous substances or materials hazardous to one's health, does not apply as under the intended use.

As to the materials being used in the furnace, it must be known whether they are aggressive to or even could lead to destruction of the insulation or of the heating rods.

It is prohibited to remove safety devices, bypass them or put them out of operation.

Unauthorised modifications/alterations of the machine shall also void the warranty.

The set-up instructions and the safety regulations must be followed, as otherwise the furnace is considered as not being used as intended for, and any and all claims against the manufacturer shall become void.

4.1.3 CE Declaration of Conformity

This product was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety.

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

EC Directives

2006/95/EC	2004/108/EC
------------	-------------

Applied harmonised standards

DIN EN ISO 12100	DIN EN ISO 13732-1
DIN EN 61010-1	DIN EN 61000-6-1, DIN EN 61000-6-3

Tab. 1

The product's Declaration of Conformity can also be viewed on the manufacturer's website.

4.1.4 Components and interfaces

EN


Fig.1 Machine overview, front

- | | |
|--|---------------------|
| 1 Exhaust air pipe | 6 Insertion plate |
| 2 Furnace lift door | 7 Shroud insulation |
| 3 Thermoelement | 8 SD card slot |
| 4 Heating rod | 9 USB port |
| 5 Indicator for flow of protective gas | 10 Controller |


Fig. 2 Machine overview, rear

- 1 Protective gas connection
- 2 Mains switch
- 3 Power supply connection
- 4 Fuse
- 5 Switchgear fan
- 6 Service unit
- 7 Compressed-air connection

Display, operating level


Fig. 3 Display, operating level

- P1...P4** Fast keys, program 1 ... 4
- Program start key
- Program stop key
- Key for calling up the configuration level
- 1 Temperature curve
- 2 Current temperature
- 3 Current program status
- 4 Current program
- 5 Indication of the program run-time
- 6 Remaining time until program ends

Display, configuration level


Fig. 4 Display, configuration level

- "Program" key
- "Function" key
- "Selection" key
- "Info" key
- "Save" key

- 1** Function indication
- 2** Selected program
- 3** Selected segment of program
- 4** Corresponding parameter

Handling the programs:

- ▷ To select a program: Press the key and the respective number key.
- ▷ To take over the program: Press the key. The display switches to the operating level.
- ▷ To display the segments of the program: Press the key.

Handling the functions:

- ▷ To select a function: Press the key and the respective number key.
 - Key 1: Function "Gas1"
 - Key 2: Function "Gas2"
 - Key 3: Function "Air"
 - Key 4: Function "Preset time for program"
- ▷ To start or stop function 1 - 3: Press the key.

To program function 4:

- ▷ Press the key.
- ▷ Press the key.
- ▷ Using the number keys, enter the desired time [hh:mm] until the program starts.

After entering the last digit, the cursor will jump back to the first numerical location, so that possibly required corrections can still be made.

- ▷ Press the key.
- ▷ Press the key.
- ▷ Press the number key for the respective program.
- ▷ Press the key. The program is taken over with the time preset. The display switches to the operating level. After the program starts, the preset time runs off first, then the program starts.

4.2 Program overview, shortcut keys

The Argotherm 2 has a permanently stored program for sintering Sintron.

Program description

- _ Program 1 [P1]
 - Standard program for Sintron

5 Installation

5.1 Setting up

NOTE:

Due to the weight of 30 kg, it is recommended that two persons set up the unit.

▷ When setting up, observe the valid country-specific accident prevention regulations.

- _ The machine is intended exclusively for use within dry, closed rooms.
- _ The machine must be positioned below an exhaust pipe or an exhaust hood (as used for pre-heating furnaces).

NOTE:

Malfunction!

▷ Do not install an exhaust pipe directly near the exhaust air pipe of the furnace!

- _ Sufficient ventilation must be ensured at the set-up location, in order to draw off the heat and possibly developing exhaust gas.
- _ A clearance of at least 50 cm is to be observed sideways, to the rear and upwards of the furnace. The lateral clearance to non-flammable materials can be reduced to 20 cm. Space requirements for the Argotherm 2 (incl. connections):
 - Width: 1.4 m (0.8 m for non-flammable materials)
 - Depth: 0.9 m
 - Height: 1.0 m
- _ The set-up surface must be level.
- _ The empty weight of the machine is 30 kg. The set-up surface must withstand an appropriate load.

- _ During operation, the room temperature must be between 5 °C and 40 °C; avoid large variations in temperature.
- _ The floor covering must be made of non-flammable material, so that hot material falling out of the furnace can not inflame the covering.
- ▷ Set up the machine on a sturdy work table or work bench (on a non-flammable surface).
- ▷ Remove all packaging material, including from the furnace compartment.

When setting up, please also take into consideration:

- _ Assembly and repair instructions, Argovent 2 (34416-FB-NTR)
- _ Assembly and repair instructions, installation of heating elements (34219-FB-NTR)

DANGER:

For carrying out the further setting-up steps, a qualified professional in the area of electrical safety is mandatory. The installation and repair of heating elements must be performed by trained and certified personnel.

5.2 Protective gas connection

- ▷ Connect the gas regulator/gauge to the protective-gas cylinder.
- ▷ At the bottom left of the furnace rear wall, insert the protective gas hose into the gas connection. The resistance of the seal must be overcome, in order to ensure tightness against leaks.


Fig. 5

- ▷ Connect the other end of the gas hose to the gas regulator/gauges of the protective-gas cylinder. For this, insert the protective gas hose into the push-in connection by turning and pushing in, until the resistance of the seal is overcome.


Fig. 6

- ▷ Completely open the right-hand valve on the pressure reducer.

To check the gas hose's tightness against leakage:

- ▷ Open the gas cylinder and gas regulator/gauge.
- ▷ Widely pivot/move about the gas hose at the gas regulator/gauge.

There should not be any hissing noise whatsoever.

Requirements for the protection-gas supply

For operation, the protective gas must fulfil the following minimum requirements:

- _ Argon 4.6 (or higher) from commercially available cylinders (preferably 20 l)
- _ Cylinder pressure 200 bar/2900 psi

 When the left-hand manometer pointer points to 0 or below, the cylinder is empty and must be exchanged. The cylinder then only holds a residual pressure of approx. 5 bar/72.5 psi, which may no longer be used for the sintering. The furnace is equipped with a compressed-air and protective-gas monitoring feature. When there is no pressure on the gas side, the furnace interrupts the sintering procedure or the program will not start.

5.3 Compressed-air connection

To provide operation with air cooling, the furnace must be connected to the compressed-air supply.

Requirements for the compressed-air supply

For operation, the compressed air must fulfil the following minimum requirements:

- _ Flow rate: 100 l/min (3.54 CFM)
- _ Pressure: 6 bar (87 psi)

Requirements on air purity:

- _ Solid contaminants: Class 3; better than 5 µm for solids
- _ Water content: Class 4; Max. pressure dew point +3 °C
- _ Total oil content: Class 2; Max. oil content 0.1 mg/m³

Establishing the compressed-air connection

- ▷ Insert the compressed-air hose to the stop into the service unit.


Fig. 7

- ▷ Connect the other end of the compressed-air hose with the compressed-air connection in the laboratory.

5.4 Electrical connection

The electrical connection is to be made via a socket outlet with earthing contact providing 220-240 V and 50-60 Hz. This socket outlet must be located close to the furnace and be easily accessible.

When using an extension cord or a multiple socket outlet, do not exceed the respective max. electrical load rating.

**WARNING:**

The mains cable may not be damaged!

- ▷ Do not place objects on the mains cable.
- ▷ Lay the mains cable in such a manner that nobody can step on or trip over it.

**NOTE:**

Do not insert the mains plug before the heating elements have been inserted.

- ▷ Set the mains switch to the **O** position.
- ▷ Connect the enclosed mains cable to the socket intended for this purpose on the rear of the furnace.


Fig. 8

- ▷ Connect the other end of the cable to the socket outlet with earthing contact.

The mains cable may only be exchanged by an approved identical one.

6 Application and Operation

Before starting-up the furnace, it should acclimate at the set-up location for 24 hours. Swift heating up in cold condition can cause damage to the heating rods (for example, burning-through).

During the first heating up, there may be odour nuisance due to binding agent fumes escaping from the insulation material.

- ▷ Provide for thorough ventilation of the furnace set-up location when heating up for the first time.

**DANGER:**

Hot air at the outlet opening!

- ▷ Do not reach with your hand at or above the outlet opening.
- ▷ Do not place objects on the furnace!

**WARNING:**

Danger of fire and explosion!

The furnace surfaces become very warm.

- ▷ Do not place objects on the furnace!

**NOTE:**

Insufficient sintering results due to faulty temperature control!

- ▷ Prior to every process, check the thermoelement for damage.

6.1 Function test

The individual functions, such as cooling and gas flow, can be manually set and checked via the controller.

A function test must be carried out before using the furnace for production:

- ▷ Check the furnace's protective gas function (see 6.1.1).
- ▷ Check the furnace's compressed-air cooling (see 6.1.2).

6.1.1 Checking the furnace function for the protective gas

The correct quantity of protective gas must be ensured prior to the first program run. It is manually set and checked with the *F1* and *F2* functions.

- ▷ Open the valve of the gas cylinder.
- ▷ Access the configuration level with the  key.
- ▷ Press the **F** key.
- ▷ Press the **1** and/or **2** key.
- ▷ Press the  key.

The function starts.

- ▷ Upon completion of the test, press the  key to end the function.

6.1.2 Checking the furnace function for the compressed-air cooling

	WARNING:
Danger of burning!	
▷ Check the air flow only when the furnace is cold.	

The furnace is factory-pre-set when shipped. When the air flow is active, a pressure of 1 bar is automatically given.

The air flow must be checked prior to the first program run.

- ▷ Access the configuration level with the  key.
- ▷ Press the **F** key.
- ▷ Press the **3** key.
- ▷ Press the  key.

Air flows into the sintering compartment.

- ▷ Upon completion of the test, press the  key to end the function.

6.2 Starting the sintering process

6.2.1 Preparations for the sintering

- ▷ Equip the heating furnace with Argovent (see separate Assembly and repair instructions, Argovent 2)
- ▷ Connect compressed air.
- ▷ Connect protective gas supply and adjust.

6.2.2 Sintering

- ▷ Shut the furnace door.
- ▷ Open the valve of the gas cylinder.
- ▷ Provide for the compressed-air supply.
- ▷ Switch the furnace on.
- ▷ Select program 1 on the furnace keypad by pressing the **P1** key.

- ▷ Start the program with the  key.

The furnace carries out an automatic pressure test for a few seconds. The pressure reading of 1.5 bar is indicated on the display.

Fig. 9 Indication of pressure test

The protective gas monitoring checks, if the gas cylinder is open and if protective gas is present for the sintering. If this is not the case, the furnace will not start and will output an error message. The gas pressure is continuously monitored in the running program.

6.2.3 Post-preparations for the Sintering

When the sintering program is over, a temperature of 50 °C is indicated and the compressed-air cooling is stopped.

- ▷ Switch the furnace off.
- ▷ Using the bowl/cover tongs, carefully remove the protective-gas retort out of the furnace compartment and place it down on a non-flammable surface.

NOTE:

Tilting off or falling down of the sintering bowl, due to its high weight.

- ▷ Apply bowl/cover tongs ensuring a firm and good grip.
 - ▷ Using the bowl/cover tongs, carefully remove the sintering bowl out of the furnace compartment and place it down on a non-flammable surface.
 - ▷ Remove the sintering cover and place it down on a non-flammable surface.
 - ▷ Remove the sintered frame.
- Before continuing:
- ▷ Blast off the sintered frame with commercially available aluminium corundum (50 - 110 µm).
 - ▷ Evaluate the sintering result based on the Sinteron sintering chart (see User Manual Argovent 2).

6.3 Finalizing the Sintering Program

When the sintering program is over, a temperature of 50 °C is indicated and the compressed-air cooling is stopped.

- ▷ Switch the furnace off.
- ▷ Shut the main valve of the gas cylinder.
- ▷ When not using for longer periods, leave the Argovent components in the furnace.

7 Cleaning and Maintenance

7.1 Cleaning

DANGER:

Risk of electric shock!
 ▷ Prior to any cleaning, pull the mains plug.

NOTE:

Possible damage to the furnace!
 ▷ Do not use a high-pressure cleaner for the cleaning!
 ▷ Use only commercially available, water-based or non-flammable, solvent-free cleaning agents for the cleaning.

Prior to cleaning:

- ▷ Allow the furnace to cool down to the room temperature.
- ▷ Clean the various components according to Table 2.

Component	Cleaning
Housing (exterior)	Wipe off with moist, lint-free cloth
Interior	Carefully vacuum out with vacuum cleaner
Insulation materials	Carefully vacuum out with vacuum cleaner
Door gasket	Wipe off with moist, lint-free cloth
Instrument panel	Wipe off with moist, lint-free cloth

Tab. 2

7.2 Daily checks

- ▷ Check the system for any exterior damage.
- ▷ Check the function of the safety devices (e.g. the mains switch).
- ▷ Check the tightness and proper seating of the gas line.
- ▷ Check the function of the fan.
- ▷ Visually check the following components:
 - Sealing surfaces on door and furnace shroud
 - Visible part of the heating rods in the furnace compartment
 - Visible part of the thermoelement in the furnace compartment

7.3 Monthly maintenance

- ▷ Carefully vacuum out and clean the furnace compartment, exhaust openings and exhaust pipe.
- ▷ Check the furnace compartment, exhaust openings and exhaust pipe for damage.

8 Malfunctions, Repairs and Warranty

Cracks in the insulation

The furnace insulation consists of a high-grade, fire-resistant material. Due to the thermal expansion, cracks develop in the insulation after only a few heating cycles. However, these do not have an influence on the function or quality of the furnace.

8.1 Malfunctions

▷ In case of malfunctions, contact the manufacturer's Technical Service.

8.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

Replacing a fuse

- ▷ Switch the furnace off and pull the mains plug.
- ▷ If required, wait until the furnace has cooled to room temperature.
- ▷ On the rear side of the furnace, turn the fuse holder in anticlockwise direction and pull it out (bayonet lock).
- ▷ Replace the defective fuse against an identical one.
- ▷ Insert fuse holder with new fuse in furnace rear wall and lock it by twisting in clockwise direction.

8.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

In addition, it also applies that warranty and liability claims are excluded for personal and property damage, when these can be attributed to one or more of the following causes:

- _ Each person involved with operating, mounting, maintaining, repairing the system must have read and understood the User Manual. No liability shall be assumed whatsoever for damages and operating malfunctions resulting from non-observance of the User Manual.
- _ Use of the system not in accordance with the intended use.
- _ Improper mounting, starting-up, operation and maintenance of the system.
- _ Operation of the system with defective safety devices or improperly mounted or inoperable safety and protective devices.
- _ Non-observance of the notes in the User Manual in terms of transport, storage, assembly, starting-up, operation, maintenance and equipping the system.
- _ Unauthorised modifications to the system.
- _ Unauthorised modification of the operating parameters.
- _ Unauthorised changing of parameters and settings as well as program changes.
- _ Original parts and accessories are designed especially for Argotherm 2. When exchanging components, use only original parts. Otherwise, the warranty shall become void. For damages caused by the use of non-original parts, the manufacturer shall exclude any and all liability.
- _ Disaster cases caused through foreign bodies and force majeure.

9 Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

10 Technical Data and Accessories/Spare Parts

Subject to changes.

Technical data

	Unit	Value
Art. No.	–	178740
Dimensions (D × W × H) incl. maintenance unit	mm	500 × 461 × 480
Weight	kg	32
Electrical connections	V/Hz	220-240/50-60
Max. output	kW	3.5
Fuse (time-delayed)	A	15
Degree of protection	–	IP20
Thermal protection class according to DIN EN 60519-2	–	Class 0
Sound pressure level	dB(A)	< 70
Ambient conditions:		
– Temperature	°C	+5 - +40
– Humidity	%	80

Tab. 3

Accessories/Spare parts

Art. No.	Designation
178741	Argotherm 2 heating rod
178702	Argotherm thermoelement
178703	Cylinder pressure reducer, EU
178704	Pressure reducer, cylinder DK
178705	Pressure reducer, cylinder FR
178706	Pressure reducer, cylinder GB/SG
178707	Pressure reducer, cylinder US
178708	Cylinder pressure reducer, AU
178719	Cylinder pressure reducer, CN
178745	Argovent 2
178746	Sintering base Argovent 2
178747	Sintering retort Argovent 2
178748	Sintering bowl Argovent 2
178749	Sintering cover Argovent 2
178750	Retort weight Argovent 2
178751	Separation disc Argovent 2
178752	AV2 base plate for Argotherm 2
178753	Connector Argovent 2
178754	Pliers, Ceramill Argovent 2
178755	Sintering pearls Sintron Argovent 2 150 g
178729	Auxiliary powder, Sintron

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Amann Girrbach AG