dental Scanner

DE Betriebsanleitung	3 - 26
EN User Manual	27 - 49
FR Mode d'emploi	50 - 74
IT Istruzioni d'uso	75 - 99
ES Modo de empleo	100 - 124
PT Instruções de operação	125 - 149

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Tab. 3

- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	28	6	Application Tips	36
2	General Safety Instructions	29	6.1	Scanning of Splitex models	37
3	Suitable Personnel	29	6.2	Scanning with the universal plate .	38
4	Machine Specifications	30	6.2.1	Universal plate	39
4.1	Delivery scope of the complete system	30	6.2.2	Positioning and affixing single models	40
4.2	Intended Use	30	6.2.3	Positioning and affixing model pairs	40
4.3	CE Declaration of Conformity	30	6.3	Articulator scan	42
4.4	Components and Interfaces	31	6.3.1	Positioning articulators	43
4.4.1	Desktop scanner	31	6.3.2	Locating surface for Artex CR	44
4.4.2	Scanner accessories	32	6.4	M-Die	44
4.4.3	PC	32	6.4.1	Product Description	44
5	Installation	33	6.4.2	Usage	45
5.1	Transport	33	6.5	All-in-bar	46
5.2	Setting up	33	6.5.1	Product Description	46
5.3	Installation of the PC	34	6.5.2	Usage	47
5.4	Installation of the software on the PC	34	7	Cleaning and Maintenance	48
5.5	Installation of the scanner	34	7.1	Cleaning	48
5.6	Activating additional modules	34	7.2	Maintenance	48
5.7	Calibration of axes	35	8	Malfunctions, Repairs and Warranty ..	49
			8.1	Malfunctions	49
			8.2	Repairs	49
			8.3	Warranty	49
			9	Environmental Protection	49
			10	Technical Data	49

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

_ **NOTE** means that property damage can occur.

_ **CAUTION** means that light to fairly serious personal injury can occur.

_ **WARNING** means that serious personal injury can occur.

_ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Fuse

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

CAUTION:

A defective machine can lead to malfunctions! On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

CAUTION:

When scanning, there is danger of crush injuries between the scanner arm and housing, as well as the risk of injury from clothing, jewellery or hair being drawn in.

- ▷ Do not reach into the scan chamber during scanning.
- ▷ Observe warning stickers on the unit.

CAUTION:

Possible health hazard through strip light! The scanner operates with stripe light. Continuous visual contact with stripe light can trigger epileptic seizures, migraines, etc.

- ▷ Do not operate the scanner in the presence of persons with such a health predisposition.

NOTE:

Possible damage to the scanner and models due to missing or incorrect fastening!

The movement of the axes in the scanner can cause unsecured or incorrectly secured models to fall down during scanning.

- ▷ Do not place objects on or inside the scanner.
- ▷ Do not place models in the scanner without fastening them.
- ▷ Use only the provided object holders or approved accessories for fastening the models.
- ▷ Affix occlusion models only with new, sturdy rubber bands. Thin or porous rubber bands can tear.

NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

3 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 Machine Specifications

4.1 Delivery scope of the complete system

- _ Desktop scanner
 - USB cable
 - Mains cable
 - Operating Instructions
 - USB stick with calibration data
 - Calibration model
 - _ Additional accessories:
 - M-Die
 - All-in bar
 - Universal plate
 - Fixation spring lever
 - _ PC
 - Windows operating system
 - Mains cable (1 x)
 - Monitor (1 x display port cable, 1 x mains cable)
 - Key board, mouse, headset
 - _ CAD software package
 - CAD dongle
 - Quickguide Installation
 - CAD software CD
- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

4.2 Intended Use

The desktop scanner is a PC-controlled 3D scanner for scanning dental tooth models.

NOTE:

The scan result depends on the materials being scanned. Metal surfaces, for example, can lead to artefacts. Scan spray can improve the scan capability of materials.

The measuring accuracy of the scanner is only guaranteed if the scanner is calibrated. A calibration model and the corresponding standard values are included in the delivery scope.

For improper use of the system, the manufacturer shall not assume any liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

Failure to follow the instructions in this Manual may impair the function or safety. Therefore, please read the Manual before starting operation.

4.3 CE Declaration of Conformity

This product was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety.

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

EC Directives	Applied harmonised standards
2014/35/EG	
2014/30/EG	DIN EN 61010-1:2010 DIN EN 61326-1:2013

Tab. 1

The product's declaration of conformity is available upon request from the manufacturer.

4.4 Components and Interfaces

4.4.1 Desktop scanner


Fig. 1 Machine overview, front

- 1 Base plate
- 2 Anti-slip mat
- 3 Orientation lines
- 4 Splitex retainer
- 5 Operation button
- 6 Height lines
- 7 Locating surface for Artex CR


Fig. 2 Machine overview, rear

- 1 Main switch
- 2 Power supply connection
- 3 USB port

4.4.2 Scanner accessories


Fig. 3

- 1 Fixation spring lever
- 2 Fixation pin
- 3 Universal plate
- 4 M-die
- 5 Compartments for abutments
- 6 Holding pin
- 7 All-in-bar

4.4.3 PC


Fig. 4 Icon image

- 1 Network connection, Internet
- 2 USB connections for keyboard, mouse and dongle
- 3 Digital connection, monitor
- 4 Power supply connection, PC

5 Installation

5.1 Transport

NOTE:

Possible damage of the scanner due to improper lifting.
 Too much load on the scanner base can cause it to warp.
 ▷ Lift and carry the scanner only in the area of the feet.


Fig. 5

▷ For transport, have two persons lift the unit by the floor plate in the area of the feet.

5.2 Setting up

NOTE:

Possible measuring errors due to reflections on the measuring object!
 Strong ambient light leads to undesired reflections on the measuring object. This impairs the accuracy of the measurements.
 ▷ Set up the scanner in such a manner that it faces away from windows or strong artificial lighting.

- _ The machine is intended exclusively for use within dry, closed rooms.
- _ Minimum space requirements for the desktop scanner (incl. connections):
 - Width: 424 mm
 - Depth: 415 mm (excl. connection cable)
 - Height: 469 mm
- _ The set-up surface must withstand a load of 25 kg.
- _ During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls; this is important due to natural oscillation).
- ▷ The opening of the scan chamber must face away from windows and strong artificial lighting.
- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.

5.3 Installation of the PC

- ▷ Connect the mains cable to the power supply connection of the PC as well as to a mains socket outlet.
- ▷ Connect a network cable to the network connection of the PC and the Internet router/modem.
- ▷ Connect the monitor and PC using the provided display port cable.
- ▷ Connect the mouse and the keyboard to the PC.

5.4 Installation of the software on the PC

- ▷ Execute installation routine of the provided software as administrator.

5.5 Installation of the scanner

 As each scanner is individually factory-calibrated, the data on the USB stick are only valid for the respective scanner.

Prior to installation:

- ▷ Check if the serial number of the scanner on the type plate is identical with the serial number stated on the USB stick!


Prior to installation:

- ▷ Make sure that the main switch of the scanner is positioned on **0**.
- ▷ Connect the provided USB cable to the scanner and PC (USB 2.0).
- ▷ Insert the CAD Dongle in a free USB port on the rear side of the PC.
- ▷ Connect the mains cable to the power supply connection of the scanner as well as to a mains socket outlet.

5.6 Activating additional modules

Additional modules must be activated in the scanner software before they can be used.

- ▷ Click on *Options > Access > Add Activation Code...* to enter the received scanner-specific activation code.

- ▷ Restart the software to activate the module.

5.7 Calibration of axes

NOTE:

Possible inaccurate measurements without calibration or calibration with a damaged calibration model!

The measuring accuracy of the scanner is only guaranteed if the scanner is calibrated. For calibration, a calibration model and the corresponding standard values are included in the delivery scope.

- ▷ Calibrate the unit
 - _ after putting into operation
 - _ during operation whenever the software prompts this action.
- ▷ Use the calibration model only when in proper condition.
- ▷ Check if the calibration model is damaged at any central location.
- ▷ Do not start the calibration until the values entered in the software correspond with the values on the calibration model.

Following the set-up and connection to the PC, the desktop scanner must be calibrated using the supplied calibration model.


Fig. 6


Fig. 7

- ▷ Insert the calibration model
- ▷ Create a new patient case and start the scanner software from the data base.
- ▷ Cancel the subsequent command prompt for the 2D scan.
- ▷ Start the axes calibration under *Scanner Options > Axis calibration*.


Fig. 8

- ▷ Follow the software instructions.

If required:

- ▷ Re-enter the calibration model values already entered during the installation by clicking on *Application > Register Calibration Object*.

6 Application Tips

With the scanner's comprehensive object holder system, individual jaw models, tooth abutments, quadrant models and occlusion models can be securely attached and scanned.


Fig. 9

- 1 Base plate
- 2 Anti-slip mat
- 3 Orientation lines
- 4 Splitex retainer
- 5 Operation button
- 6 Height lines
- 7 Locating surface for Artex CR

Element	Function
Base plate	Holds the articulator
Anti-slip mat	Secure placement of the articulator
Orientation lines	Orientation aid for secure placement of articulators on the anti-slip mat
Splitex retainer	Retains models with Splitex plates and for exact positioning of scanner accessories in the unit.
Operation button	Starts the next scan step (after command prompt through the software)
Height lines	Determination of the model height; input of the height in software possible
Locating surface	Facilitates the alignment of the Artex CR

Tab. 2

6.1 Scanning of Splitex models


Fig. 10 Model on Splitex plate in scanner

No additional attachment is required for jaw models on a Splitex plate. The model can be placed directly on the scanner's Splitex retainer. Through the integrated magnet, the model is sufficiently affixed without the requirement of additional fastening elements.

► To remove out of the scanner: Grasp the model at the rear or by the side and lift it off.

Determining the model height

In order to measure models, these must be located in the measuring field. The scanner has an automatically moving axis for alignment of the models in the measuring field. In conjunction with the scanner software, determining the model height is necessary for some scanning processes.

The model height must be determined, in order for the unit to align the model at the correct height. The determined height must be specified in the scanner software before starting the scanning process.


Fig. 11 Height lines

- 1 Single line
- 2 Double line

After inserting a model into the scanner:

- ▷ Check the height of the line (single or double line) at which the denture is located.


Fig. 12 High model (single line)


Fig. 13 Low model (double line)

- ▷ In the software, select the line at which the height of the denture is located.


Fig. 14

- 1 Button for single line (high model)
- 2 Button for double line (low model)

The scanner now automatically aligns the model to the required height in the measuring field.

- ▷ Start the scan procedure.

6.2 Scanning with the universal plate


Fig. 15 Model fixation

- 1 Universal plate (delivery scope)
- 2 Modelling wax
- 3 Fixation pin incl. fixation spring lever
- 4 Splitex retainer

The Splitex retainer is in the desktop scanner. The universal plate or alternatively models with Splitex counter plate are placed on it, depending on requirement.

6.2.1 Universal plate

The universal plate is used for transferring models without Splitex retainer/plate into the scanner. It can be used both to affix individual models as well as model pairs in occlusion.


Fig. 16 Universal plate

- 1 Top side
- 2 Metal pins
- 3 Holes for fixation pins

_ The top side of the universal plate [1] is level; it is used to apply modelling wax onto which the dental models are affixed.

_ Two metal pins are located at the rear side of the level surface [2]. They represent the rear limit for placement of the models.

_ The universal holder comprises 16 holes [3]. These are used for placement of the fixation pins including the fixation spring levers.

Placement in the scanner

▷ Place the universal plate directly on the Splitex retainer in the scanner.

Through the integrated magnet, the universal plate is sufficiently fixed on the Splitex holder without the requirement of additional fastening elements.

▷ To remove the universal plate from the scanner, grasp it at the rear or at the sides and lift it out.

Fixation pins incl. fixation spring levers

The fixation pins incl. the fixation spring levers are used to affix model pairs in occlusion onto the universal plate.

▷ For this, insert the fixation pins into optimally matching holes until the metal ring faces against the universal holder.

The fixation pins are placed loosely into the holes and achieve their stability through proper use with the fixation spring levers.

NOTE:

Only insert the fixation pins into the universal plate when being used for model fixation. Otherwise, they may fall out during the scanning process and damage the scanner.

Function method

The fixation via the fixation spring levers takes place by means of canting/chocking the spring levers with the fixation pins.

▷ Position and align the fixation spring lever with respect to the object.

▷ Using one finger, centrally press/push down the on the spring lever.

This causes the fixation spring lever to deform. After pressing/pushing down the fixation spring lever, it cants/chocks with the fixation pins and thus firmly clamps the object against the universal holder.


Fig. 17 Use of fixation spring levers

6.2.2 Positioning and affixing single models

Always position models in such a manner that they are within the scan focus.

The placement of the model on the universal plate should be performed outside of the scanner.

- ▷ Apply modelling wax uniformly on the universal plate.
- ▷ Place the model centrally on the universal plate and press it firmly against the modelling wax.

To ensure that the model being scanned is sufficiently affixed to the universal plate and does not become loose and fall off during the scan process, the following must be observed:

- ▷ Press the model with both hands firmly against the modelling wax.
- ▷ Pay attention that the complete bottom surface of the model base has contact with the modelling wax.

NOTE:

The model must not project out beyond the rear limit of the universal plate (metal pins).

- ▷ Afterwards, place the universal plate together with the model on the scanner's Splitex retainer and press it on.

NOTE:

When using modelling wax not supplied by the manufacturer, it can occur that the model becomes loose/detaches during the scan process.

Insufficient scan results or damage of mechanical components and/or the model are possible.

- ▷ Make sure that the modelling wax provides sufficient holding force during the scan process.

6.2.3 Positioning and affixing model pairs

Always position model pairs in such a manner that they are within the scan focus.

Model pairs can be affixed either with a rubber band or with use of the fixation spring levers.

The placement of model pairs on the universal plate should be performed outside of the scanner.

Placement of model pairs with integrated articulation elements

For model pairs with integrated articulation elements (e.g. triple tray), the jaw assignment is already specified. These model pairs can be applied directly to the universal plate.

- ▷ Apply modelling wax uniformly on the universal plate.
- ▷ Place the model pair centrally on the universal plate and press it firmly against the modelling wax.

Placement of model pairs without integrated articulation elements

For model pairs without integrated articulation elements, it is recommended to firstly fasten only the lower jaw to the universal plate.

- ▷ Apply modelling wax uniformly on the universal plate.
- ▷ Place the lower jaw model centrally on the universal plate and press it firmly against the modelling wax.
- ▷ Afterwards, bring the upper jaw model in occlusion to the affixed lower jaw model.

Fixation with rubber band

For the fixation of model pairs with a rubber band, a suitable and appropriate household rubber band can be used.

The contact surface of the universal plate to the Splitex retainer has an opening that provides sufficient space for a rubber band.

- ▷ Place the model pair onto the universal plate.
- ▷ Fasten/tension suitable rubber band around model pair and universal plate.

NOTE:

The rubber band must be tensioned tight enough, to prevent the models from loosening or moving during the scan process.

The contact surfaces of the Splitex plate may not be enclosed by the rubber band, as this would influence the hold or grip of the universal plate in the scanner.


Fig. 18 Fixation with rubber band

NOTE:

Tearing rubber bands during scanning can lead to damage of the scanner and the models. Especially old rubber bands can become brittle and easily tear.

- ▷ It is recommended to use several strong, new rubber bands.

Fixation with spring levers

After placing the model pairs on the universal plate, mount the fixation spring levers.

NOTE:

It is recommended to use at least two fixation spring levers to achieve a uniform distribution of the hold-down pressure and to prevent the model pairs from tilting.

The position of the fixation pins and the fixation spring levers can be adjusted.

- ▷ Insert the respective fixation pin into a hole of the universal holder that is suitable for fixation of the model.


Fig. 19 Use of fixation spring levers

6.3 Articulator scan

The articulator is used for affixing model pairs that are in occlusion. By positioning the articulator in the scanner, the model pair in occlusion is transferred to the software.

The articulator must be placed in the scanner in the specified preferential direction. With the Artex CR, a condyloid-related transmission is possible.

For vestibular scans, any articulator can be used without further aids. In this case, the assignment of upper and lower jaws is based on average value, and not condyloid-related as with the Artex CR.

Positioning the articulator in the scanner


Fig. 20 Occlusion model in the articulator

- 1 Upper supporting pin
- 2 Front side
- 3 Gripping points

- ▷ Ensure that the occlusion model is correctly articulated.
- ▷ Unscrew upper supporting pin [1] from the articulator.

Inserting an articulator without the corresponding software command can cause the articulator to fall down. This can cause damage to the scanner, the articulator or the models.

- ▷ Wait until the software prompts you to insert the articulator.
- Only then is it ensured that the axis movements are minimal during the following scanning process.

6.3.1 Positioning articulators

The scanner's base plate is moved into the horizontal position and additionally rotated.

The anti-slip mat has 2 lines to aid the positioning of articulators.


Fig. 21

- 1 Orientation line for the rear side model placement in the articulator
- 2 Orientation line for articulator center axis
- 3 Incisal pin

NOTE:

Possibility of insufficient scanning results and damage to the articulator!

If the articulator is positioned incorrectly in the scanner, it can be moved or fall out of the unit during the scanning process.

- ▷ Align the articulator in such a manner that the **front** part of the articulator (left of line [1]) is inside the base plate and does not protrude beyond the edge of the base plate. The **rear** parts of the articulator (right of line [1]) may protrude beyond the edge of the base plate.

- ▷ Align the longitudinal axis of the articulator on line [2]. In this, the incisal pin [3] must point to the left in the direction of the sensor optics.
- ▷ Align rear side of model holder in articulator along line [1].

The feet of the articulator should not be placed directly on the line [1], but usually to the right of the line. They must still be positioned on the anti-slip mat.

- ▷ When the articulator is correctly aligned: Start the scan procedure.

When scanning is complete:

- ▷ Remove articulator from scanner with both hands.
- ▷ Follow the next steps in the software.

6.3.2 Locating surface for Artex CR

For easier positioning of an Artex CR, a locating surface is integrated in the axis cover. It serves as an additional orientation aid for the positioning of an Artex CR. When the Artex CR is correctly positioned, its right side faces against the locating surface.


Fig. 22

1 Locating surface for Artex CR

6.4 M-Die

6.4.1 Product Description


Fig. 23 Delivery scope

- 1 M-Die
- 2 Compartments for abutments
- 3 Holding pin
- 4 Fixing compound

M-Die is an additional component for the “Intelligent Multi-Die” and “Multi-Cap” scanning processes:

- _ With “Intelligent Multi Die”, up to twelve abutments of a patient can be captured in a single scan.
- _ With “Multi-Cap”, up to twelve jaw and patient-independent single abutments can be captured.

The indication “Multi-cap” is limited to offset caps.

 M-Die can only be used in conjunction with the universal plate.

6.4.2 Usage

The usage is described in the “Intelligent Multi Die” mode.

 A gap is provided for between the M-Die and the universal plate. This gap allows for placement of the M-Die onto the universal plate without having to firstly remove the modelling wax from the plate.

Global Scan

 The placement of the model on the universal plate should be performed outside of the scanner.

▷ Fasten the model to be scanned (upper or lower jaw) centrally on the universal plate.

 Always position models in such a manner that they are within the scan focus.

▷ Position the universal plate into the scanner.
 ▷ Scan the models.

Positioning abutments

▷ Mount the M-die [4] onto the universal plate [5].
 ▷ Fill fixing compound [1] into the M-Die compartments.

▷ At the position indicated in the software [3], work the abutments to be scanned [2] centrally into the M-Die (here compartment 7).


Fig. 24 Positioning abutments

- 1 Fixing compound
- 2 Abutment
- 3 Compartment 7
- 4 M-Die
- 5 Universal plate

▷ Press the abutments into the fixing compound in such a manner that they face flat and level against the surface of the M-Die ([A]).

- Ensure that both fixing compound and abutment do not project out above the surface of the M-Die/compartment [B].
- Make sure that there is no gap between abutment and compound [C].
- Do not press abutment too deep into compound [D].


Fig. 25 Press-in-depth of the abutments

Scanning

- ▷ Depending on scanner, place the universal plate with M-Die directly onto the Splitex retainer or model plate.

The integrated magnet provides for sufficient fixing of the universal plate.


Fig. 26 Universal plate with M-die in the scanner

- ▷ Scan the abutments in the “Intelligent Multi Die” mode.

 The scanned abutments are automatically assigned to the correct position in the dental arch when scanning in “Intelligent Multi Die” mode.

- ▷ To remove the universal plate from the scanner: Grasp the plate at the rear or by the side and lift it off.

6.5 All-in-bar

6.5.1 Product Description


Fig. 27 Delivery scope

- 1 All-In Bar
- 2 Holding pin
- 3 Modelling wax

The All-In Bar is used for positioning of the single elements of a quadrant model (triple-tray model) for the scanning process. With the All-In Bar, single abutments as well as upper and lower jaw segments are affixed to the universal plate. In this, abutments are positioned higher on the All-In Bar in order to prevent shadings during the scanning. Up to three abutments can be placed on the All-In Bar.

 The All-In Bar can only be used in conjunction with the universal plate.

 There is a gap between the All-In Bar and the universal plate (see Fig. 30 on page 25). This gap allows for placement of the All-In Bar onto the universal plate without having to firstly remove the modelling wax remainders from the plate.

6.5.2 Usage

- The models are scanned in two sequences:
- _ Scanning of the upper and lower jaw models in occlusion.
 - _ Scanning of the adjoining and freely-positioned upper and lower jaw models in the All-In Bar.

Scanning the models in occlusion

 The placement of the model on the universal plate should be performed outside of the scanner.

- ▷ Centrally affix the model to be scanned with the fixing spring lever (alternatively using a rubber band) on the universal plate.


Fig. 28 Model in occlusion

- 1 Fixation spring lever
- 2 Triple-tray model
- 3 Abutment
- 4 Modelling wax
- 5 Universal plate

 Always position models in such a manner that they are within the scan focus.

- ▷ Position the universal plate into the scanner.
- ▷ Scan the models.

Scanning the freely-positioned models

 The placement of the model on the universal plate should be performed outside of the scanner.

- ▷ Insert the All-In Bar via the holding pins in the respective holes of the universal plate.


Fig. 29 Mounting All-In Bar on universal plate

- 1 All-In Bar
- 2 Universal plate

- ▷ Fill compartments of the All-In Bar with modelling wax.
- ▷ Disassemble quadrant model in upper and lower jaw.
- ▷ Remove abutment to be scanned.
- ▷ Affix abutment within the modelling-wax-filled compartment.

- ▷ Affix upper and lower jaw with to model base at the left and right of the All-In Bar on the universal plate using modelling wax.
 - Pay attention that the complete bottom surface of the model base has contact with the modelling wax.
 - Bring model segments into contact with the side surfaces of the All-In Bar, so that the model segments are within the scan field during the scanning process.
 - Press the model segments with both hands firmly against the modelling wax.


Fig. 30 Model with freely-positioned abutment

- 1 Model segment
- 2 Abutment
- 3 Modelling wax
- 4 All-In Bar
- 5 Model base
- 6 Universal plate

NOTE:

The model must not project out beyond the rear limit of the universal plate (metal pins).

- ▷ Position the universal plate into the scanner.
- ▷ Scan the models.

7 Cleaning and Maintenance

7.1 Cleaning

The scanner must be cleaned regularly.

- ▷ Switch the scanner off.

NOTE:

Possible damage to the 3D sensor:

- ▷ Do **not** clean the optics of the 3D sensor in the top of the scanner!

- ▷ Carefully vacuum out the interior of the scanner.
- ▷ Wipe off the housing with a moist cloth. Do not use aggressive cleaning agents.

7.2 Maintenance

Maintenance may be carried out only by trained, specialized personnel.

8 Malfunctions, Repairs and Warranty

8.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the scanner.
- ▷ Restart the PC.
- ▷ Check the cable connections; disconnect and then re-connect, as required.

8.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

8.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

10 Technical Data

Subject to changes.

	Unit	Desktop scanner
Art. No.	–	179560N
Dimensions (D × W × H)	mm	415 × 424 × 469
Weight	kg	25
Electrical connections		
_ Voltage	V/Hz	100-240/50-60
_ Fuse	A	2 × 1.6
Power output	W	60
Accuracy	µm	< 4
Allowable temperature range (operation)	°C	18 - 30
Interfaces	–	USB 2.0

Tab. 3

9 Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the scanner via the public disposal systems in an environmentally-friendly manner.

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Amann Girrbach AG