# therm 3

<b>DE</b> Betriebsanleitung	3 - 19
<b>EN</b> User Manual	20 - 35
<b>FR</b> Mode d'emploi	36 - 51
<b>IT</b> Istruzioni d'uso	52 - 68
<b>ES</b> Modo de empleo	69 - 84
<b>PT</b> Instruções de operação	85 - 101

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Tab. 4

## Table of Contents

<b>1</b>	<b>Explanation of Symbols</b> .....	<b>21</b>	<b>9</b>	<b>Environmental Protection</b> .....	<b>35</b>
<b>2</b>	<b>General Safety Instructions</b> .....	<b>22</b>	<b>10</b>	<b>Technical Data and Accessories/Spare Parts</b> .....	<b>35</b>
<b>3</b>	<b>Suitable Personnel</b> .....	<b>22</b>			
<b>4</b>	<b>Machine Specifications</b> .....	<b>23</b>			
	4.1 Delivery Scope .....	23			
	4.2 Intended Use .....	23			
	4.3 CE Declaration of Conformity .....	23			
	4.4 Components and Interfaces .....	24			
	4.5 Program overview, shortcut keys ..	26			
<b>5</b>	<b>Installation</b> .....	<b>28</b>			
	5.1 Setting up .....	28			
	5.2 Inserting a floor plate .....	29			
	5.3 Electrical connection .....	29			
<b>6</b>	<b>Application and Operation</b> .....	<b>30</b>			
	6.1 Selecting the sintering program ..	30			
	6.2 Programming the sintering program .....	31			
	6.3 Functions .....	32			
	6.4 Operation .....	32			
<b>7</b>	<b>Cleaning and Maintenance</b> .....	<b>33</b>			
	7.1 Cleaning .....	33			
	7.2 Daily checks .....	33			
	7.3 Monthly maintenance .....	33			
<b>8</b>	<b>Malfunctions, Repairs and Warranty</b> ..	<b>34</b>			
	8.1 Malfunctions .....	34			
	8.2 Repairs .....	34			
	8.3 Warranty .....	34			

## 1 Explanation of Symbols

### Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

When there is a risk of burns from hot surfaces, the exclamation mark in the warning triangle is substituted by a heat-alert symbol.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that life-threatening personal injury can occur.

### Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

### Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

## 2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

### DANGER:

This product does not comply with the ATEX Directive.

- ▷ Do not use in explosive atmospheres!
- ▷ Do not operate with explosive gases or mixtures!

### WARNING:

Fire hazard!

- ▷ In case of unexpected processes in the furnace (for example, strong smoke development or odour nuisance), switch the furnace off immediately and allow the furnace to cool down to room temperature.

### CAUTION:

A defective machine can lead to malfunctions! On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

### WARNING:

Operation with power sources, products, consumables, auxiliary aids, etc., which are subject to the Ordinance on Hazardous Substances or in any way have an impact on the health of the operating personnel, is not permitted.

### NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

### NOTE:

- Increased wear of door gasket, insulation, heating rods and furnace housing!
- ▷ Avoid opening the furnace when in hot condition (above 200 °C)!

## 3 Suitable Personnel

### NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

## 4 Machine Specifications

### 4.1 Delivery Scope

- \_ High-temperature furnace
- \_ Heating rod (4 x)
- \_ Exhaust air pipe
- \_ Floor plate
- \_ Sintering bowl, stackable
- \_ Sintering pearls
- \_ Tongs for sintering bowl
- \_ Mains cable
- \_ Allen key

▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

### 4.2 Intended Use

High-temperature furnace for the final sintering of dental oxide ceramics, particularly for zirconium oxide blanks.

The sintering furnace is pre-programmed by the manufacturer with 3 temperature programs, which are matched to the zirconium oxide frame material of Amann Girrbach AG.

247 further programming positions are freely programmable.

The stackable sintering bowl, filled with sintering pearls, is used for sintering the zirconium oxide frames. Any other or further-going use, such as the processing of other products than those intended for, as well as the handling of hazardous substances or materials hazardous to one's health, does not apply as under the intended use.

As to the materials being used in the furnace, it must be known whether they are aggressive to or even could lead to destruction of the insulation or of the heating rods.

It is prohibited to remove safety devices, bypass them or put them out of operation.

Unauthorised modifications/alterations of the machine shall void the warranty.

The set-up instructions and the safety regulations must be followed, as otherwise the furnace is considered as not being used as intended for, and any and all claims against the manufacturer shall become void.

### 4.3 CE Declaration of Conformity

This product was designed and manufactured based on careful selection of the harmonised standards to be observed, as well as additional technical specifications. It thus corresponds with the state-of-the-art and ensures maximum safety.

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

#### EC Guidelines

2006/95/EC	2004/108/EC
------------	-------------

#### Applied harmonised standards

DIN EN ISO 12100	DIN EN ISO 13732-1
DIN EN 61010-1	DIN EN 61000-6-1, DIN EN 61000-6-3

Tab. 1

The product's Declaration of Conformity can also be viewed on the manufacturer's website.

## 4.4 Components and Interfaces

EN


Fig. 1 Machine overview, front

- |                     |                     |
|---------------------|---------------------|
| 1 Exhaust air pipe  | 6 Shroud insulation |
| 2 Furnace lift door | 7 SD card slot      |
| 3 Thermoelement     | 8 USB port          |
| 4 Heating rod       | 9 Controller        |
| 5 Floor plate       |                     |


Fig. 2 Machine overview, rear

- 1 Mains switch
- 2 Power supply connection
- 3 Fuse
- 4 Switchgear fan

## Display, operating level


Fig. 3 Display, operating level

- P1...P4** Fast keys, program 1 ... 4
- Program start key
- Program stop key
- Key for calling up the configuration level
- 1 Temperature curve
- 2 Current temperature
- 3 Current program status
- 4 Current program
- 5 Indication of the program run-time
- 6 Remaining time until program ends


Fig. 4 Display, configuration level

**Legend to Fig. 4, Fig. 5 and Fig. 6:**

-  "Program" key
-  "Function" key
-  "Selection" key
-  "Info" key
-  "Save" key
-  "Move to list up" key
-  "Move to list down" key
-  "Delete character to left" key
-  "Move on digit to right" key
-  "Move on digit to left" key
-  "Delete program name" key


Fig. 5 Display, list menu

**4.5 Program overview, shortcut keys**

The high-temperature furnace has 3 permanently set programming positions for sintering zirconium oxide frames of varying indication, as well as 247 freely programmable programming positions for sintering dental oxide ceramics of other manufacturers (observe manufacturer's information).

**Program description**

- \_ Program 1 [P1] - standard
  - Standard program for zirconium oxides
  - Universal sintering program for all common indications




Fig. 6 Display, entering a program name


\_ Program 2 [P2] - fast

- Quick-sintering program for zirconium oxides - single units.
- Quick-sintering program only for single units - Max. furnace filling load is one sintering bowl; otherwise, sufficient through-sintering of the objects cannot be achieved.


\_ Program 4 [P4] - default

- Parameter freely assignable
- The program is preset to 1500 °C, 2 hours dwell time and a heating rate of 10 °C/min (not material-specific).

\_ Program 3 [P3] - slow

- Sintering program with long-term cooling for zirconium oxides.
- Sintering program with dwell phase at 900 °C and defined cooling. Particularly for large and massive frames on which the connected remaining blank is also sintered.

When the maximal temperature of 1530 °C is set, a maximal heating rate of 8 °C/min is recommended, in order to prevent ageing of the heating rods. Continuous operation of the furnace at 1530 °C leads to premature ageing and glazing of the heating rods!

## 5 Installation

### 5.1 Setting up

#### NOTE:

Due to the unit's weight of 30 kg, two persons are recommended for setting up the unit.

► When setting up, observe the valid country-specific accident prevention regulations.

- \_ The machine is intended exclusively for use within dry, closed rooms.
- \_ The machine must be positioned below an exhaust pipe or an exhaust hood (as used for pre-heating furnaces).

#### NOTE:

Malfunction!

► Do not install an exhaust pipe directly near the exhaust air pipe of the furnace!

- \_ Sufficient ventilation must be ensured at the set-up location, in order to draw off the heat and possibly developing exhaust gas.
- \_ A clearance of at least 50 cm is to be observed sideways, to the rear and upwards of the furnace. The lateral clearance to non-flammable materials can be reduced to 20 cm. Space requirements for the high-temperature furnace (incl. connections):
  - Width: 1.4 m (0.8 m for non-flammable materials)
  - Depth: 0.9 m
  - Height: 1.0 m
- \_ The set-up surface must be level.
- \_ The empty weight of the machine is 30 kg. The set-up surface must withstand an appropriate load.

- \_ During operation, the room temperature must be between 5 °C and 40 °C; avoid large variations in temperature.
  - \_ The maximal operating elevation (standard elevation zero) is 2000 m.
  - \_ The floor covering must be made of non-flammable material, so that hot material falling out of the furnace can not inflame the covering.
- Set up the machine on a sturdy work table or work bench (on a non-flammable surface).
- Remove all packaging material, including from the furnace compartment.

#### DANGER:

For carrying out the further installation steps, a qualified professional in the area of electrical safety is mandatory.

The installation and repair of heating elements must be performed by trained and certified personnel.

### 5.2 Inserting a floor plate

▷ Insert the floor plate.


Fig. 7

### 5.3 Electrical connection

The electrical connection is to be made via a socket outlet with earthing contact providing 220 - 240 V and 50 - 60 Hz. This socket outlet must be located close to the furnace and be easily accessible.

When using an extension cord or a multiple socket outlet, do not exceed the respective max. electrical load rating.

**WARNING:**

The mains cable may not be damaged!

- ▷ Do not place objects on the mains cable.
- ▷ Lay the mains cable in such a manner that nobody can step on or trip over it.
- ▷ Pay attention that the cable does not have contact with hot objects or surfaces.
- ▷ In case of damage to the insulation: Replace the mains cable.

**NOTE:**

Do not insert the mains plug before the heating elements have been inserted.

- ▷ Set the mains switch to the **0** position.
- ▷ Connect the enclosed mains cable to the socket intended for this purpose on the rear of the furnace.


Fig. 8

- ▷ Connect the other end of the cable to the socket outlet with earthing contact.

The mains cable may only be replaced against an original mains cable from the manufacturer.

## 6 Application and Operation

Before starting-up the furnace, it should acclimate at the set-up location for 24 hours. Swift heating up in cold condition can cause damage to the heating elements (for example, burning-through).

During the first heating up, there may be odour nuisance due to binding agent fumes escaping from the insulation material.

- ▷ Provide for thorough ventilation of the furnace set-up location when heating up for the first time.

### WARNING:

Danger of fire and explosion!

The furnace surfaces become very warm.

- ▷ Do not place objects on the furnace!

### NOTE:

Insufficient sintering results due to faulty temperature control!

- ▷ Prior to every process, visibly check the thermoelement for damage.

### DANGER:

Hot air at the outlet opening!

- ▷ Do not reach with your hand at or above the outlet opening.
- ▷ Do not place objects on the furnace!
- ▷ Do not hang objects directly above the furnace!

### 6.1 Selecting the sintering program

- ▷ To select a program: Press the **P** key.  
The list menu is indicated.

Fig. 9 Display, list menu

- ▷ In the list, select the desired program with the keys

-or-

- ▷ directly select the program by pressing the **P** key and the respective number key.
- ▷ To take over the program: Press the key.  
The display switches to the operating level.
- ▷ To display the segments of the program: Press the key.

## 6.2 Programming the sintering program

Information/tips on programming / controller behaviour:

- \_ The input options of the heating rate is limited from 60 °C/h (min.) to 900 °C/h (max.).
- \_ The temperature input is limited from 25 °C (min.) to 1530 °C (max.).
- \_ In case of incomplete entry of parameters within a segment (e.g. heating rate 0 °C/h), either the complete program will be discarded and not taken over and the standard program P1 is indicated, or, in case of an incomplete programmed sequence segment, this part of the program will always be discarded (all following parameters are set to 0.)
- \_ The remaining-time indication is put out rounded exactly to 10 minutes.
- \_ Depending on the furnace filling, the indication of the program status (red dot) can lightly deviate from the shown curve shape in the range of the cooling curve.

A sintering program consists of up to 8 segments. Each segment contains 3 parameters:

- \_ Heating rate in °C/h (max. 60 °C/h)
- \_ Temperature in °C (max. 1530 °C/h)
- \_ Dwell time in minutes

In order to set a program (4 to 250):

- ▷ Press the **P** key.
- ▷ Press the number key for the respective program.

**-or-**

- ▷ In the list, select the desired free program location with the **▲▼** keys.
- ▷ To enter the segment definition: Press the **→** key.
- ▷ Enter the first parameter of the segment (heating rate in °C/h) via the die number key.

- ▷ Press the **→** key.
- ▷ Enter the second parameter of the segment (temperature in °C) via the die number key.
- ▷ Press the **→** key.
- ▷ Enter the third parameter of the segment (dwell time in minutes) via the die number key.
- ▷ Press the **→** key.
- ▷ Enter further segments in the same manner.
- ▷ To save/store the program: Press the **☐** key.
- ▷ Assign a name to the program using the alphanumeric keys.
  - Navigate in the title with the **◀▶** keys.
  - delete individual characters with the **←** key
  - delete the complete name with the **CLR** key


Fig. 10 Display, entering a program name

- ▷ To save/store the program: Press the **☐** key. The program is taken over, the display changes to the operator level.

### 6.3 Functions

#### Handling

- ▷ To select a function: Press the **F** key and the respective number key.
  - Key 4: Function “Preset time for program”

#### Programming

To program function 4:

- ▷ Press the **F** key.
- ▷ Press the **4** key.
- ▷ Using the number keys, enter the desired time [hh:mm] until the program starts. The time is rounded exactly to 10 minutes.

After entering the last digit, the cursor will jump back to the first numerical location, so that possibly required corrections can still be made.

- ▷ Press the **→** key.
- ▷ Press the **P** key.
- ▷ In the list, select the desired program with the **▲** **▼** keys

**-or-**

- ▷ directly select the program by pressing the **P** key and the respective number key.
- ▷ Press the **⏏** key.
 

The program is taken over with the time preset. The display switches to the operating level. After the program starts, the preset time runs off first, then the program starts.

### 6.4 Operation

#### NOTE:

Damage to the heating rods!

- ▷ When loading and unloading the furnace, do not bump against or come in contact with the heating rods!

- ▷ When sintering zirconium frames, use a sintering bowl filled with sintering pearls.
- ▷ Using the sintering bowl tongs, carefully place the sintering bowl into the furnace (up to 3 bowls are stackable).

The placement of zirconium frames in the sintering bowl is described in the instructions of the zirconium-oxide blanks.

- ▷ Load the furnace with the units.
- ▷ Shut the furnace door.
- ▷ Switch the furnace on with the mains switch.
- ▷ Select the sintering program (see page 30).
- ▷ Define the sintering program, if required (see page 31).
- ▷ Start the sintering program with the **▶** key.
 

After the program has elapsed, the heating cycle automatically switches off.

#### WARNING:

Danger of burning!

- ▷ To remove the sintering bowl, use the sintering-bowl prongs or other suitable pliers.

## 7 Cleaning and Maintenance

### 7.1 Cleaning

**DANGER:**

 Risk of electric shock!

▷ Prior to any cleaning, pull the mains plug.

**NOTE:**

 Possible damage to the furnace!

▷ Do not use a high-pressure cleaner or compressed air for the cleaning!

▷ Use only commercially available, water-based or non-flammable, solvent-free cleaning agents for the cleaning.

Prior to cleaning:

- ▷ Allow the furnace to cool down to the room temperature.
- ▷ Clean the various components according to Table 2.

Component	Cleaning
Housing (exterior)	Wipe off with moist, lint-free cloth
Interior	Carefully vacuum out with vacuum cleaner
Insulation materials	Carefully vacuum out with vacuum cleaner
Door gasket	Wipe off with moist, lint-free cloth
Instrument panel	Wipe off with moist, lint-free cloth

Tab. 2

### 7.2 Daily checks

- ▷ Check the furnace for any exterior damage.
- ▷ Check the function of the safety devices (e.g. the mains switch).
- ▷ Check the function of the fan (check for noise).
- ▷ Visually check the following components:
  - Sealing surfaces on door and furnace shroud
  - Visible part of the heating rods in the furnace compartment
  - Visible part of the thermoelement in the furnace compartment

### 7.3 Monthly maintenance

- ▷ Carefully vacuum out and clean the furnace compartment, exhaust openings and exhaust pipe.
- ▷ Check the furnace compartment, exhaust openings and exhaust pipe for damage.

## 8 Malfunctions, Repairs and Warranty

### Cracks in the insulation

The furnace insulation consists of a high-grade, fire-resistant material. Due to the thermal expansion, cracks develop in the insulation after only a few heating cycles. However, these do not have an influence on the function or quality of the furnace.

### 8.1 Malfunctions

▷ In case of malfunctions, contact the manufacturer's Technical Service.

### 8.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

#### Replacing a fuse

- ▷ Switch the furnace off and pull the mains plug.
- ▷ Wait until the furnace has cooled to room temperature.
- ▷ On the rear side of the furnace, turn the fuse holder in anticlockwise direction and pull it out (bayonet lock).
- ▷ Replace the defective fuse against an identical one (see Page 35).
- ▷ Insert fuse holder with new fuse in furnace rear wall and lock it by twisting in clockwise direction.

### 8.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

In addition, it also applies that warranty and liability claims are excluded for personal and property damage, when these can be attributed to one or more of the following causes:

- \_ Each person involved with operating, mounting, maintaining, repairing the system must have read and understood the User Manual. No liability shall be assumed whatsoever for damages and operating malfunctions resulting from non-observance of the User Manual.
- \_ Use of the system not in accordance with the intended use.
- \_ Improper mounting, starting-up, operation and maintenance of the system.
- \_ Operation of the system with defective safety devices or improperly mounted or inoperable safety and protective devices.
- \_ Non-observance of the notes in the User Manual in terms of transport, storage, assembly, starting-up, operation, maintenance and equipping the system.
- \_ Unauthorised modifications to the system.
- \_ Unauthorised modification of the operating parameters.
- \_ Original parts and accessories are designed especially for this furnace. When exchanging components, use only original parts. Otherwise, the warranty shall become void. For damages caused by the use of non-original parts, the manufacturer shall exclude any and all liability.
- \_ Disaster cases caused through foreign bodies and force majeure.

 Each furnace is subject to a factory test run on behalf of the manufacturer. The test run may cause discolouration to the furnace compartment or sheet-metal parts.

## 10 Technical Data and Accessories/Spare Parts

 Subject to changes.

## 9 Environmental Protection

### Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

### Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

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