 **ceramill<sup>®</sup> zirconia**

 **ceramill<sup>®</sup> zolid fx**

 **ceramill<sup>®</sup> zolid fx**  
multilayer

<b>DE</b> Gebrauchsanweisung	3 - 26
<b>EN</b> Instruction Manual	27 - 49
<b>FR</b> Instructions d'utilisation	50 - 75
<b>IT</b> Istruzioni per l'uso	76 - 100
<b>ES</b> Instrucciones de uso	101 - 123

**AMANN GIRRBACH**

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

## Table of Contents

Explanation of Symbols .....	28
General Safety Instructions .....	31
Suitable Personnel .....	31
Properties .....	32
Warranty/Exclusion of Liability .....	34
Testing the Blanks .....	34
Mounting to the Worktable .....	35
Application .....	35
Process steps in the dental laboratory .....	36
Process steps at the dentist's .....	46
Environmental Protection .....	49
Download information .....	49

EN

## Explanation of Symbols

### Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- \_ **NOTE** means that property damage can occur.
- \_ **CAUTION** means that light to fairly serious personal injury can occur.
- \_ **WARNING** means that serious personal injury can occur.
- \_ **DANGER** means that life-threatening personal injury can occur.

### Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.



**Other symbols in the Manual**

<b>Symbol</b>	<b>Meaning</b>
▷	Item of an operation description
–	Item of a list
▪	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

## Other symbols on the product

Symbol	Meaning
	Medical device
	Unique Device Identifier
	Catalog number
	Batch code
	Manufacturer
	Country of manufacture
	Consult instructions for use
	Use by YYYY-MM-DD or YYYY-MM
<b>Rx only</b>	Symbol that may be used in place of statement "CAUTION". US Federal law restricts this device to sale by or on the order of a dentist.
	Serial number
	Store in a dry place

## General Safety Instructions

### CAUTION:

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

### NOTE:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Suitable Personnel

### NOTE:

The product may only be processed by professionally trained dental technicians and dentists.

## Properties

### Intended Use

Zirconium-oxide blanks for permanent and partly removable dental prosthetics.

### Product description

The blanks are made of zirconium oxide ( $ZrO_2$ ) for type II, class 4 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing permanent and partly removable prosthetic restorations (e.g. crowns and three-unit bridges to the molar region as well as implant-supported restorations) using CNC milling machines (e.g. Ceramill milling machines).

After completion of the specified end-sintering, the final product meets the requirements of DIN EN ISO 6872.

Intended patient group: Suitable for patients of all ages and gender.

Clinical benefit of the product:

- \_ Restoration of missing tooth and gingiva portions
- \_ Restoration of chewing function and aesthetics

### Safety data sheet/Declaration of conformity/SSCP

The safety data sheet and the SSCP of the products (according to the requirements of the Medical Device Regulation (MDR)) can be downloaded from the Amann Girrbach website under [Services > Downloads > Additional documents](#). The declaration of conformity of the products are available upon request from the manufacturer.

**Technical data**

	Unit	Ceramill Zolid FX Ceramill Zolid FX ML
Bending strength	MPa	≥ 500
E-module	GPa	≥ 200
Thermal expansion coefficient (CTE) (25 - 500°C)	10 <sup>-6</sup> /K	10.1 ± 0.5
Chemical solubility	μg/cm <sup>2</sup>	< 100
Vickers hardness	HV10	1300 ± 200

The delivery scope can vary country-specific.

**Chemical composition**

Oxide	Mass percentage	
	Ceramill Zolid FX	Ceramill Zolid FX ML
ZrO <sub>2</sub> + HfO <sub>2</sub> + Y <sub>2</sub> O <sub>3</sub>	≥ 99.0	
Y <sub>2</sub> O <sub>3</sub>	9.15 - 9.55	8.5 - 9.5
HfO <sub>2</sub>	≤ 5	
Al <sub>2</sub> O <sub>3</sub>	≤ 0.5	
Other oxides	≤ 1	

**Shelf life of the blanks**

When appropriately stored, the blanks can be used 5 years after the manufacturing date.

**Storage**

Store the blanks in the original packaging in a dry location.

**Frequency of use**

As an intermediate product, the respective product is not intended for single use. The final product (dental prosthesis) is patient-specific.

**Warranty/Exclusion of Liability**

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development and are tested according to the state of the art and to legal requirements. We reserve the right to make any resulting changes in handling and composition.

The respectively valid instruction-manual version can be found under [www.amangirrbach.com/instruction-manuals](http://www.amangirrbach.com/instruction-manuals). This version replaces all previous versions.

**Testing the Blanks**

Upon receipt, it is imperative to visually check the proper condition of the blanks. After using a (transport-)damaged blank, complaint claims are no longer possible.

## Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

## Application

### Indication ranges

- \_ Anatomically reduced and fully anatomical (monolithic) crowns in the anterior and posterior tooth range (e.g. single-tooth crowns, inlays, onlays, veneers).
- \_ **Up to three-unit** anatomically reduced and fully anatomical (monolithic) bridges with a maximum of one pontic in the anterior and posterior tooth range.
- \_ Implant-supported hybrid dentures (e.g. hybrid abutments and hybrid abutment crowns).

### Contraindications

- \_ Insufficient tooth-structure availability
- \_ Insufficient preparation
- \_ Insufficient oral hygiene
- \_ More than one bridge unit
- \_ Known incompatibilities with respect to the components
- \_ Heavily discoloured hard tooth structure
- \_ Provisional insertion
- \_ All indications that are not listed under "Indications"

### Undesirable side effects/residual risks

- \_ Mechanical failure (fracture of the restoration, chipping of the ceramic veneer or detachment of the restoration/loss of retention) with low risk of possible ingestion or reversible mucosal injuries.
- \_ Biological complications (plaque accumulation, intolerance reactions, marginal/secondary caries, loosening of the supporting teeth) resulting in possible loss of the restoration.

## Process steps in the dental laboratory

### Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating zirconium oxide frames in high-density sintered condition:

Minimum frame thickness in mm	Connector cross-section in mm <sup>2</sup>		Maximum number of connected bridge units	
	Anterior	Posterior	Anterior	Posterior
0.5	≥ 12	≥ 12	1	1

## Detailed minimum wall thicknesses and connector cross-sections

Indication	Pattern	Total number of units	Number of connected bridge units	Wall thickness in mm		Connector cross-section in mm <sup>2</sup>
				Incisal/occlusal	Circular	
Primary components/Double crowns	–	1	–	0.7	0.5	–
Single crown	–	1	–	0.5	0.5	–
Front-tooth bridge	OXO	3	1	0.5	0.5	≥ 12
Posterior bridge	OXO	3	1	0.7	0.5	≥ 12

O Abutment crown  
X Bridge unit


## Crown border parameters in the CAD software

Pay attention to sufficient crown border thicknesses Recommended values:

- \_ 1. Crown border thickness (Horizontal): 0.2 mm
- \_ 2. Angled: 0.3 mm
- \_ 3. Angle: 60°
- \_ 4. Height (Vertical): 0 mm

For tangential preparations, increase the “crown border thickness” and “Angled” parameters, e.g.:

- \_ 1. Crown border thickness (Horizontal): 0.25 - 0.3 mm
- \_ 2. Angled: 0.35 - 0.4 mm


Observe zirconium-oxide-efficient design forms and do not create pointed edges.


## Positioning in the CAM software

- ▷ Apply enough sufficiently thick holding ridges  
(at least 3 ridges for a single crown; thickness  $> 1.7$  mm).


- ▷ Position holding ridges approx. 1 - 2 mm away from the crown margin, especially for long crown margins.

This provides support to the crown margin.

- ▷ Position the crown in the blank in such a manner that the clearance from the crown margin to the blank surface is at least 1 mm.

This ensures that the crown has no contact to the sinter skin of the blank.

## Procedure when fabricating dental framework

For maximum fitting accuracy, the specific expansion/shrinkage value is given on the blanks.

When machining the blanks in milling machines, it is recommended to only use cutters with sufficient cutting edge quality.

- ▷ Depending on requirement in the respective CAM software, enter the appropriate value.  
In case of uncertainties, contact your CAM supplier.
- ▷ Mill the blanks with the milling machine applying the milling strategies provided for the material.

## Removal from the blank

Utmost care must be taken when removing the frames from the blank:

- ▷ Separate the frames by grinding them out using a turbine grinder (without water-cooling) and a pointed diamond bit (do not use cutting discs).

For trimming and dressing the connectors and smoothing out the frames, we recommend using a polishing set.

- ▷ Clean the frames from residual zirconia dust using a brush or compressed air.

### **Colouration of the white zirconium oxide**

- ▷ Depending on requirement, partially or completely stain the white zirconium oxide frames prior to the sintering with dye solution (e.g. Ceramill Liquid FX).

### **Colouration of the pre-stained polychrome zirconium oxide**

Ceramill Zolid FX Multilayer is available in 16 A-D Vita shades including bleach shades. The polychrome coloring with integrated color gradient integrated colour gradient gives the blanks a natural appearance immediately after sintering.


- ▷ Depending on requirement, partially stain the pre-colored polychrome zirconia restorations (e.g. Ceramill Liquid FX).
- ▷ After successful final sintering: color with stains or glaze materials or use as a framework for the veneering technique.

### **Shade selection and placement of the restoration in the blank**

The tooth shade of the blank covers off the imprinted area (e. g. 0/A1).

When a higher blank is selected, the respective lighter or darker shade can be achieved more significantly by vertically moving in the nesting. The closer the work is positioned to the stained side, the larger the brighter shade component.


*Fig. 3 Placement of the restoration in the blank*

**VS1** Tooth shade 1 (lighter)

**VS2** Tooth shade 2 (darker)

*Fig. 4 Placement of the restoration in the blank (block form)*

**VS1** Tooth shade 1 (lighter)

**VS2** Tooth shade 2 (darker)

When using the block form, the shade concept can deviate for larger front and posterior teeth, as these cannot be moved up or down in the zirconium block. Depending on the desired shade, the next lighter or darker zirconium block must thus be selected.

By applying the stain, further tooth shades can be achieved.

## Final sintering

A high-temperature sintering furnace, e.g. Ceramill Therm or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

Highly translucent zirconium oxide and opaque zirconium oxide can be sintered together in one sintering process.

### NOTE:

When highly translucent zirconium oxide and opaque zirconium oxide are sintered together in the same sintering bowl and with the same sintering pearls, then this can lead to changes in the translucency. Therefore, it is advisable for highly translucent zirconium oxide and opaque zirconium oxide, to use individual sintering pearls in separate sintering bowls.

- ▷ Place the frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

### NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

- ▷ Perform the final sintering of the frames according to the following program:
- Heat-up phase: Room temperature to final temperature 1450 °C; heating rate 5 - 10 K/min
  - Dwell time at final temperature: 2 hrs
  - Cooling phase: Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

## Frames for veneering

### Post-processing

If required, the zirconium oxide frames can be reworked after the final sintering without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 µm).

### Frame preparation for veneering

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 50 µm, pressure ≤ 2 bar) and subsequent steaming off.

### Veneer ceramics

For information on the thermal-expansion coefficient of the blanks, please refer to the Technical data (see page 33).

## Monolithic dentures

### Post-processing

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40  $\mu\text{m}$ ).
- ▷ Polish the contact surfaces to the antagonist and to the neighbouring teeth to a high gloss when in high-density sintered condition.

### Preparation for staining and glazing

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 50  $\mu\text{m}$ , pressure  $\leq$  2 bar) and subsequent steaming off.

The exterior crown surfaces are not corundum blasted.

### Staining and glazing

- ▷ Stain and glaze the units afterwards.

For bridges with massive intermediate units, we recommend step-by-step long-term cooling to 500 °C.

## Process steps at the dentist's

### Processing

All prosthetic components must be cleaned, disinfected and/or sterilised prior to each use (see "Amann Girrbach guideline for cleaning, disinfecting and sterilising prosthetic components").

### Fixation

Due to the high translucency of this zirconium oxide, the selection of cement shade and possible discolourations of the tooth stump are to be observed, as these may affect the final colour result.

Ceramill Zolid zirconium oxide restorations can be cemented conventionally with zinc phosphate cement or glass ionomer cement (e.g. GC Fuji Plus, GC) or adhesively using cement composites (Panavia<sup>®</sup> 21 or Panavia<sup>®</sup> F 2.0, Kuraray Noritake).

▷ With conventional cementing, ensure sufficient retention and a minimum stump height of 4 mm!

Temporary attachment is not recommended as scaffolding may be damaged upon removal.

Multilink Hybrid Abutment (Ivoclar Vivadent) is recommended for the extraoral attachment of a two-part abutment made of Ceramill Zolid zirconium oxide on a titanium adhesive base (hybrid abutment or a hybrid abutment crown). For attachment, follow the instructions of the manufacturer of the titanium adhesive base and the manufacturer of the attachment material.

### Post-processing at the dentist's

When grinding is required by the dentist while installing the work into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding points (recommended grain approx. 40  $\mu\text{m}$ ). For polishing, we recommend using diamond polishing bits.

#### NOTE:

When the work is not sufficiently polished, the antagonist may become damaged due to abrasion!

### Follow-up examination at the dentist's

It is urgently recommended to check the dentures in the patient's mouth once a year. In this, the remaining dentition, the antagonists and the soft tissue are also to be inspected. Corrective measures are to be undertaken as required. Here, it is also important that the objects are polished afterwards to a high gloss.

## Information to be passed on to the patient

(e.g. by healthcare professionals).

Possible complications and residual risks for zirconium-oxide-based restorations in general and not specific to Amann Girrbach restorative materials are:

- \_ Mechanical failure (fracture of the restoration, chipping of the ceramic veneer or detachment of the restoration/loss of retention) with low risk of possible ingestion or reversible mucosal injuries.
- \_ Biological complications (plaque accumulation, intolerance reactions, marginal/secondary caries, loosening of the supporting teeth) resulting in possible loss of the restoration.

To reduce possible complications/residual risks, it is recommended to have the denture checked by the dentist during the regular dental check-up (once a year).

## Environmental Protection

### Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

### Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

## Download information

Further instructions/manuals can be downloaded under [www.amanngirrbach.com/instruction-manuals](http://www.amanngirrbach.com/instruction-manuals).

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**Amann Girrbach AG**
**Amann Girrbach GmbH**