 **giroinvest super**

Gebrauchsanweisung
User Manual
Mode d'emploi
Istruzioni d'uso
Modo de empleo


AMANNGIRRBACH

Deutsch	04 - 09
English	10 - 15
Française	16 - 21
Italiano	22 - 27
Español	28 - 33

Thermische Expansion Giroinvest Super

Thermal Expansion Giroinvest Super

Expansion thermique Giroinvest Super

Espansione termica Giroinvest Super

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GIROINVEST SUPER

Giroinvest Super is a phosphate-bonded universal investment for crown, bridge and inlay work and pressable porcelain. Its very broad expansion range (up to 4.0% overall expansion) allows all precious and non-precious metal alloys (e.g. Girobond NB 781600) to be cast accurately by precisely regulating the expansion. Giroinvest Super produces very smooth casting surfaces in the speed cast technique. Giroinvest Super is suitable to be cast either with casting ring or ringless.

Giroinvest Super:

Order no.	
781670	Pack contents 40 x 150 g bag (= 6 kg)
781685	Pack contents 2 x 4 kg bag (= 8 kg)
781680	Pack contents 50 x 100 g bag (= 5 kg)

Giroinvest Super Liquid:

Order no.	
724090	Pack contents 1l

Technical Data:

Setting expansion	2.8%
Thermal expansion	1.2%
Overall expansion	4.0%

Figures relate to 100% Giroinvest Super Liquid

Mixing ratio

150 g Giroinvest Super:	37 ml liquid
100 g Giroinvest Super:	25 ml liquid

Working time	5 min
--------------	-------

Mixing:

30 s manual premix,
15 s evacuate without mixing,
60 s mix under vacuum

Setting time	30-35 min
--------------	-----------

Times relate to 20°C working temperature

Material data in accordance with

DIN EN ISO 9694:

Initial setting	9-11 min
Compression strength	
after 2 h	4-8 MPa

WARNINGS AND SAFETY INSTRUCTIONS:

- _ This investment contains quartz and cristobalite. May cause damage to the lung through prolonged or repeated exposure. Route of exposure: Inhalative. Do not breathe dust. In case of inadequate ventilation wear respiratory protection. Get medical advice / attention if you feel unwell.


Recommendations:

- _ Use a protective face mask, type FFP 2-EN 149:2001.
- _ Open the foil bag with scissors and avoid creating dust when filling the mixing bowl.
- _ Rinse the empty foil bag with water before disposal.
- _ Always use a damp cloth to remove investment dust or dried mixing liquid from the working area.
- _ Soak the cooled casting ring in water before devesting the casting.
- _ Do not open the preheat furnace during the heating stage, as the wax vapours given off can ignite with air! Furnace gases released should be extracted and directed outside.
- _ Speed cast technique: Place all casting rings quickly into the furnace and do not open the furnace for 15 min!

Warning

INSTRUCTIONS FOR PRECIOUS AND NON-PRECIOUS METAL ALLOYS

Preparation:

- _ Position the wax pattern at least 5 mm away from the ring and 10 mm below the upper rim of the casting ring.
- _ Line metal rings with a 1mm thick ring liner without overlapping the edges. Ensure the liner is flush with the upper rim of the casting ring.

Casting ring size 1 and 3: 1 ring liner

Casting ring size 6: 2 ring liners

Mixing (Mixing ratio 150 g Giroinvest

Super: 37 ml liquid):

- _ A higher ratio of Giroinvest Super Liquid increases the overall expansion.
- _ Always use distilled water for diluting Giroinvest Super Liquid mixing concentrate.

Liquid concentration:

	Ratio of Giroinvest Super Liquid / Ratio of distilled water			
	Inlays/Onlays Partial crowns	Crowns Telescope crowns (Primary parts)	Bridges	Secondary units GC Pattern Resin
High-gold-content yellow gold (gold content > 70%)	50% / 50% 18,5 ml / 18,5 ml	55% / 45% 20 ml / 17 ml	65% / 35% 24 ml / 13 ml	80% / 20% 30 ml / 7 ml
Semi-precious yellow gold (gold content < 55%)	53 / 47 20 / 17	57 / 43 21 / 16	68 / 32 25 / 12	83 / 17 31 / 6
High-gold-content, silver-coloured bonding alloys (gold content > 70%)	63 / 37 23 / 14	70 / 30 26 / 11	75 / 25 28 / 9	85 / 15 31 / 6
High-gold-content yellow gold bonding alloys (gold content > 80%)	67 / 33 25 / 12	78 / 22 29 / 8	80 / 20 30 / 7	87 / 13 32 / 5
Semi-precious bonding alloys	57 / 43 21 / 16	68 / 32 25 / 12	73 / 27 27 / 10	87 / 13 32 / 5
Pd-based alloys	55 / 45 20 / 17	63 / 37 23 / 14	73 / 27 27 / 10	85 / 15 31 / 6
Non-precious metal alloys	70 / 30 26 / 11	85-95 / 15-5 31-35 / 6-2	70-80 / 30-20 26-30 / 11-7	0° 95-100 / 5-0 35-37 / 2-0
				2°-6° 80-90 / 20-10 30-33 / 7-4

All ratios relate to 150 g Giroinvest Super: 37 ml liquid at 20°C working temperature. The ratios are approximate. Different techniques, accessories and equipment can alter the results. In case of the overnight technique, use a 10% lower liquid concentration.

Working time:

5 min (working temperature 20°C)

- _ Rinse a clean mixing bowl with water and wipe it before filling with investment. Dirty or dry mixing bowls reduce the liquid content of the investment.
- _ Pour in the liquid and then sprinkle in the powder.
- _ Manually mix the powder and liquid thoroughly with a spatula for 30s.
- _ Use a vacuum mixer for further mixing. (Set vacuum mixers with a variable motor speed to 400 rpm). Evacuate the mixture for 15s to degas fully before activating the mixing paddles. Mix under vacuum for 60s. Reducing the mixing time causes micropearls, a poor surface and an inaccurate fit.

Investing:

- _ Fill the casting ring using minimum vibration. Once the pattern is covered with investment, continue filling the casting ring without vibration.
- _ Recommendation: invest without pressure
- _ Adhere strictly to the setting time and temperature for placing the casting ring in the furnace with the speed cast technique.

Preheating:

- _ Grind the top of the mould after setting or roughen with knife.
- _ Always place the mould upright on a grooved furnace tray with the sprue facing downwards. Ensure adequate space (approx. 1cm) between the casting ring and furnace tray in furnaces with floor heating.

Setting time:

30-35 min (working temperature 20 °C)

	Preheating speed cast	Preheating conventional
Setting time	exactly 30-35 min after filling the mold	minimum 30-35 min after filling the mold
Heat rate	-	3-5 °C / min with linear furnace control
Insertion temperature	850 °C in case of different final/casting temperature, heat up or cool down to the final temperature after having placed the invested item into the furnace	Room temperature
1 st hold time at 290 °C	-	30-45 min depending on size and number of casting rings
2 nd hold time at 580 °C	-	30-45 min depending on size and number of casting rings
Hold time at final temperature (maximum 1080 °C)	35-60 min depending on size and number of casting rings	35-60 min depending on size and number of casting rings

- _ Adding casting rings later reduces the furnace temperature and increases the preheat time!
- _ Speed cast technique: Place all casting rings quickly into the furnace and do not open the furnace for 15 min!

GB

Preheating overnight:

- _ As overnight heating produces the full setting expansion, the Giroinvest Super Liquid concentration should be reduced by 10 % depending on the type of alloy and restoration.

After casting:

- _ Allow the casting ring to bench cool to hand heat after casting.
- _ Do not quench in water!
- _ The investment contains quartz and cristobalite. Do not inhale the dust, as there is a risk of damage to the lungs (silicosis, lung cancer)! (See Warnings and Safety Instructions)

INSTRUCTIONS FOR PRESSABLE PORCELAIN

Preparation:

- _ Coat the sprue former thinly with Vaseline so that it separates easily from the investment after setting.
- _ Always use materials that burn out without residue for contouring the patterns.
- _ Units per casting ring: 100g Casting ring - maximum 3 (similar) wax units; 200g Casting ring - maximum 5 (similar) wax units

Mixing (Mixing ratio 100g Giroinvest Super: 25 ml liquid):

- _ A higher ratio of Giroinvest Super Liquid increases the overall expansion.
- _ Always use distilled water for diluting Giroinvest Super Liquid mixing concentrate.

Liquid concentration:

Pattern	Giroinvest Super Liquid / distilled water
Crowns	65 % / 35 % 16 ml / 9 ml
Inlay 1 surface	40 / 60 10 / 15
Inlay 2-3 surfaces	60 / 40 15 / 10

All ratios relate to 100g Giroinvest Super: 25ml liquid at 20°C working temperature. The ratios are approximate. Different techniques, accessories and equipment can alter the results.

- _ Pour in the liquid and then sprinkle in the powder.
- _ Manually mix the powder and liquid thoroughly with a spatula for 30s.

_ Use a vacuum mixer for further mixing. (Set vacuum mixers with a variable motor speed to 400 rpm). Evacuate the mixture for 15 s to degas fully before activating the mixing paddles. Mix under vacuum for 60s. Reducing the mixing time causes micropearls, a poor surface and an inaccurate fit.

Investing:

- _ Fill the casting ring using minimum vibration. Once the pattern is covered with investment, continue filling the casting ring without vibration.
- _ 100g Casting rings: 16min setting time without pressure / 200 g Casting rings: 18min setting time without pressure.
- _ After setting, remove the sprue former, stabilising ring and paper sleeve and place the mould into the furnace preheated to 850 °C.

Preheating:

- _ Always place the mould upright on a grooved furnace tray with the press sprue facing downwards. Ensure adequate space (approx. 1cm) between the casting ring and furnace tray in furnaces with floor heating.

	Preheating speed procedure	
Casting ring size	100g Casting ring	200g Casting ring
Setting time	16 min after filling the mold	18 min after filling the mold
Hold time at final temperature	45 min	60 min

_ Adhere to the pressable porcelain manufacturer's instructions for use for further procedures.

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Amann Gırrbach AG
Amann Gırrbach GmbH