



User Guide

Set up the watch

ZEUS 6 PRO is suitable for Android system 5.1 and iOS 8.2 or above. Through the following QR code, you can first download and install FitCloudPro.



Android: Search for **FitCloudPro** in the app store to download.

IOS: Search **FitCloudPro** in App Store to

How To Charge Your Watch:

1) Please Charge The Watch For 1-2 Hours Before Use; **2)** If The Charger Connection Is Unstable, Please Reconnect It; **3)** Unable To

Turn On Or Charge, Please Use 5V/2A Or Higher To Charge For More Than 20 Minutes.

Operation Introduction

Swipe Left: Sports data, heart rate, sleep, customization; **Swipe Right:**




Shortcut menu; **Slide Up:** Show information content; **Slide Down:**

Brightness, Do Not Disturb, Settings and other modes

Bluetooth connection operation steps

App Searches For Bluetooth Device "ZEUS 6 PRO"

Function Introduction



Dial Mode: Press And Hold On The Watch Face Interface For Three Seconds To Switch Between Different Watch Faces.

Call Function: Make Calls, Talk Hands-Free, Hang Up, Adjust Volume.

Menu: 5 different UI styles, switch at will.

Information Push: The Watch Will Vibrate To Remind You Of Incoming Calls, Text Messages, Qq, Wechat, Twitter, Facebook, Whatsapp, And Skype Messages.





Weather: App Connection Displays Real-Time Weather Temperature And Type.

Music: Click To Control Music Playback On Your Phone.

Sports: Walking, Running, Aerobics And Other Modes.

Compass: can tell the direction



Heart Rate: Hold The Watch Close To Your Wrist And Wait Patiently

For 30 To 60 Seconds. The Normal Value Is 50 To 100 Times/Minute.

Blood Pressure: Hold The Watch Close To Your Wrist And Wait Patiently For 30-60 Seconds. Adult Diastolic Blood Pressure: 60-80 (MmHg), Adult Systolic Blood Pressure: 100-120 (MmHg).

Sleep: You Can View Sleep Duration, Deep Sleep, Light Sleep And Other Data.

Blood Oxygen: Hold The Watch Close To Your Wrist And Wait Patiently For 30 To 60 Seconds. The Normal Value Range Is Generally Between 95% And 100%.

Common Problem

Watch Cannot Be Turned On

Please Press And Hold The Power Button For More Than 3 Seconds Or The Battery May Be Too Low. Please Charge It In Time.

Bluetooth Is Not Connected Or Cannot Connect

- 1) Please Try To Restart The Watch And Reconnect.
- 2) Please Try To Restart Your Phone' s Bluetooth And Connect Again.
- 3) Do Not Connect Your Phone To Other Bluetooth Devices At The Same Time.

Measure Heart Rate/Blood Pressure/Blood Oxygen

- 1) Generally Caused By Poor Contact Between The Watch' s Sensor

And The Human Body During Measurement.

2) When Measuring, Please Make Sure That The Sensor Is In Full Contact With Your Wrist.

3) For People With Darker Skin And More Arm Hair, Please Turn On Enhanced Measurement In App [Device] > [Enhanced Measurement].

Sleep Data Is Not Accurate Enough

1) Monitoring Is To Simulate a Person's Natural Sleep And Wake-Up Time And Needs To Be Worn Normally.

2) If You Fall Asleep Too Late Or Wear It When You Fall Asleep, Errors May Occur.

3) Data Is Not Monitored During The Day, And The Default Monitoring Time Is From 9:30 Pm To 12:00 Noon The Next Day.

For More Frequently Asked Questions, Please View App [My] > [Faq].