INSTALLATION
BEDIENUNG
WARTUNG

ARTIKEL 8004/8005

NOBU
NOTRUFBUTTON

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INSTALLATION
OPERATION
MAINTENANCE

ARTICLE 8004/8005

NOBU
EMERGENCY CALL
BUTTON

Dear users,

You have chosen the IQcare NOBU emergency call button so that you are informed via the call system if your residents are in an emergency situation and need your help.

If handled properly, you can rely on this quality product at all times. Therefore, please observe the information in these installation and operating instructions.

Your IQcare team

FOR QUESTIONS:
JUST CALL OR EMAIL!

FON 0 23 54 / 944 99 69

E-MAIL INFO@IQFY.DE

NO ELEKTROSMOG

All IQcare sensors and switches transmit only very briefly when they are activated. This means that there is no electrosmog!

AREAS OF APPLICATION

Indoors and outdoors, e.g. room, bathroom, balcony.
As a stand-alone alarm system, wherever the residents are where the residents are at the moment.

SCOPE OF DELIVERY (depending on order)

NOBU

NOBU
Emergency
call button

NAVI
Radio
receiver

NAVI GO
Mobile
Radio receiver

PREPARATION


Please decide now which situation you have and which products you want to install in addition:

- 1 Preparations for an independent emergency call on a call system** *Continue on page 19*
- 2 NOBU as an independent emergency call on a call system (working current)** *Continue on page 20*
- 3 NOBU as stand-alone emergency call on a call system (closed-circuit current)** *Continue on page 21*
- 4 NOBU to existing radio receiver NAVI** *Continue on page 22*
- 5 NOBU on mobile receiver NAVI GO** *Continue on page 23*
- 6 Error handling** 26
- 7 Maintenance** 27
- 8 Desinfection** 27
- 9 Technical Data** 27

1 PREPARATION

If there is an adapter cable in the scope of delivery that has not yet been connected, please do so.

Important note: This is a connector, please do not twist the plug onto the connector. You will find a white line on both the radio receiver and the adapter cable you will find a white line as a mark - these must fit over each other.

TIP

If you have several IQcare wireless solutions in your home, number these units using our helpful identification stickers.

2 NOBU as an independent emergency call on a call system (working current)

- Find a socket that is close to the connection of your alarm system so that you can plug in the free end of the adapter cable at **the end of the teach-in process.**

Plug the radio receiver into this the radio receiver. The CLR-button now flashes green, the radio receiver is in the teach-in mode.

Press 3x the NOBU emergency call button.

The LRN-button now lights up permanently red. NOBU and NAVI now form one unit. The connection is confirmed by an acoustic signal.

Connect the adapter cable with your call system.

3 NOBU as an independent emergency call on a call system (closed-circuit current)

- Find a socket that is close to the connection of your alarm system so that you can plug in the free end of the adapter cable **at the end of the teach-in process**. Plug the radio receiver into this the radio receiver. The CLR-button now flashes green, the radio receiver is in the teach-in mode.
- Press once the LRN-button, the receiver is now in standby mode. The CLR-button flashes red.
- Press 3x the NOBU emergency button.
- The LRN-button now lights up permanently red. NOBU and NAVI now form one unit. The connection is confirmed by an acoustic signal.
- Connect the adapter cable with your call system.

4 Teach-in to an existing system

- You already use an IQcare system such as the SESY bed escape system - accordingly the radio receiver NAVI is already present.

Then proceed as follows:

- Unplug the Navi from the socket, wait 3 seconds (until the light goes out) and plug it back in.
- Press the LRN button for approx. 10 sec. until the CLR-button flashes: green for working current red for quiescent current
- Press the NOBU 3 x briefly.
- The receiver signals the connection with an acoustic signal.
- Finally, press the CLR-button once.

5 NOBU with mobile receiver NAVI GO

- Plug the radio receiver NAVI GO into a socket.
- The CLR-button flashes green.
- Press 3 x briefly in succession on the emergency call button.

3 X

- The CLR-button goes off, the NAVI GO lights up red and is ready for use!
- As soon as the NOBU is pressed an acoustic signal sounds and the NAVI GO flashes red. The signal must be confirmed by pressing the LRN / Off-button to confirm the signal.

Teaching of an additional emergency call button NOBU

- Unplug the NAVI GO mobile radio receiver from the power socket and wait 6 seconds.
- Plug the radio receiver NAVI GO into a socket.

- NAVI GO now lights up red.

- Press the LRN-button for approx. 10 seconds.

10 SEC.

- Now the CLR-button flashes and the NAVI GO is in its basic colour.

- Press the NOBU button 3 x.
- NOBU to be taught in additionally.

3 X

- Press 1 x on the CLR-button.

- NAVI GO now lights up permanently red and is ready for use again.

- As soon as one of the NOBUs is pressed, an acoustic signal sounds and the NAVI GO flashes red. The signal must be confirmed by pressing the LRN / Off-button.

- The NOBUs and the NAVI GO are now ready for use and can be used as desired. If you change the location of the NAVI GO receiver or the NOBUs, the affiliation remains the same.

ERROR HANDLING

- In rare cases, the installation may not work as described.
- **In this case, it is recommended to completely reset the radio receiver to factory settings. To do this, please proceed as follows:**
 1. by unplugging the radio receiver and plugging it back into the socket, the setting can be deleted within 30 seconds.
 2. Press the LRN-button for approx. 10 seconds, until the CLR-button flashes.
 - 3a. **NAVI**

Press the CLR-button for approx. 3 seconds until the light briefly goes out before it flashes green again. The radio receiver is now back in the delivery state and can be reset according to the operating instructions.
 - 3b. **NAVI GO**

Press the CLR-button for approx. 3 seconds. The Navi Go lights up red again briefly before it is back in the delivery state. The CLR-button flashes green again!

Attention: Please note that you must then reprogram all the previously connected components must then be reprogrammed!

ATTENTION:

If you have tuned in several accessory components to a Navi / Navi Go radio receiver and would like to remove some of them, please call our technical support and you will receive instructions on how to proceed by telephone!

MAINTENANCE AND CARE

- IQcare products can be cleaned carefully with a soft cloth and normal household cleaning agents.

DESINFECTATION

- Use commercially available disinfectants. A list with recommended disinfectants can be found on the homepage under the downloads ...
- ... or you can download it directly via the QR code on the right

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IQfy GmbH