

MEDICAL BREAKTHROUGH 9™

User Manual



HOW-TO VIDEOS

To make things easier for you we provided these easy-to-follow instructional videos that will help you set up and use your new Medical Breakthrough™ chair.

Note: On the next page, we included a list of Frequently Asked Questions. Please read them because they will answer most questions you may have.

▶ HOW TO INSTALL THE BREAKTHROUGH 9

We understand most people do not like reading manuals, so we made a video that will show you how to install the Breakthrough 9. Please go here to watch the video that will show you step-by-step how to install the Breakthrough 9.

www.medicalbreakthrough.org/mb9install.php

▶ HOW TO USE THE BREAKTHROUGH 9

The Breakthrough 9 is very easy to use. There are only a few buttons on the remote, but the first time you use it might be overwhelming. So we made a video that shows what each button does. With this YouTube video, you can rewind it and play it back if you miss anything. You may have a lot of questions the first few times you use the Breakthrough 9, but if you use this video as a reference guide for a few weeks, you will know everything about the chair.

www.medicalbreakthrough.org/mb9demo.php

**The videos and FAQ's will answer most questions. But, if you have further questions, feel free to contact us directly at www.medicalbreakthrough.org/contact-us.php*

***We are open 10AM-5PM PST (California Time). Thank you again for purchasing the Breakthrough 9.*



FREQUENTLY ASKED QUESTIONS

These are the most commonly asked questions we get from our customers. Please read these questions and you'll have a better understanding of the Breakthrough 9 Massage Chair in no time.

1. **What do I do with all the pads included with the chair?**

A lot of customers ask us how to use all the pads we give and why we include so many pads. Many massage chair owners buy 3 different massage chairs. They start with a light chair, move on to the medium chair, and when their body gets used to it, move on to an extremely strong chair. We don't want you to buy 3 different chairs. One chair is good enough. That's why we include a lot of comfort pads. You can place these pads wherever you want a lighter massage. They're there to help you.

Most massage chairs are very rough and does not include any pads. We include a lot of pads for you to make the chair as comfortable as possible. These are provided completely free. You can use them anywhere on your body. Everyone is different -- some people have extra sensitivity in their back, others experience it around their feet. For that reason, we included all of the pads.

Mix and match these pads to your liking. It's a lot of work for us to make each one, and although

you may not need it, we provide it for you just in case a part of your body is extra sensitive. Use them on various parts of your body and in an assortment of positions. Although these pads aren't necessary for your enjoyment, they're available for you to add an extra level of luxury to your massages.

2. **If I have questions and need some help who do I call?**

Our customer support and tech team is open Monday to Friday 10AM to 5PM PST (California Time). The chair is very easy to use, the first time you use it, you might have a lot of questions, but after a few weeks you will know the chair very well. Just play around with the settings on the chair and after a while you will understand every button. If you have questions, we are here to help. Please call (661) 775-4866 Monday to Friday. Remember, we are in California, so please call when we are in, which is 10AM to 5PM PST.



▶ Frequently Asked Questions (FAQs)

- 3. Why is the back heating not working after I turn it on?**

It does not work immediately. The back infrared heating system takes 3 minutes to warm up before you will feel any kind of warm sensation.
- 4. What does the width button do?**

This button allows you to adjust the massage hands closer to the center of your back or further out during "Tapping" or "Shiatsu" massage functions. There are 5 levels to choose from.
- 5. Why are there 2 speakers on each side of the upper arm?**

Only 1 of those is a speaker on each side. The other is where the negative oxygen ions are emitted.
- 6. How do I adjust the speed of the rollers on my feet?**

You can adjust the speed by pressing the Sole Roller button. There are 3 different speed levels to choose from.
- 7. Can I make my massages longer?**

Yes, you can. By pressing the button that looks like a clock, this is the time button. Each time it is pressed it will add 5 more minutes to your massage for up to 40 minutes.
- 8. How do I make my back massage lighter?**

The back massage intensity can be adjusted by pressing the Upper-Body Intensity buttons. Press the down arrow to make the massage lighter. With 5 levels, this allows you to find the perfect back massage for your needs
- 9. The back massagers are too fast so how can I change this?**

Press the Speed button on the right side of the remote just below the Upper-Body Intensity button. There are 5 different speeds available so press it a couple times to find the one that suits you the best.
- 10. The massage rollers are too high on my neck, how do I move them down?**

In the circle group of buttons just below the screen, the top up and down arrows are labeled "Shoulder Position." By pressing the down arrow, you make the rollers move further down on your neck.
- 11. Can I turn off the arms if I only want a back massage?**

Yes, you can. Press the Arm Air Pressure button 3 times to turn it off; it should disappear on your remote screen.
- 12. How can I save a massage?**

If you found a massage setting that you really enjoy, then press the power button. The screen will display "SA & MEM" within 10 seconds. Then, choose the M1 or M2 button on the bottom of the remote, this will save your massage. The next time you get in the chair, turn on the power and press the button you saved it under to begin your personal massage.
- 13. What does the 3D button do?**

This button moves the massage hands further out of the chair for a deeper more humanistic massage. The speed and the intensity of this function are adjustable, so play around with it some until it fits your needs.
- 14. How can I make the footrest come up without the backrest going back more?**

In the buttons grouped in the circle below the screen, the middle round one is the Zero-Gravity button. The one just to the upper left of it will make your footrest come up without moving the backrest. You will have to press and hold this button until it is where you need it and then release the button to stop the movement.
- 15. Can I move the rollers during my Back Stretch Massage?**

Yes, you can. On the second panel of buttons, in the 3rd row of buttons just below the hands buttons, you will see an up arrow, a back outline button, and a down arrow button. In the fixed or partial back stretch position, you can move the rollers either up or down.



► Frequently Asked Questions (FAQs)

16. Can I turn the feet rollers off?

Yes, you can. By pressing the Sole Roller button 3 times, it will turn off.

17. How do I turn on the Zero-Gravity Mode?

The button for this function is in the very middle of the round set of buttons just below the screen. By pressing the middle button that looks like a squiggle line, your chair will recline back into Zero-Gravity Mode.

18. How do I make the Arm Massage lighter?

On the bottom right-hand side of the remote, the very last blue button states "Arm Air Pressure" below it. By pressing this button, you can adjust the intensity of your arm massage. There are 3 levels of intensity to choose from.

19. How do I get the Shiatsu massage?

On the second panel, you will see buttons with hands on them in the 1st and 2nd row. In the second row, you can see the hand that looks like it's doing an upside down "Thumbs Up" motion. Press this button and it will turn on the Shiatsu massage mode 1 and press it again to go into mode 2.

20. Can I adjust the strength of the Auto massages?

Yes, you can by pressing either the Upper Body Intensity buttons or the blue intensity button on the second panel.

21. How do I make the backrest come up further?

In the round set of buttons, you will see the middle squiggle button. The button to the upper-right of the squiggle allows the backrest to come up. Hold this button and the backrest will come up slowly and then release the button when you've reached the desired position.

22. How do I know which automatic massage program is running?

On the top of the screen, you will see the time and next to the time you will see: FO Auto, F1 Auto, F2 Auto, F3 Auto, F5 Auto, F6 Auto, F7 Auto, or F8 Auto. FO is the Demo massage, F1 is the Comfort massage, F2 is the Fatigue Recovery massage, F3 is

the Ache Relieve massage, F5 is the Full Air massage, F6 is the Waist Stretch massage, F7 is the Synchronous Music massage, and F8 is the Relax massage.

23. How do I make the Shoulder Grasp faster?

Once in this massage method, all you have to do is press the speed button. There are 5-speed levels to choose from so finding the right one should be easy.

24. Why isn't my chair going down after my massage ended?

Once the chair has run through the timed massage, the chair will not return to its sitting position just yet. This screen will be asking you to save the massage. If you choose not to save it just press the power on/off button to make the chair return to a normal sitting position.

25. How do I turn on the Negative Oxygen Ion function?

Below the built-in speaker there is another thing that looks like a speaker on the upper arm headrest. This is where the negative oxygen ions come out. Next to this are 2 buttons, the top button turns on this function and the button below it adjusts the speed.

26. How do I make the massage stronger on my feet?

On the second panel of buttons, you will see the set of 4 buttons in blue. The second blue button states "Intensity" below it. Press this button to adjust the strength of the massage.

27. What does the Arthrosis Traction button do?

This button activates our Calf Stretch function. Pressing this button will turn it off if it is currently on or turn it on if it is currently off.



► Frequently Asked Questions (FAQs)

28. During the Full Air Massage can I turn off the Calf, Feet, and Seat functions?

Yes. During this massage, if you want to turn off the calf, feet, and seat, all you have to do is press the 1st blue button on the second panel. It's at the very bottom left-hand side. It will state "Lower Body" below it. Press it 3 times to turn it off completely.

29. Which massage is the best to fully relax me?

To fully relax your body, you want to choose the Comfort massage. Your screen will display F1 Auto on top to let you know you've chosen the right one. This massage is designed to relax the muscles and bones of your body.

30. Which functions can I adjust the speed in?

The Speed can be adjusted in: Shoulder Grasp massage, Kneading massage, Tapping massage, Kneading & Tapping massage, Shiatsu massage, 3D massage, and the Sole Roller massage.

31. How do I get the chair to stretch me?

To do this, look on the 1st panel of buttons. In the 2nd curved row press, the middle button that states "Waist Stretch" below it. The chair will squeeze your calves and pull them down gently while the backrest reclines to fully stretch you.

32. How do I make the backrest recline more?

In the set of round buttons, you will see the middle button that looks like the squiggle which is your Zero Gravity button. To the bottom right of this button, you will see a squiggle with an arrow pointing at it. Hold this button to slowly recline backward and then release it when you've gotten to where you'd like it to be at.

33. What does the Fatigue Recovery button do?

This function improves your fatigue and recovers your body's energy. It goes through the different types of hand massages on your back while also squeezing your calves, massaging your buttocks and providing a sole rolling massage.

34. What does the Relax button do?

This button provides different massages that will help to promote blood circulation throughout your body. It will go through a combination of backhand massages as well as swinging your hip from side-to-side and using the airbags throughout the chair to stimulate the blood flow in your body.

35. What does the Ache Relief button do?

This function is designed to relax your normal aching areas. It provides you with a deeper massage on your body's acupressure points by going through a mixture of the various backhand massages. The other functions, like the arms lower body and sole rollers, continue to work as well.

36. Why isn't my chair massaging after I've turned it on?

You have to press a massage function first before the chair will begin. After you've pressed one of the automatic massages or something like the kneading function, the chair will begin its body scan and then beep 7 times before the massage will actually start.



Thank you for purchasing the Medical Breakthrough 9™ Massage Chair. Please read this manual carefully for instructions on how to use it correctly. Please pay close attention to the safety and maintenance sections and keep this manual handy for future reference.

Note: Dr. Recommendation, USA reserves the right to revise the design and description of the Medical Breakthrough 9™ & the manual without any further notice. Note that the color of your massage chair may be different than the color shown in this manual.



▶ Table of Contents

Maintenance and Cleaning	Page 02
Safety Precaution	Page 03
Name and Functions of Components	Page 05
Usage and Functions	Page 06
Method of Usage	Page 07
Product Specifications	Page 16
Accessories	Page 17
Responsibilities of Others	Page 19
Error	Page 20



▶ Maintenance and Cleaning of Medical Breakthrough's High Quality Synthetic Leather

Medical Breakthrough's Synthetic Leather is made from a polyurethane resin after coating, drying, and curing. Compared with genuine leather, Medical Breakthrough's High Quality Synthetic Leather is more resistant to acid, alkali, heat and sunshine. However, because of its own physical and chemical properties, it should not be put in a place where there is dust, high temperature, high humidity, low temperature, strong sunshine, or an environment containing acid and alkali agents or solutions. Otherwise it will easily crack, age quickly and become hydrolyzed.

■ Cleaning Methods and Precautions for Medical Breakthrough's High Quality Synthetic Leather:

1. Organic agents, grease solutions and compounds, strong solvents (e.g. alcohol, turpentine, and thinner) cannot be used for cleaning the leather surface, otherwise the leather surface will become corroded.
2. Harsh or powdered detergents cannot be used for cleaning the leather surface, otherwise there will be cracks
3. Environment-friendly detergent which contains no abrading and with a PH value of 5-7 should be used for cleaning the leather surface.
4. If the leather surface becomes overly dirty, use non-alcoholic baby wipes to wipe the dirt. Brushes should not be used for cleaning the surface of the chair.

■ Maintenance Methods and Precautions of Medical Breakthrough's High Quality Synthetic Leather:

1. Good room ventilation is important because an environment that is too dry or too humid will accelerate the aging of the leather.
2. Avoid exposing the chair to direct sunshine or air conditioner output, otherwise the leather surface will harden. Direct exposure to sunshine may also cause the color of the leather to fade.
3. During overly hot days, the leather may absorb the sweat produced by the body. The organic matter in the sweat and the leather will have a chemical reaction during high temperatures and humid environments, and there will be a peculiar smell. For this, a soft cloth should be used for wiping frequently.
4. While wiping the leather, an alkali cleaner should not be used because there is an acid treatment during the leather making and the alkali will decrease the pliability of the leather and there will be Cracks if the alkali cleaner is used.



▶ Safety Precautions

Important Safety Warnings:

1. To prevent injuries, children and or pets should not play around the chair during operation.
2. Users must stop and turn the power off immediately upon feeling ill.
3. If the chair is functioning abnormally, Turn off the power switch and unplug the power cord and call our tech support department immediately. Do not attempt to disassemble or repair the chair yourself.
4. Keep all body parts away from the internal mechanisms.
5. Stop using the chair and unplug the power cord if any part of the fabric cover is torn or worn out (exposure to mechanics.) Please contact tech support immediately.
6. Users with recent, chemically treated hair, or those using any topical hair products, should place a towel between their hair and the upholstery to protect the material.
7. It is not recommended for those under the age of 16 to use the chair.
8. The warranty will void for users who weigh over 300lbs.
9. Do not use the chair if the back pad is flipped over, and the rollers are exposed.
10. Please use a well-grounding power supply that is suitable for this product.
11. Please unplug after using or before cleaning to avoid injury or damage.
12. Operate your Medical Breakthrough 9 according to this instruction manual.
13. Do not use accessories which are not recommended.
14. Do not use this product outdoors.
15. Please read this instruction manual carefully before operating your Medical Breakthrough 9.
16. The recommended time for each session in a Medical Breakthrough Chair is 15 minutes.
17. Please do not fall asleep while using this product.
18. Please do not use this product when intoxicated or ill.
19. Please do not use the product within one hour after having a meal.
20. Be cautious when adjusting the intensity of the massage in order to avoid injury.

Environment Usage

1. Please do not use the chair under high temperatures or in extremely humid environments.
2. Please do not use the chair in a transition between hot and cold environments. Let it

return to room temperature for at least one hour before use.

3. Please do not use the chair in a dusty environment.
4. Please allow enough space for the chair to perform during use.
5. Please use the chair in properly ventilated rooms.

Medical Conditions

Please consult your physician prior to use if you have:

- Osteoporosis
- Heart Disease
- Electronic medical equipment such as a heart pacemaker
- Fever
- Pregnancy or are menstruating
- Skin injuries or skin disease
- Mental conditions
- Heat sensitivity
- Irregular curvature of the spine or other spine abnormalities

Additional Safety Information:

- Do not use the chair if you are wet
- The massage chair has a heated surface. If you are sensitive to heat, please use caution when using the chair



▶ Safety Precautions

Safety:

1. Check the voltage to see if it is suitable to the specifications of this product.
2. Don't allow water to go into the chair to avoid electric shocks or damage to the product.
3. Do not damage the wires or change the circuit of this product.
4. Do not use a wet cloth to clean the electric parts such as the switch and plug.
5. Stay away from this product under a power cut state to avoid injury because the power could come back on suddenly.
6. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
7. Children should be supervised to ensure that they do not play with the appliance.
8. If the supply cord is damaged, it must be replaced by the manufacturer, its service

agent or similarly qualified persons in order to avoid a hazard.

Maintenance:

1. Damage caused by the user attempting to disassemble the chair or replacing spare parts themselves will void the warranty.
2. This product should be maintained by certain appointed agents. DO NOT attempt to disassemble or repair the chair yourself.
3. Please do not forget to cut the power after usage.
4. Do not use this product if the socket loosens.
5. If this product will be left unused for a long time, please curl up the wires and store this product in a dry and dustless environment.
6. Do not store this product under high temperatures or near fire and avoid longtime direct sunshine to it.
7. If the detachable cord supply Cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
8. Please clean this product with a dry cloth. Do not use thinner, benzene or alcohol.
9. The mechanical components of this product are specially

designed and made, so it is not necessary for it to be specially maintained.


10. Do not roll or pull this product above uneven ground; it should be lifted before moving.
11. Please use it intermittently; do not make this product work continuously for too long.

Solutions for Ordinary Malfunctions:

1. It is a natural movement sound to hear the motor sounds while is use.
2. If the controller cannot work normally, please check the plug and socket to see if they are connected strongly and to make sure the switch is turned on.
3. Once the set massage time is over, the chair will automatically power off
4. If this product works continuously for too long, the temperature protector will make the product power off. The chair will be usable after a thirty minute rest period.



▶ Name and Functions of Components



- 1. Pillow Pad
- 2. 30 Digital Audio
- 3. Arm Air Bag Assembly
- 4. Seat Cushion
- 5. Calf Massage Assembly
- 6. Foot Massage Assembly
- 7. Negative Oxygen Ion Function
- a. Upper Arm Air Bag Assembly
- 9. Back Cushion
- 10. Armrest
- 11. Led Colorful Lights
- 12. Controller Pocket
- 13. Back Cover
- 14. Controller
- 15. Driving Box Cover
- 16. Safety switch
- 17. Caster Wheels
- 18. Power Cord and Plug
- 19. Power socket
- 20. Fuse Box
- 21. Power Switch
- 22. Controller Cord Socket



▶ Usage and Functions



Functions:

1. Medical Breakthrough 9™ is designed with a set of smart 3D mechanical hands, which move up and down, stretch back and forth, all in a sleek, muted design.
2. The shoulder function is designed with an automatic detection and micro adjustment function. This has a body curve and massage point automatic detection that can automatically adjust to the distance of the massage hands, moving back and forth according to the detected body curve and massage points in order to make the massage more humanistic and scientific.
3. The massage functions are as follow: Comfort, Ease, Ache Relieve, Demo (fast experience), Relax, Music Sync, Waist Stretch, Full Air Pressure and other specific massage functions.
4. The chair can be set with 2 kinds of memory massage functions: M1 and M2 buttons are the saved massages.
5. The Upper-body massage has 3 massage positions: Overall, Partial and Fixed position. There are 6 different massage methods: Shoulder Grasping, Kneading, Tapping, combination of Kneading & Tapping, Shiatsu, and 3D with 5 massage speeds available. Under the Tapping and the Shiatsu massage status, the width of the massage hands can be adjusted to 5 separate levels.
6. The Arm air massage and Upper Arm air massage functions have 3 intensity adjustments.
7. The Lower-body air pressure massage function, Back air pressure massage function, Buttock air pressure massage, and Calf Rest air pressure have 3 intensity adjustments.
8. The upper arm has two grades for the width adjustment function. When you adjust to the last grade, push into the arm for your hand and arm device to automatically reset.
9. The chair's back heating function uses carbon fiber as an infrared heating source.
10. The footrest can be adjusted to suit various heights.
11. The chair has the ability to slide at the base in order to utilize extra space.
12. LED colorful lights on armrest with negative oxygen ionizer on the upper arm.
13. The chair has an MP3 player function with 3-D digital audio built into the upper arms.
14. Connect the smart device via Bluetooth, to play music.
15. There are 8 natural energy magnets placed in the chair's air bags.




▶ Method of Usage

LCD Manual Controller Controls The Massage Chair







▶ Method of Usage


Switch The Main Power Supply ON



To Start a Massage

- Press the red power button to start the massage, the massage chair will automatically recline back.
- Press an auto massage button; the massage hands proceed with an auto detection of your shoulders and a body scan. Then the calf rest roller begins to work at the same time. (Choose the auto program of rapid massage to start the massage program immediately to save half the time). If you made no choice of the auto massage program after 20 minutes, the massage chair will automatically switch off.
- Upon completion of scan, a beeping sound will begin emitting from the controller and the massage heads will gently roll on your shoulder region to demonstrate where the chair has sensed the position of your shoulders. During this period, if you would like to further adjust the detected shoulder region, you can do so by pressing the up or down shoulder position buttons seen below.
- If you have not made any adjustments to the shoulder position or selected a manual massage mode during the time after scanning, the chair will automatically start the fine massage program for you.

Buttons	Description	On Display
 Shoulder Position Adjustment (Upper Button)	Press to adjust the massage heads higher	
 Shoulder Position Adjustment (Lower Button)	Press to adjust the massage heads lower	



▶ Method of Usage

Choosing Massage Methods

Memory Massage: When the shoulder position detection begins, you can use the M1 or M2 memory buttons to start massage. There is no need to wait until the detection is finished, it will go back to the stored massage models to massage.



Characterization of Automatic Massages

Buttons	Description	On Display
	A quick massage that goes through all of the different massage functions	<i>F0 Auto</i>
	Massages the stiffness from muscles and bones, leaving you feeling refreshed	<i>F1 Auto</i>
	Helps the body recover after a day of stress on the body	<i>F2 Auto</i>
	Deep massage that targets pressure points throughout the body	<i>F3 Auto</i>
	Uses multiple air bags to gently squeeze your arms, shoulders, waist, and legs	<i>F5 Auto</i>
	The backrest reclines while your legs are pulled down to fully stretch your body, then the chair returns to its original position, then repeats the stretch	<i>F6 Auto</i>
	Plays soothing music that syncs to the massage	<i>F7 Auto</i>
	A full-body massage designed to completely relax your body and relieve stress	<i>F8 Auto</i>

























APP Program Auto Calling

Buttons	Description	On Display
	APP automatic function selection: 3 auto modes	<i>A1~A3</i>
	Confirm Selection: If the current APP mode is invalid and you press the confirm button, the horn sound will remind you it's an invalid selection	<i>A1~A3</i> ^{APP} Auto



▶ Method of Usage













Upper Body Manual Operation, Lower Body Manual Operation

Buttons	Description	On Display
 SHOULDER GRASP	A deep massage targeting your shoulders. Speed and intensity can be adjusted.	
 SWEDISH KNEAD	A gentle, kneading massage. Speed and intensity can be adjusted.	
 THAI TAP	Thai Tap: Rhythmic tapping along the body. 2 massage modes with adjustable speed, intensity, and width.	
 WORKOUT	A combination of knead and tap massages. Speed and intensity can be adjusted.	
 JAPANESE SHIATSU	The rollers steadily move in, then out to target specific pressure points. 2 massage modes with adjustable speed, intensity, and width.	
 3D BURMESE 3D	4 different deep tissue massages with adjustable speed and intensity.	
 SPEED	Speed: 5 Adjustable Levels	 Speed
 WIDTH	Width: 5 Adjustable Levels	 Width
	Fixed position massage	 H0
	Back and forth massage on partial position	 H4
	Back and forth massage on upper body	 H1
	Back and forth massage on upper back	 H2
	Back and forth massage on lower back	 H3
 LOWER ADJUST	During fixed or partial position massage mode, the massage hands can adjust downward	Message hand mov downward
 UPPER ADJUST	During fixed or partial position massage mode, the massage hands can adjust upward	Message hand mov upward
		Message position adjustment (under fixed and partial position massage mode, it is adjustable)



▶ Method of Usage


Air Compression Massage Functions

Buttons	Description	On Display
 FOOT REFLEXOLOGY/ GLUTE MASSAGE	Lower body air compression with 3 modes: Lower body, back and waist, and calf	Feet  Back  Buttocks 
 INTENSITY	Lower body, back and waist, and calf air compression with 3 adjustable intensities	Feet  Back  Buttocks 
 ARTHROSIS TRACTION	Lower body and leg air compression mode, also uses the Calf Stretch function	
 4D ARM REFLEXOLOGY	Arm air compression mode with 3 adjustable intensities	Arm 



Other Functions

Buttons	Description	On Display
 UPPER BODY INTENSITY WEAKEN	Decreases the intensity of the massage hands	 Upper
 UPPER BODY INTENSITY STRENGTHEN	Increases the intensity of the massage hands	 Upper
 SOLE ROLLER	Adjusts the speed of the sole rollers with 3 adjustable levels	
 HEATER	Uses infrared heating therapy to heat your back after 3 minutes	

Massage Time Adjustment

Buttons	Description	On Display
 TIME	The default time for a massage is 20 minutes. Pressing this button will add 5 minutes to your massage. The massage time can range from 5 to 40 minutes.	20 Min






Negative Oxygen ion (controller button is in the upper arm)

Buttons	Description	On Display
NEGATIVE OXYGEN ION FUNCTION	Cleanses the air and body of toxins and impurities to help you physically and mentally feel better.	
SPEED	2 adjustable speeds available	



▶ Method of Usage


The Angle of backrest and Calf Rest

Buttons	Description
	Zero Gravity Mode: Fully reclines the chair so you are laying down, then continues with a selected massage.
	Calf Rest Lifting Button: Hold this button and the calf rest lifts slowly, then stops when you release it.
	Calf Rest Decline Button: Hold this button and the calf rest declines slowly, then stops when you release it.
	Backrest Lifting Button: Hold this button and the backrest lifts slowly, then stops when you release it.
	Backrest Rest Decline Button: Hold this button and the backrest declines slowly, then stops when you release it.


Shut off Power and Finish Massage

Pressing the power button during a massage will stop all massage functions immediately. The chair can also power off when the selected massage is finished. The backrest and footrest won't return to its original position yet.

- Press M1 or M2 button within 10 seconds to store the massage mode.
- Or press the power button within 10 seconds. The massage hands, backrest and footrest will return to their original position automatically.
- Or if you wait for 10 seconds, the massage hands return to their original position.



Turn Off Power (Cut off the power of the machine as shown in figure)



▶ Method of Usage

Wireless Bluetooth Music Player Operation

- Audio devices with a Bluetooth function can pair with this chair. Any audio transferred wirelessly will be played through the chair's speaker.

Arm Assembly Installation

1. In the side of the chair (where the armrest will go), find the electrical wires and air tube and gently bring them toward you. **(Figure 1)**
2. Take a look at the armrest and notice the similar wires and air tube in the inner part of the armrest. Align the wires and tube with the corresponding ones on the chair itself. **(Figure 2)**




Figure 1




Figure 2

3. Once you know each wire and tube match between the armrest and the chair, connect them accordingly. Be sure that you are installing the armrests in the correct position. **(Figure 3)**
4. After connecting the wires and air tube, gently insert them into the bottom of the seat frame. Be careful not to fold the air pipe when fully installing the armrest. **(Figure 4)**




Figure 3






Figure 4



▶ Method of Usage

- Now that the wires are connected, notice the hooks on the armrest and how they correspond to grooves in the steel frame of the chair. **(Figure 5)**
- Mount the armrest onto the correct side of the chair and make sure the armrest hooks to the corresponding parts of the chair's steel frame. **(Figure 6, 7)**



Note:
The connect wire must be inserted into the right place, avoid to get stuck.

- Secure the armrest to the steel frame using screws and a hex wrench. (Figure 8, 9)
- Follow these instructions for the other armrest and ensure both armrests are securely fastened.




Figure: Hex head screw fixed in front part of armrest




Figure: hex head screw fixed in back part of armrest



▶ Method of Usage

Head Cushion and Backrest

1. Using the head cushion can reduce the intensity of neck and shoulder massages. It is optional so it depends on whether you want the extra cushioning. The head cushion is easily removable and the backrest is removable using a zipper. **(Figure 1)**
2. The AB flannelette is the portion of the head cushion you will insert through the belt of the backrest. **(Figure 2)**
3. Lay the head cushion on the backrest and adjust it to where it feels best. **(Figure 3)**




Figure 1




Figure 2






Figure 3



Floor Protection


Placing a heavy massage chair down on a wooden floor for a long time may damage the floor. Please place a carpet or other object on the floor to avoid damage.



Note: Do not move the chair with its wheels on the wooden floor, rough ground, or in a narrow space. Two people are required to lift the chair in this case. Carry from the calf front sealing plate and back cover (please don't lift the armrest).

Moving The Chair

Make sure all the wires are far above the ground, tilt the backrest backward until the chair is resting on the caster wheels, and push the chair forward or backward. Finally, return the chair to the normal position in a slow and gentle manner.



Note: You must turn off the power switch, then unplug the power cord and remote control wire before moving the chair.



▶ Product specification

Model: Medical Breakthrough 9™ Massage Chair

Description: Massage Chair

Rated Voltage: 110-120V-60Hz
220-240V -50Hz/60Hz

Rated Power Input: 240W

Rated Time: 20 Min

Noise: ≤ 58db

Standby power: < 0.5W

Safety Design: Class I

Material introduction: PVC, PA, steel parts, electric, and electronics parts

MATERIAL	PU	LEATHER	CLOTH	WOOD



▶ Accessories For The Breakthrough 9

If you are interested in any of these Medical Breakthrough Accessories, or if you want to see any special deals we may have, please visit www.MedicalBreakthrough.org



Medical Breakthrough Leather Protector

To protect your chair, we handcrafted a durable cover that is designed to perfectly fit your chair. This Leather Protector will help your chair look its best throughout its entire lifespan without the cracks or tears you might find in other massage chairs.



Podiatrist by Medical Breakthrough

250 million people suffer from foot pain, that's why we worked with a team of podiatrists to engineer the Podiatrist by Medical Breakthrough, a compact foot massage system designed to help heal your feet and calves on-the-go.



▶ Accessories For The Breakthrough 9



Facial Relaxation Eye Massager

Eye and head massages have been proven to relieve stress and reduce migraines, so we developed this revolution in relaxation technology. Facial Relaxation is an eye massager designed to aid you against persistent headaches, lack of sleep, and fatigued eyes.



Full Body Pain Reliever

Designed to deliver an innovative, versatile massage, the Full Body Pain Reliever will loosen tight muscles throughout your body, easing any type of pain or soreness. It's also portable so you can bring instant pain reliever wherever you may need it.



▶ Responsibilities of Others


- Unless otherwise expressly agreed in writing, In the event of a problem identified prior to installation, the unit must not be installed, and Company is not responsible for failures or damage that could have been discovered, repaired, or avoided by proper inspection prior to installation.
- Damage occurring in transit is the responsibility of the carrier.
- The user or installer MUST open the crate and inspect the unit for damage when it is delivered.
- If damage is discovered, it must be reported immediately to the seller and the carrier.
- Failure of the carrier to respond should be reported to the seller and the carrier.
- Your freight claims should be filed promptly thereafter.
- Damage occurring to the unit during installation is the responsibility of the installer, contractor or user, and damage occurring thereafter is the responsibility of the user. Failure of any non-factory installed components, parts and accessories is the sole responsibility of the user.
- Any replacement parts shall be covered only by the original equipment manufacturer warranty, if any. The Company is not responsible for costs to modify any product to obtain any code approval, such as city, county, or state building codes.
- Failure to abide by the instruction explicitly stated in the manual is the sole responsibility of the user.
- Damage to the installed equipment, or injury to any persons, as a result of, but not limited to, misuse, modifications, and negligence, is the sole responsibility of the user.



▶ Error

01. Foot Kneading Unit Signal Lost
02. Back Roller Kneading Function
03. Back Roller Up/Down Function
04. N/A > doesn't appear in current generation chairs anymore
05. Back roller moves but with limitations
06. N/A > doesn't appear in current generation chairs anymore
07. N/A > doesn't appear in current generation chairs anymore
08. Leg Wheels don't work at all
09. Leg Up/Down (Telescopic) function limited
10. Leg Up/Down (Telescopic) function failure
11. Calf Motor Failure
12. Bevel Motor Failure (Back Actuator #2 = Deeper Actuator furthest away from board)
13. Back Motor Failure (Back Actuator #1 = Closest to main board)





Manual Version 7.20