


**ECG/Blood glucose/Blood pressure Smart Watch**

User Manual

**Device's structure**




- 1.83 inch HD big screen
- Zinc alloy
- Electrode plate
- stainless steel clasp
- Silicone strap
- Charging port
- Electrode plate
- Switch on/off
- Sensor
- PTT hole
- Electrode plate


**Disassembly and wear**

**1. Disassembly**  
Align the metal pin on left end of the watch strap with the hole at the connection point of the bracelet, flip the metal ball on the right, and align the hole on the right to complete the installation and disassembly.

**2. Wear**  
Better wear the watch at 1-2CM behind the ulnar neck of the wrist. When wearing, make sure that the bottom of the watch is close to the skin to avoid light leakage of the sensor and affect the accuracy of the measurement.




**Switch on/off:** Long press the screen to power on.




**Install APP and Bluetooth connection**

1. Turn on the device, turn on the Bluetooth of the mobile phone, scan the QR code below or search "H Band" in the App store / Google play store to download and install.




2. Register and log in


Please register with an email or mobile number if you don't have a H Band account, or you can choose "Login without an account" to enter.



3. Open the APP. To the connection interface please tap "Scan Device", turn on the Bluetooth, then find the corresponding device to connect. Also can go to the data panel to get connect.



**Operation guide of the device function interface**



**Shortcut set**  
Swipe down  
Device interface  
Swipe left and right  
Notification  
Swipe up

Step  
Heart rate  
Blood pressure  
ECG  
Blood glucose  
Blood oxygen  
Body temperature  
Sleep  
Weather

**Device function introduction**

**1. Time, date & week**  
For the first time, please connect to the APP to sync the correct time, date, and week.

**2. Dial switch**  
Long press the screen to switch the dial style.  
Note: To customize the watch face, please go to APP Settings (Mine -> EPG -> Dial Settings).


**3. Shortcut set**  
Swipe down on the main time interface to enter the drop-down shortcut interface. The specific functions are as follows:  
Find phone/Do not disturb/Brightness/About/Menu style/Settings

**4. Workouts**  
Press the button on the right side of the watch to find the sport. As follows: Outdoor running, outdoor walking, indoor running, indoor walking, hiking, stepper, outdoor cycling, indoor cycling, elliptical, machine, rowing machine.


**5. Pedometer**  
Steps, calorie burnt and distance will be shown on this interface.

**6. Heart rate monitoring**  
Switch to this interface, Heart rate automatic measurement, The icon starts to flash, then shows the result.


**7. Blood pressure monitoring**  
Switch to this interface, blood pressure Manual measurement, The icon starts to flash, then shows the result.




**8. ECG test**  
Switch to the interface of ECG. Click the button to the ECG test state, and put your finger against the left electrode pad. (Please make sure the device is worn on your wrist properly and the skin is moist till the result is out). The result by this manual testing on bracelet will not be saved on App. If you want to save the data, please start the ECG testing on App. The data will be shown on a graph.



9. PTT  
9.1. Put chest sticker on the connecting wire, and then install the connecting wire on watch.  
9.2. After the watch connected to the APP, find out PTT in the side button function of the watch, and then tap the button to test.  
9.3. You can see the dynamic ECG on the APP.



9.4. When the PTT running normally, there is a small ECG icon appear in the upper right corner of the APP data panel, tap the small icon to enter and start recording.



**10. Blood glucose monitoring**  
Switch to this interface, Tap the triangle button below, The icon starts to flash, then shows the result.

**11. Blood oxygen monitoring**  
Switch to this interface, Blood oxygen automatic measurement, The icon starts to flash, then shows the result.

**12. Body temperature monitoring**  
Switch to this interface, Tap the triangle button below, The icon starts to flash, then shows the result.

**13. Sleep monitoring**  
Switch to this interface, Check the duration of deep sleep, light sleep, duration of total sleep and sleep quality.

**14. Weather**  
When the device connected to the APP, the weather will automatically be sync.

**15. Notification**  
Synchronize incoming calls, text messages, QQ, WeChat and other application software messages, never missing any important messages anytime, anywhere.

**16. Alarm**  
Switch to this function interface, To find clock, when the alarm not set, it is shown on the right side.

APP find "M", Tap "EP03" to find "Alarm Settings", tap on the upper right corner "", after setting completed will show on the device.

**17. Stopwatch**  
Tap to enter the stopwatch interface, Tap the button to start or pause, To recalculate, press the reset button, and press the exit button to return.

**18. Timer**  
Find out timer in the function list, adjust the time to start timing.

**19. Music**  
When the device connected to the APP, the music player software has been installed on the phone just tap play to realize the synchronization of the device. Tap the middle play/pause button to control the phone to play music.


**20. Find phone**  
After Bluetooth connected, tap on the device find phone function, the phone will vibrate and ringing.

**21. Female**  
Save personal information as women, This function can be set according to personal circumstances.

**22. Calculator**


**23. Camera**  
After the Bluetooth connection, enter the watch camera interface, then you can shake the device to take pictures.

**24. Settings**  
Find out set in the function list, As follows: Screen/Display/Health/Monitor/Switch/Date&time/Menu style/Notifications/Language/H Band/System.




**APP functions and instructions**

**1. Dashboard**  
① Steps and time; ② More (Heart rate / Blood pressure and more); ③ Steps, distance and calories data interface; ④ Sleep data interface; ⑤ Heart rate monitoring data interface; ⑥ Blood pressure monitoring data interface; ⑦ ECG monitoring data interface; ⑧ Blood oxygen monitoring data interface; ⑨ HRV data interface; ⑩ Body temperature data interface; ⑪ Women's menstrual cycle data interface; ⑫ Blood glucose data interface.




**2. Movement**  
① Select the GPS Movement button in the App.  
② Select the map software installed on the phone.  
③ Select Start sport button.  
④ Go to the open area and select Pause button if the GPS signal is weak.




When GPS exercise completed, select movement of the device, Slide to pause exercise, while the exercise data will be automatically saved in the APP.

**3. Mine**  
**Personal information:** Tap the avatar and nickname to enter. The content are as follows: nickname, gender, date of birth, weight, height, BMI, skin color.  
**My device:** the content are as follows:  
**Message notification:** receive/close the application message reminder.  
**Alarm settings:** check or create an alarm.  
**Sedentary setting:** Set a sedentary reminder time.  
**Heart rate alert:** Set the upper limit of resting heart rate, and automatically alarm when it reaches or exceeds it.  
**Find device:** find the device.  
**Pair device:** binding the device.  
**BP Private Mode:** the private mode to obtain more accurate data, for individually, normally use general mode.  
**Blood glucose Private Mode:** turn on the switch. To the blood sugar private mode settings, type your blood sugar figure before meals and save it. After the blood sugar private mode turned on, switch on the automatic blood sugar monitoring, wear on the watch, make the watch automatically monitor for two days, and then check the data trend in the APP.  
**Screen-on Time:** Set screen time.  
**Weather setting:** Set temperature unit/switch on weather push function.  
**Switch Setting:** H: Automatic monitoring; BP: Automatic monitoring; Blood Oxygen: Nighttime monitoring; Disconnect Alert: Scientific sleep/Automatic temperature monitoring.  
**Take photos:** after connecting Bluetooth and clicking app to take photos, the watch can control the mobile phone to take photos.  
**Reset device password:** reset the device password.  
**Dial settings:** support choose to use the interested dial of the device,



or choose the "photo dial", (customize to upload pictures as you like) also can choose opposite watch face in "More Watch Faces".

**Firmware update:** current firmware information and update firmware.  
**Clear data:** Clear all data on the device and shut down.  
**Reset:** Disconnect device.



**Help:** check answers of "Cannot connect to the device", "Easy to disconnect the device", and "Message cannot be pushed".  
**Permission management:** Obtain part of mobile phone permissions during APP use.  
**Steps goal:** Set target steps.  
**Sleep:** Set target sleep time.

**Unit setting:** Set imperial or metric units.  
**APP style:** Set the APP background style.  
**About us:** check current APP version information.  
**Clear data:** Clear all data on the device and shut down.  
**Log out:** Log out of the current account.

**Parameter**  
Chip: GR5515  
Display: 1.83-inch HD color screen  
Material: Tempered glass/Zinc alloy  
Weight: 47g  
Charge time: 2hour  
Support OS: IOS 8.0 or above  
Android 4.4 or above  
Battery: 220mAh  
Waterproof: IP67  
Standby: 15days  
Bluetooth version: 5.1

**Note**  
1. Before measuring, please make sure that you have saved correct personal information on the APP, such as height, age, gender, weight, etc., and synchronize data;  
2. Sit for 5-10 minutes before measurement to keep your heartbeat calm and breathing stable;  
3. Please keep relaxed during measurement, panting, talking walking and exercise will affect the accuracy of the measurement results;  
4. Keep a correct sitting posture, make the hands flat and keep it parallel with heart.  
The data collected by this device is for reference only, not for medical purpose

