



Proshape™ One

Reciprocating Files

INSTRUCTIONS FOR USE

COMPOSITION

The instrument is made of a nickel-titanium alloy.

Proshape™ One Indications for Use

These files are used during root canal therapy for canal shaping prior to obturation. ProShape is to be used in a reciprocating motion.

Warnings

- A rubber dam system should always be used during root canal therapy
- The Proshape™ One files are non-sterile and must be sterilized before patient use. Steam autoclave is the preferred method.
- Proshape™ One files are intended as single use instruments and should be disposed of after use.

Precautions for Use

1. Straight-line access is critical for proper rotary file use.
2. Use electric torque-controlled handpiece.
3. Irrigate and lubricate the canal frequently throughout the procedure.
4. Take each rotary file to length only one time to reduce the risk of canal transportation.
5. Inspect Instruments frequently for deformaties.

Adverse Reactions

- File Separation

Sterilization

Files are non-sterile and must be sterilized prior to use. Files are single use and should not be sterilized and reused after a procedure.

- Used files should be discarded after use.

Electric Handpiece

The Proshape™ One files can only be used in an electric handpiece and motor with torque control.

Proshape™ One Technique

- Prepare access for straight line access.
- In canals where a number 10 file is snug, begin with the red(primary) Proshape to begin shaping process.
- In larger canals use the green(medium) Proshape.
- Use file until resistance is met and move to smaller file after you have reconfirmed patency.

Obturation of Canal Systems

- Use preferred obturation method to fill the canals.

Speed and Torque

- Use an electric, torque controlled handpiece with the below parameters.

Speed	Torque
300-500rpm	175-400g cm