# IR-8M and IR-16M

User manual

Art. nr: 4112 4136

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## Introduction

### Product description

IR-8M/16M 24V/120/230V is an 8/16-channel programmable IR receiver with 8/16 alternating relay switches. The receiver can be programmed with all GewaLink channels, 4096 codes and Infra-code. The receiver relays can independently be set for bistable (latching) or monostable (non-latching) function.

### Intended purpose

IR-8M/16M is intended to support people with motor, cognitive or visual impairments to remotely control electronic equipment at home. It can receive IR and GEWA radio signals from remote controls and use relay outputs to control electronic equipment. The device is intended for indoor use.

### Detector Connection

**CAUTION:** Be careful when connecting the leadings. Making the wrong connection may damage the equipment.

The wires are colour-coded. Connect the wires as in the figure. If older types of detectors without yellow leadings are connected, a resistor at each detector input must be removed (contact Abilia for further information). It is possible to connect up to two detectors.

**NOTE:** The range is reduced if the detector eye is covered or in any other way obscured by curtains or furniture for example.

### Relay Connection

The receiver has 8/16 voltage-free alternating relay switches which are connected according to the figure. See technical data for max relay load. When connecting inductive loads, motors and relays for instance, an interference suppression device should be used. To meet the isolation requirements when connecting both extra low and low voltage, there must be a distance between these leadings that corresponds to one relay output (3 terminal blocks). Alternatively, connect extra low and low voltage to different circuit boards (applies to IR-16M). The leadings for extra low voltage shall have the same isolation level as the low voltage leadings. Extra low voltage leadings can also be protected by additional isolation.

### Fuse

The in- and output on the IR-receiver must be externally protected by a 10 A fuse in a fuse box.


### Supply Voltage

The supply voltage of the receiver should be 24 V AC/DC (+/-20%), or 120/ 230 V AC. See instructions near the power supply terminal blocks.

**NOTE:** The power must be switched off when changing fuses or in any other way working on the receiver.

## Bistable Function

There is a switch for monostable or bistable relay function. Use position "On" for bistable function and "OFF" for monostable function.

# Program IR-8M/16M

When using programming or quick programming mode all 8 switches to set channel groups **must be** in "off" position on upper and lower board. Switch the programming/permanent position switch to "on". Decide whether the relay is to have a bistable or a monostable function by using the bistable/monostable function switch. Activate the desired channel on the IR transmitter at the same time as you press the programming button for relay 1. The receiver relay is activated and the indicator lamp flashes twice thus confirming the programming. Relay 1 is now programmed. Do the same thing for the second relay of the receiver but this time using programming button for Relay 2 and another transmitter channel and so on. If you are not satisfied with the selected channel repeat the process.

**NOTE:** The corresponding relay is activated during programming. GewaLink channels 0-127 (channel 62 has no function), 4096 codes and Infra-code can be programmed.









## Quick Programming

Switch the programming/permanent position switch to "On". Decide whether the relay is to have a bistable or a monostable function by using the bistable/monostable function switch. Activate the first channel of the IR transmitter at the same time as you press two optional programming buttons. All the receiver relays are activated and the indicator lamps flashes twice thus confirming the programming. It is possible to reprogram each relay afterwards from monostable to bistable function or vice versa.

## Permanent Channel Setting

The channel switches 1 - 8 can also be permanently set without programming. Set the programming/permanent setting switch to "Off". Set the channel groups according to the enclosed table. Choose monostable or bistable function using the monostable/bistable switch.

**NOTE:** All relays will work according to the chosen setting, i. e. they will have monostable or bistable function.

Switches	Channel groups	Switches	Channel groups
On 1 8 	00-07	On 1 8 	32-39
On 1 8 	08-15	On 1 8 	40-47
On 1 8 	16-23	On 1 8 	48-55
On 1 8 	24-31	On 1 8 	56-63 (Channel 62 has no function)
Filled point ● shows switch position			

## Cleaning

Clean the case of the receiver using some washing-up liquid. Use only a damp rag, not a wet one.

## Product specifications

Supply Voltage	IR-8M/16M 24V: 24V AC/DC (+/-20%) 50-60 Hz
	IR-8M/16M 120V: 120V AC, 60 Hz
	IR-8M/16M 230V: 230V AC, 50 Hz
Max connected voltage:	250V AC
Fuse:	300 mA slow
Max Current Consumption	IR-8M/16M 24V: 300 mA
	IR-8M/16M 120V: 150 mA
	IR-8M/16M 230V: 75 mA
Max Load each Relay:	7 (2) A load at 250V AC
	1 A resistive load at 24V DC
	1 A inductive load at 24V DC
	(according to standard 947-5, DC13)
Total Max Load all Relays:	2000 VA
Max Detectors:	2
Programming Channels:	GewaLink channels 0-127 (channel 62 has no function), 4096 codes and Infra-code. Infra-code applies from serial number: 382763.
Ambient Temperature Range:	-20° - +35°C, ( -4° - +95°F)
User environment:	Indoor (IP-20)
Estimated lifetime:	5 years

## Storage and disposal

### Storage

Keep the product indoors at room temperature and away from moist.

### Recycling Waste Electrical and Electronic Equipment

The product contains electronic components that must be recycled. When the product is no longer of use, hand it in to a recycling centre for waste electrical and electronic equipment.

# Important information

Before you use the product, read and understand the contents of this user manual.

## Safety definitions

 **WARNING:** Risk of injury or death if the instructions are not followed.

 **CAUTION:** Risk of damage to the product if the instructions are not followed.

**NOTE:** Information that is necessary in a given situation.

## Serious incidents

If a serious incident related to the product occurs, report this to the manufacturer and the competent authority of your country. A serious incident directly or indirectly led, might have led or might lead to the death of a patient, user or other person; the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health threat.

## Conformance with requirements for medical devices

IR-8M/16M is marketed as a technical aid for people with disabilities. It complies with all necessary requirements, regulations and directives for medical devices.

## Product liability

The product is developed and risk-assessed according to ISO 14971. The manual and the technical specifications are made according to the safety assessments in the risk analysis. Always consider the safety when using the product. Abilia AB does not take responsibility for any consequences following incorrect use or installation of the product. Any form of damage or tampering with the product invalidates the warranty and disclaims the manufacturer's liability. The estimated lifetime of the product applies when used in accordance with its intended use and maintenance instructions.

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Abilia AB, Råsundavägen 6, 169 67 Solna, Sweden
Abilia UK Ltd, Unit 12 Buckingham Business Park, Swavesey,
Abilia AB, Råsundavägen 6, Solna, Sweden