

CLEANING & CARE INSTRUMENTS



Use circular motion

The fiber optic light transmitting surfaces on both ends of the fiber optic clip or wand should be wiped clean with isopropyl alcohol before autoclaving or chemiclaving.



In the event that a fiber optic clip becomes loose, *gentle* pressure can be applied with pliers to tighten clips.

STERILIZATION

Your fiber optic instrument may be autoclaved or chemiclaved. Do not exceed 275°F (135 C). Please note that exposure to heat will reduce the life of your fiber optic instrument. The average life expectancy of fiber optics is 18 to 24 months with autoclaving or chemiclaving. This time could be longer or shorter depending on frequency and temperature of autoclave or chemiclave.