





Baseplate Placement Tips

Placing an Adhesive Tracheostoma Baseplate

Step 1 Clean

Gently clean the skin around the tracheostoma with mild soap and warm water. Rinse and completely dry the skin.

Step 2 Remove

Remove adhesive baseplate from backing

Step 3
Find your centre

In most cases, try to centre the housing directly over the stoma and as close to the stoma as possible with placement.

Step 4
Press

Working from the inside of the baseplate / stoma to the outside, gently press the adhesive down around the stoma, using fingers to massage out any folds, creases or bubbles.



Step 5 Remove To remove baseplate, consider using a Remove Wipe. Wipe the Remove wipe around the outside edges of the baseplate and slowly peel the edges of the baseplate back until housing is entirely off.



- The closer the baseplate to the stoma, the better
- With deep-set stomas, consider "folding" the baseplate on
- Massage the baseplate from inside -> out
- Warm the baseplate with body heat prior to application
- Patient must be able to replicate. Have patient, or have the family take a photo to help with independent placement.



Wipes and Extras

- Extra Stick: Silicone Glue & Adhesive tape discs
- Skin Irritation: Skin prep / barrier wipes
- Oily Skin: Alcohol Wipes
- Difficulty removing: Remove Wipes



