# Blom-Singer®

voice restoration systems

## Blom-Singer® Gel Cap Insertion System

USER Instructions For Use

**INHEALTH**®  
TECHNOLOGIES

**R1**

37575-01E

37575-01E Effective March 2014 / Gültig ab März 2014 / Efectivo a partir de marzo de 2014 / Date d'entrée en vigueur : mars 2014 / Valido da marzo 2014 / Met ingang van maart 2014 / Vigente em de março de 2014 / V platnosti od března 2014 / Wchodzi w życie w marcu 2014 r.

Blom-Singer is a registered trademark of Hansa Medical Products, Inc.  
InHealth Technologies is a registered trademark of Helix Medical, LLC.  
© 2014 Helix Medical, LLC (InHealth Technologies)

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Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)

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### **BLOM-SINGER® GEL CAP INSERTION SYSTEM**

for Blom-Singer Low Pressure Voice Prosthesis

#### **Introduction**

*Please refer to the diagrams located at the back of this instruction manual. A glossary of terms appears at the end of this chapter.*

The Gel Cap Insertion System is designed to provide an easier and more comfortable method for inserting a 16 Fr. or 20 Fr. Blom-Singer Low Pressure Voice Prosthesis. Using the specially designed gel cap loading tool, the user places a capsule (gel cap) over the tip of the voice prosthesis to enable atraumatic insertion into the tracheoesophageal puncture.

#### **INDICATIONS**

The Blom-Singer Gel Cap Insertion System is designed for use by individuals who have had a laryngectomy and use a Blom-Singer 16 Fr. or 20 Fr. Low Pressure Voice Prosthesis.

#### **CONTRAINDICATIONS**

The Blom-Singer Gel Cap Insertion System is a medical product that should be used only by individuals with experience and training in its use and care.

#### **HOW SUPPLIED**

The Gel Cap Insertion System is supplied in a dual 16 Fr./20 Fr. component kit for use with a 16 Fr. or 20 Fr. Blom-Singer Low Pressure voice prosthesis. Store in a cool, dry place and out of direct sunlight.

#### **PRODUCT DESCRIPTION**

The Gel Cap Insertion System (diagram 1) is a dual 16 Fr./20 Fr. component kit that contains all the necessary items to insert a 16 Fr. or 20 Fr. Blom-Singer Low Pressure voice prosthesis.

#### **Dual Gel Cap Loading Tool**

The dual gel cap loading tool is designed to facilitate the loading of the voice prosthesis into the gel cap. At one end of the tool

is a 16 Fr. hexagon insertion hole and at the opposite end is a 20 Fr. round insertion hole (diagram 1a). The user selects the appropriate hole that corresponds to the voice prosthesis diameter being used.

## **Gel Cap**

The gel cap used for insertion of the prosthesis is made of a vegetable base gelatin. Consult with your physician as individual reactions to these materials may occur.

The gel cap (capsule) is supplied in a 16 Fr. and 20 Fr. size for each kit (diagram 1b). The user selects the appropriate gel cap that corresponds to the voice prosthesis diameter being used.

The function of the gel cap is to reduce the trauma and force that may be associated with low pressure voice prosthesis insertion by providing a smooth, rounded shape to the tip end of the voice prosthesis enabling easy entry into the tracheoesophageal puncture. The gel cap decreases the profile of the tip of the voice prosthesis by enclosing the esophageal retention flange (diagram 2a) inside the gel cap so that the prosthesis tip is fully encapsulated. This eliminates the presence of the esophageal retention flange during insertion. Replacement gel caps are available separately.

## **Pushrod**

The pushrod helps to guide the voice prosthesis through the hole of the dual gel cap loading tool and into the gel cap. Each kit is supplied with a 16 Fr. and 20 Fr. pushrod (diagram 1c). The user selects the appropriate pushrod that corresponds to the voice prosthesis diameter being used.

## **INSTRUCTIONS FOR USE**

The Blom-Singer Gel Cap Insertion System is designed for use only with Blom-Singer Low Pressure Voice Prostheses. Do not use this product with any other type devices.

**Note:** The Gel Cap Insertion System does not function with special order Low Pressure voice prostheses.

---

## Loading the Gel Cap

*Please refer to the diagrams at the back of this instruction manual.*

*The following procedure is provided by Eric D. Blom, Ph.D. as instructions for loading the voice prosthesis into the gel cap.*

The Gel Cap Insertion System, voice prosthesis, and user's hands should **always be thoroughly clean** to help avoid introducing contaminants into the tracheostoma, puncture, or esophagus.

**Note:** Prior to use of the Low Pressure voice prosthesis, verify that the flapper valve mechanism is intact and working properly. The valve should close flat inside the prosthesis.

**Important:** Be certain the tip of the voice prosthesis, gel cap loading tool, pushrod, and your hands are completely dry before beginning the following procedure. Failure to do so may cause premature dissolving of the gel cap.

Before beginning the following procedure, remove the items from the kit that correspond in size to your voice prosthesis, e.g., a 16 Fr. prosthesis uses the 16 Fr. gel caps and 16 Fr. pushrod.

1. Insert the neck strap of the voice prosthesis (diagram 2c) into the appropriate size hole of the dual gel cap loading tool on the side marked "top" (diagram 3). Be sure to use the correct hole size that corresponds to the diameter of your voice prosthesis, e.g., a 16 Fr. diameter prosthesis fits into the 16 Fr. hexagon hole and a 20 Fr. diameter prosthesis fits into the 20 Fr. round hole.
2. Holding the neck strap, gently pull the prosthesis down and through the hole (diagram 4) until the esophageal retention flange is positioned over the hole (diagram 5).
3. Grasp the tubular portion (body, diagram 2b) of the voice prosthesis at the base of the neck strap and very slowly and gently begin to pull the prosthesis down further into the hole (diagram 6) until the esophageal flange on the voice prosthesis begins to fold forward inside the hole as shown in diagram 7. Do not pull the prosthesis beyond this point as over-pulling will cause the voice prosthesis to be pulled completely through the gel cap loading tool. This will require starting again with step 1 above.
4. From the gel cap container, remove the 16 Fr. or 20 Fr. gel cap that corresponds to the diameter of your voice prosthesis. (**Note:**

Discard the printed end and use unprinted cap.) Place the open end of the gel cap in the hole of the loading tool with the folded flange of the prosthesis (diagram 8). The gel cap should fit into the groove and cover the prosthesis tip. Verify that the gel cap is securely in place.

5. Fully insert the pushrod, of the appropriate size, into the bottom hole of the voice prosthesis at the base of the neck strap. Place a finger on the tip of the gel cap to hold it in place while simultaneously pushing the voice prosthesis from the bottom (diagram 9). Push gently until you can see that the folded esophageal retention flange is fully positioned in the gel cap.

6. Remove your fingertip from the gel cap and push gently with the pushrod in the direction indicated by the arrow in diagram 10.

7. Continue to push the prosthesis until it is completely out of the gel cap loading tool (diagram 11). Hold the prosthesis at the base of the neck strap and remove the the pushrod from the prosthesis.

**Note: Never use the pushrod as an inserter. Always use the inserter provided with your prosthesis.**

The voice prosthesis is now fully loaded in the gel cap and ready for insertion in the usual manner.

## **Voice Prosthesis Insertion**

**Caution:** Care should be exercised to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the prosthesis. In the unlikely event that this should occur, the user should attempt to cough the voice prosthesis out of the windpipe and immediately seek medical attention if unsuccessful.

**Insertion of the prosthesis in the puncture should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma.**

1. Handle the voice prosthesis by the base of the neck strap to avoid touching the portion that inserts into the puncture. Place the tip of the inserter into the open end of the prosthesis (diagram 12).

2. Attach the strap on the prosthesis securely over the safety

---

peg on the inserter as shown in diagram 13. This should prevent inadvertent dislodgment of the prosthesis off the inserter during insertion into the puncture.

**3.** Dip the gel-cap tip of the voice prosthesis in a generous coating of an oil-free, water-soluble lubricant (diagram 14).

**Note:** Do not use petroleum-based products such as Vaseline® to lubricate the prosthesis as these products can damage silicone.

**4.** Refer to the user information data sheet packaged with the voice prosthesis for detailed instructions on how to properly and safely insert the voice prosthesis into the tracheoesophageal puncture.

Be certain that the gel cap is securely installed on the tip of the voice prosthesis prior to insertion into the tracheoesophageal puncture.

**5.** Hold the prosthesis in a fully inserted position for at least two minutes. The gel cap will dissolve inside the esophagus within 2—3 minutes after inserting the prosthesis. The esophageal retention flange, on an appropriately sized prosthesis that is positioned correctly, will unfold and the prosthesis will be in its proper position in the esophagus.

**6.** Remove the inserter from the voice prosthesis per instructions noted in the user information data sheet packaged with the prosthesis.

## **CLEANING AND CARE INSTRUCTIONS**

Wash the dual gel cap loading tool and the pushrod thoroughly with a mild detergent and rinse under a strong flush of tap water. Dry the devices with sterile gauze. Never use facial or toilet tissue to dry the dual gel cap loading tool, inserter, or voice prosthesis, as particles of lint or fabric could be inhaled into the airway upon insertion of the prosthesis.

Once the dual gel cap loading tool and the pushrod have been cleaned and thoroughly dried, they may be stored in a clean, resealable plastic bag.

If there are tears, cracks, or structural damage to the Gel Cap Insertion System or prosthesis, discontinue use.

**Caution:** Do not use solvents or petroleum-based products

to clean or lubricate the voice prosthesis. These materials may damage silicone.

For complete cleaning and care instructions for the voice prosthesis, please refer to the user information data sheet packaged with each prosthesis.

## **WARNINGS AND PRECAUTIONS**

The Gel Cap Insertion System is a medically prescribed device. Instruction in its use must be provided by a qualified medical professional trained in the use of this particular device.

Foreign objects should never be inserted into the voice prosthesis by the patient. Inserting objects other than the supplied device may cause dislodgment of the voice prosthesis or its components.

The Blom-Singer Gel Cap Insertion System is designed for use with Blom-Singer Low Pressure Voice Prostheses only. Do not use this product with any other type devices.

The Gel Cap Insertion System does not function with special order Low Pressure voice prostheses.

The tip of the Blom-Singer Low Pressure Voice Prosthesis, the gel cap loading tool, the pushrod, and the user's hands must be completely dry prior to using the Gel Cap Insertion System in order to prevent premature dissolving of the gel cap.

Never use the pushrod as a voice prosthesis inserter. Always use the inserter provided with the prosthesis.

When using the Blom-Singer Gel Cap Insertion System, always wait at least two minutes following insertion of the voice prosthesis into the tracheoesophageal puncture before detaching the strap from the safety peg and withdrawing the inserter.

Always inspect the voice prosthesis prior to insertion. Do not use a voice prosthesis that shows obvious signs of damage.

Do not use solvents or petroleum-based products to clean or lubricate the voice prosthesis. These materials may damage silicone.

Care should be exercised when inserting or removing the voice prosthesis to avoid displacement of the prosthesis which could result in aspiration (inhalation) of the voice prosthesis.

---

In the unlikely event that this should occur, the patient should attempt to cough the voice prosthesis out of the windpipe and immediately seek medical attention if unsuccessful.

If the voice prosthesis is accidentally dislodged from the puncture, immediately place a Blom-Singer tracheoesophageal puncture dilator, of the appropriate diameter, in the puncture to keep the puncture from closing.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician or a licensed speech pathologist.

## **COMPLICATIONS**

Although rare, the following complications have been identified to occur with the use of the Blom-Singer Gel Cap Insertion System or the Blom-Singer Low Pressure Voice Prosthesis. They include: tearing or other damage to the voice prosthesis from improper use of the gel cap loading tool; discomfort upon insertion of the voice prosthesis; stoma contamination or sepsis, which may require removal of the prosthesis and/or appropriate antibiotics; accidental aspiration (inhalation) of the voice prosthesis into the airway, which may require removal by a physician; aspiration (inhalation) of liquids through the prosthesis into the airway, which may cause the patient to cough; occasional extrusion of the voice prosthesis, requiring replacement after dilation of the tracheoesophageal puncture and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the tracheoesophageal puncture; intractable leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; dysphagia, or difficulty swallowing; tearing or other damage to the voice prosthesis from improper use; yeast deposits causing leakage through the prosthesis or valve incompetence; accidental ingestion of the prosthesis into the esophagus and/or GI tract.

## **ORDERING INFORMATION**

### **USA**

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

### **International**

Please contact our customer service department for a distributor referral.

### **Consumer Affairs**

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: [productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)

### **RETURNED GOODS POLICY**

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer's tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 14 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

### **LIMITED WARRANTY, LIMITED REMEDY AND LIMITATION OF LIABILITY**

InHealth Technologies, a Division of Helix Medical, LLC, warrants that the Product will be free from defects in workmanship and materials on the date of purchase. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES OR REPRESENTATIONS, EXPRESS OR IMPLIED, ARISING BY LAW OR CUSTOM, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SAID IMPLIED WARRANTIES BEING HEREBY EXPRESSLY

---

DISCLAIMED. Purchaser's sole and exclusive remedy for a valid warranty claim is replacement of the Product. This remedy does not include the cost of installation, removal, dismantling, delivery or reinstallation. InHealth shall have no liability to the extent Products are or have been modified or subject to misuse, abuse or improper storage, installation or maintenance or damaged in transit. No returns will be accepted without prior approval and authorization. INHEALTH SHALL NOT BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT, CONSEQUENTIAL, PUNITIVE, OR EXEMPLARY DAMAGES OF ANY KIND.

## **GLOSSARY**

### **Airway**

The air passage by which air reaches the lungs.

### **Candida albicans**

A yeast fungus that is a commonly occurring condition that can infect the mouth or the skin and may multiply and overgrow on to the voice prosthesis.

### **Clinician**

A health professional, such as a physician, licensed speech-language pathologist, or nurse who is involved in clinical practice.

### **Contraindications**

A reason that makes it inadvisable to prescribe a particular procedure or treatment.

### **Esophagus**

The muscular membranous tube for the passage of food from the pharynx to the stomach; the "food pipe."

### **Fistula; TE-fistula**

An abnormal duct or passage that connects an abscess or cavity to the body surface; a small opening between the trachea and the esophagus.

### **Indications**

A reason to prescribe a drug or perform a procedure.

### **Respiration**

The act or process of inhaling and exhaling; breathing; moving air into and out of the lungs. Also called ventilation.

**Respiratory Tract (laryngectomy)**

The air passage from the tracheostoma to the air sacs of the lungs, including the trachea, and bronchi through which air passes during breathing.

**Tracheoesophageal puncture (TEP)**

An opening made by a surgeon between the trachea and the esophagus. The opening is made at the time of initial surgery or later so that a voice prosthesis can be placed into this passage.

**Trachea**

The airway through which respiratory air passes to and from the lungs (windpipe).

**Tracheoesophageal (TE)**

Pertaining to the trachea and esophagus.

**Tracheostoma**

An opening through the neck and into the trachea. It is through this opening that the laryngectomye breathes.

**Voice prosthesis**

A one-way speech valve that is inserted into the TEP to allow air to pass from the trachea through the prosthesis and into the esophagus to produce sound.

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











---

## BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAFIA / BIBLIOGRAFIE

Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação. / Další odkazy na literaturu jsou k dispozici na požádání. / Dodatkowe informacje dostępne na żądanie.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.

- 
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
  11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
  12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
  13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
  14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
  15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
  16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

**Helix Medical, LLC**

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

CE 0344

**EMERGO EUROPE**

Molenstraat 15  
2513 BH, The Hague  
The Netherlands