Blom-Singer®

voice restoration systems

Blom-Singer® Indwelling Valved Insert and Plug Insert

USER Instructions For Use

INHEALTH®
TECHNOLOGIES

R1

37582-01E

37582-01E Effective March 2014 / Gültig ab März 2014 / Efectivo a partir de marzo de 2014 / Date d'entrée en vigueur : mars 2014 / Valido da marzo 2014 / Met ingang van maart 2014 / Vigente em de março de 2014 / V platnosti od března 2014 / Wchodzi w życie w marcu 2014 r.

Blom-Singer is a registered trademark of Hansa Medical Products, Inc.
InHealth Technologies is a registered trademark of Helix Medical, LLC.
© 2014 Helix Medical, LLC (InHealth Technologies)

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BLOM-SINGER® INDWELLING VALVED INSERT AND PLUG INSERT

Introduction

Please refer to the diagrams located at the back of this instruction manual. A glossary of terms appears at the end of this chapter.

The Blom-Singer Indwelling Valved Insert and Plug Insert are optional accessory devices that are used in conjunction with a Blom-Singer Indwelling voice prosthesis. These Indwelling accessory devices are designed for temporary use only to temporarily stop leakage of the voice prosthesis until medical attention can be obtained.

These Indwelling accessory devices should only be used by individuals who have received proper training in the effective management and care of Blom-Singer voice prostheses. Management includes, but is not limited to, handling of the device, dislodgement of the device, and cleaning of the device.

Your doctor will discuss with you the risks and benefits of using this device. The information in this manual is not intended to be used in place of consultation or instructions of a medical professional.

After placement of this Indwelling accessory device, contact your physician/clinician as soon as possible to have the voice prosthesis performance evaluated.

These products are intended for single patient use only. Not for use between patients. See Warnings, Precautions and Complications.

INDICATIONS

The Blom-Singer Plug Insert may be used on a temporary basis in conjunction with a Blom-Singer Indwelling voice prosthesis to prevent esophageal leakage of fluids through the voice prosthesis into the trachea during swallowing.

The Blom-Singer Indwelling Valved Insert Accessory Device

may be used in conjunction with the Blom-Singer Indwelling voice prosthesis to prevent leakage of fluids through the voice prosthesis into the trachea during swallowing and to maintain voice restoration on a temporary basis.

CONTRAINDICATIONS

The Plug Insert is a medical product and should be used only by individuals with experience and training in its use and care.

The Plug Insert is not intended to permanently replace a properly functioning voice prosthesis. It is not intended to be placed directly in a surgically created fistula between the esophagus and trachea. It is not intended to be used in combination with any other voice prosthesis or device other than the Blom-Singer ADVANTAGE or CLASSIC Indwelling voice prostheses.

The Blom-Singer Indwelling Valved Insert accessory device is a medical product and should be used only by individuals with experience and training in its use and care.

The Blom-Singer Indwelling Valved Insert accessory device is intended to be temporarily used only as an optional accessory device with the Blom-Singer Indwelling Voice Prosthesis. It is not intended to be placed in a surgically created fistula between the esophagus and trachea. It is not intended to be used in combination with any other voice prosthesis or device.

The Indwelling Valved Insert accessory device is not intended for use with the Indwelling voice prosthesis size 6mm due to device size constraints. The Indwelling Valved Insert will not fit properly into the Indwelling size 6mm.

Use of the Blom-Singer Adjustable Tracheostoma Valve (ATSV) is contraindicated (must not be used) when the Indwelling Valved Insert is placed within the Indwelling voice prosthesis. Use of the ATSV may result in accidental displacement of the Indwelling Valved Insert or Indwelling voice prosthesis which could result in aspiration (inhaling) of the valved insert or voice prosthesis.

HOW SUPPLIED

The Blom-Singer Indwelling Valved Insert and Plug Insert are supplied separately as accessory devices for use with a Blom-Singer Indwelling voice prosthesis. An inserter tool is included

with each device. An indwelling insert holder is included with the Indwelling Valved Insert.

PRODUCT DESCRIPTION

The Blom-Singer Indwelling Valved Insert and Plug Insert are made of silicone, a biocompatible material used in many medical products. They are available in two models as shown in the chart below.

Description	Size		Allows Voicing	Use with Blom-Singer Prosthesis
	16 Fr.	20 Fr.		
Indwelling Valved Insert		•	•	20 Fr. CLASSIC 7mm–25mm
Plug Insert	•	•		CLASSIC, ADVANTAGE 4–25mm

Indwelling Valved Insert

The Indwelling Valved Insert (diagram 1) is designed to prevent esophageal leakage through the voice prosthesis and to enable voice production until such time that the voice prosthesis and TEP is evaluated by a clinician. The Indwelling Valved Insert consists of a body that holds the valve assembly (diagram 1a), a neck strap with safety hole (diagram 1b), and an indwelling insert holder (diagram 1c). The Indwelling Valved Insert may be inserted by the user.

The Indwelling Valved Insert must not be used in conjunction with the Blom-Singer Adjustable Tracheostoma Valve (ATSV).

Plug Insert

The Plug Insert (diagram 2) is designed to prevent esophageal leakage through the voice prosthesis until such time that the voice prosthesis and TEP is evaluated by a clinician; it does not allow voicing. The Plug Insert consists of a plug tip (diagram 2a) and neck strap with safety hole (diagram 2b). The Plug Insert may be inserted by the user.

The Plug Insert must not be used in conjunction with the Blom-Singer Adjustable Tracheostoma Valve (ATSV).

WARNINGS AND PRECAUTIONS

InHealth relies on the qualified, trained medical professional to instruct patients on the use and care of this device and to provide patients with the Blom-Singer instructions for use supplied with the product.

Should the device show signs of malfunction, contact your physician/clinician as soon as possible to have the device performance evaluated.

If there are tears, cracks, or structural damage to the device, discontinue use.

Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP be evaluated by your clinician at regular intervals.

These Indwelling accessory devices are designed for temporary use only and are not intended to replace the Indwelling voice prosthesis. After placement of the device, contact your physician/clinician as soon as possible.

Do not detach or in any way modify the neck strap on the Indwelling Valved Insert or Plug Insert.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or cause the device to not work properly.

Never insert any cleaning device while the Indwelling Valved Insert or Plug Insert is in position in the neck. Always clean the device after it has been removed.

Never use facial or toilet tissue to dry this Indwelling accessory device as particles of lint or fabric may be deposited onto the device that can be inhaled into the airway upon reinsertion.

The Blom-Singer Adjustable Tracheostoma Valve (ATSV) must not be used when the Indwelling Valved Insert or Plug Insert is placed within the voice prosthesis. Use of the ATSV in conjunction with these Indwelling accessory devices may result in accidental displacement of the device or voice prosthesis, which could result

in inhalation of the device or voice prosthesis.

Never force the Indwelling Valved Insert or Plug Insert into the voice prosthesis. If the device does not insert easily on the first attempt, do not continue to try to insert it. Should you experience repeated difficulty in managing the Indwelling Valved Insert or Plug Insert, contact your physician/clinician as soon as possible.

The Plug Insert will not fully insert into a 4 mm size voice prosthesis due to the shorter length of the voice prosthesis shaft. A small portion of the Plug will extend from the entrance of the voice prosthesis. It is very important to not insert the Plug beyond this point of palpable resistance. Doing so could result in pushing the voice prosthesis into the esophagus.

Voice Prosthesis Dislodgement

Care should be exercised when inserting or removing the Indwelling Valved Insert, Plug Insert, or inserter tool to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the prosthesis. In the unlikely event that this should occur, you should attempt to cough the voice prosthesis out of the windpipe. Immediately seek medical attention if unsuccessful.

The neck strap should not, under any circumstances, be cut or removed from the Indwelling Valved Insert or Plug Insert. The neck strap is intended to help prevent accidental aspiration (inhaling) of the device should the device become dislodged.

If the voice prosthesis is accidentally dislodged from the TE puncture, immediately place a Blom-Singer Puncture Dilator, or a suitable device as recommended by your doctor, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. Contact your doctor immediately. The voice prosthesis should be reinserted within 24 hours.

Never insert foreign objects into the voice prosthesis. Inserting objects other than the supplied device may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or inhaled.

Never insert any cleaning device while the Indwelling Valved Insert or Plug Insert is in position in the neck. Always clean the device after it has been removed.

Voice Prosthesis Removal

Never attempt to remove the Indwelling voice prosthesis by yourself or allow anyone other than a qualified, trained medical professional to remove the voice prosthesis.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or a licensed speech pathologist.

COMPLICATIONS

Although rare, the following complications have been identified to occur with silicone prostheses of the Blom-Singer type. They include:

- stoma contamination or bacterial infection, which may require removal of the prosthesis and/or appropriate antibiotics;
- accidental inhalation of the voice prosthesis into the airway, which may require removal by a physician;
- inhalation of liquids through the prosthesis into the airway, which may cause the patient to cough;
- occasional movement out of position by the voice prosthesis, requiring replacement after enlargement of the TE puncture and additional supervision of the stoma care regimen;
- puncture enlargement resulting in leakage of fluids around the prosthesis;
- inflammatory reaction around the puncture site and formation of scar tissue;
- dislodgment of the voice prosthesis and subsequent closure of the TE puncture;
- chronic and persistent leakage around the voice prosthesis, requiring surgical revision or closure of the puncture;
- difficulty speaking or swallowing;
- tearing or other damage to the voice prosthesis from improper use;
- yeast deposits causing leakage through the prosthesis or inability of the valve to function properly;
- accidental swallowing of the prosthesis into the esophagus and/or gastrointestinal (GI) tract.

INSTRUCTIONS FOR USE

Your clinician will instruct you on the routine use and care of this Blom-Singer device.

Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP be evaluated by your clinician at regular intervals.

If there are tears, cracks, or structural damage to the device, discontinue use.

The Indwelling accessory device and user's hands **should always be thoroughly** clean to prevent contaminants from entering into the tracheostoma, puncture, or esophagus.

Warning: Care should be exercised when inserting or removing the Indwelling Valved Insert, Plug Insert, or inserter tool to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the prosthesis. In the unlikely event that this should occur, you should attempt to cough the voice prosthesis out of the windpipe. **Immediately** seek medical attention if unsuccessful.

The following procedural instructions are provided by Eric D. Blom, Ph.D. Please refer to the diagrams located at the back of this instruction manual.

Insertion of the Indwelling Valved Insert / Plug Insert

Note: Indwelling Valved Insert: Prior to use of the Indwelling Valved Insert, check the valve mechanism to be sure it is intact and working properly. The valve should close flat inside the device.

Insert the Indwelling Valved/Plug Insert only while using a close-up mirror and a bright light focused directly into the stoma. Be certain that your hands are clean.

1. Place the tip of the inserter tool into the open end of the Indwelling Valved/Plug Insert (diagram 3). (Avoid squeezing or damaging the valve in the tip of the Indwelling Valved Insert.) Handle the device by the base of the strap to avoid touching the portion that inserts into the voice prosthesis.

2. Attach the safety hole in the neck strap on the Indwelling Valved/Plug Insert securely over the safety peg on the inserter tool as shown in diagram 4. This should prevent inadvertent dislodgment of the Indwelling Valved/Plug Insert off the inserter tool during insertion into the voice prosthesis. (Indwelling Valved Insert: Place the indwelling insert holder with attached Indwelling Valved Insert around the neck.)

3. Hold the end of the inserter tool as shown in diagrams 5A and 5B, and align the top of the Indwelling Valved/Plug Insert partially inside the voice prosthesis with the neck strap oriented upwards, diagrams 6A and 6B.

4. Carefully and gently insert the Indwelling Valved/Plug Insert until it is fully inserted and positioned flush with the voice prosthesis (diagram 7).

Never force the Indwelling Valved Insert or Plug Insert into the voice prosthesis. **Caution:** If the device does not insert easily on the first attempt, do not continue to try to insert it. Should you experience repeated difficulty in managing the Indwelling Valved Insert or Plug Insert, contact your physician/clinician as soon as possible.

Caution: Plug Insert: The Plug Insert will not fully insert into a 4mm size voice prosthesis due to the shorter length of the voice prosthesis shaft. A small portion of the Plug will extend from the entrance of the voice prosthesis as shown in diagram 8. It is very important to not insert the Plug beyond this point of noticeable resistance. Doing so could result in pushing the voice prosthesis into the esophagus. See Warning in step 6 below.

5. After confirming the Indwelling Valved/Plug Insert is securely in place and fully inserted, the inserter tool should be carefully removed.

6. Hold the inserter tool and the neck strap together firmly in one hand. Carefully lift the neck strap off the safety peg on the inserter tool. Hold the neck strap against your neck while slowly and gently withdrawing the inserter tool with a twisting motion (diagram 9).

Warning: Caution should be exercised when removing the inserter tool to avoid accidental displacement of the Indwelling

Valved, Plug Insert or voice prosthesis which could result in aspiration (inhaling) of the device or voice prosthesis. Should aspiration occur, attempt to cough the device and/or voice prosthesis out of the windpipe. Seek medical attention **immediately** if unsuccessful.

7. Apply a strip of hypo-allergenic, medical adhesive tape over the neck strap to prevent strap movement as shown in diagram 10. **Warning:** Do not detach or in any way modify the neck strap on the Indwelling Valved/Plug Insert. **The neck strap should not, under any circumstances, be cut or removed from the Indwelling Valved/Plug Insert.** The neck strap is intended to help prevent accidental aspiration (inhaling) of the device should the device become dislodged.

Removal of the Indwelling Valved Insert / Plug Insert

Removal of the Indwelling Valved/Plug Insert should only be done while using a close-up mirror and a bright light focused directly into the stoma so that the open end of the prosthesis is clearly visible. Be certain that your hands are clean.

1. Stabilize the voice prosthesis by holding it carefully and firmly with tweezers to avoid displacement (diagram 11).
2. Grasp the neck strap of the Indwelling Valved/Plug Insert and gently pull it out from the voice prosthesis.

Cleaning and Care Instructions

Your clinician will instruct you on the routine use and care of this Blom-Singer Indwelling accessory device.

The Indwelling accessory device and user's hands **should always be thoroughly** clean to prevent contaminants from entering into the tracheostoma, puncture, or esophagus.

Always use **gentle** handling and pressure when cleaning the Indwelling accessory device to avoid deterioration of the device and damaging internal components which could cause leakage.

Do not insert foreign objects into the Indwelling voice prosthesis, Indwelling Valved Insert, or Plug Insert.

Precaution: Never insert any cleaning device while the Indwelling accessory device is in position in the neck. Always clean the device **after** it has been removed.

If there are tears, cracks, or structural damage to the device, discontinue use.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or cause the device to not work properly.

The Indwelling accessory device should be cleaned, rinsed and air-dried after each use. Place the device in 3% hydrogen peroxide after each use. Shake to remove any residual liquid. Allow the Indwelling Insert to dry in a well-lit, aerated area. If desired, these accessories can be washed in soap and water, rinsed and air-dried.

Never use facial or toilet tissue to dry an Indwelling accessory device as particles of lint or fabric may be deposited onto the device that can be inhaled into the airway upon reinsertion.

Once the Indwelling accessory device has been cleaned and dried completely, it should be stored in a clean, resealable plastic bag. Always handle with clean hands and avoid touching portions of the device that will enter the voice prosthesis.

Indwelling Valved Insert

Inspect the Indwelling Valved Insert routinely for structural damage. Warning: When the valve fails to close completely a few drops of fluid may pass back through the valve from the esophagus to the trachea, which may cause coughing. The valve closure can be verified by looking at it with a mirror, while it is still positioned in the puncture, to see if fluids leak through it during swallowing.

If it is determined that the user has a Candida condition, routine cleaning of the Indwelling Valved Insert will retard the development of yeast deposits. After removing the Indwelling Valved Insert, and after following the cleaning instructions above, place the Indwelling Valved Insert in a small container of fresh 3% hydrogen peroxide overnight. After soaking, the Indwelling Valved Insert should be thoroughly rinsed again with tap water, dried completely, and then placed in a clean, resealable plastic bag until its next use. Note: Soaking in hydrogen peroxide will not remove yeast that is already present on the device. It is a recommended cleaning measure only. If yeast is present on the device and causing it to leak, the Indwelling Valved Insert should be replaced.

Please refer to the instructions for use supplied with the Indwelling voice prosthesis for complete instructions on proper care and maintenance of the voice prosthesis.

SPECIAL ORDER PRODUCTS

If this instruction manual accompanies a Special Order Product, there may be differences in the physical characteristics between the product and the product descriptions described in this instruction manual. These differences will not affect the safety or desired beneficial effect of the special order product. Special Order products are nonreturnable.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ON-LINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

International

Please contact our customer service department for a distributor referral.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: productcomplaints@inhealth.com

RETURNED GOODS POLICY

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer's tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 14 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

LIMITED WARRANTY, LIMITED REMEDY AND LIMITATION OF LIABILITY

InHealth Technologies, a Division of Helix Medical, LLC, warrants that the Product will be free from defects in workmanship and materials on the date of purchase. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES OR REPRESENTATIONS, EXPRESS OR IMPLIED, ARISING BY LAW OR CUSTOM, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SAID IMPLIED WARRANTIES BEING HEREBY EXPRESSLY DISCLAIMED. Purchaser's sole and exclusive remedy for a valid warranty claim is replacement of the Product. This remedy does not include the cost of installation, removal, dismantling, delivery or reinstallation. InHealth shall have no liability to the extent Products are or have been modified or subject to misuse, abuse or improper storage, installation or maintenance or damaged in transit. No returns will be accepted without prior approval and authorization. INHEALTH SHALL NOT BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT, CONSEQUENTIAL, PUNITIVE, OR EXEMPLARY DAMAGES OF ANY KIND.

GLOSSARY

Airway

The air passage by which air reaches the lungs.

Candida albicans

A yeast fungus that is a commonly occurring condition that can infect the mouth or the skin and may multiply and overgrow on to the voice prosthesis.

Clinician

A health professional, such as a physician, licensed speech-language pathologist, or nurse who is involved in clinical practice.

Contraindications

A reason that makes it inadvisable to prescribe a particular procedure or treatment.

Esophagus

The muscular membranous tube for the passage of food from the pharynx to the stomach; the "food pipe."

Fistula; TE-fistula

An abnormal duct or passage that connects an abscess or cavity to the body surface; a small opening between the trachea and the esophagus.

Indications

A reason to prescribe a drug or perform a procedure.

Occlusion (Occlude)

To block, such as the stoma with a finger, thumb, or ATSV in order to produce tracheoesophageal voice.

Respiration

The act or process of inhaling and exhaling; breathing; moving air into and out of the lungs. Also called ventilation.

Respiratory Tract (laryngectomy)

The air passage from the tracheostoma to the air sacs of the lungs, including the trachea, and bronchi through which air passes during breathing.

Tracheoesophageal puncture (TEP)

An opening made by a surgeon between the trachea and the esophagus. The opening is made at the time of initial surgery or later so that a voice prosthesis can be placed into this passage.

Trachea

The airway through which respiratory air passes to and from the lungs (windpipe).

Tracheoesophageal (TE)

Pertaining to the trachea and esophagus.

Tracheostoma

An opening through the neck and into the trachea. It is through this opening that the laryngectomy breathes.

Voice prosthesis

A one-way speech valve that is inserted into the TEP to allow air to pass from the trachea through the prosthesis and into the esophagus to produce sound.

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










BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAFIA / BIBLIOGRAFIE

Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação. / Další odkazy na literaturu jsou k dispozici na požádání. / Dodatkowe informacje dostępne na żądanie.

1. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
2. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
3. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
4. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
5. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.

Helix Medical, LLC

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE

Molenstraat 15
2513 BH, The Hague
The Netherlands